1 Tyvek externo  
External Tyvek

2 Tyvek interno  
Internal Tyvek

3 Esferas Retrorreflexivas  
Retroreflective Spheres  
Esferas Retrorreflectantes

4 Blister interno  
Internal blister

5 Blister externo  
External blister

## Português - PT

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## English - EN

## Presentation form

Orbis kits - Retro-reflective Spheres are supplied in double blister packaging sealed with Tyvek. The commercial models, their codes, and specific quantities are shown in the table below:

Code	Sphere Diameter	Sphere Fitting Type	Quantity
ART-OR-1-01	Ø 11mm	Snap Fit	1
ART-OR-1-02	Ø 11mm	Snap Fit	2
ART-OR-1-03	Ø 11mm	Snap Fit	3
ART-OR-1-04	Ø 11mm	Snap Fit	4
ART-OR-1-05	Ø 11mm	Snap Fit	5
ART-OR-1-06	Ø 11mm	Snap Fit	6
ART-OR-1-07	Ø 11mm	Snap Fit	7
ART-OR-1-08	Ø 11mm	Snap Fit	8
ART-OR-1-09	Ø 11mm	Snap Fit	9
ART-OR-2-01	Ø 13mm	Threaded Fit	1
ART-OR-2-02	Ø 13mm	Threaded Fit	2
ART-OR-2-03	Ø 13mm	Threaded Fit	3
ART-OR-2-04	Ø 13mm	Threaded Fit	4
ART-OR-2-05	Ø 13mm	Threaded Fit	5
ART-OR-2-06	Ø 13mm	Threaded Fit	6
ART-OR-2-07	Ø 13mm	Threaded Fit	7
ART-OR-2-08	Ø 13mm	Threaded Fit	8
ART-OR-2-09	Ø 13mm	Threaded Fit	9

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


Artis | Instructions for use: Orbis: Retro-reflective Spheres

### Product Identification

Technical name: Marker for Surgical Navigation  
Commercial Name: Orbis - Retroreflective Spheres  
ANVISA Registration: 80463390009

### Product description

Orbis – Retro-reflective Spheres is a sterile, single-use disposable product, developed specifically to be used with passive Surgical Navigation systems. Its purpose is to reflect the infrared light emitted by the stereovision system, allowing the instruments in the surgical field to track.

### Product composition

Each Retroreflective Sphere unit is composed of a plastic body, coated with retroreflective material. Models with code ART-OR-1 are made up of balls with a diameter of Ø 11mm and press fit. Models with code ART-OR-2 are made up of balls with a diameter of Ø 13mm and a threaded fitting.

### Indications and use

Orbis kits – Retroreflective Spheres are indicated to be used in surgical navigation procedures performed with passive Surgical Navigation systems.

### Important information to the doctor

Orbis Kits – Retroreflective Balls are supplied sterile and must not be reprocessed, repackaged, or resterilized. The correct performance of retroreflective spheres depends on their surface integrity. Keep the surface of the retroreflective spheres dry and clean throughout the surgical procedure.

### Contraindications

There is no contraindication so far regarding the use of this product.

### Adverse reactions

There are no adverse reactions related to the use of this material.

### Expiration date

The expiration date is 3 years. Refer to the product label.

### Warnings

Orbis kits – Retroreflective Spheres should not be used if their packaging is violated or the surface of the balls is scratched or worn. During the surgical navigation procedure, if the surface of the spheres is compromised by blood, fluids, or any other substance that interferes with the retro-reflective property of the material, it must be immediately replaced. Do not attempt to clean a contaminated sphere. Mechanical cleaning can damage the surface by altering the material's retro-reflective properties.

### Precautions

Strict compliance with correct and aseptic surgical techniques must be observed. Check if there was the correct adaptation to the post and make sure that the balls have descended to the maximum limit of the fitting. The sterility of the piece is only guaranteed if the primary packaging (blister sealed with Tyvek) is intact.

### Storage and handling conditions

Orbis kits - Retroreflective Spheres must be stored protected from direct light and rain. Recommended storage temperature ranges between 5°C and 40°C. The maximum recommended stacking is 10 boxes. Once the Kit is opened, the spheres must be handled with care to avoid direct abrasion on its surface, as this can influence the retro-reflection capacity.

### Disposal Information

The disposal of products must be carried out in accordance with the precepts established by RDC resolution no. 222/2018 of march 28, 2018, considering that the material is classified as A1 Group. It is also noteworthy that before disposal, the material must be decharacterized. For that, it is suggested that the pieces be scratched, cut, perforated, or broken.

## Instructions for use

– Open the secondary packaging (outer blister sealed with tyvek) peeling off the tyvek. Place the primary packaging (inner blister sealed with tyvek) on the sterile auxiliary table;

– Open the primary packaging (inner blister sealed with tyvek) to access the spheres;


– Position each sphere on its respective post on the instrument to be used;

– Check that all spheres have descended correctly occupying the entire length of the post;

– At the end of the procedure, the spheres must be removed from the instruments and discarded;



– To remove the snap fit balls, place your fingers near the base of the ball and pull upwards. Balls with threaded fitting must be unthreaded by turning counterclockwise;

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