
1

**1. ART-TC-2 - DOMUS - KIT TEMPLATE PARA CRANIOPLASTIA**

ART-TC-2 - DOMUS - TEMPLATE KIT FOR CRANIOPLASTY
ART-TC-2 - DOMUS - KIT MODELO PARA CRANIOPLASTIA

1.1. ART-TC-2-1 - DOMUS - TEMPLATE

ART-TC-2-1 - DOMUS - TEMPLATE
ART-TC-2-1 - DOMUS - MODELO

1.2. ART-TC-2-4 - DOMUS - RÉPLICA DA FALHA ÓSSEA

ART-TC-2-4 - DOMUS - REPLICA OF THE BONE DEFECT
ART-TC-2-4 - DOMUS - RÉPLICA DEL DEFECTO ÓSEO

1.3. ART-TC-2-3 - DOMUS - TAMPA DO MOLDE

ART-TC-2-3 - DOMUS - MOLD TOP
ART-TC-2-3 - DOMUS - TAPA DEL MOLDE

1.4. ART-TC-2-2 - DOMUS - BASE DO MOLDE

ART-TC-2-2 - DOMUS - MOLD BOTTOM
ART-TC-2-2 - DOMUS - BASE DEL MOLDE

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English (EN)**Precautions**

Strict compliance with correct and aseptic surgical techniques must be observed. The Domus - Template Kit for Cranioplasty comes into direct contact with the patient, therefore, the maintenance of the entire aseptic chain during handling must be observed. The sterility of the piece is only guaranteed if the primary packaging (blister sealed with Tyvek) is intact.

Disposal Information

It is recommended to follow local legal regulations for the disposal of potentially contaminating products

Storage precautions

The Domus - Template Kit for Cranioplasty must be stored in its original unopened packaging, inside a storage room in a dry, clean environment, protected from light, with humidity between 5% and 70% and temperature not lower than 5 °C and not higher than 40 °C.

The product must not be stored directly on the floor and must maintain a minimum height of 20 cm. It should not be placed in very high places, close to light bulbs, as this can cause drying of the packaging or damage to the labeling, and it should also not be stored in places where contaminants are stored.

Transport precautionse

The transport of the Domus - Template Kit for Cranioplasty must be carried out with the piece properly stored in its original packaging, identified and treated as a fragile item, maintaining a temperature between 5°C and 45°C, protected from light. Transport the product carefully to avoid drops or friction in order not to damage the packaging or the product components, thereby ensuring the sterility of the material.

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English (EN)

Product description

The **Domus - Template Kit for Cranioplasty** It is a patient-specific product that uses the patient’s CT imaging to adjust a modeling shape, create a replica of the bone defect, and produce a template that simulates the reconstruction of the material to be shaped accordingly.

Template Kit for Cranioplasty Composition

The Product consists of four parts (**Domus - Template ART-TC-2-1 with a Domus - Mold Bottom ART-TC-2-2, Domus - Mold Top ART-TC-2-3 and Domus - Replica of the Bone Defect ART-TC-2-4**) and its function is to assist in cranioplasty procedures. Each part is individually packaged through blisters sealed with Tyvek paper, which makes up the primary packaging of the product. The three parts that make up the product are packed in a cardboard box that makes up the outer packaging of the product.

COMPOSITION - ART-TC-2 - DOMUS - TEMPLATE KIT FOR CRANIOPLASTY			
CODE	DESCRIPTION	QTY.	MATERIAL
ART-TC-2-1	Domus - Template	01	Polylactic Acid (PLA); Polyurethane Rubber
ART-TC-2-2	Domus - Mold Bottom	01	Polyurethane Acid (PLA); Polyurethane Rubber
ART-TC-2-3	Domus - Mold Top	01	Polylactic Acid (PLA); Polyurethane Rubber
ART-TC-2-4	Domus - Replica of the Bone Defect	01	Polylactic Acid (PLA)

Indications and use

The **Domus - Template Kit for Cranioplasty** It is indicated for modeling acrylic bone cement in skull reconstruction surgeries or cranioplasty procedures.

The **Domus - Template ART-TC-2-1** it is used to assess whether there is any tissue interference along the bone edge that could prevent the correct fitting of the modeled piece to the surgical bed. The **Domus - Mold Bottom ART-T C-2-2** and the **Domus - Mold Top ART-TC-2-3** are responsible for modeling the cranioplasty material, giving it the exact adaptation provided by the **Domus - Template ART-TC-2-1**.

The use of this product brings the following benefits to the patient:

- Accurate marginal adaptation of the final rebuilt piece;
- Eliminates patient exposure to the exothermic polymerization stage of the acrylic bone cement to be modeled;
- Reduction in the time required for in loco material modeling, reducing the total surgery time;
- Aesthetic gain when compared to manual modeling of acrylic bone cement;

General note

The polymerization time for the acrylic bone cement to be modeled may vary depending on the acrylic bone cement manufacturer chosen and the temperature conditions that can alter the total working time of the material. Carefully observe the specific instructions for use of the bone cement chosen for the procedure.

Warnings

For the correct use of the Domus - Template Kit for Cranioplasty, the technician or surgeon responsible for handling it must carefully read this instruction for use, and be familiar with the surgical technique regarding the use of acrylic bone cement in cranioplasty cases. It is strongly recommended to use an auxiliary press to apply uniform pressure to the mold top and thereby obtain the best result in the molding process.

Instructions and use

1. After exposing the surgical bed, position the **Domus - Template ART-TC-2-1** over the bone defect and check if the margins are free from interference, which ensures a good marginal adaptation of the modeled reconstruction piece;

2. Mix the polymer and monomer of the acrylic bone cement according to the manufacturer’s instructions until the beginning of the plastic phase. The plastic phase is characterized as the stage in which handling the material does not produce filaments;

3. To handle the acrylic bone cement in the plastic phase, a sterile glove moistened with saline solution should be used to prevent the cement from sticking to the glove;

4. Moisten the inside of the **Domus - Molde Bottom ART-TC-2-2** with saline solution to minimize the chance of acrylic bone cement adherence;

5. The acrylic bone cement should be placed into the **Domus - Mold Base ART-TC-2-2** at the beginning of the plastic phase and manually distributed evenly and consistently across the entire inner surface;

6. The **Domus - Mold Top ART-TC-2-3** must be correctly positioned, observing the indication of the X and Y orientations, presented on the corners of the cover and the base, which must coincide;

7. Having correctly positioned the cover in relation to the base, it is necessary to apply uniform pressure on the lid so that the material flows properly and uniformly, ensuring the correct modeling of the acrylic bone cement;

8. Note that the edges of the lid must meet the edges of the base to ensure that the correct planned thickness is properly reproduced;

9. Wait for the material to set completely, according to the time indicated by the acrylic bone cement manufacturer, then open the form;

10. Remove excess material respecting the line that delimits the contour of the piece;

11. Check the correct fit of the modeled piece using the **Domus - Replica of the Bone Defect ART-TC-2-4**. To ensure proper orientation, refer to the letters printed on the bone defect, where the letter "A" stands for Anterior and the letter "S" stands for Superior.

12. Before placing the piece onto the surgical bed, rinse it thoroughly with saline solution, then position it to be fixed to the patient’s bone edge.

Handling precautions

The opening of the sealed package must be done manually, in an appropriate place (surgical environment), and by a health care professional. The opening of the blister sealed with Tyvek should only be done during surgery to avoid contamination of the product. The product must not be used if the packaging has been violated

Contraindications

This product should not be used for any other purpose or to model any material not specified in this instruction for use. Only acrylic bone cements for cranioplasty procedures, with regular registration with ANVISA, may be used in conjunction with Domus - Cranioplasty Template Kit.

Adverse reactions

There are no adverse reactions related to the use of the Domus - Template Kit for Cranioplasty.

Important information for the Physician

The Domus - Cranioplasty Template Kit was developed to model acrylic bone cement according to the planning based on the patient’s CT images. It is essential that the entire contour of the bone defect is properly exposed and free of any soft tissue interference to ensure the correct adaptation of the piece.

Validity

The validity of the product is 1 year. Refer to the product label.

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