EASY-TO-USE TOOLS FOR CLEARED TISSUE LIGHT SHEET IMAGE ANALYSIS





JACK G. ZEITOUN¹, ANTHONY REKSOATMODJO¹, CHASE REDD¹, RICARDO AZEVEDO^{1,2}, SUNIL P. GANDHI^{1,2}, DAMIAN G. WHEELER^{1,3}, PHILIP CHEUNG¹


Translucence Biosystems, Inc. (Irvine, CA, USA) ²Univ. of California (Irvine, CA, USA) ³ Activity Signaling, LLC (San Diego, CA, USA)

049.03


Introduction Conventional 2D histology methods are limited to thin sections and regions of interest, restricting insight into brain-wide structure and function. Recent advances in tissue clearing and light sheet imaging have enabled the evaluation of whole-brains at cellular-resolution but have introduced massive data volumes that present a challenge for analysis and reproducibility. High-resolution 3D light sheet imaging requires efficient workflows for stitching, registration, and quantification across terabyte-scale datasets. Translucence Biosystems has developed a suite of Al-powered tools—Stitchy, BrainSeeker, Maverick, and Voxels—to streamline every stage of the 3D image analysis pipeline, making whole-brain analysis accessible, reliable, and reproducible. Pipeline For the Generation of 3D Anatomics Data Tissue CLEARING LIGHT SHEET IMAGING PREPRICESSING QUANTIFICATION Addicable Data













This figure shows two tools under development: 1) object-based co-localization on multi-channel imaged data and 2) our Franslucence reference brain. For our colocalization tool, segmented objects from two fluorescence channels: cFos(+) neurons are shown in Red and neurons expressing Npas4 are shown in green. Our workflow determined the closeness of the centroids of the segmented objects and the yellow dots indicate the location of neurons expressing both proteins. These calculation are done across the entire brain, but a simple cropped 2D plane is shown. The Translucence Reference Brain was generated in partnership with UNC and Kitware, using 64 light sheet autofluorescence image volumes registered together using a neural-networkbased registration method called ConstrICON. Unlike the Allen Reference Brain, generated with serial two-photon tomography, the Translucence Reference Brain originates from light sheet data allowing for superior registration of light

Summary

- The new Translucence software suite integrates image stitching, registration, and quantification into a unified, scalable workflow for 3D light sheet datasets.
- Stitchy enables automated tiling and z-plane correction across fluorescence channels, supporting all major LSFM file formats and large-batch processing.
- BrainSeeker provides rapid and accurate brain registration to common atlas spaces and advanced Al-based alignment tools.
- Maverick and Voxels empower users to design, visualize, and execute machine-learning-based analysis pipelines through a scientist-friendly, no-code user interface.
- Together, these tools accelerate data processing, improve reproducibility, and transform terabyte-scale 3D light sheet imaging datasets into actionable scientific insights.