


# WHOLE BRAIN MEASUREMENT OF PLAQUES, MICROGLIA, $\alpha$ -SYNUCLEIN AND TYROSINE HYDROXYLASE IN AD AND PD MODEL MICE

NATHANIEL GUANZON<sup>1</sup>, YESSENIA GALLEGOS<sup>1</sup>, CHASE REDD<sup>1</sup>, ERIK CASTILLO<sup>1</sup>, ERIC BLAES<sup>1</sup>, RICARDO AZEVEDO<sup>1,2</sup>, SUNIL GANDHI<sup>1,2</sup>, DAMIAN C. WHEELER<sup>1,3</sup>  
<sup>1</sup>Translucence Biosystems Inc, Irvine, CA, USA, <sup>2</sup>Center for the Neurobiology of Learning and Memory, Department of Neurobiology and Behavior, University of California, Irvine; <sup>3</sup>Activity Signaling, LLC


## Introduction

- Alzheimer's disease and Parkinson's disease involve complex, brain-wide pathology, highlighting the need for more comprehensive tools to evaluate disease-related changes across the intact brain.
- Traditional 2D histology is constrained to thin sections and select regions, offering only a partial and potentially biased assessment of disease-associated biology.
- 3D tissue clearing combined with light-sheet imaging enables whole-brain visualization at cellular resolution across intact tissues.
- When paired with AI-powered analysis, this approach enables unbiased quantification of amyloid plaques, microglial activation, TH-positive staining, and  $\alpha$ -synuclein cell bodies and fibrils across hundreds of brain regions.
- In this poster, we demonstrate how this pipeline has been adapted for AD and PD mouse models. In AD (5xFAD) mice, we segment microglia, plaque-associated microglia, and amyloid plaques. In PD models (MPTP, and PFF-injected mice), we use TH and  $\alpha$ -synuclein immunostaining to map dopaminergic degeneration and  $\alpha$ -synuclein pathology across the intact brain.


## Whole-Brain Immunostaining of Iba1(+) Microglia




## Brain-Wide Regional Quantification of $\beta$ -Amyloid Plaques in 5xFAD Mouse Brains




## Brain-Wide Quantification of Plaque-Associated Microglia



## Whole-Brain Quantification of Tyrosine Hydroxylase in a PD Model




## Whole-Brain Quantification of $\alpha$ -Synuclein Pathology in a PFF Model



## Colocalization of Microglia and Amyloid Plaques in 5xFAD Mice

In WT mice, individual homeostatic microglia are spaced out and have a ramified morphology with processes surveying the surrounding region. In diseased and inflamed states, activated microglia take on a condensed morphology. In 5xFAD Tg mice, microglia move to colocalize with  $\beta$ -amyloid plaques. Our automated methods can identify and count plaques and plaque-associated microglia (PAM) throughout the brain.



## Summary

- In 5xFAD mice, whole-brain immunostaining and quantification revealed widespread  $\beta$ -amyloid plaque pathology across anatomically distinct brain regions.
- Iba1 quantification showed that microglia colocalize with amyloid plaques and adopt an activated morphology, enabling selective detection and quantification of plaque-associated microglia (PAMs).
- In PD models, whole-brain TH immunostaining revealed diverse regional patterns of TH-positive cell bodies and neurites and demonstrated reduced TH signal in injured animals relative to controls.
- Whole-brain  $\alpha$ -synuclein immunolabeling in PFF-treated mice combined with AI-powered quantification of  $\alpha$ -synuclein-positive cell bodies and fibrils, highlighted the regional distribution of PD pathology.
- Together, these data demonstrate the value of our integrated tissue clearing, light-sheet imaging, and AI-powered analysis pipeline for unbiased whole-brain measurement of amyloid pathology, microglial activation, dopaminergic degeneration, and  $\alpha$ -synuclein aggregation.