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Typ IIR, BINDEBÄNDER

DE

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



TYPE IIR, TIE STRAPS

EN

Product: Surgical Face Mask type IIR **Size:** 17,5x9,5 cm **Model:** Tie straps, planar **Range of Application:** The surgical face mask type IIR is a fluid resistant, disposable medical device that creates a physical barrier between the mouth and nose of the wearer and the immediate environment. It is designed to prevent spreading of infections by catching pathogens shed in liquid droplets and aerosols from the wearer's mouth and nose and protect the wearer from potentially contaminated fluid splash. **Structure and Composition:** This product consists of a mask body (made of non-woven fabrics and polypropylene melt-blown fabrics), a nose clip of galvanized steel and tie straps. The product is free from latex, PVC and nickel. **Product performance:** This product fulfills the performance requirements of ISO EN 14683:2019 type IIR.

Instructions: (1) Clean and disinfect hands before use. (2) Open the package and take one mask out. The colored side is the external side. Avoid contact with the internal side of the mask. (3) With the nose clip facing upwards, attach the mask to your nose and stretch the mask around your chin. (4,5,6) Then attach the tie straps around your head and neck. (7) Use the index finger of both hands to pinch the nose clip sideways from the middle so it takes the shape of your nose bridge. Make sure the mask fits your face-shape. When removing the mask, avoid touching the external side of the mask (the external side of the mask is contaminated). Only touch the tie straps and discard of the mask into a waste container with a lid. If the mask has been used in a medical environment it should be discarded according to medical waste regulations. Clean and disinfect hands after use.

Caution: Check whether the package is intact before use. Confirm the label, production date and use-before-date. Use within the validity period. Do not use if the package is damaged. The product should be used as soon as possible after the package is opened. This product is disposable and should never be reused. After use, discard of the mask into a waste container with a lid. If the mask has been used in a medical environment it should be discarded according to medical waste regulations.

 Date of manufacture - Indicates the date when the medical device was manufactured	 Manufacturer - Indicates the manufacturer of the medical device
 Expiration Date - Indicates the date after which the medical device should not be used	 List number - Indicates the manufacturer's list number so that the medical device can be identified
 Do not expose to direct sunlight - Indicates a medical device that should not be exposed to direct sunlight	 Batch code - Indicates the manufacturer's batch code so that the batch can be identified
 Store dry - Indicates that the medical device should be protected against moisture	 Medical device - Indicates that it is a medical device
 Do not use if damaged packaging - Indicates a medical device that should not be used if the packaging has been damaged or opened	 CE marking - Indicates the products sold in the EEA have been assessed to meet high safety, health, and environmental protection requirements
 Do not reuse - Indicates that the medical device is intended for single use or for an individual patient during a single procedure	 Important - Indicates that it is necessary for the user to read the user instructions for important safety information about e.g. warnings and precautions , which for various reasons, cannot be indicated on the medical device
 Read the instruction manual - Indicates that it is necessary for the user to read the instruction manual	 Authorized Representative - Indicates the Authorized Representative of the European Community

Storage requirements: Store protected from sunlight in a dry, clean and well-ventilated environment.

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