

COVID-19 Antigen Detection Kit

Nasal swab | Self-test

DE COVID-19 Antigen-Test | Nasenabstricht
NL COVID-19 Antigentest | Neusswab
FR Test d'antigène COVID-19 | Écouvillon nasal
SV COVID-19 Antigentest | Nasal provtagning


EN Instructions for use
DE Gebrauchsanweisung
NL Gebruiksaanwijzing
FR Mode d'emploi
SV Bruksanvisning


EN Sampling swab **DE** Abstrichtupfer **NL** Bemonsteringsstaaf **FR** Écouvillon de prélèvement **SV** Testpinne

EN Tube cap **DE** Röhrchenkappe **NL** Buisdop **FR** Bouchon du tube **SV** Kork till provrör

EN Test tube with buffer liquid **DE** Probenentnahmehöhrchen mit Pufferlösung **NL** Monsterextractiebuis met buffervloeistof **FR** Tube de test avec liquide tampon **SV** Provrör med buffertvåtska

EN Foil pouch and test cartridge **DE** Folienbeutel und Testkassette **NL** Aluminium zakje en testkaart **FR** Pochette en aluminium et cartouche de test **SV** Förpackning med testkassett

EN Desiccant (discard, do not open) **DE** Trockenmittel (entsorgen, nicht öffnen) **NL** Droogmiddel (weggooien, niet openen) **FR** Déshydratant (jeter, ne pas ouvrir) **SV** Torkmedel (släng, öppna ej)

Reference Number: COVID-19-NG21 (New Gene)
Article Number: REF 5301, 5305, 5325 (Gibson Medical)

Version: GM1-v05-NS-HT
Date: 2021-12

PRODUCT INFO

IVDD Directive
(98/79/EC)

EN COVID-19 Antigen Detection Kit Nasal swab

SUMMARY: COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection. Asymptomatic virus carriers can also be infectious sources. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The most common symptoms include fever, dry cough and tiredness. Some cases also report aches and pains, sore throat, diarrhea, conjunctivitis, headache, loss of taste or smell, a rash on skin or discoloration of fingers or toes. The serious symptoms include difficulty breathing or shortness of breath, chest pain or pressure, loss of speech or movement. Without timely medical intervention, severe complications or death may occur on COVID-19 cases.

INTENDED USE: This is a in-vitro diagnostic rapid self-test for the qualitative detection of SARS-CoV-2 antigens (Ag). Suitable for people with COVID-19 symptoms to assist the early diagnosis of SARS-CoV-2 infection. Can also be used by people without COVID-19 symptoms to regularly monitor their health status.

PRINCIPLE: The COVID-19 Antigen Detection Kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein from SARS-CoV-2. The test strip is composed of the following parts: nsample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal-gold conjugated with the monoclonal antibody against the nucleocapsid protein of SARS-CoV-2; the reaction membrane contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates absorbed in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen is present in the sample, the complex of the anti-SARS-CoV-2 conjugate and the virus will be captured by the specific anti-SARS-CoV-2 monoclonal antibodies coated on the test line region (T). Absence of the T line suggests a negative result. To serve as a procedural control, a red line will always appear in the control line region (C) indicating that proper volume of sample has been added and membrane wicking effect has occurred.

WARNINGS AND PRECAUTIONS: This product is applicable for nasal sampling, other test samples may result in incorrect or invalid test results. Do not use components from other sources. Keep out of reach of children. Children under the age of 15, or people unable to perform the test themselves, should be tested by their parents, guardians or under adult supervision. For people who have recent nasal trauma or surgery, or have severe coagulopathy, gentle operation is required for nasal swab collection to avoid injuries to the nose. Do not force the swab into the nose. Contact a health care provider if your nose is damaged during the procedure. Carry out the analysis within two hours after sampling. Outdated and stale samples can cause erroneous results. Check the expiration date on the foil pouch containing the test cartridge. Do not use expired tests. Make sure the correct amount of sample is used for testing; too much or too little can cause incorrect results. Test result is displayed 15-30 minutes after the liquid is expelled into the well. Read the result in a place with sufficient light to be able to interpret the result accurately. Use a phone or watch as a timer. After 30 minutes, the result is invalid. If the test line or control line is outside the test window, do not use the test cartridge as the test result is invalid. Retest the sample with another test cartridge. This product is for single use only. Do not recycle used components. Disinfect used products with household bleach spray or a 70-75% alcohol spray, or place in a sealed waste bag before discarding in household waste. Wash your hands thoroughly after the test.

LIMITATIONS: This product is intended for self-diagnosis of COVID-19. The final diagnosis should be determined by a physician after evaluation of clinical signs and results of other examinations. A negative result indicates that there is no virus in the sample, or that the virus load is below the detection limit for this product. The possibility of viral infection cannot be completely ruled out. Incorrect sampling or low viral load may cause a false negative result. A positive result indicates that the tested sample has a viral load that is higher than the detection limit of this product. The color intensity of the test line does not correspond to the severity of infection or disease progression in the person. The amount of viral antigens in the sample will decrease with the duration of the disease. Samples taken one week after the onset of symptoms are more likely to show false negative results. The use of a liquid other than the supplied extraction buffer solution will invalidate the test. In case of further suspicion of infection, repeat the test after 1-2 days.

STORAGE AND STABILITY: Store the test kit at a temperature between 2-30°C. Avoid exposure to sunlight. The test kit is stable until the expiration date on the test cartridge packaging. Once the test cartridge package has been opened, the test cartridge should be used within one hour (1h). Prolonged exposure to hot or humid environments can give incorrect results.

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
DEAnweisungen

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SVInstruktioner

TAKING THE SAMPLE: **1** Open the package containing the sampling swab. Do not touch the swab tip. Gently insert the tip of the swab into one nostril, 2-3 cm. Press it against the nasal wall and rotate for 20 seconds to collect nasal secretion. Repeat in the other nostril. **2** Remove the sample extraction tube seal and insert the swab. Press it against the sides and rotate for 10 seconds. **3** Squeeze the sides of the tube to expel sample solution from the swab. Place the tube cap firmly onto the tube. Avoid touching the top of the tube cap. Shake the tube, or gently press the tube, to mix the sample solution. Let it rest for one minute to release viral antigens.

TEST PROCEDURE: **4** Peel to open the foil pouch and remove the test cartridge from its package. Place it on a flat surface. Avoid touching the well and window of the test cartridge. **5** Turn and tap the test tube to mix the liquid again. Avoid touching the top of the tube cap. Hold the tube upside down and squeeze to expel 3 drops of sample solution into the loading well. Testresult is displayed 15-30 minutes after the liquid is expelled into the well. After 30 minutes, the result is invalid.

RESULTS: **POSITIVE:** Red line at C and T, regardless of how faint. White line at T is considered negative. **NEGATIVE:** Red line at C only, see Limitations. **INVALID:** Without line at C, retest with another test cartridge.

AFTER TESTING: Contact a health care provider in case of a positive result. Additional testing, such as PCR test, may be required to confirm viral infection or to initiate infection tracing. Always follow local and governmental recommendations. Disinfect used products and samples, or place in a sealed waste bag, before discarding in household waste.

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