


# Surgical Face Mask

EN 14683:2019 + AC:2019 **TYPE IIR** REF: 1097 | Eye Shield + Anti-Fog



GIBSONMEDICAL.SE

## SURGICAL FACE MASK














TYPE IIR, TIE STRAPS, EYE SHIELD + ANTI-FOG

EN

**PRODUCT:** Surgical Face Mask Type IIR **SIZE:** 17.5 x 9.5 cm **MODEL:** Tie Straps, Eye Shield + Anti-Fog **PRODUCT PERFORMANCE:** EN 14683:2019+AC:2019 Type IIR. **RANGE OF APPLICATION:** Fluid resistant disposable face mask that creates a physical barrier between the mouth and nose of the wearer and the immediate environment. Designed to prevent spreading of infections by catching pathogens shed in liquid droplets and particles from the wearer's mouth and nose. Also protects the wearer from potentially contaminated fluid splash. **STRUCTURE AND COMPOSITION:** Mask body made of three layer polypropylene non-woven fabric (spunbond-melt-blown-spunbond) with anti-fog foam. Spunbond tie straps and aluminium nose clip. Eye shield is made of PET (polyethylene terephthalate). Free from latex, PVC and nickel.

**INSTRUCTIONS:** Clean and disinfect hands before use. Open the package and take one mask out. Avoid touching the internal side of the mask. The coloured side is the external side. With the nose clip facing upwards, attach the mask to your nose and stretch the mask around your chin. Tie the upper straps around the back of the head and lower straps around the neck. Use the index finger of both hands to pinch the nose clip so it takes the shape of your nose bridge. Make sure the mask fits your face shape and that it is comfortable. Do not touch the mask when in use. Only touch the tie straps when removing it. Discard into a waste container with a lid, or according to medical waste regulations in medical environments. Clean and disinfect hands after use.

**CAUTION:** Check whether the product is intact before use. Do not use a damaged product. Confirm the label, production date and use-before-date. Use within the validity period. The product is disposable and should not be reused.

 <b>Manufacturer</b> - Indicates the manufacturer of the medical device	 <b>Read the instruction manual</b> - Indicates that it is necessary for the user to read the instruction manual
 <b>Date of manufacture</b> - Indicates the date when the medical device was manufactured	 <b>List number</b> - Indicates the manufacturer's list number so that the medical device can be identified
 <b>Expiration date</b> - Indicates the date after which the medical device should not be used	 <b>Batch code</b> - Indicates the manufacturer's batch code so that the batch can be identified
 <b>Do not expose to direct sunlight</b> - Indicates a medical device that should not be exposed to direct sunlight	 <b>CE marking</b> - Indicates that products sold in the EEA have been assessed to meet high health-, safety- and environmental protection requirements
 <b>Store dry</b> - Indicates that the medical device should be protected against moisture	 <b>Authorized Representative</b> - Indicates the Authorized Representative of the European Community
 <b>Do not reuse</b> - Indicates that the medical device is intended for single use or for an individual patient during a single procedure	 <b>Financial contribution towards recovery of packaging</b> - Indicates that a financial contribution has been made towards recovery of packaging in Europe
 <b>Medical device</b> - Indicates that it is a medical device	

**STORAGE REQUIREMENTS:** Store protected from sunlight in a dry, clean and well-ventilated environment.

**RECYCLING:** Packaging is sorted as paper. Face mask as residual waste.