
EN	NEUTRAL ELECTRODES
DE	NEUTRALELEKTRODEN
FR	ELECTRODES NEUTRES
IT	ELETTRIDI NEUTRI
BG	НЕУТРАЛЕН ЕЛЕКТРОД
CS	NEUTRÁLNÍ ELEKTRODY
DA	NEUTRALELEKTRODER
EL	ΟΥΔΕΤΕΡΑ ΗΛΕΚΤΡΟΔΙΑ
ES	ELECTRODOS NEUTROS
ET	NEUTRAALELEKTROODID
FI	NEUTRAALIELEKTRODIT
HR	NEUTRALNE ELEKTRODE
HU	NEUTRÁLIS ELEKTRÓDÁK
LT	NEUTRALŪS ELEKTRODAI
LV	NEITRĀLIE ELEKTRODI
NL	NEUTRALE ELEKTRODEN
PL	ELEKTRODY NEUTRALNE
PT <sub>EU</sub>	ELÉTRODOS NEUTROS
PT <sub>BRASIL</sub>	PLACAS TERRA
RO	ELECTROZI NEUTRI
SK	NEUTRÁLNE ELEKTRÓDY
SL	NEVTRALNE ELEKTRODE
SV	NEUTRALELEKTRODER
KO	중립 전극
NO	NØYTRALE ELEKTRODER
RU	НЕЙТРАЛЬНЫЕ ЭЛЕКТРОДЫ
TH	ອົເລີກໂທຣດທີ່ເປັນກລາງ
TR	NÖTRAL ELEKTROLAR
ZH	中性电极
AR	الإلكترودات المحايدة

Leonhard Lang GmbH, Archenweg 56, 6020 Innsbruck, Austria;  
Fax +43 512 392210, medical@LeonhardLang.at

## NEUTRAL ELECTRODES FOR NEONATES

Please read this document and keep it safe. Ensure that any individual using this product is familiar with and understands the information contained in this document.

## INTENDED USE / INTENDED PURPOSE

Neutral electrodes are intended as accessories for electrosurgical devices. Neutral electrodes are non-sterile; they are intended only for single use on intact (uninjured) skin and lead energy from the patient to the electrosurgical device.

### Indications:

- Neutral electrodes return high-frequency current from the patient back to an electrosurgical device and are indicated for all monopolar surgical procedures in which high-frequency electrical current is used to cut or coagulate tissue.

### Contraindications:

- Non-conventional electrosurgical procedures (e.g., high-current operation mode).
- Damage to or soiling of the neutral electrode that could lead to increased electrical resistance, inadequate contact, or short-circuit.
- Incompatibility of the neutral electrode with the electrosurgical device that could lead to a malfunction or a hazard for the patient or user.
- If the surgeon considers that there are additional contraindications or that situations could arise based on the current specialist literature that may endanger the patient or user.

## WARNING

Improper use of neutral electrodes can cause damage to tissue. These instructions for use serve to ensure patient safety. **NOT FOLLOWING THESE INSTRUCTIONS MAY LEAD TO BURNS, PRESSURE NECROSSES, OR OTHER SKIN TRAUMA DURING USE.**

0.45 - 4.99 kg		✓ 1 - 11 lbs

- PRODUCT LIMITATIONS:** This neutral electrode for neonates is only intended for use on neonates with a body weight of 0.45 – 4.99 kg (1 - 11 lbs). For neonates that weigh more than this, use a neutral electrode for children.

- PRODUCT LIMITATIONS:** Only use this neutral electrode for neonates for currents of up to 350 mA.

- PRODUCT LIMITATIONS:** These neutral electrodes for neonates are designed for use during conventional monopolar electrosurgical procedures. Restrict the duration of activation to a maximum of 30 seconds within each 2 minute period. Exceeding this limit can inflict burns on the patient due to overloading, even if the neutral electrode is correctly and fully applied and used together with an activated contact quality monitoring system.

- Do not use the neutral electrode if the packaging is open or damaged.
- Do not use the neutral electrode if it is damaged, modified, or expired.
- If an electrosurgical unit is equipped with a neutral electrode contact quality monitoring system, always use a split electrode. Never deactivate the acoustic alarm of this monitoring system during surgery.
- Use the lowest effective safe power setting to achieve the desired surgical effect.

- Do not operate the electrodes in the immediate presence of flammable gases (including concentrated oxygen) and solvents to avoid possible explosion or fire hazard. Avoid contact between the product and solvent-based liquids.

- Determine the maximum power settings and prepare the HF generator before the first use of the neutral electrodes for infants according to the guideline given below.

## COMPATIBILITY

For any compatibility questions, particularly regarding compatibility with a specific electrode contact quality monitoring system, please contact your local distributor. In any case, the split neutral electrodes for

neonates are compatible with the dynamic contact quality monitoring systems REM™, NESSY®, and ARM™.

### CAUTION

- Surgeons should be familiar with the effects of electrosurgery on small patients, and should first consider the use of bipolar electrosurgery, which does not require a neutral electrode!
- Take particular care when using this neutral electrode for neonates on premature infants born before 34 weeks of gestation. If possible, consult a neonatologist to determine whether the neonate's skin has sufficiently developed for contact with adhesive materials.
- Do not trim down the neutral electrode to reduce its size! Do not use additional gel!
- Do not reposition the neutral electrode!
- If the adhesive border surrounding the gel is not used, or is only partially used, the neutral electrode must be consistently monitored for the ingress of fluid! It should not be placed at points at which fluids may be present!
- Before using the electrosurgical device, familiarize yourself with its instructions for use, in particular with respect to power setting restrictions and the maximum uninterrupted activation period.
- Should the coagulating or cutting effect diminish or fail during surgery, or a higher than normal power setting be required, immediately make sure that the neutral electrode is appropriately positioned and in full contact with the skin. Inspect all connections of the neutral electrode, such as clamp, cable, and connector, plus the active accessories and the generator before turning up the power setting!
- If you reposition the patient, make sure the entire surface of the neutral electrode still sticks well to the skin and check all cable connections.
- Do not reuse the neutral electrode! If re-used, the adhesive strength and the electrical properties may be insufficient, which could lead to a patient injury. There is also a risk of cross-infection from one patient to another.
- The electrosurgical devices used must be type BF (Body Floating) or CF (Cardiac Floating).
- The approved length of accessories defined for the specific electrosurgical device must be compared to the length of the neutral electrode cable. The length of the neutral electrode cable (3 m / 9.8 ft or 5 m / 16.4 ft) can be found on the packaging information.

### USER GROUP

Neutral electrodes must be applied on the patient by a physician or by appropriately trained medical professionals.

## INSTRUCTIONS FOR USE

### Determining the Maximum Power Settings and Preparing the HF Generator:

- Before the neutral electrode for neonates is used for the first time, the HF generator must be calibrated by a qualified technician to determine which power setting corresponds to a current of 350 mA for each mode. When doing so, use a commercially available HF tester (e.g. from Fluke or mtk biomed). Set the load resistance to 125 ohms. Then, the corresponding power in watts must be determined for each mode at a current of 350 mA as a maximum power setting (coagulation, electrotomy, combined forms). The thereby determined maximum power settings should be indicated clearly and visibly on the HF generator for use by the surgeon. They must under no circumstances be exceeded during use of the neutral electrode for neonates.

### Selecting a Neutral Electrode:

- Whenever the size of the patient allows, use a neutral electrode for infants in order to reduce the likelihood of skin burns.
- The neutral electrode for neonates is only intended for use on neonates with a body weight of 0.45 – 4.99 kg (1 – 11 lbs). For neonates that weigh more than this, use a neutral electrode for children.

### Selecting and Preparing the Application Site:

- The recommended application site on neonates weighing between 0.45 and 4.99 kg (1 – 11 lbs) is on the back, below the shoulder blades, and above the sacrum. Special clinical circumstances may require alternative application sites. If alternative application sites are required, ensure the neutral electrode is in contact with the skin of the neonate over the maximum possible surface area and reduce the power settings.
- Make sure that the selected site is not thermally insulated or subject to additional heat from heat sources during the operation.
- For surgeries during which the HF current may flow through parts of the body with a relatively small cross-sectional area, bipolar methods may be more suitable in order to avoid unintentional damage to tissue.
- Monitoring electrodes or other devices that may provide an alternative earthing route for the high-frequency current should be applied as far from the surgical site as possible. The exclusive use of ECG and other monitoring systems, cables, and wiring equipped with HF current limiting devices (e.g. high frequency filters or chokes) is recommended. If this is not possible, the neutral electrode must be applied closer to the surgical site than these electrodes or devices. The use of needle monitoring electrodes is not recommended.
- Thoroughly clean the selected area of skin (e.g. to remove creams or vermix caseosa). Carefully dry the area, especially if skin cleaning fluids have been used. Avoid using flammable skin prepping agents or disinfectants (e.g. acetone degreasers). **Bear in mind that failure to prepare the area can lead to skin burns.**
- Avoid skin-to-skin contact (for example between the patient's arms and body) by inserting insulation between the relevant areas, for example dry gauze.

- Position the neutral electrode cable so that the cable is not pulling the neutral electrode away from the patient's skin. The neutral electrode cable should also be positioned to avoid contact with the patient and other cables, and to ensure there are no loops. In particular, do not wrap the neutral electrode cable around one of the patient's extremities or other grounded objects, in order to prevent capacitive coupling and resulting burns.

- Check that the neutral electrode is securely adhered to the skin across its entire surface, and that the clamp is securely attached to the connection lug of the electrode. Make sure that the clamp is not exerting any unnecessary pressure on the patient's skin. If the adhesive border surrounding the gel is not used, or is only partially used, the neutral electrode must be consistently monitored to prevent the ingress of fluid.
- Before using the electrosurgical device, familiarize yourself with the instructions for use, especially with respect to the instructions for using neutral electrodes in operations on neonates.
- Connect the neutral electrode cable with the electrosurgical device in accordance with the instructions for use for the generator.

### Removing and Disposing of the Neutral Electrode:

- After use, carefully remove the neutral electrode using one hand while supporting the skin beneath it with the other. Lift a corner of the electrode and slowly detach it from the bottom upwards (not using the neutral electrode cable). Ripping off the electrode or peeling it off quickly may damage the skin. Take extra care when the skin is exceptionally sensitive, especially with newborns and premature infants. Do not use solvents to remove adhesive residue. The slow and careful removal of any adhesive residue with warm water and cotton balls is preferred.
- To open the cable clamp, lift the locking lever.
- Dispose of the product in accordance with local regulations or the guidelines of the hospital/health care facility.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Environmental conditions: +10°C to +40°C; 20% to 85% RH  
The conditions for transport and storage are stated on the primary and secondary packaging.

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LEONHARD LANG GmbH