Integrated MALDI-MSI and Glycopeptide LC-MS/MS Reveal Site- and Region-Specific Glycosylation Changes in Tumor Tissue



Hongxia Bai¹, Grace Grimsle¹, Richard Drake², Klaus Lindpaintner¹

¹Bruker-Daltonics, Life-Science Mass Spectrometry, Billerica, MA, USA; ²

Department of Pharmacology, Medical University of South Carolina

Introduction

Glycosylation plays a central role in maintaining protein structure and function, particularly in cell-cell communication and immune signaling. In this study, we integrate high-resolution mass spectrometry imaging (MSI) of released glycans with intact glycopeptide analysis to uncover spatial and molecular differences in N-glycosylation between tumor and normal prostate tissue, and to connect and resolve findings from glycomics spatial biology imaging with commensurate observations on the level of individual glycoprotein analysis. MALDI-MSI revealed prominent changes in glycan spatial distribution, including enrichment of high-mannose and bisected glycans in tumor tissue, while complex-type fucosylated and sialylated glycans predominated in normal regions. LC-MS/MS analysis of glycopeptides validated these trends and highlighted specific glycoproteins involved. This dual-modality approach provides mechanistic insights into glycosylation remodeling during malignant transformation and strongly indicates the potential relevance of glycosylation-specific biomarkers

Methods




Figure1. Adjacent normal and Pca tissue sections were processed by MSI for glycomics and by LC-MS/MS for glycoproteomics.

Results




Figure 2. Histology and MALDI-MS glycan profiling of normal (upper image and spectrum) and Pca (bottom image and spectrum) tissue.




Figure 3. Spatial distribution of representative N-glycan classes across normal and PCa prostate tissues.





Figure 4. Chord diagram and bar plots summarizing glycopeptide diversity and distribution in Normal (Upper) and Pca (Bottom) tissue.



A L G F E n A T Q A L G R

Figure 5. A representative glycopeptide spectrum match with peptide backbone ALGFENATQALGR and the glycan composition H3N2F1. This peptide is from Galectin-3-binding protein (LG3BP), which is known to play roles in cancer progression, immune modulation, and metastasis, and has been reported as a potential biomarker in prostate and other cancers.




Figure 6. Glycopeptide-protein mapping of peptides bearing the same glycan motif. Each colored cell represents a glycopeptide carrying the shared glycan composition H3N2F1 (HexNAc(2)Hex(3)Fuc(1)), identified on a corresponding protein (y-axis). Orange indicates detection in tumor tissue, blue indicates normal tissue, and red marks peptides shared by both conditions. This map highlights the diversity of peptide backbones and proteins modified with the same glycan motif across tissue types

Summary

- Dual approach: MALDI-MSI + LC-MS/MS glycoproteomics
- PCa:
 \(\) High-mannose, biantennary glycansNormal:
 \(\) Fucosylated, sialylated glycans
- Normal and PCa showed different glycoproteins profiles
- Shared glycan (H3N2F1) found on diverse proteins/peptides
- These results underscore the value of combining imaging and glycoproteomics to uncover biologically meaningful glycosylation changes in cancer.

Disclaimer: H.B,; G.G.; and K.L. are employer of Bruker;

MUSC is a Bruker Excellence Center for Glycomics.

Technology