



## Introduction

Detergent-, salt-, and contaminant-free peptide samples are essential for bottom-up proteomics. The PreOmics® Phoenix Peptide Cleanup Kit is designed to assist researchers achieving best results with few sample preparation steps and little hands-on time. For sample-specific protocols and optimization visit [www.preomics.com/resources](http://www.preomics.com/resources) or contact [info@preomics.com](mailto:info@preomics.com).

## Kit Contents

The kit contains all chemicals and plasticware to perform an efficient peptide cleanup, removing detergents, fatty acids, sugars, salts, and other contaminants.

Component	Cap	Quantity	Buffer Properties			Description	Storage
			Organic	Acidic	Basic		
STOP	●	1x 15 mL	●	●		For sample acidification to promote efficient peptide binding.	RT
WASH X	●	3x 25 mL	●	●		For removing hydrophobic contaminants.	RT
WASH 1	●	2x 25 mL	●	●		For removing hydrophobic contaminants.	RT
WASH 2	●	1x 25 mL		●		For removing hydrophilic contaminants.	RT
ELUTE	●	1x 25 mL	●		●	For eluting peptides from the cartridge.	RT
LC-LOAD	○	1x 25 mL		●		For loading peptides on reversed-phase LC-MS column.	RT
CARTRIDGE		96x				Cartridges with SPE sorbent for peptide purification from 1–100 µg protein starting material. Racked in adapter plate and sealed with silicone mat.	RT
WASTE PLATE		1x				Deep well plate for collecting waste after washes.	RT
COLLECTION PLATE		1x				LoBind® plate for collecting peptides after elution. The plate has a max. working volume of 150 µL, and can be used after validating the workflow for elution with lower buffer volumes. Alternatively, use ELUTE PLATE (see Pre-Requisites section below).	RT
ADAPTER PLATE		1x				Enables cartridges to be placed on top of 96-well plates.	RT
ADAPTER		8x				Enables a cartridge to be placed into a tube.	RT

## Pre-Requisites

Common lab equipment is required for the sample preparation.

Equipment	Quantity and Description
PIPETTE	Careful sample handling and pipetting reduces contaminations and improves quantification.
SAMPLE	Peptide solution (20–100 µL) generated from tryptic digest of 1–100 µg protein starting material. For other sample types, contact PreOmics for adapted protocols.
96 WELL PLATES	96 deep well & 96 well skirted plates to balance WASTE & COLLECTION PLATES in centrifuge.
CENTRIFUGE	Swing-bucket centrifuges are required for peptide loading, washing, and elution.
ELUTE PLATE	96 deep well plate with >250 µL for elution of peptides from cartridge (e.g., Eppendorf Deepwell plate 96/500 µL Protein LoBind®, catalogue number EU: 0030504100, US: 951032107).
SEALING MAT	Prevents sample contamination and evaporation (e.g., Eppendorf Sealing Mat, catalogue number EU&US: 0030127978).
VACUUM EVAPORATOR	To evaporate volatile buffers from the eluate before LC-MS.
ULTRASONIC BATH	Optional: can be used to resuspend peptides.

## Procedure



## Method



### 1. LOAD

- 1.1. Mix 20–100 µL SAMPLE with equal volume of **STOP** ●. *\*NOTE1\**
- 1.2. Place **ADAPTER PLATE** with **CARTRIDGE** on **WASTE PLATE**. Label plate and wells, and transfer the sample to the **CARTRIDGE**. *\*NOTE2\**

### 2. PURIFY

- 2.1. Spin **CARTRIDGE** in a **CENTRIFUGE** (2,250 rcf; 1–3 min). If needed, adjust time to ensure complete flow-through.
- 2.2. Add 200 µL **WASH X** ○ to **CARTRIDGE**, repeat step 2.1, discard flow-through.
- 2.3. Repeat step 2.2 twice.
- 2.4. Add 200 µL **WASH 1** ● to **CARTRIDGE**, repeat step 2.1, discard flow-through.
- 2.5. Repeat step 2.4 once.
- 2.6. Add 200 µL **WASH 2** ● to **CARTRIDGE**, repeat step 2.1, discard flow-through. *\*SP\**

### 3. ELUTE

- 3.1. Discard **WASTE PLATE**. Use **ADAPTER PLATE** to place **CARTRIDGE** on top of the ELUTE PLATE. Label plate and wells.  
Alternatively, use the supplied **COLLECTION PLATE** after validating the workflow for elution with lower buffer volumes.  
See "Kit Contents" for details.
- 3.2. Add 100 µL **ELUTE**  to **CARTRIDGE**, spin **CARTRIDGE** in a **CENTRIFUGE** (2,250 rcf; 1–3 min).  
Keep flow-through in ELUTE PLATE.
- 3.3. Repeat step 3.2, keep flow-through in the same ELUTE PLATE.
- 3.4. Discard **CARTRIDGE** and place and place ELUTE PLATE in a **VACUUM EVAPORATOR** (45°C; until completely dry). **\*SP\***
- 3.5. Reconstitute peptides by adding **LC-LOAD**  to ELUTE PLATE. Adjust the volume according to specific requirements.  
For example, add 50 µL **LC-LOAD** to 100 µg protein starting material.
- 3.6. Sonicate ELUTE PLATE in an **ULTRASONIC BATH** (5 min) or shake (RT; 500 rpm; 5 min).
- 3.7. Spin plate in a **CENTRIFUGE** as follows:
  - User-provided ELUTE PLATE: RT; maximum rcf recommended by manufacturer; 5–15 min.
  - **COLLECTION PLATE**: RT; 2,250 rcf; 15 min.  
Transfer the supernatant to a clean plate and avoid touching the bottom of the well during transfer. **\*NOTE3\***

#### **\*NOTE1\***

The **SAMPLE-STOP** mix is expected to be acidic (pH<3.0) when the peptide sample was stored in a buffer commonly used for tryptic digestion. If the pH is above 3.0, contact [info@preomics.com](mailto:info@preomics.com) for an optimized protocol. The pH of the mixture can be measured with pH paper (the **STOP** pH itself cannot be determined using a pH measurement device or pH paper due to the high concentration of organic solvent).

#### **\*NOTE2\***

Alternatively, use **ADAPTERS** to place single **CARTRIDGES** on top of 2 mL tubes.  
Refer to the corresponding Phoenix 4x kit protocol for subsequent steps.

#### **\*NOTE3\***

At this point, peptide concentration can be measured, or the sample can be directly injected for LC-MS analysis. Visit our FAQ for recommendations on peptide quantitation assays. The silicone mat provided with the kit is incompatible with autosamplers.

#### **\*SP\* - Storage Point:**

At this point, close the peptide containing plate or **CARTRIDGE** using the silicone mat.  
Peptides can be frozen at -20°C for two weeks. Dried peptides, prior to reconstitution in **LC-LOAD**, can also be stored long-term at -80°C.

For answers to frequently asked questions, please visit our FAQ page at [www.preomics.com/faq](http://www.preomics.com/faq).

For our General Terms and Conditions, refer to [www.preomics.com](http://www.preomics.com).

For trademark information, visit [www.preomics.com/legal/trademarks](http://www.preomics.com/legal/trademarks).