



















Introduction

Sample preparation is one of the essential steps of bottom-up proteomics. The PreOmics® iST sample preparation kit is designed to streamline this process, enabling researchers to achieve optimal results with reduced steps and hands-on time. For specific recommendations and optimization, visit www.preomics.com/resources or contact info@preomics.com.

Kit Contents

The kit includes all essential components for proteomic sample preparation: denaturing, reducing, and alkylating agents, enzymes, cartridges, and wash buffers for peptide cleanup.

Component	Cap	Quantity	Buffer Properties				Description	Storage
iST 8x (P.O.00001)			Organic	Acidic	Basic	Neutral		
DIGEST		2x					Trypsin/LysC mix to digest proteins.	-20°C
RESUSPEND		1x 2 mL					For reconstituting lyophilized proteolytic enzymes.	RT
LYSE		1x 1 mL					For denaturing, reducing, and alkylating proteins.	RT
STOP		1x 1 mL					For stopping enzymatic activity.	RT
WASH 1		1x 2 mL					For removing hydrophobic contaminants.	RT
WASH 2		1x 2 mL					For removing hydrophilic contaminants.	RT
ELUTE		1x 2 mL					For eluting peptides from the cartridge.	RT
LC-LOAD		1x 1 mL					For loading peptides on reversed-phase LC-MS column.	RT
CARTRIDGE		8x					Cartridges with SPE sorbent for peptide purification from 1–100 µg protein starting material.	RT
WASTE TUBE		8x					2 mL tube for collecting waste after washing steps.	RT
COLLECTION TUBE		8x					1.5 mL tube for collecting peptides after elution.	RT
ADAPTER		8x					Enables a cartridge to be placed into a tube.	RT

Pre-Requisites

Common lab equipment is required for the sample preparation.

Equipment	Quantity and Description
PIPETTE	Careful sample handling and pipetting reduces contaminations and improves quantification.
SAMPLE	Mammalian tissue, 1–2 mg wet weight. For most mammalian tissue types, 1–2 mg of wet-weight tissue sample is equivalent to 100–200 µg protein starting amount. Visit our FAQ for more information on tissue starting amounts.
REACTION TUBE	1.5 mL microcentrifuge low protein binding tubes are recommended (e.g., Eppendorf Protein LoBind® Tubes, catalogue number EU: 0030108116, US: 022431081)
SONICATOR	For tissue homogenization using bead-based sonication (e.g., Diagenode Bioruptor®). Alternatively, use BeatBox® for tissue homogenization. Consult the BeatBox protocols at www.preomics.com/resources for details.
GLASS BEADS	Protein extraction glass beads (Diagenode #C20000021; Ø<1 mm) to facilitate tissue lysis.
HEATING SHAKER	Two heating shakers for tubes are recommended for protein denaturation and digestion.
CENTRIFUGE	Benchtop centrifuge for 1.5 or 2 mL tubes is required for loading, washing, and elution.
VACUUM EVAPORATOR	To evaporate volatile buffers from the eluate before LC-MS.
ULTRASONIC BATH	Optional: can be used to resuspend peptides.

Procedure



Method

1. LYSE

- 1.1. Place **SAMPLE** in a clean **REACTION TUBE**. Add 40–50 mg **GLASS BEADS** and 100 µL **LYSE** to sample. ***NOTE1***
- 1.2. Shear sample in a **SONICATOR** (10 cycles; 30 sec ON/OFF).
- 1.3. Place sample in a pre-heated **HEATING SHAKER** (95°C; 1,000 rpm; 10 min). For harder tissue like heart or muscle, repeat steps 1.2 and 1.3 once.
- 1.4. Optional: Spin down droplets (RT; max. 300 rcf; 10 sec).

2. DIGEST

- 2.1. Determine protein concentration of homogenate. Visit our **FAQ** for recommendations on protein quantitation assays. Proceed with digestion of up to 100 µg of extracted protein in a maximum total volume of 50 µL (if the sample volume is less than 50 µL, adjust to 50 µL with **LYSE**).
- 2.2. Optional: Spin down lyophilized enzyme mix in the **DIGEST** tube (RT; max. 300 rcf; 10 sec).
- 2.3. Add 210 µL **RESUSPEND** (yellow circle) to **DIGEST** (red circle) (1 tube for 4 reactions), shake (RT; 500 rpm; 10 min) and pipette up/down. ***NOTE2***
- 2.4. Add 50 µL **DIGEST** (red circle) to **REACTION TUBE** and place it in a pre-heated **HEATING SHAKER** (37°C; 500 rpm; 3 h).
- 2.5. Add 100 µL **STOP** (black circle) to **REACTION TUBE** (precipitation may occur), shake (RT; 500 rpm; 1 min) and pipette up/down. ***SP***
- 2.6. Spin sample in **CENTRIFUGE** (16,000 rcf; 1 min).
- 2.7. Assemble **ADAPTER** and **CARTRIDGE** to the top of **WASTE TUBE**. Label all tubes and transfer supernatant to the **CARTRIDGE**.

3. PURIFY

- 3.1. Centrifuge the **CARTRIDGE** (RT; 3,800 rcf; 1–3 min). If needed, adjust time to ensure complete flow-through.
- 3.2. Add 200 µL **WASH 1** (blue circle) to **CARTRIDGE**, repeat step 3.1.
- 3.3. Add 200 µL **WASH 2** (green circle) to **CARTRIDGE**, repeat step 3.1. ***SP***
- 3.4. Transfer assembled **ADAPTER** and **CARTRIDGE** onto clean **COLLECTION TUBE**. Discard **WASTE TUBE**. Label all tubes.
- 3.5. Add 100 µL **ELUTE** (pink circle) to **CARTRIDGE**, repeat step 3.1, keep flow-through in **COLLECTION TUBE**.
- 3.6. Repeat step 3.5, keep flow-through in the same **COLLECTION TUBE**.
- 3.7. Discard **CARTRIDGE** and place **COLLECTION TUBE** in a **VACUUM EVAPORATOR** (45°C; until completely dry). ***SP***
- 3.8. Reconstitute peptides by adding **LC-LOAD** (white circle) to **COLLECTION TUBE**. Adjust the volume according to specific requirements. For example, add 50 µL **LC-LOAD** to 100 µg protein starting material.
- 3.9. Sonicate **COLLECTION TUBE** in an **ULTRASONIC BATH** (5 min) or shake (RT; 500 rpm; 5 min).
- 3.10. Spin **COLLECTION TUBE** in a **CENTRIFUGE** (RT; 16,000 rcf; 5 min). Transfer the supernatant to a clean vial and avoid touching the bottom of the **COLLECTION TUBE** during transfer. ***NOTE3***

- *NOTE1*** Alternatively, use BeatBox® for tissue homogenization and lysis. Consult the BeatBox protocols at www.preomics.com/resources for details.
- *NOTE 2*** Resuspended DIGEST can be stored for up to two weeks at 4°C. For longer storage periods, visit our FAQ.
- *NOTE 3*** At this point, peptide concentration can be measured, or the sample can be directly injected for LC-MS analysis. Visit our FAQ for recommendations on peptide quantitation assays.
- *SP* - Storage Point:** At this point, close the peptide containing **COLLECTION TUBE** or **CARTRIDGE** using the silicone lid. Peptides can be frozen at -20°C for two weeks. Dried peptides, prior to reconstitution in **LC-LOAD**, can also be stored long-term at -80°C.

Data analysis

Consider the following as fixed modifications in your database search:

MODIFICATION	DESCRIPTION	COMPOSITION	SPECIFICITY	MASS	UNIMOD #
ALKYLATION	Carbamidomethyl on cysteine	C ₂ H ₃ NO	[C]	+57Da	4

For answers to frequently asked questions, please visit our FAQ page at www.preomics.com/faq.

For our General Terms and Conditions, refer to www.preomics.com.

For trademark information, visit www.preomics.com/legal/trademarks.