

# Quantitative and precise plasma proteome analysis with enhanced coverage

Maura Grey<sup>1</sup>, Zehan Hu<sup>2</sup>, Katharina Limm<sup>2</sup>, Sandra Schär<sup>3</sup>, Arthur Viode<sup>3</sup>, Roland Bruderer<sup>3</sup>

<sup>1</sup>PreOmics Inc., Billerica, US | <sup>2</sup>PreOmics GmbH, Martinsried, Germany | <sup>3</sup>Biognosys AG, Zurich, Switzerland

## INTRODUCTION

Plasma is a valuable source for the discovery of biomarkers, as it contains proteins such as cytokines and tissue leakage markers that reflect the health status of patients. However, mass spectrometry (MS)-based plasma proteomics is challenged by the vast dynamic range of protein concentrations and signal suppression from high-abundance proteins, which can obscure low-abundance analytes. Enrichment strategies improve coverage but often distort the original protein concentrations, making it difficult to interpret the quantification. The new combined P2 paramagnetic nano-particle and iST technology further extends proteome coverage while preserving quantitative accuracy. In this study, we evaluated its performance using spike-in experiments with different plasma sources.

## MATERIALS & METHODS

**Samples:** A species dilution experiment was performed by mixing K2EDTA bovine plasma (Neo Biotech) with pooled K2EDTA human plasma (Proteogenex) at seven defined volume ratios (1:0 to 0:1), each in triplicate.

For protein-specific spike-in experiments, purified IL-6 (Absea, cat. PP-496), IL-1 $\beta$  (Absea, cat. PC-490), and MAPT (Absea, cat. PC-006) were added to 100  $\mu$ L K2EDTA human plasma to generate a five-fold serial dilution series from 1250 ng to 0.4 ng, each in triplicate.

**Sample preparation:** Neat plasma (2  $\mu$ L) was processed using the iST-BCT protocol. P2-iST Plasma samples were prepared from 100  $\mu$ L plasma on a KingFisher™ Flex, including enrichment, lysis, digestion, and peptide cleanup.

**LC-MS analysis:** 800 ng peptides were analyzed using the Evosep™ One system with Evotip Pure™ and an Aurora Elite CSI 75  $\mu$ m  $\times$  150 mm C18 column, coupled to a timsTOF HT mass spectrometer operating in diaPASEF® mode (Whisper™ Zoom 40 SPD).

**Data processing:** Raw files were analyzed in Spectronaut® 20 using Pulsar-generated spectral libraries; fold changes were derived across the dilution series. Specific normalization and noise filtration strategies were applied.

## KEY TAKEAWAYS

### High precision at greater depth

P2-iST Plasma markedly expands proteome coverage while maintaining low variability, enabling precise quantification for many more proteins than neat plasma.


### Accurate relative quantification

Fold-change accuracy is preserved across multispecies dilutions and spike-ins, with quantitative performance comparable to neat plasma at the proteome level.


### Linear signal-concentration relationship

P2-iST Plasma provides robust linearity and improved LOD/LOQ across a broad concentration range, supporting both discovery and targeted proteomics.


## RESULTS




**Experimental design of (1) species dilution and (2) protein-specific spike-ins for benchmarking the quantitative performance of the P2-iST Plasma workflow.**



**Enhanced proteome coverage and precision using P2-iST Plasma in a species dilution study.** (A) P2-iST Plasma enables exceptional plasma dynamic range (7–8 orders of magnitude) by enriching low-abundance proteins, increasing the number of precisely quantified proteins (CV <20%) by 5–7-fold compared with neat plasma, depending on species and mixing ratio. These results demonstrate the high repeatability of the P2-iST Plasma workflow relative to neat plasma.



**High relative quantitation accuracy demonstrated by P2-iST Plasma.** (B) Both P2-iST Plasma and neat plasma samples show high fold-change (FC) accuracy and strong quantitative concordance at the full-proteome level. (C) Fifteen commonly detected FDA-approved biomarker bovine orthologues exhibit accurate FC trends in both workflows, with a single outlier (LDHB, Q5E9B1) in the P2-iST Plasma dataset, potentially attributable to carryover from randomized acquisition order or high sequence conservation between bovine and human proteins. (D) Protein-specific spike-in responses for IL1B, IL6, and MAPT show consistent and accurate FCs derived from pairwise dilutions, confirming that P2-iST Plasma supports robust quantitative measurement of cytokines and key biomarkers relevant to neurodegenerative diseases.



**High linearity delivered by P2-iST Plasma.** (E) Matrix-matched, calibration-based LOD/LOQ analysis demonstrates that P2-iST Plasma achieves higher sensitivity and improved linearity compared with neat plasma, enabling detection and quantification of ~3.5-fold more proteins, particularly at intermediate dilution levels. (F) Spiked-in cytokines IL1B and IL6, as well as the biomarker MAPT, are detectable at tens of pg/ $\mu$ L due to low background noise and are quantitatively measured from hundreds of pg/ $\mu$ L up to at least 12.5 ng/ $\mu$ L with a highly linear response, highlighting the suitability of P2-iST for targeted proteomics applications.

## CONTACT & MORE

Contact:  
Maura.Grey@preomics.com

**Conflict of Interest Disclosure**  
Grey, M. is employed by PreOmics Inc.  
Hu, Z., Limm, K. are employed by PreOmics GmbH.  
Schär, S., Arthur, V., Bruderer, R. are employed by Biognosys AG.

**References:**  
• Pernemalm, M. et al. eLife, 8:e41608 (2019).  
• Soni, R.K. Clin Proteom, 21, 43 (2024).  
• Huang, T. et al. bioRxiv (2023). <https://doi.org/10.1101/2023.08.28.555225>  
• Pino, L.K. et al. J. Proteome Res. 19 (3), 1147 (2020).

