

## nAspro & nAspro BBY





(Nano Aspheric Precision Optics)

Hydrophilic Monofocal Aspheric Intraocular Lens  
with Delivery System

### DESCRIPTION OF THE DEVICE

**NASPRO** Aspheric lens is an optical Implant, designed to be surgically implanted in the human eye as a replacement of the crystalline lens for visual correction of aphakia.

### MATERIAL AND PRODUCT CHARACTERISTICS

Material Characteristics		
Brand	Naspro	Naspro BBY
Model	NAS207	NASY207
Material	Copolymer of Hydroxy Ethyl Methacrylate and Methyl methacrylate	
Colour	Clear	Yellow
Water Content	26%	
Refractive Index	1.46	
Minimum Resolution Efficiency	70%	
YAG Laser Compatible	Yes	
Light Transmittance	NAS207	
	NASY207	
Product Characteristics		
Optics Design	Bi Convex / Convex-Concave	Bi Convex / Convex-Concave
Haptic Design	Dual Haptic	Dual Haptic
Drawing		
Optic Diameter	6.00mm	6.00mm
Overall Length	12.50 mm	12.50 mm
Power Range	-15.00 D to +53.00 D (In steps of 0.5D)	-15.00 D to +53.00 D (In steps of 0.5D)

### INTENDED USE AND MODE OF ACTION

**NASPRO** Aspheric IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This position allows the lens to function as a refractive medium in the correction of aphakia.

Throughout life, the crystalline lens protects the retina from hazardous UV light. Cataract extraction removes this protective barrier and may accelerate retinal photo toxicity to control the amount of radiation transmitted to the retina after cataract surgery, IOLs were created with the capability of blocking UV radiation. All IOLs today block short wavelength UV radiation up to 400nm.

Yellow IOLs approximate the UV and blue-light attenuating properties of the human crystalline lens. This IOL contains covalently bound chromophores, that providing increased protection against harmful UV light while allowing optimal blue light transmission, improving contrast sensitivity in low light conditions especially in eyes at risk of age-related macular degeneration.

Blue-absorbing chromophores added to IOL material to expand protection and decrease the amount of UV-blue energy (400 nm to 500 nm range) reaching the retina. It is thought that adding a blue-blocking chromophores increased protection from photic retinopathy, because it simulated the properties of the natural crystalline lens.

### SPHERICAL ABERRATION

When parallel rays of incident light pass through an ideal lens they are brought to a common focus. When such rays pass through a spherical lens of positive power, rays near the lens periphery are brought to a focus closer to

the posterior lens surface than those rays passing through, or near, the optical centre. This condition is termed positive spherical aberration. In the young eye the naturally occurring positive spherical aberration of the cornea is compensated by the naturally occurring negative aberration of the crystalline lens. As age advances the crystalline lens becomes more positive and does not fully compensate for the spherical aberration of the cornea.

The anterior surface of the NASPRO Aspheric IOL is modified to have a negative aspheric surface. This refracts all parallel rays incident on the anterior lens surface in such a way that they are brought to a common focus. NASPRO Aspheric IOLs can therefore improve contrast sensitivity and functional visual activity in the pseudophakic patient in comparison to a standard spherical lens.

### INDICATIONS

**NASPRO** Aspheric IOL indicated for the replacement of the human lens to achieve visual correction of aphakia following cataract ( congenital or traumatic or senile cataract) surgery. These lenses are intended for placement in the capsular bag.

### CONTRAINDICATIONS

Apart from non-specific contraindications related to any form of ocular surgery, the following specific contraindications must be respected:

- Microphthalmia
- Active ocular diseases (e.g chronic severe uveitis, proliferative diabetic retinopathy, chronic glaucoma not responsive to medication)
- Corneal de-compensation or corneal endothelial cell insufficiency
- Corneal Opacification

### PRECAUTIONS

As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: Temporary Corneal Edema, Secondary Cataract Formation, Vitreous Herniation into the Formation, Mal-positioned Lens, Aphakic Glaucoma, Temporary Flat Anterior Chamber, Retinal Detachment, Lens Implants Loop Amputation, Endothelial Conical Dystrophy Anterior Chamber, Ocular Infection, Pupillary Block, Corneal Dystrophy and IOL Opacification

### WARNINGS

1. Do Not Resterilize. Resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.
2. Do not store below 5° C and above 40° C to avoid shocks & fragile.
3. Lens should not be altered in any manner.
4. Lens should not be repackaged by anyone.
5. Single use IOLs and single use injectors cannot be reused, as they are not designed to perform as intended after the first and only usage to avoid infection.
6. Do not use after the expiry date.
7. Non-toothed, polished instruments must be used if handling the IOL.
8. Do not allow the IOL to contact substances that are unsterile or ocular-incompatible prior to placement into the eye.
9. Do not use unsterile surgical instruments or instruments that may carry a risk of contamination.
10. Do not allow the IOL to dehydrate during the procedure.
11. Once closed, do not reopen the flaps of the injector.
12. Following implantation, irrigate/aspirate to eliminate any viscoelastic residues from the bag, especially between the IOL and posterior capsule.

### CALCULATION OF LENS POWER

Accurate Keratometry and biometry is essential to successful visual outcomes. Preoperative calculation of the required lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. The A-constant listed on the outer label is presented as a guideline and is a starting point for implant power calculations. This A-constant has been theoretically derived. Lens constants must be "Personalized" to compensate for the differences in instrumentation, measurement technique, and IOL power calculation methods. A convenient initial estimate can be obtained by referencing to the personalized lens constant for a similar lens model.

### DISCLAIMER OF LIABILITY

AI Optics (P) Ltd., accepts no liability for improper model selection by the physician, for improper handling, use, surgical technique applied or for any other iatrogenic error caused by the implanting surgeon. AI Optics (P) Ltd., makes no expressed or implied warranties in connection with sale of the IOL.

### STERILIZATION

Sterilization is done By Moist Heat Sterilization Method. Validity of sterilization applies as long as the sealed inner peel-pouch is not disturbed or damaged. Any damage to the peel-pouch or any accidental opening of the peel-pouch is to be declared as "NOT STERILE".

### PATIENT IDENTITY CARD

The packaging contains product identification stickers for maintaining a record of the IOL implantation. Surgeons are requested to give the "Patient ID Card" to the patient after implantation and advise them to carry the card at all times.

### DISPOSAL OF DEVICE

There are NO specific guidelines for disposing of Non sterile or contaminated medical devices. Follow standard procedure for discarding outdated or contaminated products.

### RETURN GOODS POLICY

AI Optics (P) Ltd accepts returned lenses for exchanges only. NO cash refunds will be issued. To return lenses, you must first obtain a return authorization number from customer services department. No returned goods will be accepted without proper authorization number. Returned lenses should be shipped by traceable method. No credit will be given to lost or damaged lenses in shipment. Lenses will be replaced as long as they are returned within six months of their original invoice date.

### CUSTOMER FEEDBACK AND ADVERSE EVENTS REPORTING

Please report any feedback, adverse events/complaints to AI OPTICS (P) LTD., 84/1, Aavin Dairy Farm Road, SIDCO Industrial Estate, Ambattur, Chennai-600098, INDIA. Email: info@appasamy.com.

### HOW SUPPLIED

AI OPTICS (P) Ltd, IOLs are supplied with Delivery system depends on customer requirement.

The Package Contains

One Sterile Intra Ocular Lens, One sterile Lens Delivery system, Patient ID Card, Labels and Information for use (IFU).


AI OPTICS (P) Ltd recommends using the qualified IOL Delivery System for patient's safety. The use of an unqualified Delivery system may cause damage to the lens and potential complications during the implantation process.

If use AI OPTICS (P) Ltd qualified Lens Delivery system, please follow the "Instruction for Use of Delivery System" (see below).

**NOTE:** To avoid dehydration, leave lens immersed in saline until ready to fold and implant. The lens should be folded and implanted within 3 minutes after removal from the saline.

### INSTRUCTION FOR OPENING OF PRIMARY PACK

After opening the outer box, examine the device package for any damage, and verify that information (lens model, power, serial number, Expiry Date, etc.) Provided on the device primary pack is same as the information provided on the outer box label. Peel the Primary Pack and hold the PP cup in sealing free side by one hand. Slowly peel off the aluminium foil by other hand.




### INSTRUCTION FOR USE OF AI OPTICS DELIVERY SYSTEM

Sterilization of Delivery system is done By EO Sterilization Method.


#### STEP: 1

Open the cartridge using thumb and index finger. It is most convenient to load the lens.




#### STEP: 2

Apply a thin layer of viscoelastic in the lens loading chamber and in the barrel to serve as a lubricant for the free movement of lens through the cartridge during injection. Please Refer "IOL loading Method" for place the lens into the chamber of the cartridge.




### IOL LOADING METHOD

OPTIC DESIGN	HAPTIC DESIGN	METHOD
Bi-Convex (More than 7.0 D) & Convex-Concave (Less than 7.0 D)	Dual Haptic	The lens should be placed into the cartridge barrel as Reverse 'E'. The anterior side of the lens should be like this 


#### STEP: 3

Use flat forceps to gently press the lens and at the same time, partially close the cartridge. Then gently push the leading and trailing haptic over the optic and close cartridge.




#### STEP: 4

If it is unable to close the cartridge, and then open the cartridge, remove the lens and reload the lens following steps 2, 3 and 4. Make sure neither the haptic nor optic is caught between the shutters.




#### STEP: 5

Ensure that there is no gap between the shutters before loading the cartridge into the injector. Load the cartridge into the injector through the groove provided.


















#### STEP: 6

Gently push the plunger until the lens moves towards the cartridge tip. The lens will completely come out of the cartridge.




### SYMBOLS USED ON LABELING

<b>D</b>	Diopter (Spherical equivalent)
<b>IOL</b>	Intraocular lens
	Sterilized using Steam
	Do not Re-use the product
	Do not Re-sterilize the product
<b>SN</b>	Serial number
	Do not use if package is damaged
	Store between 5°C to 40°C
	Consult Instructions for use
	Keep away from sunlight
	Keep dry
	Sterilized using Ethylene oxide
	Manufacturer
	Date of Manufacture
	Expiry Date
<b>LOT</b>	LOT Number
<b>REF</b>	Reference (Model) Number
<b>ØB</b>	Optic Diameter of IOL
<b>ØT</b>	Overall length of IOL
	Single sterile barrier system
	Double sterile barrier system
	Caution
<b>CE ###</b>	CE conformity Marking and the Notified body number
<b>EC REP</b>	Authorized Representative in the European Community/ European Union

**EC REP Amstermed B.V.**  
Saturnusstraat 46-62, Unit 032,  
2132 HB Hoofddorp, The Netherlands  
Email: info@amstermed.nl



Manufacturing Licence No. MFG/MD/2025/000107  
Manufactured by:

 **AI OPTICS PRIVATE LIMITED,**  
84/1, Aavin Dairy Farm Road,  
SIDCO Industrial Estate, Ambattur,  
Chennai – 600098, INDIA.