English

Paragonix SherpaPak® Cardiac Transport System Part#106,L-621 ver. 4

2026-03-18 **MD**

INSTRUCTIONS FOR USE

1. INDICATIONS FOR USE

The Paragonix SherpaPak® Cardiac Transport System is intended to be used for the static hypothermic preservation of hearts during transportation and eventual transplantation into a recipient using cold storage solutions indicated for use with the heart.

The intended organ storage time for the Paragonix SherpaPak® Cardiac Transport System is up to 4 hours.

Donor hearts exceeding clinically accepted static hypothermic preservation times should be evaluated by the transplant surgeon to determine transplantability in accordance with accepted clinical guidelines and in the best medical interest of the intended recipient.

1.1 CONTRAINDICATIONS

The Paragonix SherpaPak® Cardiac Transport System should not be used in the following circumstances:

1. When safe connection with the aorta cannot be made, e.g., due to lack of enough length of the aortic root.

1.2 WARRANTY STATEMENT

Paragonix Technologies, Inc., (Paragonix®) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling and storage of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Paragonix's control directly affect the device and the results obtained from its use.

Paragonix's obligation under this warranty is limited to the repair or replacement of this device and Paragonix shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device.

Paragonix does not authorize any other person to assume for it, any other or additional liability or responsibility in connection with this device. Paragonix assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.

2. SAFETY REQUIREMENTS

⚠ Caution: Donor hearts exceeding 4 hours static hypothermic preservation time will require transplant surgeon evaluation to determine transplantability.

Should any serious incident occur in relation to the device, please contact Paragonix and the

appropriate regulatory authority in the county of use where the incident occurred. Paragonix can be reached via your local sales rep or contacted via:

email: support@paragonixtechnologies.com
telephone: +1.781.428.4828

2.1 IMPORTANT INFORMATION





It is important that all personnel who will operate the Paragonix SherpaPak® Cardiac Transport System (Paragonix SherpaPak® CTS) read and understand these instructions for use before operating the device.

All personnel should follow all warnings and precautions outlined below, for their safety and the safety of those around them.

The Paragonix SherpaPak® Cardiac Transport System is a device for use by hospitals that perform cardiac recovery procedures to store and transport donor hearts to the transplant facility. The Paragonix SherpaPak® Cardiac Transport System is for use within acute care facilities that have an existing agreement with an Organ Procurement Organization (OPO) requiring the hospital to notify the OPO or third party designated by the OPO in a timely manner about all deaths and imminent deaths that occur in the hospital. Such hospitals must also have documented protocols and procedures for determining death or imminent death, obtaining family consent, and trained personnel with fully equipped surgical recovery operating rooms to perform the recovery procedure. The facility must also be able to provide Advanced Cardiac Life Support, should that be required.

2.2 SYMBOL DEFINITIONS

Table 1: Symbols used in the labeling of Paragonix SherpaPak® CTS and their definitions.	
Symbols	Definitions
	Use by YYYY- MM-DD
	Do not reuse
	Do not use if package is damaged or open
	Indicates the need to consult the instructions for use
	Batch code
	Caution, consult accompanying documents
	Manufacturer
	Medical Device
	Device or device component is non-sterile
	Sterile using irradiation
	Fill Port
	Vent Port
	Fragile, handle with care
	Keep Dry

	Temperature limits
	Pressure limits
	Humidity limits
	Unique Device Identifier
	Serial Number
	Model Number
	Country of Manufacture
	Date of Manufacture
	Authorized representative in the EU
	Authorized representative in Switzerland
	European Importer

2.3 WARNINGS ⚠

- **Caution: Federal (US) law restricts this device to sale by or on the order of a physician.**
- **Use aseptic technique as appropriate.** Certain components of the Paragonix SherpaPak® CTS are provided sterile. Maintain aseptic technique when handling sterile components of the Paragonix SherpaPak® CTS.
- **Chilled (4°C) preservation solution, cleared for use with the heart, must be available.**

⚠ Do NOT reuse any component of Paragonix SherpaPak® CTS. SherpaPak® CTS is intended for single use only. DO NOT RE-USE. Re-use of the device may result in infection and other complications due to the loss of sterility. Certain components of the Paragonix SherpaPak® CTS are sterile as supplied (sterilization method is gamma irradiation). SherpaPak® should be disposed of in accordance with local guidelines for biomedical waste.

- **Use institution-specific precautions with the donor heart and preservation solution when operating the Paragonix SherpaPak® CTS.** The heart and preservation solution may carry undetected pathogens from the donor. Use universal precautions for bloodborne pathogens in handling the heart, and in handling and disposing of the Paragonix SherpaPak® CTS and preservation solution to prevent the possible transmission of pathogens to personnel. As appropriate, this may include use of personal protective equipment (e.g. gloves, masks, gowns, goggles or equivalent eye protection) and disposal of materials as potentially infectious biohazard waste.
- **Prior to use, inspect all components of the Paragonix SherpaPak® CTS. Do not use if any component is loose, broken, or damaged.**
- **Do not open the Paragonix SherpaPak® CTS during organ transport.**

- **No modification of the Paragonix SherpaPak® CTS is allowed.**

2.4 OPERATING CONDITIONS ⚠

Operating Temperature: 22°C
Recommended Operating Room temperature is 22°C. Direct sunlight and outdoor temperature extremes (high and low) can affect Paragonix SherpaPak® CTS internal temperature (4 to 8°C). During exposure to temperature extremes, Paragonix SherpaPak® CTS temperature must be frequently monitored.

Operating Pressure: Sea Level to 8000 ft.

The operating pressure range conforms to commercial airliner transport. Extreme pressure levels can impact performance. During exposure to extreme pressures (altitudes), Paragonix SherpaPak® CTS operation must be frequently monitored.

Operating Humidity: 40-60% Relative Humidity

Extreme humidity levels can impact performance. During exposure to extreme humidity levels, Paragonix SherpaPak® CTS operation must be frequently monitored.

2.5 ACCURACY ⚠


Temperature Accuracy
The Paragonix SherpaPak® CTS comes with a pre-installed datalogger capable of reporting both the ambient temperature and the temperature within the Organ Assembly.

Temperature Accuracy is ±0.5°C from 0° to 50°C.

Time Accuracy
The pre-installed datalogger logs temperature against time with an accuracy of ±1 minute per month.

2.6 PRECAUTIONS ⚠

Do not wet temperature probe connectors on Paragonix SherpaPak® CTS
A series of three (3) key areas are identified on the Paragonix SherpaPak® CTS as being regions to "Keep Dry". These regions are called out by the following label:

These labels are affixed in the below-indicated locations, and the warning must be followed in order to ensure proper device function.

Table 2: Identification of temperature probe connector locations on Paragonix SherpaPak® CTS to be protected from water.

Shipper (1)	
Organ Lid	
Shell Base	

During Preparation

Should the user notice that during preparation of the system any of the connectors become wet, the Organ Canister and Shell Assembly should be discarded and replaced with a new Organ Canister and Shell Assembly.

Note: It is not required to re-attach the heart to a new Heart Connector.

Replacement of the Organ Canister and Shell Assembly should be performed in accordance with Sections 4.7 and 4.8 below.

After Preparation

If after complete packaging of the donor heart into the shipper and activation of the datalogger, the temperature reading is above 15°C, this indicates that the connectors have been wetted during the preparation steps. If this occurs, the user is to remove the Organ Canister and Shell Assembly from the Shipper and replace with a new Organ Canister and Shell Assembly.

Note: It is not required to re-attach the heart to a new Heart Connector. Replacement of the Organ Canister and Shell Assembly should be performed in accordance with Sections 4.7 and 4.8 below.

Keep the Paragonix SherpaPak® CTS primarily upright during transportation.

The Paragonix SherpaPak® CTS is designed to be transported upright. Temporary tilting ±45° from horizontal in any direction is acceptable.

Avoid direct sunlight and hot or cold temperature extremes

The Paragonix SherpaPak® CTS is designed to be transported under the same environmental conditions as is appropriate for people. Avoid extended exposure to outdoor conditions (sunlight, heat, or cold).

Ensure adequate connection between aorta and Heart Connector.

For proper heart preservation, the aorta must be connected correctly (Section 4.7).

Use Caution when lifting the Paragonix SherpaPak® CTS

A fully loaded Paragonix SherpaPak® CTS weighs ~30 lbs. Use proper lifting practices.

Exercise Caution when using in presence of other electronic devices

While no significant risk of reciprocal interference has been identified, Paragonix SherpaPak® CTS should always be closely monitored when in use.

2.7 SIDE EFFECTS

All surgical procedures and medical devices have potential risks. The potential surgical risks of a transplant with a donor heart are similar for the Paragonix SherpaPak Cardiac Transport System and traditional ice storage preservation. There is a risk of receiving a heart that does not function properly after transplant. There is also a risk that the donor heart may be damaged during preservation.

2.8 SIDE EFFECTS ASSOCIATED WITH PARAGONIX SHERPAPAK CTS

- It is possible that after preservation using the Paragonix SherpaPak CTS, the transplant doctor may decide that the donor heart is not suitable for transplantation.
- The Paragonix SherpaPak CTS is prepared and operated by trained medical professionals. Like with many medical technologies, there are inherent risks including injury to the aortic valve leaflets and coronary orifices, vascular perforation, infection or delayed organ function.

2.9 CLINICAL BENEFITS (EU)

The Paragonix SherpaPak Cardiac Transport System is designed to improve organ preservation during transportation and storage from organ donor to recipient. Improved organ preservation can enhance donor organ quality, leading to reduced primary graft dysfunction (PGD) and improved longterm outcomes.

3. PRE-ASSEMBLY CHECKLIST

Table 3: Pre-Assembly Checklist for Paragonix SherpaPak® CTS for identifying component counts, locations, and sterility disposition.

Box	Item (Quantity)
LARGE Non Sterile	Shipper(1)
SMALL Non Sterile	SherpaCool® Pouch (1) ⚠ CAUTION: PRECONDITION AT -20°C OR BELOW FOR AT LEAST 48 HOURS PRIOR TO USE.
	SherpaCool® Ribbon (3) ⚠ RIBBONS CAUTION: PRECONDITION AT -20°C OR BELOW FOR AT LEAST 48 HOURS PRIOR TO USE.

MEDIUM Sterile	Organ Canister (1)
	Organ Lid (1)
	Shell Base (1)
	Shell Top (1)
	Heart Connector (1) 10 mm
	Heart Connector (1) 18 mm
Heart Connector (1) 25.4 mm	
Heart Connector (1) 32 mm	

3.1 PARAGONIX SHERPAPAK® CTS SET-UP AND PREPARATION

- Inspect all parts on arrival for any signs of damage that may have occurred during transport.
- Report any damage or concerns about the condition of the Paragonix SherpaPak® CTS immediately to Paragonix Technologies Inc.®
- When in transit, the SherpaCool® Pouches and Ribbons may be stored in the shipper for up to 18 hours.

Note: If SherpaCool® and Ribbons are set up in the shipper, place the ribbons in their intended locations as shown in section 4.6. Place the pouch in the center of the shipper.












Figure 14: Preservation Solution
⚠ CAUTION: FOLLOW DIRECTIONS IN THE ENCLOSED IFU REGARDING STORAGE AND PREPARATION FOR USE

4. OPERATING INSTRUCTIONS
4.1 GENERAL INFORMATION
Before using in a clinical setting, operators must be trained in the use and functional understanding of the Paragonix SherpaPak® CTS.

4.2 OVERVIEW USING THE PARAGONIX SHERPAK® CTS

Using the Paragonix SherpaPak® CTS involves performing the following procedures:

1. Removal of packaging, preconditioning of SherpaCool® Box (Section 4.3) containing SherpaCool® Pouch and SherpaCool® Ribbons at or below -20 °C for at least 48 hours.
2. Setup of datalogging software for use with Paragonix SherpaPak® CTS (Section 4.4, optional, recommended)
3. Transport of Paragonix SherpaPak® CTS to the recovery site (Section 4.5)
4. Preparing the Paragonix SherpaPak® CTS and Shipper for deployment at recovery site (Section 4.6)
5. Preparing the Paragonix SherpaPak® CTS for heart receipt at recovery site (Section 4.7)
6. Heart recovery, aortic connection, and preservation (Section 4.8)
7. Traveling with the Paragonix SherpaPak® CTS to the transplant site (Section 4.9)
8. Removing the heart from Paragonix SherpaPak® CTS for transplant (Section 4.10)

These instructions may be modified based on actual use. The instructions are designed for implementation by two operators. These instructions may be modified for a single operator, provided that proper aseptic technique is used in conjunction with procedures to be performed on the sterile field with sterile Paragonix SherpaPak® CTS components

4.3 REMOVAL OF PACKAGING AND PRECONDITIONING OF SHERPACOOL® BOX CONTAINING SHERPACOOL® RIBBONS AND SHERPACOOL® POUCHES

Paragonix SherpaPak® CTS must be maintained in a ready-to-use condition, so that it can be available to the heart recovery team at all times.

Make the following preparations:

- Do not remove SherpaCool® Ribbons and Pouches from the SherpaCool® Box.
- Write the date and time of SherpaCool® Box transition to -20°C storage on the Label affixed to the SherpaCool® Box.
- Place the SherpaCool® Box into a -20°C (or colder) freezer for a minimum of 48 hours.
- Remove SherpaCool® Box after a minimum preconditioning time of 48 hours at -20°C (or colder).
- Do not remove SherpaCool® Box until recovery team departs to donor hospital.
- Do not remove SherpaCool® Box until all other components of Paragonix SherpaPak® CTS and associated equipment and materials have been prepared for transport.
- SherpaCool® Box must be transported to the donor site on ice OR Within the shipper as described in section 3.1 above.
- Chilled (4°C) preservation solution, cleared for use With the heart, must be available.

4.4 SETUP DATALOGGING SOFTWARE FOR USE WITH PARAGONIX SHERPAK® CTS (OPTIONAL, RECOMMENDED)

In order to recover logged temperature data from Paragonix SherpaPak® CTS, the user must first download and install the latest version of the mobile Paragonix App at either the App Store or Google Play store. This requires use of a Bluetooth® -enabled iOS or Android device.

NOTE: Paragonix SherpaPak® CTS comes pre-configured with appropriate temperature monitoring settings installed in the on-board temperature monitoring system. This configuration must NEVER be modified by the user. After setting up an account for use of the Paragonix App:

- After logger is started (Section 4.7, Step 9 of this manual) it is possible to visualize the logged data by connecting the Paragonix SherpaPak® CTS to the Paragonix App.
 - Log in to the Paragonix App on your mobile device.
 - Select the Datalogger matching the serial number shown at the back of the Shipper from the Paragonix software application.
 - When prompted, enter the Datalogger PassKey provided on the back of the Shipper adjacent to the datalogger serial number.


Figure 15: Indication of datalogger Serial Number and PassKey on Shipper

- Once connected to the datalogger, logged data may be viewed locally and/or uploaded to alternative devices via the Paragonix software application.
- For assistance with data recovery, contact Paragonix directly via provided telephone number and have available the Shipper Serial Number and Data Logger Serial Number.

- Do not dispose of Shipper unit until desired data has been obtained.


4.5 TRANSPORT OF PARAGONIX SHERPAK® CTS TO THE RECOVERY SITE

Transport various supplies of Paragonix SherpaPak® CTS to the recovery site by:


- Place SherpaCool® Box onto wet ice in a cooler for transportation OR in the shipper as described in section 3.1.
- SherpaCool® Box may be transported in the same cooler as other adjunct supplies (preservation solutions, etc.)
- Maximum allowable transit time for SherpaCool® is 4 hours for wet ice transport, 18 hours within the shipper.
- Carry Paragonix SherpaPak® CTS (i.e. sterile components) in a separate bag or, optionally, place on ice in cooler While still boxed.
- Unbox Shipper prior to transportation to the recovery site.
- Wheel or carry Shipper to recovery site.

4.6 PREPARATION OF PARAGONIX SHERPAK® CTS FOR DEPLOYMENT AT RECOVERY SITE


1. Place box containing Shell Assembly with nested Organ Assembly and Heart Connectors onto a non-sterile table in the operating room (OR). Follow hospital procedures for moving equipment into the OR.
2. Using checklist (See Section 3.1), double-check for presence of all Paragonix SherpaPak® CTS components PRIOR to unpacking the components.
3. Open SherpaCool® Box, leaving SherpaCool® Ribbons and SherpaCool® Pouch on ice until installation in the shipper.
4. Place one SherpaCool® Ribbon into the cavity of the Shipper. Ensure that SherpaCool® Ribbon is dry by wiping with an absorbent cloth as necessary prior to installation.


5. Press the SherpaCool® Ribbon out against the walls of the Shipper to form a ring. Ensure that SherpaCool® Ribbon is dry by wiping with an absorbent cloth as necessary prior to installation.

6. Repeat steps 5 and 6 as above twice more, inserting the remaining two SherpaCool® Ribbons into the Shipper in concentric rings. Ensure that SherpaCool® Ribbon is dry by wiping with an absorbent cloth as necessary prior to installation. A total of three (3) SherpaCool® Ribbons should be installed in the Shipper in this fashion. Place SherpaCool® Pouch in center of SherpaCool® Ribbons.

7. Close the lid of the Shipper and leave closed until the Organ and Shell Assembly are fully loaded and prepared for installation.

4.7 PREPARING THE PARAGONIX SHERPAK® CTS FOR HEART RECEIPT AT RECOVERY SITE

1. Open shipping box containing the nested Organ Assembly and Shell Assembly in the large sterilization bag and Heart Connectors in the small sterilization bags.
2. Remove sterilization bag containing the nested Organ Assembly and Shell Assembly.

NON-STERILE FIELD


NON-STERILE FIELD

3. Remove the nested Organ Assembly and Shell Assembly from sterilization bag and, without compromising sterility, introduce into the sterile field.


NOTE: ALL COMPONENTS OF THE ORGAN ASSEMBLY AND SHELL ASSEMBLY SHOULD REMAIN STERILE THROUGHOUT THE RECOVERY PROCEDURE


4. Unclasp all four Shell Assembly latches by opening two sets of opposing latches. First **Group 1** of opposing latches, then **Group 2** of opposing latches (Closed configuration shown below)


5. Remove the Shell Top (part #4) from the Shell Base (part #3).


6. Remove the Organ Canister (part #1) and Organ Lid (part #2) from the Shell Base (part #3).

7. Unclasp all four latches on the Organ Canister and Organ Lid by opening two sets of opposing latches. First **Group 1** of opposing latches, then **Group 2** of opposing latches

8. Remove the Organ Lid (part #2).


9. Visually inspect all material for damage and discard if necessary.


STERILE FIELD

4.8 HEART RECOVERY, AORTIC CONNECTION, AND PRESERVATION


1. Fill the Organ Canister with approximately 3 L of cold (4°C) heart preservation solution per institution recommendation without over-filling, allowing for displacement of heart (solution level will be topped up prior to final disposition).

2. Remove Heart Connectors from the Paragonix SherpaPak® CTS box (medium box).

3. Identify most appropriate Heart Connector size from among 10mm, 18mm, 25.4mm, and 32mm Heart Connectors provided. Remove the selected sterile Heart Connector from its sterilization bag and introduce into sterile field.

4. The Heart Connector has a bayonet joint that fits into socket at the center of the Organ Lid (Part #2). It can be twisted in either direction to lock into the Lid.

5. Insert the Heart Connector into the aorta and tie the aorta to the connector.

NOTE: THE 18mm, 25.4mm, AND 32mm HEART CONNECTORS ARE 38.1mm IN LENGTH. THE 10mm HEART CONNECTOR IS 35mm IN LENGTH. MEDICAL PERSONNEL SHOULD ENSURE THAT THE AMOUNT OF RECOVERED ASCENDING AORTA WILL BE SUFFICIENT TO HAVE ADEQUATE AORTIC LENGTH FOR THE ANASTOMOSIS AND AVOID INJURY TO THE AORTIC VALVE LEAFLETS AND CORONARY ORIFICES.

6. Securely tie the vessel into position with umbilical tape above the slippage restraint or use a tie gun to affix a tie around the aorta. Ensure there are no leaks.


7. Attach the Heart Connector to the central socket on the underside of the Organ Lid (Part #2) by twisting the Heart Connector in either a clockwise or counterclockwise direction while supporting the heart. A tactile click indicates the Heart Connector is engaged to the Organ Lid.

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING


8. While supporting the heart, insert it into the Organ Canister (Part #1) and seat the Organ Lid (Part #2) onto the lip of the Organ Canister (Part #1).

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING


9. Lock the Organ Lid (Part #2) onto the Organ Canister (Part #1) by closing two groups of opposing latches. First Group 1 of opposing latches, then Group 2 of opposing latches


10. Open the **VENT PORT** and the **FILL PORT** on the Organ Lid (Part #2) (so labeled).

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING

11. Using an infusion set (not supplied) connected to a bag of cold (4°C) preservation solution, fill the compartment in the Organ Lid (Part #2) with preservation solution through the **FILL PORT** (so labeled) located on the top surface of the Lid. Filling is performed directly from a line attached to the bag of preservation solution. Continue filling until fluid appears at the **VENT PORT**.


12. Securely close the **VENT PORT** and the **FILL PORT** on the Organ Lid (part #2) using the attached Luer caps.


13. Insert the assembled and filled Organ Canister into the Shell Base (part #3). The Organ Canister seats into position when the latches align with the associated notches on the Shell Base.

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING

14. Place the Shell Top (Part #4) onto the upper lip of the Shell Base (Part #3).

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING.

10. Open the **VENT PORT** and the **FILL PORT** on the Organ Lid (Part #2) (so labeled).

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING

11. Using an infusion set (not supplied) connected to a bag of cold (4°C) preservation solution, fill the compartment in the Organ Lid (Part #2) with preservation solution through the **FILL PORT** (so labeled) located on the top surface of the Lid. Filling is performed directly from a line attached to the bag of preservation solution. Continue filling until fluid appears at the **VENT PORT**.

12. Securely close the **VENT PORT** and the **FILL PORT** on the Organ Lid (part #2) using the attached Luer caps.


13. Insert the assembled and filled Organ Canister into the Shell Base (part #3). The Organ Canister seats into position when the latches align with the associated notches on the Shell Base.

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING

14. Place the Shell Top (Part #4) onto the upper lip of the Shell Base (Part #3).

NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING.

15. Close the four latches of the Shell Assembly by closing two sets of opposing latches.
First **Group 1** of opposing latches, then **Group 2** of opposing latches

20. Seat and lock each of the 4 Shipper latches securely. The Paragonix SherpaPak® CTS is now ready for transport.

16. Transfer the fully loaded Shell Assembly with nested Canister Assembly to the back table and place into the pre-assembled Shipper being sure to remove SherpaCool® Pouch from Shipper. Orient so that the tabs on the Shell Assembly align with the tabs on the Shipper ring.

NON-STERILE FIELD

4.9 TRAVELING WITH THE PARAGONIX SHERPAPAK® CTS TO THE TRANSPLANT SITE


Confirm that the Shipper is securely closed with all 4 latches in place. In the event of transport by vehicle, roll the Paragonix SherpaPak® CTS to the vehicle and place the Paragonix SherpaPak® CTS in a flat location secured against shifting or tipping during transit. Secure the Paragonix SherpaPak® CTS as necessary to accomplish this. If transported by air, either on a helicopter or airplane, follow the crew's instructions and secure the Paragonix SherpaPak® CTS to prevent shifting or tipping during transit.

Upon arrival at the transplant site, follow hospital procedures for moving equipment into the transplant OR. Identify a non-sterile table in the OR for the Paragonix SherpaPak® CTS.


Upon arrival to the recipient site, continue use of the Paragonix SherpaPak® CTS until the recipient has been prepared to accept the donor heart and that the Shipper is securely closed with all 4 latches in place.

4.10 REMOVING THE HEART FROM THE PARAGONIX SHERPAPAK® CTS FOR TRANSPLANT


1. Remove the Shipper lid and Paragonix SherpaCool® Pouch from atop the Shell Assembly.

2. Pinch lock arm tabs releasing Shell Assembly


3. Using the handle lift the Shell Assembly out of the Shipper and place on table.


NON-STERILE FIELD


NOTE: ENSURE ELECTRICAL CONNECTORS ARE NOT WETTED DURING HANDLING. ENSURE THAT THE SHELL ASSEMBLY IS INTRODUCED VERTICALLY. DO NOT USE EXCESSIVE FORCE WHEN INSTALLING THE SHELL ASSEMBLY INTO THE SHIPPER UNIT.

17. Insert the SherpaCool® Pouch beneath the handle of the Shell Assembly.


18. Close the lid of the Shipper.


19. With the tip of the finger, press and hold Button 1 on the datalogger for 10 seconds to power on the unit and initiate temperature and time logging.


NOTE: FOR THE FIRST HOUR FOLLOWING ASSEMBLY OF PARAGONIX SHERPAPAK® CTS YOU MAY EXPERIENCE TRANSIENT TEMPERATURE BELOW 4°C OR ABOVE 8°C. THIS IS DUE TO SYSTEM STABILIZATION IMMEDIATELY FOLLOWING FINAL ASSEMBLY.

4. Remove the Shell Top (Part #4) from the Shell Base by unclasping two sets of opposing latches as done previously.


5. Touching ONLY the Shell Base, position the Canister Assembly such that a gownned person in the sterile field can remove the Canister Assembly. The Canister Assembly may now be slowly drawn from the Shell Base and introduced into the sterile field by a gownned individual using proper aseptic technique.


STERILE FIELD


6. Place the Canister Assembly in a basin on a flat surface in the sterile field


7. Open the FILL PORT located on the Organ Lid.


8. Unlatch the four latches on the Canister Assembly by opening two opposing sets of latches.

9. Lift the Organ Lid (part #2), detach the Heart Connector from the Organ Lid (Part #2) by rotating either clockwise or counter-clockwise.

10. Remove the heart from the Heart Connector by cutting the aorta just below the tie.

Allow for heart preparation for transplantation as per transplant facility protocol. Following use, dispose of heart preservation solution, and discard the entire Paragonix SherpaPak® CTS per facility protocol.

5 TROUBLESHOOTING

Table 4. Troubleshooting

TROUBLE: Visible bubbles or air pockets in Organ Canister
PROBABLE CAUSE: Insufficient purging
ACTION:
1. Under aseptic conditions, open fill port, connect tubing connected to preservation solution; at the same time, open purge port until fluid is leaving the lid via the vent port.
2. If needed, tip Canister Assembly slightly to the side to dislodge trapped bubbles. Repeat fill and vent process.

TROUBLE: Organ Canister Leaking preservation solution
PROBABLE CAUSE: Loose fitting between Organ Canister and Organ Lid
ACTION: Reposition Organ Lid on Organ Canister and ensure latches are closed.

TROUBLE: No temperature displayed on datalogger LCD screen
PROBABLE CAUSE: Datalogger not activated and collecting data
ACTION:
1. First attempt to start datalogger by pressing Button 1 and holding for 10 seconds (repeat step 4.8.18).
2. If the datalogger does not activate and display temperature, refer to Frequently Asked Questions (Section 5).
3. If unable to remedy and temperature still not displayed, contact Paragonix immediately via phone for further assistance.

6. STORAGE
Unopened Paragonix SherpaPak® CTS and Shippers should be stored indoors in a dry location out of direct sunlight under normal temperature and humidity conditions. Unopened SherpaCool® Ribbons and Pouches should be stored at -20°C in preparation for use.

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L-621 (DOC-3612) Ver. 4

Approved By:

(CO-26196) Additional Translations to CTS and BARO IFUs

Description

Updating L-621 CTS Multilingual IFU to include Polish translations and L-659 BARO Multilingual IFU to include Polish and Czech translations. Version numbers and issue dates were updated.

Justification

Additional countries needed per Getinge request.

Assigned To:	Initiated By:	Priority:	Impact:
Giovanna Fonseca	Giovanna Fonseca	Medium	Minor

Version History:

Author	Effective Date	CO#	Ver.	Status
Giovanna Fonseca	April 7, 2026 10:59 AM EDT	CO-26196	4	Published
Julia Bryant	July 22, 2025 11:07 AM EDT	CO-25193	3	Superseded
Grant Real	September 17, 2024 1:10 PM EDT	CO-23628	2	Superseded
Kathryn Jenkins	June 14, 2024 9:56 AM EDT	CO-23491	1	Superseded
Grant Real	May 1, 2024 2:43 PM EDT	CO-23429	0	Superseded