INSTRUCTIONS FOR USE
1. INDICATIONS FOR USE

BAROguard® System is intended to be used for the static hypothermic preservation of lungs during transportation and eventual transplantation into a recipient using cold storage solutions indicated for use with the lungs.

The intended organ storage time for BAROguard® System is up to 8 hours.

Donor lungs exceeding clinically accepted static hypothermic preservation times should be evaluated by the transplant surgeon to determine transplantability in accordance with accepted clinical guidelines and in the best medical interest of the intended recipient.

Note: Partial lungs can be transported via BAROguard® System by packaging lungs per institutional protocol and UNOS guidelines.

1.1 WARRANTY STATEMENT
Paragonix Technologies, Inc.® (Paragonix) warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling and storage of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Paragonix's control directly affect the device and the results obtained from its use. Paragonix's obligation under this warranty is limited to the repair or replacement of this device and Paragonix shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. Paragonix does not authorize any other person to assume for it any other or additional liability or responsibility in connection with this device. Paragonix assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such devices.

2. SAFETY REQUIREMENTS

⚠ Caution: Donor lungs exceeding 8 hours static hypothermic preservation time will require transplant surgeon evaluation to determine transplantability.

Should any serious incident occur in relation to the device, please contact Paragonix and the appropriate regulatory authority in the county of use where the incident occurred. Paragonix can be reached via your local sales rep or contacted via:

email: support@paragonixtechnologies.com
telephone: +1.781.428.4828

2.1 TRAINING

Schedule training prior to first use of the device. To schedule, contact Paragonix by email (support@paragonixtechnologies.com) or by phone at +1.781.428.4828. This training is for intended users to review the steps recommended by the Instructions for Use in concert with appropriate Paragonix Technologies staff.

2.2 IMPORTANT INFORMATION



It is important that all personnel who will operate the BAROguard® System read and understand these instructions for use before operating the device.

All personnel should follow all warnings and precautions outlined below, for their safety and the safety of those around them. BAROguard® System is a device for use by hospitals that perform lung recovery procedures to store and transport donor lungs to the transplant facility. BAROguard® System is for use within acute care facilities that have an existing agreement with an Organ Procurement Organization (OPO) requiring the hospital to notify the OPO or third party designated by the OPO in a timely manner about all deaths and imminent deaths that occur in the hospital. Such hospitals must also have documented protocols and procedures for determining death or imminent death, obtaining family consent, and trained personnel with fully equipped surgical recovery operating rooms to perform the recovery procedure.

2.3 SYMBOL DEFINITIONS

Table 1: Symbols used in the labeling BAROguard® System and their definitions.	
Symbols	Definitions
	Use by YYYY-MM-DD
	Do not reuse
	Do not use if package is damaged or open
	Indicates the need to consult the instructions for use
	Batch code
	Caution, consult accompanying documents
	Manufacturer
	Medical Device

	Device or device component is non-sterile
	Sterilized using irradiation
	Fragile, handle with care
	Keep Dry
	Temperature limits
	Pressure limits
	Humidity limits
	Unique Device Identifier
	Serial Number
	Model Number
	Country of Manufacture
	Date of Manufacture
	Authorized representative in the EU
	Indicates the Swiss authorized representative.
	Indicates European Importer

2.4 WARNINGS 

- **Caution: Federal (US) law restricts this device to sale by or on the order of a physician.**
- **Use aseptic technique as appropriate. Certain components of BAROguard® System are provided sterile. Maintain aseptic technique when handling sterile components of BAROguard® System.**
- **Chilled (4°C) preservation solution, cleared for use with the lungs, must be available.**

⚠ Do NOT reuse any component of the BAROguard® System. BAROguard® System is intended for single use only. DO NOT RE-USE. Re-use of the device may result in infection and other complications due to the loss of sterility. Certain components of the BAROguard® System are sterile as supplied (sterilization method is gamma irradiation). BAROguard® System should be disposed of in accordance with local guidelines for biomedical waste.

- **Use institution-specific precautions with the donor lungs and preservation solution when operating the BAROguard® System.** The lungs and preservation solution may carry undetected pathogens from the donor. Use universal precautions for bloodborne pathogens in handling the lungs, and in handling and disposing of the BAROguard® System and preservation solution to prevent the possible transmission of pathogens to personnel. As appropriate, this may include use of personal protective equipment

(e.g., gloves, masks, gowns, goggles, or equivalent eye protection) and disposal of materials as potentially infectious biohazard waste.

- **Prior to use, inspect all components of the BAROguard® System. Do not use if any component is loose, broken or damaged.**
 - **Do not open the BAROguard® System during organ transport.**
 - **No modification of the BAROguard® System is allowed.**
- WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.**
- WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the BAROguard® System device. Otherwise, degradation of the performance of this equipment could result.**
- NOTE: BAROguard® System is compliant with IEC 61000-4-2: 8KV contact, 15KV air, IEC 61000-4-3: 3V 80-2700MHz, Table 9 per the standard, IEC 61000-4-8 30A/M and additionally the following:**

Table 2. EMC/EMI Safety	
Test Standard	Test Level
TCISPR 11	Group 1 Class B, 30-1000 MHz
IEC61000-4-2	± 8kV contact, ± 2kV, ±4kV, ±8kV, ±15kV air
IEC61000-4-3	Radiated RF EM fields 3V/m, 80MHz-2.7GHz, 80% AM at 1 kHz
IEC61000-4-8	30 A/m, 50Hz or 60 Hz
IEC6100-4-39	9 kHz to 13.56 MHz

2.5 TRANSPORT AND STORAGE CONDITIONS 

- **Transport and Storage Temperature: -20°C to 40°C**
- **Transport and Storage Pressure: 625hPa – 1060hPa**
- **Transport and Storage Humidity: 5-80% Relative Humidity**

2.6 OPERATING CONDITIONS

⚠ Operating Temperature: 22°C Recommended Operating Room temperature is 22°C. Direct sunlight and outdoor temperature extremes (high and low) can affect the BAROguard® System internal temperature (4 to 8°C). During exposure to temperature extremes, the BAROguard® System temperature must be frequently monitored.

⚠ Operating Pressure: Sea Level to 8000 ft. The operating pressure range conforms to commercial airliner transport. Extreme pressure levels

can impact performance. During exposure to extreme pressures (altitudes), the BAROguard® System operation must be frequently monitored.

Operating Humidity: 40-60% Relative Humidity
Extreme humidity levels can impact performance. During exposure to extreme humidity levels, the BAROguard® System operation must be frequently monitored.

2.7 DATALOGGER COMMUNICATION AND ACCURACY

Temperature Accuracy Wireless Data Standard: Bluetooth Low Energy (Bluetooth Smart)
BAROguard® System datalogger monitors temperature and pressure and can transmit to a mobile device using optional mobile app (section 4.4) using Bluetooth® Low Energy Technology.
Datalogger Radio Power: 1 mW (0 dBm)
Datalogger Transmission Range: Approximately 30.5 m (100 ft) line-of-site

Temperature Accuracy
BAROguard® System comes with a pre-installed datalogger capable of reporting the temperature within the Organ Assembly. Temperature Accuracy is ±0.5°C from 0° to 50°C.

Pressure Accuracy
BAROguard® System comes with a pre-installed datalogger capable of reporting airway pressure of the donor lung. Pressure accuracy is ±1.0 cmH₂O.

Time Accuracy
The pre-installed datalogger logs temperature against time with an accuracy of ±1 minute per month.

2.8 PRECAUTIONS

- **Keep the BAROguard® System primarily upright during transportation.**
The BAROguard® System is designed to be transported upright. Temporary tilting ±45° from horizontal in any direction is acceptable.
- **Do not place BAROguard® System in the cargo hold of an airplane.**
BAROguard® System is designed to be transported in pressurized cabin during aircraft transportation.
- **Avoid direct sunlight and hot or cold temperature extremes.**

The BAROguard® System is designed to be transported under the same environmental conditions as is appropriate for people. Avoid extended exposure to outdoor conditions (sunlight, heat or cold).

- **Use caution when lifting the BAROguard® System.**
A fully loaded BAROguard® System weighs ~30 lbs. Use proper lifting practices.
- **Exercise Caution when using in presence of other electronic devices and electromagnetic emitters.**

BAROguard® System has been tested for radiated immunity only at selected frequencies listed in Section 2.4 of this manual. While no significant risk of reciprocal interference has been identified, BAROguard® System should always be closely monitored when in use and reoriented or relocated if

abnormal performance is observed. Electromagnetic emitters include radiofrequency identification (RFID) readers, electronic security systems (e.g., metal detectors, electronic article surveillance), near-field communications (NFC) systems, wireless power transfer (WPT), Cellular 5G, and unique medical emitters such as electrocautery, MRI, electrosurgical units, and diathermy equipment. **Do not use BAROguard® System in oxygen rich environment.** BAROguard® System is not designed or evaluated for use in an oxygen rich environment. **Do not use BAROguard® System with flammable anesthetics.** BAROguard® System is not designed or evaluated for use with flammable anesthetics.

2.9 SIDE EFFECTS
All surgical procedures and medical devices have potential risks. The potential surgical risks of a transplant with donor lungs are similar for the Paragonix BAROguard System and traditional ice storage preservation. There is a risk of receiving lungs that do not function properly after transplant. There is also a risk that the donor lungs may be damaged during preservation.

2.10 SIDE EFFECTS ASSOCIATED WITH BARO GUARD

- It is possible that after preservation using the Paragonix BAROguard® system, the transplant doctor may decide that the donor lungs are not suitable for transplantation.
- The Paragonix BAROguard® system is prepared and operated by trained medical professionals. Like with many medical technologies, there are inherent risks including injury to the donor tissue, air leaks, infection or delayed organ function.

2.11 CLINICAL BENEFITS (EU)
The Paragonix BAROguard® System is designed to improve organ preservation during transportation and storage from organ donor to recipient. Improved organ preservation can mitigate the risks of cold ischemic time on reperfusion injury, leading to reduced post-operative complications and improved long-term outcomes.

3 PRE-ASSEMBLY CHECKLIST

Table 3: Pre-Assembly Checklist for BAROguard® System identifying component counts, locations, and sterility disposition.

Included with BAROguard® system	Part Number	Item
Yes (Non-sterile)	PRGNX-7007-001	BAROguard® System Shipper (1)

Yes (Non-sterile)	PRGNX-7007-002	BAROguard® System Paragonix SherpaCool® Materials (1) CAUTION: PRECONDITION AT -20°C OR BELOW FOR AT LEAST 48 HOURS PRIOR TO USE
Yes (Sterile)	PRGNX-7007-003	BAROguard® System Sterile Components (1) 4 total Liters of fluid exclusively for packaging the donor lung(s) (e.g. 2L of chilled (4°C), sterile preservation solution for the inner bag and 2L of chilled (4°C), sterile or preservation solution for the middle bag)
No (Sterile)	Varies	Hemostat (atraumatic DeBakey clamp or similar) or other equivalent airway securement tool (1)


Figure 2: Shipper

Figure 3: SherpaCool® Materials.

CAUTION: PRECONDITION AT -20°C OR BELOW FOR AT LEAST 48 HOURS PRIOR TO USE.


Figure 4: Preservation Solution and/or Saline

CAUTION: FOLLOW DIRECTIONS IN THE ENCLOSED IFU REGARDING STORAGE AND PREPARATION FOR USE

Figure 5: Sterile Components

CAUTION: FOLLOW DIRECTIONS IN THE ENCLOSED IFU REGARDING STORAGE AND PREPARATION FOR USE

3.1 BAROguard® System Set-up and Preparation

1. Inspect all parts on arrival for any signs of damage that may have occurred during transport.
 2. Report any damage or concerns about the condition of BAROguard® System immediately to Paragonix Technologies Inc.®
 3. When in transit, the SherpaCool® Pouches may be stored in the shipper for up to 18 hours and with a total maximum of 7 liters of saline and/or preservation solution based on required practice.
1. **Note: If SherpaCool® is set up in the shipper, place the pouches in the intended locations as shown in section 4.8.**
 2. Paragonix SherpaCool® pouches should not be placed in direct contact with the temperature probe.


Figure 1: BAROguard® System, and as shipped (in boxes)

4 OPERATING INSTRUCTIONS

4.1 GENERAL INFORMATION
Before using in a clinical setting, operators must be trained in the use and functional understanding of the BAROguard® System.

4.2 OVERVIEW USING BARO GUARD® SYSTEM
Using BAROguard® System involves performing the following procedures:

1. Removal of packaging, preconditioning of SherpaCool® Box (Section 4.3) containing SherpaCool® Pouches at or below -20 C for at least 48 hours.

2. Setup of datalogging software for use with the BAROguard® System (Section 4.4, optional, recommended)
3. Transport of the BAROguard® System to the recovery site (Section 4.5)
4. Preparing the BAROguard® System and Shipper for deployment at recovery site (Section 4.6)
5. Preparing BAROguard® System for lungs receipt at recovery site (Section 4.7)
6. Donor lung(s) recovery, packaging, and preservation (Section 4.8)
7. Traveling with BAROguard® System to the transplant site (Section 4.9)
8. Removing the lung from the BAROguard® System for transplant (Section 4.10)

These instructions may be modified based on actual use. The instructions are designed for implementation by two operators. These instructions may be modified for a single operator, provided that proper aseptic technique is used in conjunction with procedures to be performed on the sterile field with sterile BAROguard® System components.

4.3 REMOVAL OF PACKAGING AND PRECONDITIONING OF SHERPACOOl® BOX CONTAINING SHERPACOOl® POUCHES

BAROguard® System must be maintained in a ready-to-use condition, so that it can be available to the lung recovery team at all times.

Make the following preparations:

- Do not remove SherpaCool® Pouches from the SherpaCool® Box.
- Write the date and time of SherpaCool® Box transition to -20°C storage on the label affixed to the SherpaCool® Box.
- Place the SherpaCool® Box into a -20°C (or colder) freezer for a minimum of 48 hours.
- Remove SherpaCool® Box after a minimum preconditioning time of 48 hours at -20°C (or colder).
- Do not remove Paragonix SherpaCool® Box until recovery team departs to donor hospital.
- Do not remove Paragonix SherpaCool® Box until all other components of BAROguard® System and associated equipment and materials have been prepared for transport.
- Paragonix SherpaCool® Box must be transported to the donor site on ice OR within the shipper as described in section 3.1 above.
- 4 total Liters of fluids exclusively for packaging the donor lung(s) (e.g. 2L of chilled (4°C), sterile preservation solution cleared for use with lungs for the inner bag and 2L of chilled (4°C), sterile saline or preservation solution for the middle bag)

4.4 SETUP OF DATALOGGING SOFTWARE FOR USE WITH BAROGUARD® SYSTEM (OPTIONAL, RECOMMENDED)

In order to recover logged temperature data from BAROguard® System, the user must first download and install the latest version of the mobile Paragonix App at either the App Store or Google Play store. This requires use of a Bluetooth®-enabled iOS or Android device.

NOTE: BAROguard® System comes pre-configured with appropriate temperature monitoring settings installed in the on-board temperature monitoring system. This configuration must NEVER be modified by the user. After setting up an account for use of the Paragonix App:

- After logger is started (Section 4.8, Step 20 of this manual) it is possible to visualize the logged data by connecting the BAROguard® System to the Paragonix App.
 - Log in to the Paragonix App on your mobile device.
 - Scan the QR code on the back of the Shipper to pair the device or select the Datalogger matching the serial number shown at the back of the Shipper from the Paragonix app.
 - If prompted enter the Datalogger PassKey provided on the back of the Shipper adjacent to the datalogger serial number.


Figure 6: Indication of datalogger Serial Number and Pass-Key on Shipper

- **Once connected to the the datalogger, logged data may be viewed locally and/or uploaded to alternative devices via the Paragonix software application.**
- **For assistance in data recovery, contact Paragonix® directly via provided telephone number and have available the Shipper Serial Number and DataLogger Serial Number.**
- **Do not dispose of Shipper unit until desired data has been obtained.**

4.5 TRANSPORT OF BAROGUARD® SYSTEM TO THE RECOVERY SITE

Transport various supplies of BAROguard® System to the recovery site by:

- Place SherpaCool® Box onto wet ice in an ice chest for transportation OR in the shipper as described in section 3.1.
- Paragonix SherpaCool® Box may be transported in the same ice chest as other adjunct supplies (preservation solutions, etc.)
- Maximum allowable transit time for SherpaCool® is 4 hours for wet ice transport, 18 hours within the shipper.
- Carry BAROguard® System Nested Bag Assembly and (sterile components) in original packaging
- Unbox Shipper prior to transportation to the recovery site.
- Wheel or carry Shipper to recovery site.

4.6 PREPARATION OF BAROGUARD® SYSTEM FOR DEPLOYMENT AT RECOVERY SITE

1. Using Pre-Assembly Checklist (See Section 3 Table 3), double-check for presence of all BAROguard® System components PRIOR to unpacking the components.
2. Open Paragonix SherpaCool® Box, leaving SherpaCool® Pouches on ice until installation in the Shipper.


NON-STERILE FIELD

4.7 PREPARING BAROGAURD™ FOR LUNGS RECEIPT AT RECOVERY SITE

1. Aseptically remove Sterile Components (nested bags, trachea connectors, cable ties, cable installation tool) and place on sterile table

Cable Installation Tool

Trachea Connectors (15mm, 18mm, 20mm, 23mm)


Nested Bag Assembly, Cables Ties, Bag Strings

2. Visually inspect all material for damage and discard if necessary.

STERILE FIELD

STERILE FIELD

3. Remove Shipper Lid and Upper Tray from Shipper

NON-STERILE FIELD

4.8 LUNG RECOVERY, PACKAGING, AND PRESERVATION

1. Prepare nested bag within back table basin.

STERILE FIELD

2. Add 2 liters of chilled (4°C), sterile preservation solution to inner bag.

STERILE FIELD

3. Prepare lungs as appropriate for insertion of trachea connector

⚠ Note: Maintain airway pressure during recovery from donor cavity and back table prep airway securement tool (e.g., hemostat, atraumatic DeBakey or similar clamp)

⚠ Note: Only package lungs that have been examined and accepted during final assessment without air leaks per institutional protocol.



STERILE FIELD

4. Insert appropriate size trachea connector into the trachea and secure with 2 cable ties

⚠ Note: Using less than 2 cable ties may lead to arterial leaks.


⚠ Note: Keep hemostat applied. Do not remove airway securement tool (e.g. hemostat, atraumatic DeBakey or similar clamp) until step 7.
NOTE: DO NOT REMOVE HEMOSTAT

STERILE FIELD

5. Remove cap from bag tubing connector while depressing the white button

⚠ Note: This may require two hands.

STERILE FIELD


6. Insert trachea connector into bag tubing connector.

⚠ Note: Keep hemostat applied. Do not remove airway securement tool (e.g. hemostat, atraumatic DeBakey or similar clamp) until step 7.

STERILE FIELD

7. Remove airway securement tool (e.g. hemostat, atraumatic DeBakey or similar clamp).

STERILE FIELD

8. Remove air from inner bag and close using institutional protocol with either woven ties or cable ties.

⚠ Note: Consider institutional protocol at recipient center for unpackaging donor lung(s) when selecting closure method.
OPTION 1: WOVEN TIES

STERILE FIELD

OPTION 2: CABLE TIES

STERILE FIELD

9. Add 2 liters of chilled (4°C), sterile saline or chilled (4°C) sterile preservation solution to the middle bag.

⚠ Note: Do not use sterile slush or ice.

STERILE FIELD

10. Remove air from middle bag and close using institutional protocol with either woven ties or cable ties.

⚠ Note: Consider institutional protocol at recipient center for unpackaging donor lung(s) when selecting closure method.

STERILE FIELD

11. Remove air from outer bag and close using institutional protocol with either woven ties or cable ties.

⚠ Note: Consider institutional protocol at recipient center for unpackaging donor lung(s) when selecting closure method.

STERILE FIELD

12. Turn on pump by depressing the indicated button.

⚠ Note: Turn on pump prior to connecting bags to the shipper to ensure lungs remain inflated. Illuminated blue LED switch and audible tone will confirm pump operation.

NON-STERILE FIELD

13. Lift only the front edge of the Bottom Tray towards the datalogger.

⚠ Note: Use caution, temperature probe is affixed to Bottom Tray. Do not pull.

NON-STERILE FIELD

14. Place Paragonix SherpaCool® Pouch 1 underneath the Bottom Paragonix SherpaCool® Tray and lower the tray.

NON-STERILE FIELD

15. Place the packaged donor lungs within the device and connect the tubing to the airline port within the shipper.

⚠ NOTE: Ensure the pump is ON prior to connecting the tubing within the device to prevent potential partial lung deflation.

ENSURE PUMP IS TURNED ON

NON-STERILE FIELD

16. Load the 2 Paragonix SherpaCool® Pouch "2" materials into the Upper Paragonix SherpaCool® Tray.

NON-STERILE FIELD

17. Replace the Upper Tray in the Shipper.

NON-STERILE FIELD

18. Close the lid of the device.

NON-STERILE FIELD

19. Seat and lock each of the 4 Shipper latches securely.

NON-STERILE FIELD

20. Press and hold Button 1 on the datalogger for 10 seconds to power on the datalogger and initiate the logging of temperature, airway pressure, and time.

NON-STERILE FIELD

Press Button 1 to toggle between airway pressure, lung cavity probe temperature, battery voltage, and internal BAROguard® System system temperature.

⚠ NOTE: For the first hour following assembly of BAROguard® System you may experience transient temperature below 4°C or above 8°C. This is due to system stabilization immediately following final assembly.

⚠ NOTE: Time to target pressure stabilization may vary depending on initial pressure maintained during procurement from donor cavity.

4.9 TRAVELING WITH BARO GUARD® SYSTEM TO THE TRANSPLANT SITE

Confirm that the device is securely closed with all 4 latches in place.

In the event of transport by vehicle, roll BAROguard® System to the vehicle and place BAROguard® System in a flat location secured against shifting or tipping during transit. Secure BAROguard® System as necessary to accomplish this. If transported by air, either on a helicopter or airplane, follow the crew's instructions and secure BAROguard® System to prevent shifting or tipping during transit.


⚠ NOTE: Use caution when loading and unloading BAROguard® System to ensure device power remains ON.

Upon arrival at the recipient transplant site, follow institutional protocol for moving equipment into the transplant OR. Identify a non-sterile surface in the OR for staging BAROguard® System. Continue use of BAROguard® System until the recipient has been prepared to accept the donor lungs and that the device is securely closed with all 4 latches in place.

4.10 REMOVING THE LUNGS FROM BARO GUARD® SYSTEM FOR TRANSPLANT

1. Remove the device lid and the Upper Tray containing Paragonix SherpaCool® materials.

⚠ Note: DO NOT TURN OFF PUMP


NON-STERILE FIELD

2. Disconnect nested bag assembly using the quick connect.

⚠ Note: DO NOT TURN OFF PUMP

⚠ Note: Disconnect tubing **PRIOR** to lifting the package donor lung(s).

NON-STERILE FIELD

3. **AFTER** disconnecting the quick connect, remove the packaged donor lung(s) from the device.


NON-STERILE FIELD

4. Touching **ONLY** the outer bag, a nonsterile operator will open it such that a sterile operator can remove the middle bag per institutional protocol.

NON-STERILE FIELD

5. Using appropriate aseptic technique, the sterile operator will disconnect the outer bag from the middle bag. Carefully remove the middle bag assembly and bring to a basin on the sterile table.

STERILE FIELD

6. Open the middle bag per institutional protocol

⚠ NOTE: Middle bag and inner bag remain connected

STERILE FIELD

7. Open the inner bag per institutional protocol.

STERILE FIELD

8. Remove first lung from inner bag for preparation of transplant per institutional protocol.

STERILE FIELD

9. As the first lung is being transplanted, repackage second lung back into shipper per institutional protocol.

STERILE FIELD

10. Unpack the second donor lung per institutional protocol once it is ready for transplant.

STERILE FIELD

11. Following use, dispose of lungs preservation solution, and discard the entire BAROguard® System per institutional protocol via the biohazard waste stream.

5. TROUBLESHOOTING

Table 4. Troubleshooting

TROUBLE:
No temperature or pressure displayed on datalogger LCD screen.

PROBABLE CAUSE:
Datalogger not activated and collecting data.

ACTION:

1. First attempt to start datalogger by pressing Button 1 and holding for 10 seconds (repeat step 4.8.20).
2. If the datalogger does not activate and display temperature, contact Paragonix immediately via the phone for further assistance.

TROUBLE:
Datalogger alarm is going off.

PROBABLE CAUSE:
Pressure in airway is less than 10 cm H₂O.

ACTION:

1. Ensure the quick connector of the nested bags is connected to the quick connect in the Shipper wall.
2. If the alarm continues for another minute disconnect the quick connect of the nested bags from the quick connect in the wall of the shipper

6. STORAGE

Unopened BAROguard® System Shippers should be stored indoors in a dry location out of direct sunlight under normal temperature and humidity conditions. Unopened SherpaCool® Pouches should be stored at -20°C in preparation for use.

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L-659 (DOC-4313) Ver. 2

Approved By:

(CO-26196) Additional Translations to CTS and BARO IFUs

Description

Updating L-621 CTS Multilingual IFU to include Polish translations and L-659 BARO Multilingual IFU to include Polish and Czech translations. Version numbers and issue dates were updated.

Justification

Additional countries needed per Getinge request.

Assigned To:	Initiated By:	Priority:	Impact:
Giovanna Fonseca	Giovanna Fonseca	Medium	Minor

Version History:

Author	Effective Date	CO#	Ver.	Status
Giovanna Fonseca	April 7, 2026 10:59 AM EDT	CO-26196	2	Published
Julia Bryant	July 22, 2025 11:06 AM EDT	CO-25311	1	Superseded
Julia Bryant	October 8, 2024 10:49 AM EDT	CO-23631	0	Superseded