SELECT LANGUAGE:

ENGLISH

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SPANISH • ESPAÑOL

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GERMAN • DEUTSCH

FINNISH • SUOMI

SWEDISH • SVENSKA

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ARM GUARD

SINGLE-ARM PROTECTION & SUPPORT

Arm Guard is a single-use, freestanding positioner designed to position the patient's arm in an extended or bent position at the patient's side or over the patient's chest for medical procedures.

To reorder the following disposable products, please contact
sales@bonefoam.com

- 925.SAS-WS-BX10 - Arm Guard - Box of 10

Note for RWS:

Items in the footer on this cover page have been resized and repositioned to accommodate the new EU info. This should be carried over to other languages.

New content that needs to be translated or carried over and inserted into non-English chapters is highlighted in magenta within the ENG chapter.

BoneFoam, Inc.
20175 County Road 50
Corcoran, MN 55340
USA

Email: customerservice@bonefoam.com
Phone: (763) 559-1830
Toll-free: 1 (877) 861-2663
Fax: (763) 559-1822


BoneFoam B.V.
P.O. Box 10110, NL-5060 GA Oisterwijk
Kantstraat 19, 5076 NP Haaren
EAR@qnetbv.eu

bonefoam.com

SAS-eIFU Rev.0 2025-07-01

INSTRUCTIONS FOR USE

1. Place arm in Arm Guard.
2. Close upper and lower wraps to secure arm with antecubital space exposed.
3. Pull strap to bend arm and Arm Guard to desired flexion.
4. Secure strap to table rail.

WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- Risk of cross-contamination if not disposed of immediately after use.
- Incorrect positioning of the patient with this product can increase the risk of:
 - Pressure injury resulting in markings or ulcers
 - Peripheral nerve injury, which can lead to temporary or permanent paralysis
 - Venous pressure, which can lead to ischemia
 - Instability, which can lead to patient falls
- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- **Ensure device is compatible with work surface, procedural instruments, and equipment.**
- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
 - Foam: Rated either UL 94 HBF or CAL 117

CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.
- Prior to using this product, read the Instructions for Use.
- Discontinue use if foam has degraded or has lost resilience.

SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

PACKAGING

The product is packaged in a clear polythene bag inside a corrugated cardboard box. Components may be compressed to minimize packaging and storage footprint.

SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
925.SAS-WS-BX10	Arm Guard Box of 10	Nonsterile, single-use positioner intended to secure the arm over the body in the supine or beach chair positions.	Indicated for use by healthcare professionals to position the patient's nonoperative arm during surgery.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
 bonetfoam.com/eIFU	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the device is intended for one single-use only
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

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