
▲BONEFOAM™


TROOP

AIRWAY MANAGEMENT
ELEVATION SYSTEM

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▲ BONEFOAM™

TROOP

AIRWAY MANAGEMENT ELEVATION SYSTEM

The Troop product line offers reusable and single-use tabletop positioning sets designed to position patients, particularly obese patients, for optimum airway management during medical procedures.

To reorder the following Troop Elevation System disposable products, please contact sales@bonefoam.com

- 10-910-00 - Troop Elevation System - Disposable - Box of 3
- 10-910-03 - Disposable Barrier Covers - Box of 30
- 10-910-02 - Head Cradle - Disposable - Box of 4

BoneFoam, Inc.
2601 49th Ave N Suite 100
Minneapolis, MN 55430
USA

Email: customerservice@bonefoam.com
Phone: (763) 559-1830
Toll-free: 1 (877) 861-2663
Fax: (763) 559-1822


BoneFoam B.V.
P.O. Box 10110, NL-5060 GA Oisterwijk
Kantstraat 19, 5076 NP Haaren
EAR@qnetbv.eu

bonefoam.com

TEP-eIFU Rev.1 2026-03-15

INSTRUCTIONS FOR USE

1. Place Troop Airway Management Elevation System at head of table and cover with barrier cover. Place head cradle on pillow.
2. Position patient on Troop Airway Management Elevation System, with shoulders positioned at shoulder bump, and head positioned on head cradle.
3. Replace armboard pads with Troop Armboard pads and position arms in channels

*The Troop Airway Management Elevation System Addition may be placed on top of the Troop Airway Management Elevation System when treating patients 450-500 lbs or patients whose BMI is >50.

CLEANING AND DISINFECTING INSTRUCTIONS

- Do not sterilize the product. This product is not intended to be sterilized.
- For coated products, wipe the product with a clean wipe and disinfectants after every use. The following disinfecting materials are known to be compatible with coated BoneFoam products:
 - Rubbing alcohol (70% isopropyl alcohol)
 - General purpose/all purpose cleaner
 - Read and follow the cleaning product's instructions.
- Make sure that the coated product is dry before you store it or use it again.

WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- For single-use products or accessories, risk of cross-contamination if not disposed of immediately after use.
- Incorrect positioning of the patient with this product can increase the risk of:
 - Pressure injury resulting in markings or ulcers
 - Peripheral nerve injury, which can lead to temporary or permanent paralysis
 - Venous pressure, which can lead to ischemia
 - Instability, which can lead to patient falls

- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- Ensure device is compatible with work surface, procedural instruments, and equipment.
- This product has a vinyl coating. While this coating is intended to enhance product durability and usability and is safe for skin contact, it may cause allergic reactions in individuals sensitive to vinyl or its derivatives.
- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
 - Coating: UL 94 HBF
 - Foam: Rated either UL 94 HBF or CAL 117

CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Allow coated products to reach room temperature after storage or transportation in cold temperatures.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid storing where sharp objects could pierce the coating of the product.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.
- Prior to using this product, read the Instructions for Use.
- Discontinue use if foam has degraded or lost resilience, or if coating is torn, cracked, or loose and wrinkles, as this may lead to acute pressure points.

SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

PACKAGING

The product is packaged in a clear polythene bag inside a corrugated cardboard box.

SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
 <small>bonefoam.com/eIFU</small>	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the device is intended for one single-use only
	Indicates the temperature limits to which the medical device can be stored, transported or used
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
REUSABLE SETS / POSITIONERS			
10-910-06-C	Complete System: Troop Elevation Pillow - Coated, Troop Elevation Pillow Addition - Coated, Head Cradle - Coated, Armboard Pads - Coated, and single-use Disposable Barrier Covers	Nonsterile, reusable positioning set, including nonsterile, single-use barrier covers, intended to support and stabilize the upper body in the head-elevated laryngoscopy position (HELP).	Indicated for use by healthcare professionals to position the head and torso of super morbidly obese patients (>400 lbs) for use during intubation and for the duration of the medical procedure.
10-910-01-C	Standard System: Troop Elevation Pillow - Coated, Head Cradle - Coated, Armboard Pads - Coated, and Disposable Barrier Covers		Indicated for use by healthcare professionals to position the head and torso of morbidly obese patients (<400 lbs) for use during intubation and for the duration of the medical procedure.
10-910-13-C-I	3-Set: Troop Elevation Pillow - Coated, Troop Elevation Pillow Addition - Coated, and Head Cradle - Coated	Nonsterile, reusable positioning set intended to support and stabilize the upper body in the head-elevated laryngoscopy position (HELP).	Indicated for use by healthcare professionals to position the head and torso of super morbidly obese patients (>400 lbs) for use during intubation and for the duration of the medical procedure.
10-910-12-C-I	2-Set: Troop Elevation Pillow - Coated, and Head Cradle - Coated		Indicated for use by healthcare professionals to position the head and torso of morbidly obese patients (<400 lbs) for use during intubation and for the duration of the medical procedure.
10-910-02-C	Head Cradle - Coated Box of 4	Nonsterile, reusable positioner intended to support and stabilize the head in the supine position.	Indicated for use by healthcare professionals to position the patient's head in the supine position during medical procedures. Can be used as a standalone positioner, or with the Troop Airway Management Elevation Complete System, Standard System, 3-Set, or 2-Set.
10-910-07-C	Armboard Pad - Coated Set of 2	Nonsterile, reusable, positioners intended to support and stabilize the arms in the in the head-elevated laryngoscopy position (HELP).	Indicated for use by healthcare professionals to position the patient's arms with the Troop Airway Management Elevation Complete System, Standard System, 3-Set, or 2-Set.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
SINGLE-USE SETS / ACCESSORIES			
10-910-00*	Troop Elevation System - Disposable Box of 3	Nonsterile, single-use positioning set intended to support and stabilize the upper body in the head-elevated laryngoscopy position (HELP).	Indicated for use by healthcare professionals to position the head and torso of super morbidly obese patients (>400 lbs) for use during intubation and for the duration of the medical procedure.
10-910-03*	Disposable Barrier Cover Box of 30	Nonsterile, single-use covers intended to provide a barrier between the patient and the positioner.	Indicated for use by healthcare professional to cover the Troop Elevation Pillow or Troop Elevation Pillow Addition.
10-910-02*	Head Cradle - Disposable Box of 4	Nonsterile, single-use positioner intended to support and stabilize the head in the supine position.	Indicated for use by healthcare professionals to position the patient's head in the supine position during medical procedures.
REPLACEMENT COMPONENTS			
10-911-06-C	Troop Elevation Pillow - Coated	Nonsterile, reusable positioner intended to ramp up and elevate the upper body in the head-elevated laryngoscopy position (HELP).	Indicated for use by healthcare professionals to position the torso of morbidly obese patients (<400 lbs) for use during intubation and for the duration of the medical procedure. Part of the Troop Airway Management Elevation Complete System, Standard System, 3-Set, or 2-Set.
10-911-07-C	Troop Elevation Pillow Addition - Coated		Indicated for use by healthcare professionals to position the torso of super morbidly obese patients (>400 lbs) for use during intubation and for the duration of the medical procedure. Part of the Troop Airway Management Elevation Complete System and 3-Set. Optional addition to be used with the Troop Airway Management Elevation Standard System or 2-Set.
10-911-08-C	Head Cradle - Coated	Nonsterile, reusable positioner intended to support and stabilize the head during medical procedures.	Indicated for use by healthcare professionals to position the patient's head in the supine position during medical procedures. Can be used as a standalone positioner, or with the Troop Airway Management Elevation Complete System, Standard System, 3-Set, or 2-Set.
10-911-09-C	Armboard Pad - Coated	Nonsterile, reusable, positioner intended to support and stabilize the arms in the in the head-elevated laryngoscopy position (HELP).	Indicated for use by healthcare professionals to position the patient's arms with the Troop Airway Management Elevation Complete System, Standard System, 3-Set, or 2-Set.

*  applies only to these products.  does not apply to these products.

▲ BONEFOAM™

TROOP

AIRWAY MANAGEMENT ELEVATION SYSTEM

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