


▲BONEFOAM™

SIDELINE

LATERAL LUMBAR
POSITIONING SYSTEM

104.SLPS-LB

SELECT LANGUAGE:

ENGLISH

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SPANISH • ESPAÑOL

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GERMAN • DEUTSCH

FINNISH • SUOMI

SWEDISH • SVENSKA

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SIDELINE

LATERAL LUMBAR POSITIONING SYSTEM

The Sideline Lateral Lumbar Positioning System is a tabletop positioning set with both reusable and single-use components designed to position a patient in the lateral decubitus position for lateral lumbar procedures.

The following disposable products are required for use with the Positioning System, but are not included:

- 104.SLPS-SUK-BX5 - Single-Use Kits - Box of 5

To reorder, please contact
sales@bonefoam.com

BoneFoam, Inc.
2601 49th Ave N Suite 100
Minneapolis, MN 55430
USA

Email: customerservice@bonefoam.com
Phone: (763) 559-1830
Toll-free: 1 (877) 861-2663
Fax: (763) 559-1822


BoneFoam B.V.
P.O. Box 10110, NL-5060 GA Oisterwijk
Kantstraat 19, 5076 NP Haaren
EAR@qnetbv.eu


bonefoam.com

SLPS-eIFU Rev.1 2026-03-15


INSTRUCTIONS FOR USE


- 1** Replace table pads with torso base and lower base and secure with table straps. Replace armboard pad with Shark base.


- 2** Attach one end of hip strap and one end of cross traction straps to lower base. Attach one end of torso strap to torso base. Secure other ends of straps to themselves so straps are hanging. Place facility draw sheet over table break between hip strap and torso strap.

- 3** Transfer patient to base pieces and into lateral position following facility transfer protocol. Adjust patient position with facility draw sheet so greater trochanter is in line with hip strap and place head guard under patient head.


- 4** Place Shark top on Shark base and place top arm in channel on Shark top.

- 5** Lift patient with facility draw sheet to position Ax Bump and Hip Bump. The Hip Bump should be positioned between greater trochanter and iliac crest. Place Leg Guard between patient's leg.

- 6** Adjust height of Head Guard to desired head position.

- 7** Secure torso strap over patient. Position hip strap over greater trochanter and secure tightly. Attach traction straps to hip strap and tighten to desired amount of traction.

CLEANING AND DISINFECTING INSTRUCTIONS

- Do not sterilize the product. This product is not intended to be sterilized.
- For coated products, wipe the product with a clean wipe and disinfectants after every use. The following disinfecting materials are known to be compatible with coated BoneFoam products:
 - Rubbing alcohol (70% isopropyl alcohol)
 - General purpose/all purpose cleaner
 - Read and follow the cleaning product's instructions.
- Make sure that the coated product is dry before you store it or use it again.

WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- For single-use accessories, risk of cross-contamination if not disposed of immediately after use.
- Incorrect positioning of the patient with this product can increase the risk of:
 - Pressure injury resulting in markings or ulcers
 - Peripheral nerve injury, which can lead to temporary or permanent paralysis
 - Venous pressure, which can lead to ischemia
 - Instability, which can lead to patient falls
- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- Ensure device is compatible with work surface, procedural instruments, and equipment.
- This product has a vinyl coating. While this coating is intended to enhance product durability and usability and is safe for skin contact, it may cause allergic reactions in individuals sensitive to vinyl or its derivatives.
- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
 - Coating: UL 94 HBF
 - Foam: Rated either UL 94 HBF or CAL 117

CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Allow coated products to reach room temperature after storage or transportation in cold temperatures.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid storing where sharp objects could pierce the coating of the product.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.

- Prior to using this product, read the Instructions for Use.
- Discontinue use if foam has degraded or lost resilience, or if coating is torn, cracked, or loose and wrinkles, as this may lead to acute pressure points.

SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

PACKAGING

The product is packaged in a clear polythene bag inside a corrugated cardboard box.

SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
REUSABLE SETS			
104.SLPS	Sideline: Torso Base, Lower Base, Ax Bump, Hip Bump, the Shark Small with Straps, and one Single-Use Kit	Nonsterile, reusable positioning system with a nonsterile, single-use positioning set intended to secure, support, and stabilize the patient in the lateral decubitus position and provide traction on the pelvis anatomy for lateral lumbar spine access.	Indicated for use by healthcare professionals to position the patient's full body for lateral lumbar procedures.
SINGLE-USE SETS			
104.SLPS-SUK-BX5*	Single-Use Kit Box of 5	Nonsterile, single-use positioning set intended to stabilize the patient with the desired amount of lateral cervical bending, pad between the patient's legs, secure the patient's torso, and provide the desired amount of permanent hip traction in the lateral decubitus position.	Indicated for use by healthcare professionals with the Sideline to position the patient's head, pad the legs, secure the arms, and manipulate the pelvis for lateral lumbar procedures.
REPLACEMENT COMPONENTS			
925.SHARK-S-WS	Shark Small	Nonsterile, reusable positioner intended to support and stabilize the arms and block the torso from rolling anteriorly in the lateral position.	Indicated for use by healthcare professionals to position the patient's arms and torso on a standard armboard during medical procedures not involving the arms.
104.SLPS-TB	Torso Base	Nonsterile, reusable base component intended to support and stabilize the upper body in the lateral position. The Torso Base pads the upper body, secures to the operating room table using table straps, connects to a single-use torso strap, and acts as a platform for the adjustable Lateral Head Support, Ax Bump, and Hip Bump.	Indicated for use by healthcare professionals to position the patient's upper body for lateral lumbar procedures. Part of the Sideline Positioning System.
104.SLPS-LB	Lower Base	Nonsterile, reusable base component intended to support the lower extremities and provide traction to move the pelvis anatomy away from the ribs in the lateral position. The Lower Base pads the legs, secures to the table using a single-use hip strap and cross-traction straps, and acts as a platform for the Sideline Hip Bump.	Indicated for use by healthcare professionals to position the patient's pelvis and lower body for lateral lumbar procedures. Part of the Sideline Positioning System.
104.SLPS-AB	Ax Bump	Nonsterile, reusable positioner intended to elevate the axilla and support and stabilize the torso in the lateral decubitus position.	Indicated for use by healthcare professionals to position the patient's axilla for medical procedures.
104.SLPS-HB	Hip Bump	Nonsterile, reusable positioner intended to elevate the pelvis and support and stabilize the trunk in the lateral decubitus position.	Indicated for use by healthcare professionals to position the patient's pelvis for medical procedures. Part of the Sideline Positioning System.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
ADDITIONAL COMPONENTS			
925.SHARK-WS	Shark Regular	Nonsterile, reusable positioner intended to support and stabilize the arms and block the torso from rolling anteriorly in the lateral position.	Indicated for use by healthcare professionals to position the patient's arms and torso on a standard armboard during medical procedures not involving the arms.

*  applies only to these products.  does not apply to these products.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
 <small>bonefoam.com/eIFU</small>	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the device is intended for one single-use only
	Indicates the temperature limits to which the medical device can be stored, transported or used
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

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