SELECT LANGUAGE:

ENGLISH

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SPANISH • ESPAÑOL

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GERMAN • DEUTSCH

FINNISH • SUOMI

SWEDISH • SVENSKA

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▲ BONEFOAM™
RAMP
 LEG ELEVATION

The Ramp product line offers reusable, freestanding positioners designed to position the patient's operative leg in an elevated position compared to the nonoperative leg for lower extremity procedures either:

- **parallel to the operative surface in the supine position**
 (applies to 925.TR, 925.TRS, 925.PTR, 925.SCTR and 925.PSCTR)
- **with the knee bent in the supine position** (applies to the 925.SER)
- **in the prone position**
 (applies to 925.TRPP2)

The following disposable products are recommended for use with the Ramps, but are not included:

- 925.TRSB-BX20 - Ramp Sterile Bags - Box of 20 (fit 925.TR, 925.TRS, 925.PTR, 925.SCTR, 925.PSCTR and 925.SER)

To order, please contact
sales@bonefoam.com

BoneFoam, Inc.
 2601 49th Ave N Suite 100
 Minneapolis, MN 55430
 USA

Email: customerservice@bonefoam.com
 Phone: (763) 559-1830
 Toll-free: 1 (877) 861-2663
 Fax: (763) 559-1822


BoneFoam B.V.
 P.O. Box 10110, NL-5060 GA Oisterwijk
 Kantstraat 19, 5076 NP Haaren
 EAR@qnetbv.eu

bonefoam.com

RAMP-eIFU Rev.1 2026-03-15

INSTRUCTIONS FOR USE

CLEANING AND DISINFECTING INSTRUCTIONS

- Do not sterilize the product. This product is not intended to be sterilized.
- For coated products, wipe the product with a clean wipe and disinfectants after every use. The following disinfecting materials are known to be compatible with coated BoneFoam products:
 - Rubbing alcohol (70% isopropyl alcohol)
 - General purpose/all purpose cleaner
 - Read and follow the cleaning product's instructions.
- Make sure that the coated product is dry before you store it or use it again.

WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- Incorrect positioning of the patient with this product can increase the risk of:
 - Pressure injury resulting in markings or ulcers
 - Peripheral nerve injury, which can lead to temporary or permanent paralysis
 - Venous pressure, which can lead to ischemia
 - Instability, which can lead to patient falls
- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- Ensure device is compatible with work surface, procedural instruments, and equipment.
- This product has a vinyl coating. While this coating is intended to enhance product durability and usability and is safe for skin contact, it may cause allergic reactions in individuals sensitive to vinyl or its derivatives.

- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
 - Coating: UL 94 HBF
 - Foam: Rated either UL 94 HBF or CAL 117

CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Allow coated products to reach room temperature after storage or transportation in cold temperatures.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid storing where sharp objects could pierce the coating of the product.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.
- Prior to using this product, read the Instructions for Use.
- Discontinue use if foam has degraded or lost resilience, or if coating is torn, cracked, or loose and wrinkles, as this may lead to acute pressure points.

SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.


PACKAGING


The product is packaged in a clear polythene bag inside a corrugated cardboard box.

SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
925.TR	Ramp	Nonsterile, reusable patient positioner intended to elevate, support, and stabilize the leg parallel to the operative surface in the supine position.	Indicated for use by healthcare professionals to position the portion of the patient's operative leg below the knee for lower extremity procedures.
925.TRS	Ramp Shorty		
925.PTR	Ramp Pediatric		
925.SCTR	Ramp Surgery Center		
925.PSCTR	Ramp Surgery Center Pediatric		
925.SER	Ramp Semi-Extended	Nonsterile, reusable patient positioner intended to elevate, support, and stabilize the leg with the knee bent in the supine position.	Indicated for use by healthcare professionals to position the patient's operative leg for lower extremity procedures.
925.TRPP2	Ramp Prone	Nonsterile, reusable patient positioner intended to elevate, support, and stabilize the legs in the prone position.	Indicated for use by healthcare professionals to position the portion of the patient's operative leg below the knee and position the foot of the nonoperative leg for lower extremity procedures.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
 bonefoam.com/eIFU	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the temperature limits to which the medical device can be stored, transported or used
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

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