
▲ BONEFOAM™


# NEXTEND

ANTERIOR CERVICAL  
POSITIONING SYSTEM

104.NPS

104.NPL

## SELECT LANGUAGE:

ENGLISH

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SPANISH • ESPAÑOL

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GERMAN • DEUTSCH

FINNISH • SUOMI

SWEDISH • SVENSKA

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Nextend | System

Nextend | Lite

The Nextend Anterior Cervical Positioning System is a tabletop positioning set with both reusable and single-use components designed to position a patient for anterior cervical procedures.

The following disposable products are required for use with the Nextend | System and Nextend | Lite:

- 104.NPS-SUK-BX5  
Nextend | System: Single-Use Kit  
Box of 5
- 104.NPL-SUK-BX5  
Nextend | Lite: Single-Use Kit  
Box of 5

To reorder, please contact  
[sales@bonefoam.com](mailto:sales@bonefoam.com)

BoneFoam, Inc.  
2601 49th Ave N Suite 100  
Minneapolis, MN 55430  
USA

Email: [customerservice@bonefoam.com](mailto:customerservice@bonefoam.com)  
Phone: (763) 559-1830  
Toll-free: 1 (877) 861-2663  
Fax: (763) 559-1822


BoneFoam B.V.  
P.O. Box 10110, NL-5060 GA Oisterwijk  
Kantstraat 19, 5076 NP Haaren  
[EAR@qnetbv.eu](mailto:EAR@qnetbv.eu)

[bonefoam.com](http://bonefoam.com)


NPS-eIFU Rev.1 2026-03-15

## INSTRUCTIONS FOR USE

104.NPS


104.NPL

1

Place positioner(s) on table and secure with table straps. Place single-use head support on scapula bolster. Attach black disposable patient straps on top two lateral strap slots on one side of torso base (applies only to 104.NPS).

104.NPS

104.NPL

2

Transfer patient to Nextend following facility transfer protocol. Position so occupant's head is in head bolster hole and shoulders overhang edge of scapula bolster. Knees should be positioned over apex of Dual Leg Positioner (applies only to 104.NPS).


3

Secure patient with black arm straps (applies only to 104.NPS).


4

Attach blue Y-shaped traction strap to torso base, with gray side facing patient's skin, and place over patient. Adjust horizontal yoke so straps are over patient's AC joints.

5

Remove layers of head support to desired amount of cervical extension.

6

If using permanent traction, secure traction anchor to the foot of the bed. Feed the Y-strap through the slot of the traction anchor, pull desired traction, and secure y-strap to itself.

## CLEANING AND DISINFECTING INSTRUCTIONS

- Do not sterilize the product. This product is not intended to be sterilized.
- For coated products, wipe the product with a clean wipe and disinfectants after every use. The following disinfecting materials are known to be compatible with coated BoneFoam products:
  - Rubbing alcohol (70% isopropyl alcohol)
  - General purpose/all purpose cleaner
    - Read and follow the cleaning product's instructions.
- Make sure that the coated product is dry before you store it or use it again.

## WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- For single-use accessories, risk of cross-contamination if not disposed of immediately after use.
- Incorrect positioning of the patient with this product can increase the risk of:
  - Pressure injury resulting in markings or ulcers
  - Peripheral nerve injury, which can lead to temporary or permanent paralysis
  - Venous pressure, which can lead to ischemia
  - Instability, which can lead to patient falls
- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- Ensure device is compatible with work surface, procedural instruments, and equipment.
- This product has a vinyl coating. While this coating is intended to enhance product durability and usability and is safe for skin contact, it may cause allergic reactions in individuals sensitive to vinyl or its derivatives.
- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
  - Coating: UL 94 HBF
  - Foam: Rated either UL 94 HBF or CAL 117

## CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Allow coated products to reach room temperature after storage or transportation in cold temperatures.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid storing where sharp objects could pierce the coating of the product.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.
- Prior to using this product, read the Instructions for Use.
- Discontinue use if foam has degraded or lost resilience, or if coating is torn, cracked, or loose and wrinkles, as this may lead to acute pressure points.

## SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

## PACKAGING

The product is packaged in a clear polythene bag inside a corrugated cardboard box.

## SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
<b>REUSABLE SETS</b>			
104.NPS	Nextend   System: System Bolster, Dual, Traction Anchor, and one Nextend Single-Use Kit	Nonsterile, reusable positioning system which includes a nonsterile, single-use positioning set intended to secure, support, and stabilize the patient in the supine position with desired amount of cervical flexion/extension and provide the desired amount of intermittent or permanent shoulder traction.	Indicated for use by healthcare professionals to position the patient's full body for anterior cervical procedures.
104.NPL	Nextend   Lite: Lite Bolster, Traction Anchor, and one Nextend Lite Single-Use Kit	Nonsterile, reusable positioning system which includes a nonsterile, single-use positioning set intended to secure, support, and stabilize the patient in the supine position with desired amount of cervical flexion/extension and provide the desired amount of intermittent or permanent shoulder traction.	Indicated for use by healthcare professionals to position the patient's full body for anterior cervical procedures.
<b>SINGLE-USE SETS</b>			
104.NPS-SUK-BX5*	Nextend   System: Single-Use Kit Box of 5	Nonsterile, single-use positioning set intended to stabilize the patient with the desired amount of cervical flexion/extension, secure the patient's arms, and provide the desired amount of intermittent or permanent shoulder traction necessary in the supine position.	Indicated for use by healthcare professionals with the Nextend Positioning System, to position the patient's head and secure their arms for anterior cervical procedures.
104.NPL-SUK-BX5*	Nextend   Lite: Single-Use Kit Box of 5	Nonsterile, single-use positioning set intended to stabilize the patient with the desired amount of cervical flexion/extension and provide the desired amount of intermittent or permanent shoulder traction necessary in the supine position.	Indicated for use by healthcare professionals with the Nextend   Lite, to position the patient's head and shoulders for anterior cervical procedures.

\*  applies only to these products.  does not apply to these products.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
<b>REPLACEMENT COMPONENTS</b>			
925.DLS-WS	Duall	Nonsterile, reusable lower extremity positioner intended to support and stabilize the legs in the supine position.	Indicated for use by healthcare professionals as a nonoperative surface to position the patient's legs for medical procedures not involving the legs.
104.NPSB	Nextend   System: Bolster	Nonsterile, reusable bolster component intended to support the patient's head, neck, and torso during surgery in the supine position. The Bolster secures to the table with a table strap, and connects to single-use patient arm and shoulder traction straps.	Indicated for use by healthcare professionals to position the patient's torso for anterior cervical procedures. To be used with the Nextend Single-Use Kit, which is part of the Nextend Positioning System.
104.NPLB	Nextend   Lite: Bolster	Nonsterile, reusable bolster component intended to support the patient's head, neck, and upper torso during surgery in the supine position. The Bolster secures to the table with a table strap, and connects to single-use shoulder traction straps.	Indicated for use by healthcare professionals to position the patient's torso for anterior cervical procedures. To be used with the Nextend   Lite Single-Use Kit, which is part of the Nextend   Lite.
104.NPS-TA	Traction Anchor	Nonsterile, reusable plate, which is secured to foot side of table using a table strap. Features a connection point for securing a patient traction strap for sustained traction.	Indicated for use by healthcare professionals to position the patient's shoulders for anterior cervical procedures. To be used with the Single-Use Kit, which is part of the Nextend Positioning System.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the device is intended for one single-use only
	Indicates the temperature limits to which the medical device can be stored, transported or used
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

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