SELECT LANGUAGE:

ENGLISH

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SPANISH • ESPAÑOL

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GERMAN • DEUTSCH

FINNISH • SUOMI

SWEDISH • SVENSKA

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UNIVERSAL PRONE

POSITIONING SYSTEM

The Universal Prone Positioning System is a two-piece, reusable tabletop positioner designed to position a patient for prone medical procedures.

The following disposable products are recommended for use with the Universal Prone Positioning System:

- 104.UPC-SUK-BX5 – Universal Prone | Cervical – Single Use Kit – Box of 5
- 925.FGC-BX10 – Face Guard Cushion – Box of 10
- 925.SG-BX10 – Shin Guard – Box of 10
- 925.UPS-DBC-BX20 – Disposable Barrier Covers – Box of 20 (Compatible with 925.UPS-S and 925.UPS-R; not compatible with 925.UPS-L)

To order, please contact
sales@bonefoam.com

BoneFoam, Inc.
2601 49th Ave N Suite 100
Minneapolis, MN 55430
USA

Email: customerservice@bonefoam.com
Phone: (763) 559-1830
Toll-free: 1 (877) 861-2663
Fax: (763) 559-1822


BoneFoam B.V.
P.O. Box 10110, NL-5060 GA Oisterwijk
Kantstraat 19, 5076 NP Haaren
EAR@qnetbv.eu

bonefoam.com

UPS-eIFU Rev.1 2026-07-10


INSTRUCTIONS FOR USE



Recommended table set up for general prone procedures.

1

Place upper and lower body pieces on table over table break.

Recommended table set up for prone spine procedures with pedestal positioned out of the way of imaging equipment areas.


2

Cover positioner with disposable barrier cover.


3

Transfer patient onto positioner following hospital transfer protocols.

4

Adjust patient to ensure breast tissue, abdomen, and genitalia are positioned in corresponding cavities.


5

Route catheters through lower body channels.


*When breaking the table at the patient's hips be aware the patient may drift or move caudal down the positioner. If intending to break the table and lower the patient's legs and/or flex their hips, position the patient higher up the positioner (more cephalad) before breaking the table.


925.UPS-L

925.UPS-R

925.UPS-S

CLEANING AND DISINFECTING INSTRUCTIONS

- Do not sterilize the product. This product is not intended to be sterilized.
- For coated products, wipe the product with a clean wipe and disinfectants after every use. The following disinfecting materials are known to be compatible with coated BoneFoam products:
 - Rubbing alcohol (70% isopropyl alcohol)
 - General purpose/all purpose cleaner
 - Read and follow the cleaning product's instructions.
- Make sure that the coated product is dry before you store it or use it again.

WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- For single-use accessories, risk of cross-contamination if not disposed of immediately after use.
- Incorrect positioning of the patient with this product can increase the risk of:
 - Pressure injury resulting in markings or ulcers
 - Peripheral nerve injury, which can lead to temporary or permanent paralysis
 - Venous pressure, which can lead to ischemia
 - Instability, which can lead to patient falls
- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- Ensure device is compatible with work surface, procedural instruments, and equipment.
- This product has a vinyl coating. While this coating is intended to enhance product durability and usability and is safe for skin contact, it may cause allergic reactions in individuals sensitive to vinyl or its derivatives.
- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
 - Coating: UL 94 HBF
 - Foam: Rated either UL 94 HBF or CAL 117

CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Allow coated products to reach room temperature after storage or transportation in cold temperatures.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid storing where sharp objects could pierce the coating of the product.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.
- Avoid using adhesive directly on the products coating as this can damage the coating over time.
- Prior to using this product, read the Instructions for Use.
- Discontinue use if foam has degraded or lost resilience, or if coating is torn, cracked, or loose and wrinkles, as this may lead to acute pressure points.

SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

PACKAGING

The product is packaged in a clear polythene bag inside a corrugated cardboard box.

SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
 bonefoam.com/eIFU	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the device is intended for one single-use only
	Indicates the temperature limits to which the medical device can be stored, transported or used
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
REUSABLE POSITIONERS			
925.UPS-L	Universal Prone Large	Nonsterile, reusable positioning set intended to support and stabilize the patient with abdominal decompression, breast and genitalia cavities, and enable hip flexion/extension adjustments in the prone position. Intended to provide the desired amount of intermittent or permanent shoulder traction in the prone position when used with a nonsterile, single-use-strap set and nonsterile, reusable traction anchor.	Indicated for use by healthcare professionals to position the torso, pelvis, and thigh of morbidly obese or larger patients during prone medical procedures on standard OR tables with table width extensions. May be used with the Posterior Cervical Strap Kit and Traction Anchor.
925.UPS-R	Universal Prone Regular		Indicated for use by healthcare professionals to position the patient’s torso, pelvis, and thigh during prone medical procedures on standard OR tables. May be used with the Universal Prone Disposable Barrier Covers and/or the Posterior Cervical Strap Kit and Traction Anchor.
925.UPS-S	Universal Prone Small		
SINGLE-USE ACCESSORIES			
925.UPS-DBC-BX20*	Disposable Barrier Cover Box of 20 (not compatible with 925.UPS-L)	Nonsterile, single-use covers intended to provide a barrier between the patient and the positioner.	Indicated for use by healthcare professionals to cover the Universal Prone Regular and Small.
REPLACEMENT COMPONENTS			
925.UPS-L-UB	Upper Body Large	Nonsterile, reusable upper body component intended to support and stabilize the patient with abdominal decompression and breast cavities in the prone position. Intended to provide the desired amount of intermittent or permanent shoulder traction in the prone position when used with a nonsterile, single-use strap set.	Indicated for use by healthcare professionals with the Universal Prone Lower Body Large (part of the Universal Prone Large), to position the torso and pelvis of morbidly obese or larger patients during prone medical procedures on standard OR tables with table width extensions.
925.UPS-L-LB	Lower Body Large	Nonsterile, reusable lower body component intended to support and stabilize the patient with genitalia cavities and enable hip flexion/extension adjustments in the prone position.	Indicated for use by healthcare professionals with the Universal Prone Upper Body Large (part of the Universal Prone Large), to position the pelvis and thighs of morbidly obese or larger patients during prone medical procedures on standard OR tables with table width extensions.

(replacement components continued on next page)

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
REPLACEMENT COMPONENTS			
925.UPS-R-UB	Upper Body Regular	Nonsterile, reusable upper body component intended to support and stabilize the patient with abdominal decompression and breast cavities in the prone position. Intended to provide the desired amount of intermittent or permanent shoulder traction in the prone position when used with a nonsterile, single-use strap set.	Indicated for use by healthcare professionals with the Universal Prone Regular Lower Body (part of the Universal Prone Positioning System Regular), to position the patient's torso and pelvis during prone medical procedures on standard OR tables.
925.UPS-R-LB	Lower Body Regular	Nonsterile, reusable lower body component intended to support and stabilize the patient with genitalia cavities and enable hip flexion/extension adjustments in the prone position.	Indicated for use by healthcare professionals with the Universal Prone Regular Upper Body (part of the Universal Prone Positioning System Regular), to position the patient's pelvis and thighs during prone medical procedures on standard OR tables.
925.UPS-S-UB	Upper Body Small	Nonsterile, reusable upper body component intended to support and stabilize the patient with abdominal decompression and breast cavities in the prone position. Intended to provide the desired amount of intermittent or permanent shoulder traction in the prone position when used with a nonsterile, single-use strap set.	Indicated for use by healthcare professionals with the Universal Prone Small Lower Body (part of the Universal Prone Positioning System Small), to position the patient's torso and pelvis during prone medical procedures on standard OR tables.
925.UPS-S-LB	Lower Body Small	Nonsterile, reusable lower body component intended to support and stabilize the patient with genitalia cavities and enable hip flexion/extension adjustments in the prone position.	Indicated for use by healthcare professionals with the Universal Prone Small Upper Body (part of the Universal Prone Positioning System Small), to position the patient's pelvis and thighs during prone medical procedures on standard OR tables.

*  applies only to these products.  does not apply to these products.

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