SELECT LANGUAGE:

ENGLISH

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SPANISH • ESPAÑOL

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GERMAN • DEUTSCH

FINNISH • SUOMI

SWEDISH • SVENSKA

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UNIVERSAL PRONE

CERVICAL | SINGLE-USE KIT

The Universal Prone Cervical Single-Use Kit is a single-use positioning strap kit designed to achieve shoulder traction and arm restraint of a patient during posterior cervical procedures.

The following product is designed to be used with the Universal Prone | Large, Universal Prone | Regular, Universal Prone | Small:

- 104.UPC-SUK-BX5 – Universal Prone | Cervical – Single Use Kit – Box of 5

To reorder this disposable product, please contact sales@bonefoam.com

BoneFoam, Inc.
2601 49th Ave N Suite 100
Minneapolis, MN 55430
USA

Email: customerservice@bonefoam.com
Phone: (763) 559-1830
Toll-free: 1 (877) 861-2663
Fax: (763) 559-1822


BoneFoam B.V.
P.O. Box 10110, NL-5060 GA Oisterwijk
Kantstraat 19, 5076 NP Haaren
EAR@qnetbv.eu

bonefoam.com


UPC-SUK-eIFU Rev.0 2026-07-10

INSTRUCTIONS FOR USE

1

Attach Y-shaped traction strap to the connection points on upper body portion of the Universal Prone so the hook-and-loop sections where the horizontal yoke portion is connected are facing away from the patient's skin. Adjust horizontal yoke so straps are over patient's AC joints.


2

Lay the H-shaped strap centered on the patient's back and align the cutout with the antecubital space. Weave strap arms underneath patient's arms.

3

Route straps back over the top of the patient's arms and secure hook portion to the middle section.

4

If using permanent traction, secure traction anchor to the foot of the bed. Feed the Y-strap through the slot of the traction anchor, pull desired traction, and secure Y-strap to itself.

CLEANING AND DISINFECTING INSTRUCTIONS

- Do not sterilize the product. This product is not intended to be sterilized.
- For coated products, wipe the product with a clean wipe and disinfectants after every use. The following disinfecting materials are known to be compatible with coated BoneFoam products:
 - Rubbing alcohol (70% isopropyl alcohol)
 - General purpose/all purpose cleaner
 - Read and follow the cleaning product's instructions.
- Make sure that the coated product is dry before you store it or use it again.

WARNINGS/PRECAUTIONS/CONTRAINDICATIONS

- Intended to contact intact skin only.
- Risk of cross-contamination if not disposed of immediately after use.
- Incorrect positioning of the patient with this product can increase the risk of:
 - Pressure injury resulting in markings or ulcers
 - Peripheral nerve injury, which can lead to temporary or permanent paralysis
 - Venous pressure, which can lead to ischemia
 - Instability, which can lead to patient falls
- Product is to be used by licensed medical professionals only.
- Not intended for patient transport.
- Ensure device is compatible with work surface, procedural instruments, and equipment.
- Possibility of artifacts in imaging applications depending on setup and imaging settings.
- Risk of fire. Below are flammability ratings of the materials:
 - Laminate Straps: Rated FMVSS-302

CARE, HANDLING AND USEFUL LIFE

- Avoid using the product if package is damaged during transportation.
- Transport and store in dry conditions and out of direct light where possible.
- To prevent patient and/or user injury and/or equipment damage, examine the product for potential damage or wear prior to use. Do not use the product if damage is visible, if parts are missing, or if it does not function as expected.
- Avoid stacking items on top of product when storing.
- Do not submerge the product in water. Equipment damage can occur.
- Prior to using this product, read the Instructions for Use.
- Discard and do not use product if laminate strap material has degraded or delaminated.

SAFE DISPOSAL INSTRUCTIONS

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

PACKAGING

The product is packaged in a clear polythene bag inside a corrugated cardboard box. Product may be compressed to minimize packaging and storage footprint.

SAFETY INFORMATION

The techniques detailed in this manual are only manufacturer's suggestions. The final responsibility for patient care with respect to this product remains with the attending physician. Any serious incident that has occurred in relation to this product should be reported to the manufacturer and regulatory authority, where required. Keep this manual available for future reference.

SKU	DESCRIPTION	INTENDED USE	INDICATIONS
104.UPC-SUK-BX5	Universal Prone Cervical – Single Use Kit – Box of 5	Nonsterile, single-use positioning set intended to secure the patient's arms and provide the desired amount of intermittent or permanent shoulder traction necessary in the prone position.	Indicated for use by healthcare professionals with the Universal Prone Positioning Systems to distract the patient's shoulders and secure their arms for posterior cervical procedures.

SYMBOL	DESCRIPTION
	Indicates the device is a medical device
	Indicates the medical device manufacturer
	Indicates the country where the medical device was manufactured and date of manufacture
	Indicates the manufacturer's lot number
	Indicates the medical device Global Trade Item Number
	Indicates the manufacturer's catalogue number
	Indicates the unique product identifier information for the device
 bonefoam.com/eIFU	Indicates the need for the user to consult the instruction for use
	Indicates the device does not contain natural rubber or dry natural rubber latex

SYMBOL	DESCRIPTION
	Indicates the device is intended for one single-use only
	Indicates a medical device should not be used if the package has been damaged or opened
	Indicates quantity in package
	Indicates the European Authorized Representative name and address in the European Community
	Signifies European technical conformity
	Signifies Great Britain technical conformity
	Indicates a UK Responsible Person on medical device information supplied by the manufacturer
	Indicates an Authorized Representative in Switzerland on medical information supplied by the manufacturer

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