

PICOHI™ 300

300 ps

Experience photoacoustic effect  
PICOHI 300 picosecond laser



For Export





## Experience, the real photoacoustic effect

# 1

Next generation pico laser **Real 300ps pulse duration**

Stable technology providing pulse duration of 300ps.

# 2

Non-photothermal effect, **Photomechanical effect treatment**

Upgraded from photothermal effect treatment, photomechanical effect treatment allows selective treatment on targeted areas.


# 3





Various handpieces for **improving of skin rejuvenation**

DOE & MLA handpieces can provide various types of treatment in multiple sizes and depths depending on the treatment area.

# Why 300ps?


- Stress Relaxation Time(SRT) of  $1\mu\text{m}$  melanosome is 300ps. The shock wave caused by the picosecond pulse acts only on the melanosome without being transmitted to the surrounding tissues.
- Thermal Relaxation Time(TRT) of 10–100nm tattoo ink particle is 100ps–10ns. 300ps is equipped with appropriate pulse duration to destroy smaller particles than other picoseconds lasers.
- Micro bubble occurs in the skin more strongly by high peak power related to shorter pulse duration. High peak power with 300ps pulse duration increases LIOB formation and skin rejuvenation effect.




| Nano Laser   | Pico Laser  |
|--|---|
| Photothermal effect  | Photomechanical effect  |
| Not precise destruction of melanin pigment<br>Damages on near by tissues   | Precise destruction of melanin pigment<br>Targeting on the pigment is possible with great damage effects  |
| <div data-bbox="377 531 1156 853">  </div> <div data-bbox="389 936 698 1186">  </div> <div data-bbox="377 1193 698 1222"> <p>[Marks after Nanosecond laser]</p> </div> <div data-bbox="769 1009 1149 1120"> <p>The heat damages melanin pigment and affects the surrounding tissue.</p> </div> | <div data-bbox="1396 531 2181 853">  </div> <div data-bbox="1409 936 1717 1186">  </div> <div data-bbox="1409 1193 1717 1222"> <p>[Marks after Picosecond laser]</p> </div> <div data-bbox="1763 1006 2219 1115"> <p>Physical destruction on melanin pigment while having no affects on the surrounding tissue.</p> </div> |

# High Power

300 picosecond pulse duration



Peak power higher than Q-switch or other pico lasers provide maximum effect on tattoo removal, pigment treatment and skin rejuvenation.






# High Effect

## LIOB effect with PICOHI VMLA handpiece (biopsy result)

Bubble within the skin helps to regenerate tissues which create new collagen for skin rejuvenation effect. MLA & DOE handpieces of PICOHI provides effective rejuvenation treatment.



## Rejuvenation



### VMLA H/P

Precise treatment at different layers from Deep dermis to Epidermis with 13mm beam size.

Depth depending on the target can be adjusted (0.5mm)



### ZMLA H/P

Cover the various lesion sizes and providing effective results.

Adjustable Beam Size depending on the treatment (4~12mm)



### DOE H/P

Synergy rejuvenation effect with MLA handpieces.

Very even energy uniformity with 10mm(7\*7 dots) spot size

## Pigmentation



### Collimated H/P

Stable and effective toning treatment with same spot size regardless of laser separation distance.

1064nm Collimated / 10mm spot size



### Zoom H/P

Various pigmentation treatment possible.

532nm · 1064nm basic zoom / 2~10mm spot size

## Tattoo Removal

### Red and Yellow

Color treatment

Zoom H/P  
532nm · 1.5~7.5mm spot size

### Black

Color treatment

Zoom H/P  
1064nm · 2~10mm spot size

## Clinical Data

### Lentigine Treatment | after 1 treatment



532nm Zoom : 3mm · 0.3J/cm<sup>2</sup> · 1pass

### Tattoo Removal Treatment | after 1 treatment



1064nm Zoom : 3mm · 1.96J/cm<sup>2</sup> · 1pass


### Scratched Scar Treatment | after 1 treatment



ZMLA : 5mm · 1.3J/cm<sup>2</sup> · Stacking




**Zoom H/P**  
1064nm, 532nm




**Collimated H/P**  
1064nm



**VMLA H/P**  
1064nm, 532nm



**ZMLA H/P**  
1064nm, 532nm



**DOE H/P**  
1064nm



**DOE H/P**  
532nm





VMLA

PICOHI™300

ZMLA

Depth depending on the target can be adjusted  
(0.5mm)

Feature

Adjustable Beam Size depending on the treatment  
(4~12mm)


Precise treatment at different layers  
from Deep dermis to Epidermis

Function

Adjustable Beam Size depending on the treatment  
for effective results





Beam mode image



[ ZMLA LASER TOP VIEW ]

|             |   |    |    |    |   |
|-------------|---|----|----|----|---|
| Beam Size   | 4 | 5  | 6  | 7  | 8 |
| Energy (mj) |   |    |    |    |   |
| Beam Size   | 9 | 10 | 11 | 12 |   |
| Energy (mj) |   |    |    |    |   |

| MLA(VMLA, ZMLA)  |   | PICO Fractional   | DOE(532nm, 1064nm)                          |  |
|--|---|-------------------|---|--|
|  | Very High                                   | Energy Efficiency | Relatively Low                              |  |
|  | Relatively Low                              | Energy Uniformity | Very Even                                   |  |
|  | Scar Curing<br>Rejuvenation<br>Pigmentation | Scope of Use      | Pigmentation<br>Rejuvenation<br>Scar Curing |  |
| Real MLA Beam  |   |                   |   | Real DOE Beam  |

Output Energy : 500mJ / 5Hz

PICOHI provides **MLA**, **DOE** for selective use for effective treatment

Comparison of Minimum Energy at Zoom HP(1064nm)

| Spot Size | Minimum Energy(mJ) |             | Minimum Fluence(J/cm²) |             |
|-----------|--------------------|-------------|------------------------|-------------|
|           | 450ps              | PICOHI™ 300 | 450ps                  | PICOHI™ 300 |
| 2mm       | 131                | 3.1         | 4.1                    | 0.1         |
| 3mm       | 134                | 7.1         | 1.9                    | 0.1         |
| 4mm       | 138                | 12.6        | 1.1                    | 0.1         |
| 5mm       | 137                | 19.6        | 0.7                    | 0.1         |
| 6mm       | 141                | 28.3        | 0.5                    | 0.1         |
| 7mm       | 153                | 38.5        | 0.4                    | 0.1         |
| 8mm       | X                  | 50.2        | X                      | 0.1         |
| 9mm       | X                  | 63.6        | X                      | 0.1         |
| 10mm      | X                  | 78.5        | X                      | 0.1         |

Comparison of Minimum Energy at Zoom HP(532nm)

| Spot Size | Minimum Energy(mJ) |             | Minimum Fluence(J/cm²) |             |
|-----------|--------------------|-------------|------------------------|-------------|
|           | 450ps              | PICOHI™ 300 | 450ps                  | PICOHI™ 300 |
| 1.5mm     | X                  | 1.8         | X                      | 0.1         |
| 2mm       | 72                 | 3.1         | 2.3                    | 0.1         |
| 3mm       | 70                 | 7.1         | 1.0                    | 0.1         |
| 4mm       | 75                 | 12.6        | 0.6                    | 0.1         |
| 5mm       | 78                 | 19.6        | 0.4                    | 0.1         |
| 6mm       | 84                 | 28.3        | 0.3                    | 0.1         |
| 6.5mm     | X                  | 33.2        | X                      | 0.1         |
| 7mm       | 76                 | 38.5        | 0.2                    | 0.1         |
| 7.5mm     | X                  | 44.2        | X                      | 0.1         |

The Fluence(J/m²) which indicates "Energy per unit area" can be adjusted from 0.1J/m² to a wide energy range in a sensitive way compared to 450ps.



# High Convenience

 **RMS system** (Remote Maintenance System)

## Patient History Record

Fast and smart treatment is possible as the parameters for previous treatments will be recoded and memorized. The saved parameters can be managed through history for future treatments which may be a similar case.

**TREATMENT INFORMATION**

Find View  
JANETE SMITH

**FEEDBACK SURVEY**

07.23.2018 12:05 ZOOM 532

Very Poor Poor Average Good Excellent

😊 😐 😊 😄 😁

CANCEL OK

BACK TREATMENT

**TREATMENT INFORMATION**

Treatment history Patients information

Hand piece: Zoom Cx1 MLJ DOE DYE RLL Mode: S32 Period: 1 year Lesion: All Feedback: All

| Time           | jet  | MODE | FLUENCE | RATE | SPOT | SPOT  | LESION  | FB |
|----------------|------|------|---------|------|------|-------|---------|----|
| 07.23.18 12:05 | Zoom | S32  | 50.0    | 10   | 10   | 10000 | Freckle | 😊  |
| 05.20.18 14:32 | Zoom | S32  | 40.0    | 8    | 8    | 10000 | Tanning | 😊  |
| 02.28.18 10:24 | Zoom | S32  | 35.0    | 6    | 6    | 8000  | Tattoo  | 😊  |
| 07.23.18 12:05 | Zoom | S32  | 50.0    | 8    | 8    | 8000  | Leflip  | 😊  |
| 05.20.18 14:32 | Zoom | S32  | 40.0    | 10   | 10   | 10000 | Pixi    | 😊  |
| 02.28.18 10:24 | Zoom | S32  | 35.0    | 7    | 7    | 10000 | Freckle | 😊  |
| 07.23.18 12:05 | Zoom | S32  | 50.0    | 6    | 6    | 10000 | Tanning | 😊  |
| 05.20.18 14:32 | Zoom | S32  | 40.0    | 6    | 6    | 8000  | Tattoo  | 😊  |


STATISTICS TREATMENT

| Real Time Device Record Check |



## Remote Management through Cloud System

- Lamp and consumable replacement notice
- Protocols saved
- Self diagnosis and inspection
- Malfunction prevention and assist functions



## Easy & Intuitive GUI

### Treatment Information






### Treatment Malfunction Prevention (coloring system)



Different colored GUI (Blue for 1064nm, Green for 532nm) for exact and safe treatment.

## Easy to Use





## PICOHI Specification

|                              |                  |                                   |
|------------------------------|------------------|-----------------------------------|
| Wavelengths                  |                  | Nd:YAG 1064nm, 532nm              |
| Pulse Duration               |                  | 300ps, 275ps                      |
| Peak Power                   |                  | 1.67GW, 0.91GW                    |
| Pulse Energy[mJ] Max         |                  | 500mJ, 250mJ                      |
| Repetition                   |                  | 1, 2, 5, 10Hz                     |
| Handpiece                    | Zoom H/P         | 1064nm (Spot size : 2 to 10mm)    |
|                              |                  | 532nm (Spot size : 1.5 to 7.5mm)  |
|                              | Collimated H/P   | 1064nm (Spot size : 10mm)         |
|                              | VMLA H/P         | 1064nm (Beam size : 13mm)         |
|                              | ZMLA H/P         | 1064nm (Beam size: 4 to 12mm)     |
|                              | DOE H/P          | 1064nm (Beam size: 10 x 10mm)     |
| 532nm (Beam size: 10 x 10mm) |                  |                                   |
| Convenience system           |                  | RMS System (WIFI)                 |
|                              |                  | Wide LCD 12.1"                    |
|                              |                  | Upper shelf display zone          |
|                              |                  | Handpiece storage                 |
|                              |                  | Hidden front handle & Back handle |
| General                      | Electrical Power | 200~240 VAC, 4.4 KVA, 50/60Hz     |
|                              | Dimensions       | 455(W) x 1040(L) x 975(H) mm      |
|                              | Weight           | 150kg                             |