ACTHRD™

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PARPi

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



ACTHRD™ 報告

ACTHRD™ Report

Identifier
Project ID:
Report No.:
Report Date:

Subject		
Identifier:	Subject ID:	
Date of Birth:	Gender:	
Diagnosis: Ovarian Cancer		
Ordering Physician		
Referral Doctor:	Tel:	
Referral Institution:		
Address:		
Specimen		
Specimen ID:	Collection Site:	Specimen Type:
Date Received:	Sample ID:	D/ID:

ABOUT ACTHRD™

The test is a next-generation sequencing (NGS) based assay to detect single nucleotide variants (SNVs), small insertions and deletions (InDels), copy number alterations (loss) (CNA-loss), large genomic rearrangements (LGRs) in the *BRCA1* and *BRCA2* genes, and analyze the genomic loss of heterozygosity (LOH) status using DNA isolated from formalin-fixed, paraffin-embedded (FFPE) tumor tissue specimens. Results of the test are used as an aid in identifying cancer patients with positive homologous recombination deficiency (HRD) status for treatment with PARP inhibitors. Additionally, the test is intended to provide tumor profiling on SNVs, InDels, and CNA-loss of 22 homologous recombination repair (HRR) genes for use by qualified health care professionals in accordance with professional guidelines, and is not conclusive or prescriptive for labeled use of any specific therapeutic products.

Report Summary for Actionable Variants/Biomarkers

Homologous Recombination Deficiency (HRD) Status		
HRD Status	Gene/Biomarker	Result
Positive	<i>BRCA1</i>	Negative
	<i>BRCA2</i>	Positive
	Loss of heterozygosity (LOH)	Positive (0.57)

* Homologous Recombination Deficiency (HRD) Positive= Either *BRCA1/2* mutant or LOH positive

S Sensitive A Resistant

Variants/Biomarkers with Clinical Significance (Target Therapy)

Genomic Alterations	Evidence Level 1, 2 (FDA-approved, NCCN guideline)	Evidence Level 3A, 3B, 4 (Others)
LOH positive	S Niraparib, Olaparib	S Rucaparib
<i>BRCA2</i> Y1655*	S Niraparib, Olaparib, Rucaparib	S Talazoparib
<i>BRCA2</i> Heterozygous deletion	S Niraparib, Olaparib, Rucaparib	S Talazoparib

Note:

- The therapeutic agents and possible effects to a given drug are based on mapping the variants/biomarkers with ACT Genomics clinical knowledge database. The mapping results only provide information for reference, but not a medical recommendation.
- Please refer to the corresponding sections for more detailed information about genomic alteration and clinical relevance listed above.

ACT Genomics' laboratory is accredited by CAP (CAP number: 9228096).
ACT Genomics only provides a technical report of the test; please consult a specialist physician to determine the appropriate clinical solution and follow the instructions of the physician. *The results are only valid for the tested sample(s).
ACT Genomics Laboratory (HQ): 3F, No. 346, Xinyi 2nd Rd., Xinyi Dist., Taipei, Taiwan 10646.
Email: service@actgenomics.com | T: +886-2-2795-3660 | F: +886-2-2795-5036

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ACTHRD™ Report

Identifier
Project ID:
Report No.:
Report Date:

Testing Results of Variants/Biomarkers with Clinical Relevance

Homologous Recombination Deficiency (HRD) Status		
Positive		

* Homologous Recombination Deficiency (HRD) Positive= Either *BRCA1/2* mutant or LOH positive

Single Nucleotide and Small InDel Variants		
Gene Alterations	Allele Frequency	Classification
<i>BRCA2</i> c.4965del (Y1655*)	68.7%	Deleterious

Copy Number Alterations			
Gene	Chromosome	Variation	Copy Number
<i>BRCA2</i>	Chr13	Heterozygous deletion	1

Large Genomic Rearrangements	
Gene	Alteration
Not detected	

Loss of Heterozygosity Status	
LOH Status	LOH Score
Positive	0.57

Note:

- CNA, LGR, and LOH status in the tumor were determined based on $\geq 30\%$ tumor purity (calculated by NGS/estimated by the pathologist).
- Samples with tumor purity lower than 40% are unable to be analyzed for the LOH status.
- Samples with tumor purity lower than 30% are unable to be analyzed for CNA and LGR.
- The homologous recombination deficiency (HRD) status is defined as deleterious or suspected deleterious alterations of *BRCA1* and *BRCA2*, and/or LOH status positive. The threshold for LOH status positive is set at a score ≥ 0.4 . Of note, the HRD definition is in accordance to the approved therapeutic product labeling in ovarian cancer and not in other cancer types.
- Only deleterious or suspected deleterious variants are listed in this section. For variants of unknown significance (VUS), please refer to Supplementary Information on Testing Results. Variants classified as benign, likely benign, and synonymous variants are not reported.

Supplementary Information for Therapeutic Implications

Targeted Therapies		
Genomic Alterations	Therapies	Evidence Level
LOH positive	S Niraparib, Olaparib	1
<i>BRCA2</i> Y1655*	S Niraparib, Olaparib, Rucaparib	1
<i>BRCA2</i> Heterozygous deletion	S Niraparib, Olaparib, Rucaparib	1
<i>BRCA2</i> Y1655*	S Talazoparib	3A
<i>BRCA2</i> Heterozygous deletion	S Talazoparib	3A
LOH positive	S Rucaparib	3B

Therapies associated with benefit or lack of benefit are based on biomarkers detected in this tumor and published evidence in professional guidelines or peer-reviewed journals.

Level	Description
1	FDA-recognized biomarkers predictive of response or resistance to FDA-approved drugs in this indication
2	Standard care biomarkers (recommended by the NCCN guideline) are predictive of response or resistance to FDA-approved drugs in this indication
3A	Biomarkers predictive of response or resistance to therapies approved by the FDA or NCCN guideline in a different cancer type
3B	Biomarkers that serve as inclusion criteria for clinical trials (minimal supportive data required)
4	Biomarkers that show plausible therapeutic significance based on small studies, few case reports, or preclinical studies

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ACT Genomics' laboratory is accredited by CAP (CAP number: 9228096).
ACT Genomics only provides a technical report of the test; please consult a specialist physician to determine the appropriate clinical solution and follow the instructions of the physician. *The results are only valid for the tested sample(s).
ACT Genomics Laboratory (HQ): 3F, No. 346, Xinyi 2nd Rd., Xinyi Dist., Taipei, Taiwan 10646.
Email: service@actgenomics.com | T: +886-2-2795-3660 | F: +886-2-2795-5036

ACTHRD™

ACT Genomics Co., Ltd.

Member of Delta Group

 www.actgenomics.com

 service@actgenomics.com

 +886-2-2795-3660

©2026 ACT Genomics Co., LTD. All Rights Reserved. It is possible that the test may return with no abnormal mutation identified for certain genes or part of the test results may not be available due to the technical limitations of the test itself and/or an individual's genetic differences and/or intra- or inter-tumoral heterogeneity even the test has been conducted under standard process. This material (including the test and the results) therefrom are intended for educational purposes only for the use of healthcare professionals and do not replace independent professional judgement. This material (including the test and the results) shall at no time be deemed a diagnostic or medical treatment recommendation to any individual. No representation, warranty, express or implied, is made as to, and no reliance should be placed on the fairness, accuracy, completeness or correctness of the information or opinions which may be contained herein. The test takers shall always consult their healthcare professionals for any enquiries of clinical interpretation of the test results.

 ACT GENOMICS® and 帕癌檢® are registered trademarks of ACT Genomics Co., LTD.

TW_DM_20260420_ACTHRD_Brochure_V01