Notice d'utilisation - Français (FR)	2
Instructions for use - English (EN)	71
Benutzerhandbuch - Deutsch (DE)	140
Manuale d'uso - Italiano (IT)	209
Manual de instrucciones - Español (ES)	278
Brugsmeddelelse - Dansk (DK)	347
Használati útmutató - Magyar (HU)	416

LIFEAZ

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INSTRUCTIONS FOR USE

INDEX

1. REGULATORY INFORMATION	6
1.1. Symbols	6
1.2. Glossary	9
2. BOX CONTENTS	10
3. CLARK DESCRIPTION	12
4. GENERAL INFORMATION	14
4.1. Intended purpose	14
4.2. Indications for use	14
4.3. Clinical benefit	14
4.4. Intended patient population	15
4.5. Intended user population	15
4.6. Contraindications and precautions	16
5. SETTING UP YOUR CLARK AED	23
5.1. Upon reception of your CLARK AED	23
5.2. Setting up your CLARK AED	23
5.3. Configuring a Wi-Fi network	27
5.4. Demonstration mode	30
6. HOW TO USE YOUR CLARK AED IN CASE OF EMERGENCY	31
6.1. In case of emergency (A victim is unconscious and not breathing)	31
6.2. If emergency medical services request the data from your Clark	36
7. AFTER USING CLARK AND AFTER THE ARRIVAL OF EMERGENCY MEDICAL SERVICES	37
8. MAINTENANCE	39

8.1. Replacing the battery	39
8.2. Change the electrodes	39
8.3. Upkeep and regular maintenance	40
8.4. Clark connection interfaces	43
8.4.1. USB connection	43
8.4.2. Wi-Fi connection when configuring a Wi-Fi network	45
9. TROUBLESHOOTING	47
9.1. Alarms	47
9.2. Alerts	48
9.3. Other problems	53
10. TECHNICAL INFORMATION	54
10.1. CLARK technical features	54
10.2. CLARK therapy characteristics	58
10.3. Electromagnetic compatibility	62
11. CONNECTIVITY	66
11.1. Mobile connectivity	66
11.2. Wi-Fi Connectivity	67
11.3. GPS	67
12. SEPARATE ACCESSORIES	68

**CLARK, YOUR AUTOMATED
EXTERNAL DEFIBRILLATOR, HAS ARRIVED.**

**READ THIS INSTRUCTIONS FOR USE CLOSELY BEFORE
PROCEEDING INSTALLING AND USING YOUR DEVICE.**

1. REGULATORY INFORMATION

1.1. SYMBOLS

Medical device

AED sticker. International symbol for automated external defibrillators (AED)

IP54

Protection level against water and particle penetration :

5 = protection from dust and other solid particles

4 = protection from water and liquid infiltration

Batch number

Catalog number

Commercial product reference

CE marking : fulfills the european requirements for medical devices. The four numbered product identification number refers to the notified body involved in the device conformity assessment.

Type BF applied part resistant to defibrillation

Remark : classification of electromedical devices according to the ISO 60601-1 standard.

Do not reuse / single-use device

(symbol located on the electrodes pouch)

Do not use if the packaging is damaged

General warning symbol

Humidity limit : Indicates the humidity range the device can be exposed to safely. The humidity limit is located next to the upper and lower horizontal lines.

Atmospheric pressure limit : Indicates the Atmospheric pressure range the device can be exposed to safely.

Unique device identifier

Do not plug in your Clark

Consult the instructions for use

Dangerous voltage

Do not expose to direct sunlight

Keep in a dry area

Serial number

Temperature limit : Indicates the temperature range the device can be exposed to safely. The upper and lower temperature limits are indicated next to the upper and lower horizontal lines. The lowest acceptable temperature is indicated on the left bottom side, the highest acceptable temperature is indicated on the right upper side.

Manufacturer name and address

Use-by date / expiration date : Accompanied by the date that indicates that the device should not be used after the indicated date.

WEEE - Waste from Electrical and Electronic Equipment
European directive concerning waste electrical and electronic equipment.

Non-sterile device (*symbol located on the electrode pouch*)

Latex-free device (*symbol located on the electrode pouch*)

A set of two electrodes is available in each electrode pouch (*symbol located on the electrode pouch*)

Temperature limit to ensure correct electrode use : For immediate use on a victim if the temperature is within 0 °C and 50 °C.

Lithium - manganese dioxide (*battery type*)

Do not deform or damage (*battery label*)

Do not discard in fire (*battery label*)

Do not open or disassemble (*battery label*)

Do not recharge the battery

Distributor name and address (*box label*)

Consult the instructions for use (*box and battery label*)

THROUGHOUT THIS INSTRUCTIONS FOR USE, THE FOLLOWING SYMBOLS ARE DEFINED AS:

WARNING

This symbol is used to indicate a warning.

It indicates a warning, risk, or danger that could lead to serious injury or death.

A precaution describes something that could cause minor bodily injuries, damage to your defibrillator, loss of data stored in CLARK, or lower the chances of successful defibrillation.

This indicates useful information, recommendation, advice, etc.

A summary of safety and clinical performance is available in the Eudamed database, via the following link:

<https://ec.europa.eu/tools/eudamed>

This summary is related to the following Basic UDI-DI:

- 3770022399CLARKXB

- 3770022399PROCEDUREPACKWL

1.2. GLOSSARY

AED	Automated external defibrillator. A device that analyzes a victim's cardiac rhythm and delivers an electric shock if a shockable cardiac rhythm is detected.
CPR	Cardiopulmonary resuscitation Chest compressions performed on a victim of sudden cardiac arrest that replace the pump function of the heart.
Defibrillation	Delivery of an electric shock to the heart in order to reverse ventricular fibrillation or ventricular tachycardia.
ECG	Electrocardiogram A record of the heart's electrical activity.
Fibrillation	Chaotic cardiac activity. It causes the heart's chambers to quiver (rapid and chaotic contractions) and prevents the heart from pumping normally.
Impedance	Resistance to the electric shock delivered by an AED.
LED	Light-emitting diode
Witness or user	Person helping a victim of sudden cardiac arrest
SCA	Sudden cardiac arrest. A sudden loss of heart function. A victim of SCA is unconscious and is not breathing.
Ventricular fibrillation	A life-threatening, abnormal cardiac rhythm caused by quivering ventricles. See the definition of fibrillation above.
Ventricular tachycardia	Rapid ventricular cardiac rhythm. Note : not all tachycardias are pathological. Only rapid ventricular tachycardia requires defibrillation.
Victim	A person suffering from sudden cardiac arrest

2. BOX CONTENTS

ELECTRODES

BATTERY PACK

FIRST AID KIT

INSTRUCTIONS
FOR USE

ELECTRODES INS-
TRUCTIONS FOR USE


The Lifeaz CLARK defibrillator is designed to be used only with the accessories provided by Lifeaz. Your CLARK may not perform correctly if it is used with accessories that have not been approved. Do not connect to another device.

See section 12 Spare accessories for more details about the accessories available for order.

3. CLARK

DESCRIPTION


- | | | | |
|---|---------------|---|-------------|
| 1 | Top surface | 5 | Front side |
| 2 | Rear side | 6 | Left side |
| 3 | Right side | 7 | Status LED |
| 4 | Inferior side | 8 | AED sticker |


Above - closed


Back


Back - Silicone stopper removed


Above - open

Inferior side - Closed battery hatch

Inferior - Opened battery hatch

- | | | |
|------------------------------|--|-----------------------------|
| 1 ON button | 6 Language selected LED | 11 Child mode button |
| 2 ON LED | 7 Language selection buttons | 12 CPR LED |
| 3 Electrode connector | 8 "Do not touch the victim" LED | 13 Status LED |
| 4 Speaker | 9 "Shock delivered" LED | 14 Silicone stopper |
| 5 Battery LED | 10 Child mode LED | 15 Micro USB port |

4. GENERAL INFORMATION

4.1. INTENDED PURPOSE

CLARK is an automated external defibrillator (AED), a medical device designed to automatically assess the patient's heart rhythm and automatically administer a defibrillation shock to victims of sudden cardiac arrest, if required. Using a CLARK defibrillator to deliver the therapeutic electric shock through the heart can restore normal heart rhythm and blood flow.

4.2. INDICATIONS FOR USE

CLARK is indicated for use on victims of sudden cardiac arrest (SCA).

Sudden cardiac arrest is a life-threatening condition that occurs when a person:

- ▶ is not breathing
- ▶ is unconscious

4.3. CLINICAL BENEFIT

CLARK's clinical benefit is the administration of an electric shock to a patient in shockable sudden cardiac arrest, to terminate the arrhythmia and achieve restoration of spontaneous circulation.

Ensure that your device is correctly installed and maintained according to the recommendations located in section 5 Setting up your CLARK AED and 8 Maintenance.

4.4. INTENDED PATIENT POPULATION

CLARK can be used on any patient victim of SCA, adult or child older than 1 year.

4.5. INTENDED USER POPULATION

CLARK is usable by any person above ten (10) years old, witness of a SCA (no need to be a health care professional to use the device, or to have been trained to CPR). Once turned on, CLARK is entirely automatic. It is equipped with a visual guide, and vocally instructs the user step by step during use on a victim of sudden cardiac arrest.

It can therefore be used by untrained users. Local regulations may differ regarding the use of an AED.

Refer to your local regulations to verify that anyone can use the device.

It is recommended that you, your colleagues and your family, learn how to use CLARK in advance. This will allow you to be ready to use it in an emergency.

4.6. CONTRAINDICATIONS AND PRECAUTIONS

There is no known contraindication for the use of an automated external defibrillator.

If you think a person is suffering from sudden cardiac arrest (meaning that the person is unconscious and not breathing), you can use an AED.

There is no risk to aggravate the situation : the defibrillator is designed to analyse a victim's cardiac rhythm and detect if a shock is necessary.

It will not deliver a shock if the victim does not need one.

However, there are certain precautions to take into account :

SECURITY RISK AND POSSIBLE DAMAGE

MRI risk : keep the defibrillator away from magnetic resonance imaging machines (MRI).

POSSIBLE ELECTROMAGNETIC INTERFERENCE WITH THE DEVICE'S PERFORMANCE

Electronic equipment kept in close proximity may emit strong electromagnetic or radio frequencies which could affect this device's performance.

Electromagnetic interference may result in the device's malfunction, a deformed ECG, inability to detect a shockable cardiac rhythm, or prematurely stop defibrillation.

Avoid using the device next to high frequency surgical equipment, mobile phones or other RF communication devices.

WARNING

RF mobile communication equipment (including peripheral telephones, antenna cables and external antennas) should not be used within 30 cm of CLARK, including the cables specified by the manufacturer.

If not, this may negatively affect the device's performance.

WARNING

Avoid using this device next to or on top of other equipment as it may function incorrectly. If such use is necessary, this device and the accompanying equipment should be monitored closely to ensure proper operation.

IMPROPER PERFORMANCE OF THE DEVICE/ ACCESSORIES POSSIBLE

Using damaged/expired equipment or accessories can cause CLARK to function improperly and/or injure the victim or user.

CLARK may not function with other manufacturers' cables or electrodes and will nullify the safety certification.

Only use the Lifeaz accessories identified in the instructions for use (electrodes, batteries).

POSSIBLE ELECTRICAL INTERFERENCE

WARNING

Using accessories other than those provided by this device's manufacturer may increase the electromagnetic emissions or decrease the electromagnetic immunity of this device, resulting in incorrect performance.

Using cables, electrodes or accessories not specified for use with this device can increase emissions or decrease resistance to electromagnetic interferences and may degrade the defibrillator's performance or the performance of other nearby equipment.

Only use the pieces and accessories specified in this instructions for use.

This device can provoke electromagnetic interference, in particular charge and energy transfers.

This interference may affect the performance of other equipment nearby.

If possible, verify the effect of the defibrillator's discharge on other equipment before using the defibrillator in an emergency.

POSSIBLE FIRE OR EXPLOSION

FLAMMABLE GAS

Do not keep the defibrillator near flammable gas or in direct contact with flammable materials. Do not use this defibrillator near flammable gas or anesthetics. Be careful when using the defibrillator near oxygen sources (like a respirator mask or ventilation tube). Turn off the gas source and move it away from the victim during defibrillation.

BATTERY

CLARK's battery is not rechargeable.

Do not try to recharge, open, damage or burn the battery, as it could catch fire or explode.

Do not connect CLARK to a power supply.

LIQUIDS / DUST

Even if Clark is protected against water and other particle infiltration (IP54, see section 1.2 Glossary), avoid spilling liquids on CLARK or its accessories. This could damage the device or result in fire or electric shock. Do not submerge any part of the defibrillator in water or other liquid. If you store Clark outdoors, protect it in a suitable case that ensures it is stored as per the specifications of this user manual. Protect CLARK from excessive dust. Do not sterilize CLARK or its accessories.

TOUCHING OR MOVING THE VICTIM

Performing CPR, touching or moving the victim during CLARK's analysis of the cardiac rhythm can delay or falsify the analysis. If CLARK announces that a shock is needed while you are touching or moving the victim, stop CPR and leave the victim as still as possible to allow CLARK to reconfirm its analysis.

MOBILE TELEPHONES

CLARK can function properly even in close proximity to two-way radios and mobile phones. Normally, using a mobile phone near the victim should not affect CLARK's performance. However, it's recommended to not keep these devices closer than necessary to the victim and CLARK.

SKIN BURNS

Do not let the electrodes touch each other or any other electrode, power cable, bandaid, medical patch, etc.

Such contact can result in an electrical arc or skin burns from the shock, and can also divert the electric current away from the victim.

During a shock, the air bubbles between the skin and the electrode can cause skin burns. To avoid air bubbles, verify that the electrodes adhere completely to the skin. Do not use dried electrodes because they will not adhere properly to the skin.

ALLERGIC REACTION / CONTACT INJURY

CLARK and its accessories shouldn't cause an allergic reaction or skin irritation. However, if you notice skin irritation after touching CLARK or its accessories, please contact Lifeaz.

ELECTRICAL SHOCK RISK

The defibrillator delivers a high energy shock. Unless it is used exactly as described in this instructions for use, this electrical energy can cause serious injury or death.

Even though no former first aid training is required to use this device, it's recommended to familiarize yourself with these instructions for use and all associated instructions, precautions, connectors and accessories before using the defibrillator.

Do not insert fingers or anything other than the battery provided into the battery connector.

Do not disassemble the device. Risk of electric shock. Only LIFEAZ authorized personnel may disassemble the device.

The device has a non-automatic discharge that can only be activated by personnel authorized by LIFEAZ and is subject to a specific internal procedure.

CHILDREN

Keep CLARK out of reach from young children to avoid the potential risk of swallowing small parts or strangulation by the electrode cables.

POSSIBLE DEVICE SHUTDOWN

Once the battery LED starts to flash, there are only a few shocks available. It's recommended to keep an additional battery available.

See section 12 Separate accessories to learn more.

POSSIBLE DEVICE MALFUNCTION WARNING

No modification to this device is authorized.

DEFIBRILLATOR'S IMPACT ON OTHER DEVICES

Before delivering a shock, it's important to remove all other electrical medical equipment from the victim, including blood flowmeters, that are not protected from defibrillation.

DEVICE USE

Using the device in any way other than described in this instructions for use can damage the defibrillator or its accessories and render the device inoperative. The device will need to be replaced, which may include additional fees. Regularly check CLARK and its accessories for damage, by following the instructions.

MAINTENANCE

Neglecting maintenance for Clark can damage the device or cause it to function incorrectly. Follow the instructions in section 8 Maintenance.

Considering the infrequent use of the device, follow the maintenance instructions closely to ensure that your device is always functional and ready to be used at any time.

Do not execute any maintenance while the device is in use.

INFORMATION RESPONSIBILITY

Regardless of the conditions of use for this device, it is the client's responsibility to provide the information contained in this instructions for use, including the general safety information, to the appropriate person or persons who may use the defibrillator.

SERIOUS INCIDENT

In the event of a serious incident involving the device, the manufacturer should be notified as well as the competent authority of the member state in which the user is established.

5. SETTING UP YOUR CLARK AED

5.1. UPON RECEPTION OF YOUR CLARK AED

OPEN THE BOX AND VERIFY THAT ALL ELEMENTS LISTED IN SECTION 2 ARE INCLUDED :

- ▶ CLARK
- ▶ 1 (one) electrode pouch
- ▶ 1 (one) battery pack
- ▶ 1 (one) instructions for use
- ▶ 1 (one) first aid kit (containing a razor, scissors, compress, gloves)
- ▶ 1 (one) electrode instructions for use

5.2. SETTING UP YOUR CLARK AED

CONNECT THE ELECTRODES

Do not open the electrode pouch.

Do not pull on the cable once they've been attached.

Do not connect the electrodes if the cable, connector or pouch is damaged.

The expiration date of your electrodes is located on the pouch next to this symbol : 

The conservation period of CLARK electrodes is around 24 months.


Upon reception, check the expiration date of the electrodes.

Clark's battery is not rechargeable. Do not connect Clark to a main power source.


You'll receive your CLARK with the electrodes already connected.
If you ever need to reconnect them, follow these instructions :

Connect the electrodes as shown in the image.
The electrodes' cable should be facing the interior (to the right)

IMPORTANT : You do not need to apply excessive force in order to correctly connect the electrodes. If you can't easily connect the electrodes, the connector might be facing the wrong direction. Turn the connector and try again according to the visual. When you are looking at your CLARK AED, the cable should be facing towards the right.

PLACE THE ELECTRODES' CABLE IN THE FASTENINGS AS SHOWN IN THE IMAGE.

INSERT THE BATTERY


- ▶ Open the cover on the bottom of your CLARK with a cross-headed screwdriver.


LOOSEN THE FOUR SCREWS BEFORE REMOVING THE COVER. DO NOT REMOVE THE SCREWS :

The screws will stay attached to the cover even after being completely loosened.


- ▶ Remove the cover

- ▶ If your battery pack is one piece, connect the battery as shown in the image.

- ▶ If your battery pack is in two parts, connect them together and then connect the whole unit to the defibrillator as shown in the image.

The connector can only be plugged in one direction. If you're unable to insert it, verify that you are plugging it in the correct direction.

REATTACH THE COVER AND THE SCREWS

Reattach the cover and screw the four screws with a cross-head screwdriver.

You don't need to use excessive force when rescrewing the cover.

TEST YOUR CLARK

Turn on the device by pressing the green button.

The device should turn on and you should hear the instructions detailed in section 6 How to use your CLARK AED in an emergency.

If a problem arises, please contact us at clients@lifeaz.co.

ON / OFF

STAND BY MODE

Put Clark in stand by mode by pressing the green button for 2 (two) seconds.

5.3. CONFIGURING A WI-FI NETWORK

The Wi-Fi network configuration mode allows you to connect CLARK to a Wi-Fi network. Once configured, CLARK will prioritize Wi-Fi over mobile network for transmitting status reports.

Refer to section 8.4.2. *Wi-Fi connection when configuring a Wi-Fi network* for specifications and limitations related to Wi-Fi use.

Notes :

- ▶ Only CLARK devices with serial numbers beginning with CLA20, CLA21, CLA22, CLA23, CLA2408, CLA2413, or CLA2435 have a built-in Wi-Fi connectivity module. On other devices, Wi-Fi network configuration mode will not be triggered.
- ▶ CLARK can store up to 3 Wi-Fi networks and will connect to the first available network among the three. If you configure a fourth network, it will replace the oldest one stored.

For the configuration to be successful, you must be near the Wi-Fi network to which you want to connect CLARK.

1. SWITCHING TO WI-FI CONFIGURATION MODE

To start Wi-Fi configuration mode, press the local language selection button and the child mode button simultaneously for at least 5 seconds.


The status LED will light up orange continuously throughout the configuration phase. Clark will then play the voice message “Wi-Fi configuration mode activated”.

Note: From this point, you have 5 minutes to complete the connection, otherwise you will need to restart the configuration mode.

2. CONNECTING TO THE TEMPORARY WI-FI NETWORK PROVIDED BY CLARK

Once Wi-Fi configuration mode is activated, a new Wi-Fi network will become available. This network has the following characteristics :

- ▶ Access point name: **LIFEAZ-CLAXXXXXXXX(RY)**, where “CLAXXXXXXXX(RY)” corresponds to the serial number written on the label on the back of the product
- ▶ Network type: public

Connect to this network with a phone, PC, or tablet.

3. ACCESSING THE CONFIGURATION PAGE

Once connected to the LIFEAZ-CLAXXXXXXXX(RX) Wi-Fi network, your phone/PC/tablet will automatically redirect you to the configuration web page.

If this does not happen, you can manually go to a search engine and type in the following URL : **<http://wifi.lifeaz.fr/>**

4. ENTERING DATA ON THE CONFIGURATION PAGE

The access settings you enter (network name and password) must match those of the Wi-Fi network (internet box) to which you want to connect CLARK.

You must copy them into the following fields (as shown in the image):

- ▶ **Network Name** : Name of the Wi-Fi network
- ▶ **Pass phrase** : Wi-Fi network password
This is hidden. If there is an error or you are unsure, delete it completely before entering it again.
- ▶ **Device Name** : Do not fill in

Once you have entered the settings, click **Connect**.

Note: The **Refresh** button displays a list of Wi-Fi networks near CLARK.

A screenshot of a web browser displaying the 'Atmel Connect to Network' configuration page. The page has a blue header with the 'Atmel' logo and the title 'Connect to Network'. Below the header, there are four input fields: 'Network Name' with a red placeholder text 'Name of your Wi-Fi box', 'Pass phrase' with a red placeholder text 'Box password', and 'Device Name'. A blue 'Connect' button is located below the 'Pass phrase' field. At the bottom of the page, there is a 'Detect Device' section with a blue 'Refresh' button. The 'Connect' button is highlighted with a red rectangle.

5. CONFIRMING THE WI-FI CONFIGURATION

After confirming the settings, CLARK will stop emitting the temporary Wi-Fi network, and your phone/PC/tablet will automatically disconnect.

CLARK will then test the connection to the Wi-Fi network you have selected.

CLARK will then:

1. Broadcast the voice message "Connection successful" or the voice message "Connection failed".
2. Turn off the status LED, indicating that CLARK is back in standby mode.

If CLARK has broadcasted the message "Connection successful" :

- ▶ This means that the Wi-Fi network connection test was successful.
- ▶ CLARK will now use Wi-Fi connectivity as a priority for sending its status reports.

If CLARK has broadcasted the message "Connection failed":

- ▶ This means that the Wi-Fi network connection test has failed.
- ▶ CLARK will therefore not be able to connect to the Wi-Fi network.
- ▶ If the cause of the error is identified, or if you have any doubts, you can repeat the configuration procedure.
- ▶ If the problem persists, contact Lifeaz services.

The list of possible causes of configuration errors are:

- ▶ Inability to connect to the target Wi-Fi network due to router incompatibility or unsupported security settings.
- ▶ Failure to transmit Wi-Fi settings due to signal loss or incorrect login credentials.
- ▶ Inability to access the configuration page from certain browsers or phones/PCs/tablets.
- ▶ Use of a target Wi-Fi network using an unsupported security mode (e.g., WPA3 or captive portal).
- ▶ Interruption of the Wi-Fi configuration phase by pressing the ON button.

5.4. DEMONSTRATION MODE

Do not open the electrode pouch during demonstration mode. These are single-use only and are intended for use in a real cardiac arrest situation.

Thanks to the demonstration mode, you can familiarise yourself with all the voice and visual instructions given by Clark during use.

To start it, press the local language and foreign language selection buttons simultaneously for at least 5 seconds.

The status LED will light up red continuously throughout the demonstration.

At the end of the sequence, Clark will automatically turn off.

This mode is intended for demonstration purposes only and does not allow for heart rhythm analysis or the delivery of an electric shock.
Therefore, do not unpack the electrodes.

6. HOW TO USE YOUR CLARK AED **IN CASE OF EMERGENCY**

6.1. IN CASE OF EMERGENCY (A VICTIM IS UNCONSCIOUS AND NOT BREATHING)

1. CALL EMERGENCY MEDICAL SERVICES

If you haven't called emergency medical services yet, call them immediately or ask a bystander to call.

If you don't remember the number for emergency medical services, it is located on CLARK's interior after opening the lid.

2. OPEN CLARK'S LID

3. PRESS THE ON BUTTON

The device will turn on and vocal instructions will start.

ON / OFF

CAREFULLY FOLLOW THE VOCAL AND VISUAL INSTRUCTIONS.

4. THE LOCAL LANGUAGE IS SELECTED BY DEFAULT

To change the default language :

CLARK includes two languages : Your local language and French.

Press the language button that corresponds to your desired language.

You can switch to French and back to your local language at any moment (the vocal instructions won't be interrupted)

5. ADULT MODE IS SELECTED BY DEFAULT


Only if you want to change the Adult mode selected by default :

Choose the Child Mode :

If you use the device on a child (between 1-8 years old), select the Child Mode by pressing the button "Child Mode". You can switch back to Adult Mode by pressing the "Child Mode" button again.

Once you've selected the Child Mode, the Child Mode LED will light up.

The light will go off if you switch back to Adult Mode.

1-8 YEARS OLD ANS

IMPORTANT

If you're not sure of the victim's age, select the Adult Mode.

6. ATTACH THE ELECTRODES

Remove the victim's clothes that cover the chest. Scissors are included with your AED.

The victim's chest must be dry.
Compresses are included with your AED in the first aid kit.

If the victim's chest is hairy, you can use a razor to remove the hair from the areas where the electrodes will be placed. A razor is included with your AED in the first aid kit.


Open the electrode pouch

Place the electrodes on the victim's chest according to the visual located on the electrodes and on the defibrillator.

If you've improperly placed the electrodes, the AED may ask you to reposition them. Follow CLARK's instructions until they're in the correct position.

WARNING : The defibrillation electrodes must be kept apart from other electrodes or metal pieces in contact with the victim. Any other device that is not protected from defibrillation must be removed from the victim during defibrillation.

7. CARDIAC ANALYSIS AND SHOCK

Once the electrodes are detected by the defibrillator, you will be vocally instructed to no longer touch the victim so the device can analyze the victim's cardiac rhythm and decide if a shock is necessary.

In this case, the following LEDs will flash :

- ▶ First, the "Do not touch the victim" LED with vocal instructions to no longer touch the victim.
- ▶ Then, the "Shock delivered" LED when a shock is delivered. A shock will only be delivered if the AED's analysis recommends a shock.

A shock will only be delivered if the AED's analysis recommends a shock.

The AED may indicate that it's unable to analyze the cardiac rhythm due to electrical noise. In this case, follow the AED's instructions to ensure that the cardiac rhythm is correctly detected.

WARNING : avoid all contact with the victim's body such as exposed skin on the head or limbs, conductive liquids including gel, blood or other bodily fluids and metal objects such as a bed frame or stretcher that could divert the electrical current of defibrillation.

8. CARDIOPULMONARY RESUSCITATION (CPR)

If a shock is delivered or if no shock is recommended after an analysis, the AED will ask you to perform CPR.

The LED shown on the visual will flash to the rhythm needed for CPR. You will also hear a rhythmic beep. Perform compressions to the rhythm of the beep.

9. ANALYSIS, SHOCK, CPR

After two minutes of CPR, the AED will analyze the cardiac rhythm again. You will receive the same instructions as steps 7 and 8 : Do not touch the victim during the AED's analysis, and start CPR afterwards for two minutes.

IF, FOR ANY REASON, THE AED STOPS WORKING, CONTINUE CPR UNTIL EMERGENCY MEDICAL SERVICES ARRIVE.

DO NOT TURN OFF THE AED OR REMOVE THE ELECTRODES FROM THE VICTIM'S CHEST. EMERGENCY MEDICAL SERVICES WILL REMOVE THEM.

10. CLARK STAND BY MODE

To put Clark in standby mode, press the ON button for at least 2 (two) seconds.

ON / OFF

DO NOT TURN OFF THE AED UNTIL EMERGENCY MEDICAL SERVICES HAVE REMOVED THE ELECTRODES FROM THE VICTIM.

6.2. IF EMERGENCY MEDICAL SERVICES REQUEST THE DATA FROM YOUR CLARK

Once the device is used in an emergency, emergency medical services may ask to collect the data stored in your CLARK AED.

IN THIS CASE PLEASE CONTACT US AT CLIENTS@LIFEAZ.CO

7. AFTER USING CLARK AND AFTER THE ARRIVAL OF EMERGENCY MEDICAL SERVICES

CONTACT LIFEAZ

You will be asked about the state of your device to ensure it wasn't damaged during use.

You may also be asked to provide the data stored in your device after use.

Follow the instructions provided by the Lifeaz customer service.

REMOVE USED ELECTRODES AND REPLACE THEM WITH A NEW PAIR

Defibrillation electrodes are single-use accessories. They must be replaced after use.

See section 5.2. Configuration of your CLARK AED to learn more.

Do not discard used electrodes in the trash. Dispose of them as electronic waste in a local electronic waste disposal center.

Do not discard your CLARK in the trash or via an electronic waste disposal system. Contact Lifeaz to discard your CLARK.

The elements of the first aid kit are single-use accessories. They must be replaced after use.

CLEAN YOUR DEVICE

It is recommended to wear gloves when cleaning after use, to avoid any risk of contamination.

Using a clean cloth, wipe all accessible surfaces with a product such as rubbing alcohol (60-90% concentration).

Dry the device with a clean cloth if necessary.

8. MAINTENANCE

Your CLARK must be ready for emergency use at all times. Maintenance should only be performed when your CLARK is not in use.

8.1. REPLACING THE BATTERY

Once the battery is depleted, your device will indicate that it needs to be changed (LED flashes RED followed by a beep, as described in section 9 Troubleshooting). Follow the instructions in section 5.2 Configuration of your CLARK AED.

**CONTACT LIFEAZ TO RECEIVE A NEW BATTERY.
PLEASE CONTACT US AT CLIENTS@LIFEAZ.CO**

Discard your battery according to local regulations.
Clark's battery is not rechargeable.
Do not connect Clark to a main power source.

8.2. CHANGE THE ELECTRODES

Electrodes are a single-use accessory and must be changed after use or upon expiration.

If the electrodes are expired, your device will indicate the need to replace them (LED flashes RED followed by a beep, as detailed in section 9 Troubleshooting).

Contact Lifeaz to receive a new pack of electrodes upon expiration.

If the electrodes were used, you must order a new pair.
Please contact clients@lifeaz.co

To change electrodes, follow the process described in section 5.2. Configuration of your CLARK AED steps 1 and 2.

Discard your electrodes according to local regulations.

8.3. UPKEEP AND REGULAR MAINTENANCE

Your CLARK must be ready for emergency use at all times.

In addition to the maintenance previously described (electrode and battery replacements), you must follow the instructions listed below to ensure proper upkeep for your device.

Do not use your CLARK or its accessories in any way other than those described in this instructions for use.
For questions, please contact Lifeaz.

WHEN ?

After every use

WHAT DO I DO ?

Follow the procedure described in section 7: "AFTER USING CLARK AND AFTER THE ARRIVAL OF EMERGENCY MEDICAL SERVICES".

HOW ?

See section 7.

WHEN ?	WHAT DO I DO ?	HOW ?
Every week	Verify that your device is ready to be used in an emergency	Verify that the LED located on the front of Clark is flashing green
Every month	Verify that the environmental conditions are optimal for storing your device.	Verify that your device is stored in an appropriate environment : Temperature : 5 ° C - 35 ° C Humidity : 20 %- 80 % Atmospheric pressure : 700 hPa- 1060 hPa
Every year	Check that Clark can be switched on and does not emit an alert message	- Open the lid - Switch on the device and make sure that the voice instructions start up normally. Check that you don't hear any alert messages (see section 9)
Every two years	Visually check that your CLARK has no mechanical damage Manually test your CLARK	- Open the lid - Check that the plastic box isn't broken, chipped or scratched. - Check that all visual instructions are legible (not scratched or blurry) - Check that the electrode pouch is sealed - Check the expiration date on the electrode pouch - Begin the manual tests by pressing the ON button for 5 seconds

WHEN ?

WHAT DO I DO ?

HOW ?

- The device will switch to Interface test mode and vocal instructions will begin : you will be instructed to perform several actions to verify the audio system, buttons and LEDs.

- Follow the vocal instructions while performing the requested actions.

Note: At some point, you may be asked to press the English language button. If so, you will have to press the French language button instead. Similarly, if you are asked to press the French language button, you will have to press the English language button.

- If you do not hear vocal instructions, do not press any button. This means that the test is over.

If the manual tests fail, see section 9 Troubleshooting for further instructions.

Regularly / If needed

Your CLARK should be clean (not covered with dust...).
You should clean it regularly.

Use a clean wipe to remove dust.
You can also use a lightly damp wipe (with water).
Even if your device is protected from water infiltration, avoid spilling on your CLARK.

8.4. CLARK CONNECTION INTERFACES

Clark is not intended to be connected to an information technology network (IT network) to operate. However, it may be connected occasionally to exchange information electronically through non-LIFEAZ devices.

There are two connection interfaces: wireless (4G LTE-M and/or Wi-Fi) and USB interfaces.

- ▶ The wireless connectivity (4G LTE-M) IT-Network end-to-end has been validated and is managed by LIFEAZ. So it is not in the scope of section 14.13 of IEC 60601-1.
- ▶ USB creates a two-nodes IT network when Clark is connected to a PC through USB cable. Details of this connection are provided in section 8.4.1. *USB connection*.
- ▶ Configuring a Wi-Fi network creates two-nodes IT network when a phone, PC, or tablet is connected to the access point emitted by CLARK. Details of this connection are provided in section 8.4.2. *Wi-Fi connection when configuring a Wi-Fi network*.
- ▶ A successful configuration enables data communication via a Wi-Fi network (WPA2) validated by LIFEAZ.

8.4.1. USB CONNECTION

1.A. PURPOSE OF CLARK CONNECTION TO AN IT-NETWORK

The USB data interface may exchange information electronically with non-LIFEAZ devices. The purpose of the USB connection is :

- ▶ ECG retrieval
- ▶ Logs retrieval
- ▶ Embedded Software Upgrade

Your Lifeaz representative can provide you with computer software tools for data exchange by USB. They will provide you with all necessary information to install and use the tools.

Note: The tools are validated by LIFEAZ so the network is considered within the control of LIFEAZ.

1.B. REQUIRED CHARACTERISTICS OF THE IT-NETWORK INCORPORATING CLARK

The required characteristics for USB communication with Clark is:

- ▶ Computer with Windows 10 or above operating system

1.C. REQUIRED CONFIGURATION OF THE IT-NETWORK INCORPORATING CLARK

The required configuration for USB communication with Clark is:

- ▶ To have Administrator rights on your computer
- ▶ To have an internet connection (recommended for optimal performance)

1.D. TECHNICAL SPECIFICATIONS OF THE NETWORK CONNECTION OF CLARK INCLUDING SECURITY SPECIFICATIONS

Exchange of data is managed through software tools that embed an encrypted key for secured connection.

1.E. INTENDED INFORMATION FLOW BETWEEN CLARK, THE IT-NETWORK AND OTHER DEVICES ON THE IT-NETWORK, AND THE INTENDED ROUTING THROUGH THE IT-NETWORK

The information flow is the following:

- ▶ ECG data from Clark to PC
- ▶ Logs data from Clark to PC
- ▶ Embedded Software Upgrade from PC to Clark

1.F. LIST OF THE HAZARDOUS SITUATIONS RESULTING FROM A FAILURE OF THE IT-NETWORK TO PROVIDE THE CHARACTERISTICS REQUIRED TO MEET THE PURPOSE OF CLARK CONNECTION TO THE IT-NETWORK

There is no unacceptable risk identified in case of failure to establish data transfer or failure during data transfer. In particular, in case of interruption of Embedded Software upgrade, Clark will roll back to its current version.

Note: Connecting Clark device to a shared IT network with other devices, or subsequent changes to the shared IT network, can lead to previously unidentified risks for patients, operators, and third parties. These risks must be identified, analyzed, evaluated, and controlled before connecting Clark into the IT network.

Changes to the IT-NETWORK include:

- ▶ changes in the IT-network configuration;
- ▶ connection of additional items to the IT- NETWORK ;
- ▶ disconnecting items from the IT- NETWORK ;
- ▶ update of equipment connected to the IT- NETWORK ; and
- ▶ upgrade of equipment connected to the IT- NETWORK.

8.4.2. WI-FI CONNECTION WHEN CONFIGURING A WI-FI NETWORK

2.A. PURPOSE OF CLARK CONNECTION TO AN IT-NETWORK

The connection to the Wi-Fi network provided by CLARK is intended exclusively for the purpose of connecting the device to a target Wi-Fi network.

This network allows the user to transfer the access settings to the target Wi-Fi network.

Note: The Wi-Fi access point provided by CLARK and the associated configuration page are validated by LIFEAZ. The network is therefore considered to be under the control of LIFEAZ.

2.B. REQUIRED CHARACTERISTICS OF THE IT-NETWORK INCORPORATING CLARK

The required characteristics of the Wi-Fi network for Clark to connect to it are :

- ▶ private Wi-Fi network (requiring a password to access it)
- ▶ WPA2 type

2.C. REQUIRED CONFIGURATION OF THE IT-NETWORK FOR CONNECTION

The prerequisites for configuring the target Wi-Fi network with CLARK are :

- ▶ Access to the Wi-Fi network connection settings (know the network name and password of the target network)
- ▶ A phone, PC, or tablet that can connect to a Wi-Fi network

2.D. TECHNICAL SPECIFICATIONS OF THE NETWORK CONNECTION OF CLARK INCLUDING SECURITY SPECIFICATIONS

The technical specifications of the Wi-Fi connection provided by CLARK are as follows:

- ▶ Supported standards: IEEE 802.11 b/g/n (2.4 GHz band)
- ▶ Operating mode: temporary Wi-Fi access point provided by CLARK
- ▶ Access: reserved for the Wi-Fi configuration phase

Data exchanges are managed via the configuration page accessible through the temporary Wi-Fi network provided by CLARK.

No patient data is transmitted or stored via the temporary Wi-Fi network provided by CLARK.

2.E. INTENDED INFORMATION FLOW BETWEEN CLARK, THE IT-NETWORK AND OTHER DEVICES ON THE IT-NETWORK, AND THE INTENDED ROUTING THROUGH THE IT-NETWORK

The flow of information is as follows:

- ▶ Settings for accessing the target Wi-Fi network from the phone/PC/tablet to CLARK.

2.F. LIST OF THE HAZARDOUS SITUATIONS RESULTING FROM A FAILURE OF THE IT-NETWORK TO PROVIDE THE CHARACTERISTICS REQUIRED TO MEET THE PURPOSE OF CLARK CONNECTION TO THE IT-NETWORK

There is no unacceptable risk identified in the event of a computer network failure during the Wi-Fi network configuration phase.

The consequences of a failure are limited to the inability to configure Wi-Fi for CLARK, or a delay in doing so.

Note: Connecting Clark device to a shared IT network with other devices, or subsequent changes to the shared IT network, can lead to previously unidentified risks for patients, operators, and third parties. These risks must be identified, analyzed, evaluated, and controlled before connecting Clark into the IT network.

Changes to the IT-NETWORK include:

- ▶ changes in the IT-network configuration;
- ▶ connection of additional items to the IT- NETWORK ;
- ▶ disconnecting items from the IT- NETWORK ;
- ▶ update of equipment connected to the IT- NETWORK ; and
- ▶ upgrade of equipment connected to the IT- NETWORK.

9. TROUBLESHOOTING

Once the battery is connected, your CLARK will auto-test regularly (daily, weekly and monthly) and will notify you in case a problem is detected. The autotests include battery and electrode tests. It also autotests every time the electrodes and battery are connected. After every autotest, if a problem is detected, you will be alerted. Below you will find all possible alerts, their cause and the recommended solution.

9.1. ALARMS

The device monitors battery status (depending on voltage, temperature, last charge duration) on a daily basis. As soon as a low battery is detected, the AED takes the following actions:

CASE 1: LOW BATTERY

- ▶ When Clark is switched ON, a "battery low" voice message is broadcasted and a "low battery" visual alarm signal blinks (battery LED)

When the low battery indication is initially given, the device can still be used (3 shocks are possible).

CASE 2: BATTERY ALMOST DEPLETED

- ▶ When Clark is switched ON, a "battery low" voice message is broadcasted and a "low battery" visual alarm signal blinks (battery LED)
- ▶ The status LED flashes red and Clark emits 5 beeps every 30 minutes

In this case, the battery capacity may not be sufficient to deliver a shock and the device operation cannot be guaranteed.

In both cases, the user should:

- ▶ Replace the battery as soon as possible
- ▶ In case of emergency, start or continue CPR and seek another defibrillator if possible until emergency services arrive

The user should check the alarm by switching ON the device, as recommended for regular maintenance.

If the "battery low" voice message is broadcasted but the battery LED is not blinking, please contact LIFEAZ.

9.2. ALERTS

01

STATUS LED FLASHES

Red

BEEP

3 beeps every 6 hours

VOCAL INSTRUCTIONS

This signifies a major device defect.
After pressing the ON button, you will hear the following messages :
"Device out of order"
"Please contact Lifeaz services"
"Device malfunction"
"Start CPR"

WHAT DOES THIS MEAN ?

Your device is in need of maintenance. Clark is non-functional and should not be used. This may be caused by a software or hardware defect. If possible, CLARK will continue to give you instructions to perform CPR.

WHAT DO I DO ?

Contact Lifeaz as soon as possible.

**STATUS LED
FLASHES**

Red
+ Battery LED flashes

BEEP

5 beeps every 30
minutes

**VOCAL
INSTRUCTIONS**

After pressing the
ON button, you will
hear the following
message :
"Battery low"

WHAT DOES THIS MEAN ?

Depleted main battery
Or internal battery.

WHAT DO I DO ?

If the battery is depleted : replace it with a new battery (process described in section 5.2. Configuration of your Clark AED). If you haven't received the replacement, see section 12 Separate accessories to order a new one.

If the problem persists after replacement, this could signify that the internal battery is low, which you can't change yourself.

Contact Lifeaz to replace the device.

**STATUS LED
FLASHES**

Red

BEEP

5 beeps every 30
minutes

**VOCAL
INSTRUCTIONS**

After pressing the
ON button, you will
hear the following
message :
"Electrodes expired"

WHAT DOES THIS MEAN ?

The electrodes are expired.

WHAT DO I DO ?

Replace the electrodes with new ones, as described in section 8 Maintenance. If you haven't received your replacement pair, see section 12 Separate accessories to order new electrodes.

**STATUS LED
FLASHES**

Red

BEEP

5 beeps every 30
minutes

**VOCAL
INSTRUCTIONS**

After pressing the ON
button, you will not
hear an alert message.
The voice instructions
will begin with:
"Call the emergency
services"

WHAT DOES THIS MEAN ?

Capacitor charging problem or battery problem

WHAT DO I DO ?

Contact Lifeaz

**STATUS LED
FLASHES**

Orange

BEEP

No vocal message

**VOCAL
INSTRUCTIONS**

After pressing the ON
button, you will hear
the following message:
"Plug the electrodes
connector"

WHAT DOES THIS MEAN ?

The electrodes are not connected.

WHAT DO I DO ?

- Connect the electrodes.
- Verify that the electrodes are properly connected (unplug/replug if needed)

06**STATUS LED
FLASHES**

Orange

BEEP

No vocal message

**VOCAL
INSTRUCTIONS**

After pressing the ON button, you will not hear an alert message. The voice instructions will begin with: "Call the emergency services"

WHAT DOES THIS MEAN ?

The wrong electrodes are connected.

WHAT DO I DO ?

Verify that the electrodes you connected are Lifeaz electrodes.

07**STATUS LED
FLASHES**

Red

BEEP

No vocal message

**VOCAL
INSTRUCTIONS**

After pressing the ON button, you will hear the following message : "Therapy disabled. Device error."

WHAT DOES THIS MEAN ?

Hardware problem

WHAT DO I DO ?

Contact Lifeaz.

**STATUS LED
FLASHES**

Orange

BEEP

5 beeps every 30
minutes

**VOCAL
INSTRUCTIONS**

After pressing the ON
button, you will not
hear an alert message.
The voice instructions
will begin with:
"Call the emergency
services"

WHAT DOES THIS MEAN ?

Hardware problem

WHAT DO I DO ?

Contact Lifeaz.

9.3. OTHER PROBLEMS

WHAT'S THE PROBLEM ?

The electrodes are expired but I haven't received the replacement electrodes.

I received an alert from the device but when I turn it on, CLARK doesn't describe the problem.

The device doesn't turn on when I connect the battery

I received new electrodes but can't connect them.

The device no longer works / you're no longer subscribed to CLARK services

WHAT CAN I DO ?

Contact Lifeaz at clients@lifeaz.co.

Take out the battery and reconnect it. If the problem isn't resolved, contact Lifeaz

Take out the battery and reconnect it. If the problem isn't resolved, contact Lifeaz

Make sure you connected the electrodes correctly, as described in the section 5.2 Configuration of your CLARK AED.

Contact Lifeaz. You will receive instructions to return your device.

Do not discard your CLARK in the trash or electronic waste system.
Return your device to Lifeaz. It will be recycled properly.

10. TECHNICAL INFORMATION

10.1. CLARK TECHNICAL FEATURES

PRODUCT NAME	CLARK	
DIMENSIONS	Length : 210 mm Width : 210 mm Height : 75 mm	
WEIGHT	1,3 kg	
PROTECTION WARRANTY (FROM LIQUIDS AND SOLIDS)	IP54 5 = Protection from dust and other microscopic particles (protection from solids) 4 = Protection from water splashing (protection from water infiltration)	
ENVIRONMENTAL CONDITIONS	Standby (nominal environment)	Temperature : 5 °C-35 °C (conditioned by the electrodes) Relative humidity : 20 %-80 % Pressure : 700 hPa-1060 hPa
	Operating (extreme environment)	For immediate use (maximum 2 hours) Temperature : 0 °C-40 °C Relative humidity : 5 %-95 % Pressure : 700 hPa-1060 hPa
	Transport and storage	Temperature : 5 °C-35 °C Relative humidity : 20 %-80 % Pressure : 700 hPa-1060 hPa

BATTERY

Battery type	Lithium - Manganese dioxide (Li/MnO ₂), 18 V
Capacity	New and fully charged : maximum number of shocks available at 20 ° C : 50 (estimated from design evaluation)
Conservation period (before insertion)	4 years
Standby life (after insertion)	5 years (standby life may be reduced in areas with poor 4G LTE-M or 2G network coverage and may be affected by product storage conditions)
Once the device is turned on	At least 1 hour with a new battery (estimated from design evaluation)
Shocks remaining after the battery LED flashes	Once the battery light turns on, 3 shocks are left. However, it is recommended to replace the battery once the battery LED turns on.

WAVEFORM

Biphasic truncated exponential with fixed energy (determined by the victim's impedance)

ELECTRODES

Surface contact

Disposable defibrillation electrode adhesives with a surface of 144 cm² each.

Expiration date

The electrodes' expiration date is indicated directly on their pouch next to the following symbol : 

The conservation period for CLARK electrodes is about 24 months.

Check the expiration date of your electrodes upon reception. Do not reuse the electrodes.

The electrodes are Applied Parts
(according to the IEC 60601-1 norm)

SHOCK CHARACTERISTICS

The maximum amount of time between the start of the cardiac rhythm analysis, with an ECG, and shock availability

<25 sec with a new battery

<35 sec with a battery after 6 shocks at 150J

The maximum amount of time between turning on the device and shock availability

<35 sec with a new battery

<45 sec with a battery after 6 shocks at 150J

Charge time of a capacitor completely discharged when CLARK is formatted with maximal energy

<8 seconds with a new battery

<8 seconds with a battery after 6 shocks at 150J

Shock cycle time

Following the recommendations of the European Resuscitation Council (ERC) 2021 : 2 minutes between two consecutive shocks.

«Shock delivered» Indicator

“Do not touch the victim” LED flashes and a vocal message

DISARM

Once charged, the defibrillator will disarm if :

- ▶ the victim's cardiac rhythm is identified as non-shockable,
- ▶ when the device is in stand by mode, by pressing on the ON button for at least two (2) seconds,
- ▶ the adhesive electrodes are removed from the victim or the electrodes are disconnected from CLARK
- ▶ the battery is depleted or removed from the device

ARTIFACT DETECTION

If electrical “noise” (artifact) is detected and interferes with the cardiac rhythm analysis, the analysis will be delayed until a clear ECG signal can be detected.

ANALYSIS PROTOCOL

Depending on the analysis's results, either a shock will be delivered or time will be allotted to perform CPR.

PROTECTION FROM DEFIBRILLATION

Your device is type BF. Meaning that the device is protected from defibrillation.

The recuperation time after your CLARK delivers a defibrillation shock is less than 1 second.

EXPECTED SERVICE LIFE

Your device has an expected service life of 5 years minimum (excluding consumables)


10.2. CLARK THERAPY CHARACTERISTICS

10.2.1. WAVEFORM

The shock delivered by your CLARK AED is a waveform Biphasic Truncated Exponential (BTE).

BTE is a state-of-the-art shock currently used in most available AEDs.

Defibrillation Waveform

Example of biphasic truncated exponential waveform (BTE)

The waveform is modified depending on the defibrillation impedance of the victim. Your CLARK delivers a shock if the impedance measured is between 25 and 175 OHM (Ω)

WAVEFORM SPECIFICATIONS :

Shock characteristics for an ADULT are the following :

VICTIM IMPEDANCE (OHM)	TOTAL PHASE DURATION (MS)	PHASE 1 DURATION (MS)	PHASE 2 DURATION (MS)	ENERGY DELIVERED (J)
25	6	3	3	140
50	9	4,5	4,5	151
75	10,8	6,5	4,3	152
100	14,1	8,5	5,6	155
125	17,5	10,5	7	157
150	20	12	8	157
175	20	12	8	152

The intermediate phase between phase 1 and phase 2 is set at 0,5 ms.

Shock characteristics for a CHILD are the following (if the child mode is selected) :

VICTIM IMPEDANCE (OHM)	TOTAL PHASE DURATION (MS)	PHASE 1 DURATION (MS)	PHASE 2 DURATION (MS)	ENERGY DELIVERED (J)
25	6	3	3	48
50	9	4,5	4,5	52
75	10,8	6,5	4,3	52
100	14,1	8,5	5,6	53
125	15	9	6	52
150	15	9	6	50
175	15	9	6	47

The intermediate phase between phase 1 and phase 2 is set at 0,5 ms.

10.2.2. ECG ANALYSIS

An AED detects and analyzes the heart's electrical signal (ECG) to determine if a shock is necessary or not.

The ECG is categorized as shockable or non-shockable.

SHOCKABLE RHYTHMS

Unless defibrillation is delivered rapidly, certain cardiac rhythms can be fatal. These rhythms include ventricular fibrillation (VF) and rapid ventricular tachycardia (VT). An unconscious victim without a pulse is always VF or almost always VT. CLARK uses multiple methods to determine if a rhythm is shockable.

REMARK : Due to security reasons for the victim, some low amplitude or low frequency rhythms will not be interpreted as VF shockable rhythms. Additionnally, certain VT rhythms generally associated with circulation will not be interpreted as shockable rhythms.

NON-SHOCKABLE RHYTHMS

Non-shockable rhythms are rhythms other than shockable rhythms

These are benign rhythms (or even normal) that should not be shocked, especially on victims with a pulse, because there is no possible advantage and may even result in a deterioration of the rhythm.

Non-shockable rhythms include the normal sinus rhythm, supraventricular tachycardia, sinus bradycardia, fibrillation and atrial flutter, heart block, idioventricular rhythms, premature ventricular contractions and other rhythms accompanied by a detectable pulse and/or a conscious victim.

In case of incorrect device use, asystole is included in this category to maximize security.

Note : The American Heart Association (AHA) established a guide that defines shockable and non-shockable rhythms. This guide equally defines the recommended sensitivity and specificity.

Guideline: "Automatic External Defibrillators for Public Access Defibrillation: Recommendations for Specifying and Reporting Arrhythmia Analysis Algorithm Performance, Incorporating New Waveforms, and Enhancing Safety, 18 march 1997"

DEVICE PERFORMANCE

Defibrillator performance is defined by two characteristics : sensitivity and specificity.

Sensitivity describes the device's ability to correctly identify shockable rhythms.

Specificity describes the device's ability to correctly identify non-shockable rhythms.

They are calculated as the following :

	RAPID VF AND VT	ALL OTHER ECG RHYTHMS
SHOCK	A	B
NO SHOCK	C	D

A true positive (A) is a correct classification of a shockable rhythm.

A true negative (D) is a correct classification of all non-shockable rhythms.

A false positive (B) is a rhythm or asystolie that is incorrectly classified as a shockable rhythm.

A false negative (C) is a VF or VT associated with cardiac arrest that is incorrectly classified as non-shockable.

The device's sensitivity to shockable rhythms is $A / (A + C)$. The true predictive value is expressed as $A / (A + B)$. The device's specificity for non-shockable rhythms is $D / (B + D)$. The rate of false positives is expressed as $B / (B+D)$.

The table below shows CLARK's performance :

RHYTHM TYPE	PERFORMANCE RESULT
Sensitivity	
Ventricular fibrillation and rapid Ventricular tachycardia (VT)	>90%
Specificity	
Normal sinus rhythm	>99%
Asystole	>95%
Other non shockable rhythms (atrial fibrillation, PVC)	>95%

The databases used to evaluate the algorithm are the following :

MITDB MIT-BIH database / CUDB database, AHA arrhythmia database.

RECORDING METHODS : these databases contain ECG recordings (Holter) acquired using a frequency of at least 250 HZ and a range of 0,05 to 125 HZ of useful data.

RHYTHM SELECTION CRITERIA : rhythms are selected depending on their classification (VF, non-shockable rhythms, asystole)

ANNOTATION METHODS AND CRITERIA : rhythm annotation is included in the evaluation database. Every ECG was categorized by arrhythmia type and entered into the algorithm. The algorithm annotates each 3-second ECG signal as shockable or non-shockable.

GLOBAL ALGORITHM RESULTS (3 SECOND SAMPLES)

	RAPID VF AND VT	ALL OTHER ECG RHYTHMS
SHOCK	A = 5299	B = 136
NO SHOCK	C = 291	D = 39301

GLOBAL SENSITIVITY $A/(A+C)$: 94.79%

GLOBAL SPECIFICITY $D/(B+D)$: 99.65%

Note :

True predictive value $A/(A+B)$: 97.50%

False positive rate $B/(B+D)$: 0.3%

Positive predictive precision $(A+D)/(A+B+C+D)$: 99.05%

10.3. ELECTROMAGNETIC COMPATIBILITY

The CLARK defibrillator must be used in an electromagnetically compatible environment specified in the tables below.

The client or person using CLARK must ensure they are in a suitable environment.

10.3.1. ELECTROMAGNETIC EMISSIONS

EMISSIONS TEST

RF
CISPR

CONFORMITY

Group I
Class B

ENVIRONMENT ELECTROMAGNETIC ADVICE

Clark's RF emissions are very weak and should not provoke interferences with electronic equipment nearby.
CLARK can be used in all establishments, including domestic establishment, and those directly connected to a low voltage public electrical power supply.

10.3.2. ELECTROMAGNETIC IMMUNITY

IMMUNITY TEST

Electrostatic
discharge (ESD)
IEC 61000-4-2

TEST LEVEL IEC60601

± 6 kV contact
 ± 8 kV air

CONFORMITY LEVEL

± 6 kV contact
 ± 8 kV air

ELECTROMAGNETIC ENVIRONMENT ADVICE

There is no particular recommendation regarding electrostatic discharge

Rapid and
transitory electrical
surges IEC 61000-
4-4

± 2 kV for electrical
source lines
 ± 1 kV for entry /
exit lines

Not applicable

There is no power
supply

IMMUNITY TEST	TEST LEVEL IEC60601	CONFORMITY LEVEL	ELECTROMAGNETIC ENVIRONMENT ADVICE
Powersurge IEC 61000-4-5	± 2 kV for electrical source lines ± 1 kV for entry / exit lines	Not applicable	There is no power supply
Voltage dip, brief outages and tension variations on power source lines.	Not applicable	Not applicable	There is no power supply
Industrial frequency magnetic fields	3 A/m	3 A/m	Industrial frequency magnetic fields must be held to the same environmental standards as commercial businesses / hospitals. There is no particular recommendation for non commercial / non hospital environments.

IMMUNITY TEST	TEST LEVEL IEC60601	CONFORMITY LEVEL	ELECTROMAGNETIC ENVIRONMENT ADVICE
---------------	------------------------	---------------------	---------------------------------------

RADIO FREQUENCY (RF)

Radiated RF electromagnetic fields RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	20 V/m	Mobile RF communication equipment should not be used next to CLARK, including cables. Interferences can occur next to equipment with the following symbol :
--	--------------------------------	--------	---

Conducted disturbances, induced by RF fields.	Not applicable	Not applicable	There is no power supply
---	----------------	----------------	--------------------------

11. CONNECTIVITY

CLARK is equipped with connectivity capabilities.

11.1. MOBILE CONNECTIVITY

Your CLARK includes a mobile connectivity module. Thanks to this feature, information regarding your device's current state are automatically sent to lifeaz, who will contact you in case a problem is detected.

FREQUENCY RANGE - CLARK with a serial number that begins with CLA20, CLA21 or CLA22

2G	EGSM900 (Class4-2W)	Tx=800MHz to 915MHz Rx=925MHz to 960MHz
	DCS1800 (Class1-1w)	Tx=1710MHz to 1785MHz Rx=1805MHz to 1880MHz

FREQUENCY RANGE - CLARK with a serial number that does not begin with CLA20, CLA21 or CLA22

4G	LTE CAT M1 (20dBm $\pm$ 2.7dB)	B8 (FDD): Tx=880MHz to 915MHz Rx=925MHz to 960MHz B20 (FDD): Tx=832MHz to 862MHz Rx=791MHz to 821MHz B28 (FDD): Tx=703MHz to 748MHz Rx=758MHz to 803MHz	ETSI EN301908 EN 62311
2G	GSM900 (Class-4 33dBm)	Tx: 880MHz to 915MHz Rx: 925MHz to 960MHz	ETSI EN301511 EN 62311

11.2. WI-FI CONNECTIVITY

Only CLARK devices with a serial number that begins with CLA20, CLA21, CLA22, CLA23, CLA2408, CLA2413 or CLA2435 include a WIFI connectivity module.

IEEE 802.11 b/g/n

Frequency range : 2400 MHz to 2483.5 MHz

Maximum transmitter power output : 19.87dBm

11.3. GPS

Only CLARK devices with a serial number that begins with CLA20, CLA21 or CLA22 include a GPS module operating in receiver mode in the following frequency bands:

- ▶ L1: 1559 MHz to 1610 MHz
- ▶ L2: 1215 MHz to 1300 MHz
- ▶ L5: 1164 MHz to 1215 MHz

The use of this feature is not available yet. You won't have to take any action to make it work.

12. SEPARATE ACCESSORIES

If you would like to order additional accessories or if you haven't received an accessory replacement, **please contact us at clients@lifeaz.co**.

REFERENCE	ITEM
LFZ00232	Electrodes
LFZ00303	Battery
LFZ00414	Clark Instructions for use
LFZ00512	First aid kit
LFZ00239	Electrodes instructions for use

Other non-medical accessories can be ordered upon request by writing to clients@lifeaz.co

Lifeaz SAS
25 Rue de Tolbiac
75013 Paris
France

Email: clients@lifeaz.co
www.lifeaz.co

REF

NAME

LFZ00499

AED CLARK 4G UK PRO

Instructions for use history

Instructions for use revision : LAB126K

Publication date : 2026-02

To download additional copies of this instructions for use, please send a mail to clients@lifeaz.co You can also download the latest applicable version of this manual from Lifeaz website: lifeaz.co.

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