


$GeneSafe^{m}$ the evolution of NIPT

A non-invasive prenatal test that screens multiple genes for mutations causing severe genetic disorders in the fetus





analyses circulating cellfree fetal DNA (cfDNA) from a maternal blood sample.

The test is performed after **10 weeks** of pregnancy.

GeneSafe works as a **complementary screen** to traditional and genomewide NIPT PrenatalSafe. It screens for several life-altering genetic disorders that are not screened with current NIPT technology, allowing a **complete picture** of the risk of a pregnancy being affected by a genetic disorder.



GeneSafe facilitates early diagnosis of single-gene disorders single-gene disorders.

It involves 3 different levels of screening:



This test screens for 5 common inherited recessive genetic disorders, such as Cystic Fibrosis, Beta-Thalassemia, Sickle cell anaemia, Deafness autosomal recessive type 1A, Deafness autosomal recessive type 1B.

Genes screened: CFTR, CX26 (GJB2), CX30 (GJB6), HBB



This test screens for 44 severe genetic disorders due to de novo mutations (a gene mutation that is not inherited) in 25 genes

Genes screened: ASXL1, BRAF, CBL, CHD7, COL1A1, COL1A2, COL2A1, FGFR2, FGFR3, HDAC8, JAG1, KRAS, MAP2K1, MAP2K2, MECP2, NIPBL, NRAS, NSD1, PTPN11, RAF1, RIT1, SETBP1, SHOC2, SIX3, SOS1



This test screens for both inherited and de novo single-gene disorders and represents a combination of the tests GeneSafe INHERITED and GeneSafe providing a complete picture of the pregnancy risk.



GENE	GENETIC DISORDER
CFTR	Cystic Fibrosis
CX26 (GJB2)	Deafness autosomal recessive type 1A
CX30 (GJB6)	Deafness autosomal recessive type 1B
НВВ	Beta-Thalassemia
НВВ	Sickle cell anemia

The inherited recessive disorders screened by GeneSafe INHERITED are the most common in the European population



identifies fetal conditions that could be **missed by traditional prenatal screening.**


GENE	SYNDROMIC DISORDERS
JAG1	Alagille syndrome
CHD7	CHARGE syndrome
HDAC8	Cornelia de Lange syndrome 5
NIPBL	Cornelia de Lange syndrome 1
MECP2	Rett syndrome
NSD1	Sotos syndrome 1
ASXL1	Bohring-Opitz syndrome
SETBP1	Schinzel-Giedion syndrome
SIX3	Holoprosencephaly
	NOONAN SYNDROMES
BRAF	Cardiofaciocutaneous syndrome 1
CBL	Noonan syndrome-like disorder with or without juvenile myelomonocytic leukemia (NSLL)
KRAS	Noonan syndrome/cancers
MAP2K1	Cardiofaciocutaneous syndrome 3
MAP2K2	Cardiofaciocutaneous syndrome 4
NRAS	Noonan syndrome 6/cancers
PTPN11	Noonan syndrome 1/ LEOPARD syndrome/cancers
PTPN11	Juvenile myelomonocytic leukemia (JMML)
RAF1	Noonan syndrome 5/LEOPARD syndrome 2
RIT1	Noonan syndrome 8
KITT	
SHOC2	Noonan syndrome-like disorder with loose anagen hair

GENE	SKELETAL DISORDERS
COL2A1	Achondrogenesis, type II or hypochondrogenesis
FGFR3	Achondroplasia
	CATSHL syndrome
	Crouzon syndrome with acanthosis nigricans
	Hypochondroplasia
	Muenke syndrome
	Thanatophoric dysplasia, type I
	Thanatophoric dysplasia, type II
	Ehlers-Danlos syndrome, classic
	Ehlers-Danlos syndrome, type VIIA
COL1A1	Osteogenesis imperfecta, type I
COLIAI	Osteogenesis imperfecta, type II
	Osteogenesis imperfecta, type III
	Osteogenesis imperfecta, type IV
COL1A2	Ehlers-Danlos syndrome, cardiac valvular form
	Ehlers-Danlos syndrome, type VIIB
	Osteogenesis imperfecta, type II
	Osteogenesis imperfecta, type III
	Osteogenesis imperfecta, type IV
	CRANIOSYNOSTOSIS SYNDROMES
FGFR2	Antley-Bixler syndrome without genital anomalies or disordered steroidogenesis
	Apert syndrome
	Crouzon syndrome
	Jackson-Weiss syndrome
	Pfeiffer syndrome type 1
	Pfeiffer syndrome type 2
	Pfeiffer syndrome type 3

GeneSafe detects de novo mutations in 25 genes causing 44 different genetic disorders. The genetic conditions screened by this innovative test are often present in the fetus in the absence of a family history of the condition. This is a paradigm shift in prenatal screening. GeneSafe screens for de novo mutations that cannot be detected by standard carrier screening, as these mutations are not present in parents' somatic cells. The genetic disorders screened by GeneSafe can cause skeletal dysplasias, cardiac defects¹⁻²⁻³, multiple congenital anomalies⁴⁻⁵, autism⁶, epilepsy⁷, and/or intellectual disability⁸⁻⁹.



GeneSafe screens for conditions common across all maternal ages



All pregnant women – **regardless of age** – are at equal risk of the genetic conditions screened by GeneSafe. Although the occurrence of each disorder is relatively rare, the cumulative rate of occurrence of these conditions (~1 in 600) is similar to that of Down Syndrome, in younger women¹⁰.

^{10.} McRae J, et al. Prevalence and architecture of de novo mutations in developmental disorders. Nature 542, 433-438



can identify conditions that may have otherwise gone undetected until after birth or into childhood



Many disorders screened with GeneSafe are not typically associated with abnormal prenatal ultrasound findings (especially in the first trimester), or may not be evident until late second/ third trimester, when confirmatory invasive testing can pose a risk of preterm birth, or baby's health after delivery.



GeneSafe[™] screens for genetic disorders associated with advanced paternal age



While traditional NIPT screens for conditions typically associated with advanced maternal age (e.g. Down Syndrome), GeneSafe DE NOVO screens for genetic disorders (e.g. Achondroplasia, Pfeiffer syndrome, Crouzen syndrome, Apert syndrome, Thanatophoric dysplasia, Osteogenesis Imperfecta, etc.), which are associated with advanced paternal age (men who are >40 years old)11 ensuring a comprehensive screen for couple of advanced age. These disorders typically are caused by **mutations** arising during spermatogenesis. As a man ages, the chance for these errors to occur substantially increases.





POSITIVE

Pathogenic / Likely Pathogenic mutation(s) detected

This result shows that the test detected one or more mutations in one or more genes.

A patient with a positive *GeneSafe* test result should be referred for genetic counselling and should always be followed-up with an invasive diagnostic test for confirmation of test results, before any medical decisions are made.

Only **known pathogenic** and **likely pathogenic** mutations are reported.




NO Pathogenic / Likely Pathogenic mutation(s) detected

This result shows the test has not detected any disease causing mutation in the targeted genes screened.

Negative screening results mean that there is a **very low risk** that the fetus has one of the disorders screened with GeneSafe although no guarantee may be given that the fetus is actually healthy.



GeneSafe[™] a groundbreaking technology allowing for a genetic analysis that is revolutionary



Advanced bioinformatic analysis to identify de novo mutations







High Resolution: coverage >550X



Exceptional performance:

Sensitivity and specificity >99%



Accurate measurement of Fetal Fraction (FF)



Low limit of detection: highly accurate at low cfDNA quantity (FF> 2%)



Low incidence of inconclusive results (<1%)





SIMPLE: a simple **blood sample (8-10 ml)** collected at **10+ weeks** of gestation is required



SAFE: it is a **non-invasive** test, **no risk** to the fetus and the mother



RELIABLE:

Sensitivity and specificity >99%



FAST:

Turnaround time of 10 days





Is intended for patients who meet any of the following criteria:

- · Advanced paternal age (men who are > 40 years old)
- Abnormal ultrasound finding(s) suggestive of monogenic disorder
- · Patients wishing to avoid an invasive diagnostic procedure
- · Patients at risk for genetic conditions screened

The test is suitable for:

- · both single and twin pregnancies.
- patients whose pregnancies have been achieved by IVF techniques, including pregnancies with egg donation or surrogacy.

SeneSafe[™] 5 easy steps



Order the test



Fill in and Sign the Test Requisition and Consent Form.



Draw a maternal blood sample (10 ml) at 10+ weeks of pregnancy



Call DHL on 0844 248 0844 to arrange collection of the sample, that is it!



Receive results in just 10 days



GeneSafe[™] includes

Free follow-up of abnormal results

Free CVS or Amniocentesis

in collaboration with reference gynecologists

Reimbursement of the test fee

for cases with inconclusive test results

Genetic Counselling for Positve Results



advanced molecular diagnostics solutions using state-of-the art technologies



Test performed in Italy (Rome or MIlan)



Fast TAT: 10 days



20 years experience in prenatal molecular diagnostics



Personalised genetic counselling with genetic counselors experts in discussing genetic test results and familial risks.



Fully **Automated work-flow,** from **cfDNA extraction** to data analysis



Test available worldwide



Over 200,000 genetic tests/year



Dedicated R&D team

Numerous peer-reviewed papers published in renowned international journals



Eurofins Clinical Diagnostics, 17 Doman Road, Camberley, GU15 3DF Laboratori e Studi Medici , Via Castel Giubileo, 11 - 00138 Roma Laboratori e Studi Medici, Via Enrico Cialdini, 16 - 20161 Milano Tel.: +39 06 8811270 | Fax: +39 0664492025

www.laboratoriogenoma.eu www.prenatalsafe.co.uk www.prenatalsafekaryo.co.uk www.genesafe.co.uk

Next Step

Please call or email to find out more UK toll free number: 0808 1691022 Email: NIPT@eurofins.co.uk

