



Synthetic Genetic Elements for Improving Recombinant Protein Expression

syngensys

www.syngensys.com
made in Sheffield



Renata Caikauskaitė, Cristina Alexandru-Crivac, Ryan Taylor,
Molly Smith, Kate Fewkes, Andy Racher, David James
Email: andy.racher@syngensys.com
SynGenSys Ltd, Sheffield, S1 2JE, UK

Abstract

High productivity cell lines are essential for improving the availability of mAb drugs. Recombinant gene transcription ultimately controls the recombinant protein production rate by mammalian cells. Synthetic promoters can significantly increase productivity and facilitate an optimal stoichiometry of recombinant polypeptide co-expression. We describe improved mAb expression in GS-KO CHO cells by using novel transcriptional hubs that utilise an optimal combination of synthetic promoters to further enhance expression without adversely affecting growth

SynGenSys's Third Generation Synthetic Promoters


- Rationally designed using rules formulated by SynGenSys
- Not variations of hCMV-MIE promoter
- Access more of CHO cell transcriptional landscape than hCMV-MIE promoter
- High level of sequence diversity
- Higher productivity, compared to hCMV-MIE, seen for a range of products when transcription driven by a CHO.SET[®] promoter

- Transcriptional Hub[™] created by combining CHO.SET[®] promoters with synthetic selectable marker promoters (SMPs) engineered to minimise overlap in TFRE species
- Productivity using Transcriptional Hub[™] higher than either CHO.SET[®] promoters or SMPs in isolation or standard hCMV-MIE/SV40E combination
- Transcriptional Hub[™] concept offers approach to increasing selection stringency without the downside of using high concentrations of selective agents

Introduction

The need to overcome limitations in patient access to biologics drives efforts to enhance productivity and reduce costs. In the bioreactor, product concentration is a function of cell specific production rate (Qp), mean viable cell concentration, and process duration. To date, most efforts have focussed upon improving the space-time yield of viable biomass with relatively few efforts focussing upon improving Qp. An assumption is that the widely used hCMV-MIE promoter, used to drive transgene transcription, cannot be improved upon. SynGenSys has formulated rules for designing mammalian promoters can be formulated and used them to build better promoters. These promoters are not variants or combinations of naturally occurring viral or endogenous promoters but true *de novo* designs leading to better Qp values compared to hCMV-MIE. Coupling CHO.SET[®] promoters with SMPs to minimise promoter interference increases productivity further.


1. The problem with standard promoter set-ups




- Low selectable marker (SM) gene expression is potentially growth limiting
- CMV competes with SV40E promoter, a strong promoter in its own right, for transcription factors causing substantial reduction in GS gene transcription
- Fixed transgene and selectable marker expression ratios when using standard CMV-CMV-SV40 set-up
- The need for high levels of selective drugs can cause selection for unwanted phenotypes

2. Transcriptional Hubs[™]: a route to higher productivities

A standard approach to obtaining more productive cell lines is to use higher concentrations of selective drugs to increase selection stringency during CLC. This risks selection for unwanted phenotypes e.g. reduced growth rate or altered cellular redox potential. Increased expression of the selectable marker gene risks damping transgene transcription as both promoters compete for the same transcription factors (TFs). Transcriptional hubs, artificial assemblies of genes whose transcription is driven by distinct promoters, create the possibility of stringent selection using synthetic promoters designed with minimal overlap of TFREs and no increased levels of selective drugs.




3. Transcriptional Hub[™]: faster recovery following transfection



- Cells transfected with Transcriptional Hub[™] vector recover faster compared to standard CMV/SV40 vector
- Cell concentration at day 10 0.5-0.7 x 10⁶ viable cells/mL for cultures using Transcriptional Hub[™] compared to 0.3 x 10⁶ viable cells/mL for CMV/SV40 vector, with culture viabilities of about 70% compared to 59% for CMV/SV40
- For Transcriptional Hub[™] vectors, cell concentration achieves 10⁹/mL between days 13-16 compared to 19-21 for CMV/SV40 vector
- Similar results seen for CHO GS-KO hosts obtained from two different vendors

4. Transcriptional Hubs[™]: higher productivities



5. Transcriptional Hubs[™]: answering the productivity question

