



900 N. 400 W. Ste. 11
North Salt Lake, UT 84054

2000 SERIES OPERATOR'S MANUAL



Produced by General Water Technologies

Thank you to our loyal customers!

Contact: 888-586-5654 M-F 8:00 AM – 5:00 PM MST (message prompt after hours for emergency response 24/7)

Service: service@genwater.com

Sales: sales@genwater.com

Orders: orders@genwater.com

Accounting: accounting@genwater.com

INDEX

Safety References	3
Specifications	4
Water Purity Conversion Charts	5
Purification Process	6
Continuous Decontamination	6
System Features	7
Flow Chart	8
DI Cartridge Replacement Procedure	9
Recommended Filter & Membrane Replacement Schedule	10
Making a connection with Push-Fit Fittings	10
Daily Checklist	11

Safety References

This water purification system complies with state-of-the-art standards in terms of technology and safety regulations. All parts have been produced and or manufactured by preferred vendors to GWT specifications, and then assembled by GWT personnel (henceforth “supplier”) according to our proprietary design.

The supplier has taken all precautions to ensure safe operation. Users must ensure that the installation is performed in a manner that the method/mode of operation will not be affected. A signature at installation verifies that this is the case. Return visits to amend this placement will be subject to additional fees.

The system is supplier factory-approved and fulfills our safety standards.

Please read the following safety precaution before start-up and comply with them during operation.

The system must be operated by trained staff. Users should read this manual prior to utilizing this equipment.

The user is allowed to carry out only that maintenance work described in this manual or provided at training by the installation team. Only use designated spare parts for maintenance work.

The system must be connected to a fused, grounded power supply.

If leaks in the system become evident, the unit must be isolated immediately from the electrical supply and the fault rectified. If required, inform your suppliers customer service department.

Carry out maintenance and inspections as recommended by the supplier at the recommended intervals to insure user-safety and optimal system performance.

The supplier does not take any responsibility for injuries to users caused by improper application of this system, including negligence or abuse.

Specifications

Basic System Info (per system)

Model	GWT 2010-10
Production Capability	5 l/h
Delivery (Peak Flow) Capability	~8~ l/h
Footprint (System face on)	24" depth x 14" wide x 23" tall
System Weight	70 lbs in operation
Door Hinge	Right, opening left to right
RO Pressure (Range)	80 - 100 PSI
Delivery Pressure	As Indicated by instrumentation, adjustable
Dispense Tap for Manual Fill	Included
Manual feed by lab in event of water outage	Not included in this model
Touchscreen Interface	Not included in this model
Live Remote Viewing	Not included in this model
Audible Leak Detection/Alarm	Not recommended
GWT 2000 DI Cartridge	Located inside door
LC 123 Carbon Filters	Inline filter
Startup Amperage	6-7 AMP
Running Amperage	3-4 AMP

Presite Requirements (per system)

Domestic Water Feed to Equipment	1/2" Supply line with 1/2" dedicated drop to unit
Connection to Domestic 3/4" Line	1/2" FNPT ball valve
Minimum Pressure	30 PSI
Pressure Gauge	Upstream side of 1/2" FNPT ball valve
Power Connection	120 VAC/20 AMP
Drain Capacity	1 gallons/minute

Water Purity Conversion Charts

Parts Per Million (PPM)	Grains Per Gallon (GPG)	Hardness Rating
0 - 75	0 – 4.3	Soft
76 - 150	4.4 – 8.7	Moderately Hard
151 - 300	8.8 – 17.5	Hard
Over 300	17.6+	Very Hard

Conductivity (Microsiemens/cm)	Resistivity (MegOhm/cm)
0.056	18
0.063	16
0.071	14
0.083	12
0.100	10
0.133	7.5
0.200	5
0.500	2
1.000	1
1.333	0.75
2.00	0.5
4.00	0.25
10.00	0.1
20.00	0.05
40.00	0.025
80.00	0.013
100.00	0.01
200.00	0.005
500.00	0.002
1000.00	0.001
2000.00	0.0005
5000.00	0.0002
10000.00	0.0001

Purification Process

The process flow with all its components is illustrated on the flow chart (see pg. 8).

When the system is started in normal operation the inlet valve opens and the reverse osmosis pump starts. The inlet valve controls the water supply into the system. The inlet valve will close, and the pump will stop when the system is shutdown, unplugged, paused, permeate tank is full, or the inlet water pressure drops.

The sediment filter (if applicable) and carbon filter absorb organic chemicals like chlorides, insecticides, pesticides and herbicides from the water and remove particles larger than 5 microns.

The reverse osmosis pump increases the input pressure to the operating pressure for the reverse osmosis process.

The reverse osmosis modules remove ~ 98%~ of the waters salt content. Moreover, nearly all bacteria, viruses, heavy metal complexes, and organic material with molecular weight exceeding 300 Daltons are rejected. The purified water product is called Permeate; the wastewater including all rejected components is called concentrate.

The concentrate is then directed the drain.

The permeate quality is controlled with a Permeate sensor. The Permeate quality is displayed on the front display.

The system is connected to a storage tank with a level sensor. The production of permeate to fill the tank is regulated by the level sensor. When the permeate reaches the maximum level (100%) the production will stop. If the permeate falls below 70% the production of the permeate will be started again and the tank will be filled to the maximum level. If the permeate level falls below 10% distribution pump will automatically shut off and illuminate "Tank Empty" light to prevent a pump failure.

Continuous Decontamination

The permeate is stored in the storage tank. The tank is equipped with a distribution / circulation pump. The pump is continually circulating the water from the storage tank through the UV disinfection (the cell walls of bacteria, virus, and protozoa are penetrated, permanently altering the DNA of the microorganisms. This effectively inactivates the microorganisms, making them unable to infect and reproduce). Then passes through a Mixed Bed DI Resin as well as a .2-micron filter. From there, water either feeds to the analyzer on a demand-basis, or back into the tank, where it will undertake the process another time, thus ensuring continuity in the decontamination process.

The quality of the pure water being delivered is constantly monitored by a Product MOhm*cm Sensor viewable on the front Display.

System Features

This water purification system is designed to produce Pure Water out of drinking water from public supply. It unifies purification technologies with Reverse- Osmosis according to the latest developments. The system is connected to a tank with level control measurement.

The production rate of pure water has reduced salt content of 98-99%. Bacteria and particles are rejected for more than 99%. Water quality and tank filling is displayed on the systems screen.

The system runs in the following modes:

RO RUNNING Mode: Pure water production and tank filling

```
RO Running
Permeate Cond    0äS
Permeate Temp 26BC
RO Runtime 00000 HRS
```

TANK FULL Mode: Production is stopped, and the system is waiting for new water demand.

```
Tank Full
Permeate Cond    0äS
(From last run)
RO Runtime 00000 HRS
```

TANK FULL DRAW DOWN Mode: Production is stopped, and the system is waiting for new water demand.

```
Tank Full--Draw Down
Permeate Cond    0äS
(From last run)
RO Runtime 00000 HRS
```


POWER OFF Mode: System has been turned off.

```
Power Off
```

All components of the RO Cabinet are assembled in an aluminum housing. The front door can easily be opened for the access of the modules. The rear panel can be accessed by removing four screws.

The pure water storage tank is 5 liters and is made of Polypropylene.

Flow Chart




DI Cartridge Replacement Procedure

Materials needed:

- 2 new GWT DI cartridges
- Filter wrench
- Disposable gloves

How to Replace DI Cartridges:

- 1.** Turn off delivery pump.
- 2.** Remove filter housing by turning counterclockwise. (Hint: hold filter wrench in your right hand and push toward left hand.)
- 3.** Remove old cartridges from filter housings.
- 4.** Remove gaskets from inside the blue caps. (DI Cartridge gaskets stick to caps)
- 5.** Drop in new DI Cartridges.
- 6.** Install Filter housing by turning housing clockwise until hand tight



- 7.** Turn delivery pump on.
- 8.** Depress the red pressure release buttons on top of the blue caps until water flows out and housings are full.
- 9.** Open Vent Tube Valve attached to in-line .2-micron filter until you get a continuous flow of water.
- 10.** Close Vent Tube Valve


Recommended Filter & Membrane Replacement Schedule

FILTER	FREQUENCY (Approx)	
DI Cartridges	Change as your own standard operating procedure dictates. GWT recommends:	<ul style="list-style-type: none"> • Type 1 water = 10.0-18.2 Megohm (replace when below 10Megohm) • Type 2 water = 1.18.2 Megohm (replace when below 1 Megohm)
Carbon Filter	1-3 Years	GWT Replaces this during yearly recertification.
RO Membrane	1-5 Years	GWT Replaces this during yearly recertification.
.2 Micron Germ Filter	1-3 Years	GWT Replaces this during yearly recertification.
Ultraviolet Light	9000 Hours	GWT Replaces this during yearly recertification.

*Replacement frequency depends on feed water quality and pure water usage.


Making a connection with Push-Fit Fittings

John Guest fittings are used throughout the system.




Step 1 – Cut Tubing

Cut the tube square ensuring it is free of score marks. Avoid damage to the internal. O-ring by removing burrs and sharp edges on the tubing. Inside the fitting are stainless steel grippers that grab the tubing.



Step 2 – Assemble the Connection

Push the pipe into the fitting to the tube stop. The collet has stainless steel grippers that hold the pipe firmly in position. The O-ring provides a permanent leak proof seal. Simply pull on the tube to check it is secure.



Step 3 – Disconnect

Depressurize the line before removing the tube. Push in the collet (ring closest to tube) squarely against the face of the fitting. The tube can be removed when the collet is pushed in. The fitting can be re-used.



Daily Checklist

Month _____ Year _____

Serial # _____

*Please contact General Water Technologies 801-294-2426 with any questions.

Date	Product Mohm*cm <i>Range 10.0 - 18.2</i>	Permeate, μS/cm <i>Range 0.0 - 50.0</i>	Runtime	CLS
1	<i>Changed?</i>			
2	<i>Changed?</i>			
3	<i>Changed?</i>			
4	<i>Changed?</i>			
5	<i>Changed?</i>			
6	<i>Changed?</i>			
7	<i>Changed?</i>			
8	<i>Changed?</i>			
9	<i>Changed?</i>			
10	<i>Changed?</i>			
11	<i>Changed?</i>			
12	<i>Changed?</i>			
13	<i>Changed?</i>			
14	<i>Changed?</i>			
15	<i>Changed?</i>			
16	<i>Changed?</i>			
17	<i>Changed?</i>			
18	<i>Changed?</i>			
19	<i>Changed?</i>			
20	<i>Changed?</i>			
21	<i>Changed?</i>			
22	<i>Changed?</i>			
23	<i>Changed?</i>			
24	<i>Changed?</i>			
25	<i>Changed?</i>			
26	<i>Changed?</i>			
27	<i>Changed?</i>			
28	<i>Changed?</i>			
29	<i>Changed?</i>			
30	<i>Changed?</i>			
31	<i>Changed?</i>			