ENGLISH

INSTRUCTIONS FOR USE - INSTRUMENTS

In accordance with the directive 93/42/EEC relative to medical devices and its amendments, this product must be handled and/or implanted by WELL-TRAINED AND QUALIFIED PERSONS, AWARE OF THESE DIRECTIVES FOR USE.

The instruments are intended for use in surgery, and should be used only for the introduction of associated Hintermann Series products.

None of the instruments should be implanted.

Only medical professionals who are thoroughly familiar with the instruments function, application, and use should use them in surgery.

Only a medical professional who can perform the orthopedic surgery required by the particular patient should use the surgical instruments.

Improper maintenance, handling, or poor cleaning procedures can render the instrument unsuitable for its intended purpose, or even dangerous to the patient or surgical staff.

These devices do not contain phthalates unless it is indicated on the label.

Products are sold either sterile or non-sterile.

1. Single use instruments:

The single use orthopedic instruments manufactured by DT MedTech, LLC must not be reused. Re-use would incur the risk of modifying the properties and performance of the instrument. DT MedTech, LLC accepts no responsibility for such re-use.

2. Packaging (produit sold STERILE only):

Instruments manufactured by DT MedTech, LLC and sold sterile have been sterilized either by gamma radiation or using ethylene oxide (ETO). The sterilization method is specified on the packaging. Components sterilized by radiation are exposed to a minimum of 25 kGy of gamma irradiation. Check packaging and labeling before use. The sterility is guaranteed as long as the packaging has not been damaged within the expiration date.

Do not use any product in which the packaging has been opened or damaged outside the operating theatre.

Instructions for re-sterilization:

Re-sterilization is only allowed for non-used products. Remove delivery packaging in compliance with current regulations to (re)sterilize non-sterile products. DT MedTech, LLC recommends to sterilize its products by the steam autoclaving procedure regularly used in the hospital. (cf. Handling and Reprocessing).

3. Re-use of the instruments:

Unless labeled for use, the instruments could be re-used. Active surgical instruments such as those used with powered equipment, impact or bone manipulation, have a limited lifespan.

4. Handling and reprocessing (NON-STERILE product or considered to be):

Check the integrity of the packaging and labeling before opening the packing.

Remove all the products from their packaging prior to sterilization.

All products should be cleaned, decontaminated, and sterilized before use.

All instruments should be cleaned, decontaminated and sterilized after use.

Re-sterilization has little effect on these products.

Preparation: Instrument assembles (e.g., internal screwdriver and associated external screwdriver) should be separated prior to cleaning.

Cleaning: Cleaning can be performed manually or ultrasonically in accordance with the specifications designated by the manufacturer of the hospital's equipment.

Manual cleaning:

Rinse the devices under cool, running tap water to remove gross soil. Manual cleaning consists of the use of a soft bristled cleaning brush (e.g., Spectrum M16) with an aldehyde-free cleaner. Special care is taken for threaded, difficult to reach, and lumened devices; fluidic devices with a large amount of the detergent solution and use a soft tumbler brush (size 4-5 mm diameter x 24") to clean the devices. Visually inspect until all areas are clean.

Note: Certain solutions such as those containing bleach or formalin may damage the devices, and they must not be used.

Use of metallic brushes or other abrasive products is also forbidden.

Steps 1 and 2

Conditions	Results	
Product used	ENZOL, or other aldehyde-free, phenol-free and QUAT - free solutions or concentrates	
Concentration	12ml/Liter	
Soak disassembled instruments	12ml/Liter	
Seeking volume used	4 Liters	
Seek bath temperature	> 35°C	
Seek time	16 minutes	
Step 3	Tap water temperature < 20°C	
Rinse instruments	1 minute per device	
Rate of water	3 Litres/minute	
Step 4, 5 and 6	Product Used	
Soak instruments	VALSULF or other aldehyde-free, phenol-free and QUAT - free solutions or concentrates	
Concentration	12ml/Liter	
Soaking volume used	4 Liters	
Seek bath temperature	< 20°C	
Seek time	15 minutes	
Ultrasound action	15 minutes	
Brushing	Yes	
Step 7	RO/DI water temperature 22-28°C	
Rinse instruments	2 minutes per device	
Rate of water	2 Liters/minute	
Step 8	Inspection	Visual clean and not slippery
Visual inspect	Sécher les instruments et les réassembler	
Step 9	Drying	Yes

Dry instruments and re-assemble

Cleaning should be immediately followed by profusely rinsing with deionized water.

Check that water flows out the cannulated parts.

Drying: Temperature should not exceed 95°C.

Packaging: No specific requirements.

Sterilization:

Stainless Steel sterilization trays *Pre-Vacuum **Pre-Vacuum
Minimum Temperature 134°C (273° F)
Exposure Time 18 minutes
Purge 2-3 minutes
Vacuum drying 20 minutes

Note: *This sterilization method is recommended for use in some countries outside of the USA from World Health Organization.

**This sterilization method is recommended for use inside the USA.

The fully loaded implant and instrument tray is recommended to be steam sterilized by the hospital using wrapping which is FDA and/or ISO compliant.

Non-Stainless Trays *Pre-Vacuum **Pre-Vacuum
Minimum Temperature 134°C (273° F)
Exposure Time 18 minutes
Purge 2-3 minutes
Vacuum drying 36 minutes

Note: *This sterilization method is recommended for use in some countries outside of the USA from World Health Organization.

**This sterilization method is recommended for use inside the USA.

The fully loaded implant and instrument tray is recommended to be steam sterilized by the hospital using wrapping which is FDA and/or ISO compliant.

The autoclave must be properly installed, maintained and calibrated. These sterilization parameters assume that all instruments have been properly cleaned and decontaminated prior to sterilization. The sterilization parameters have been validated with specific configurations as noted in the tray instructions. If other configurations are added to the tray or to the sterilizer, the recommended parameters may not be valid and new cycle parameters may need to be validated by the sterilizer. However, individuals or hospitals not using the recommended method are advised to validate the alternative method using appropriate laboratory techniques. ETO sterilization or cold sterilization techniques are not recommended.

⚠ Warning: DT MedTech, LLC cannot guarantee sterility for products that have been cleaned and/or sterilized by the purchaser or user.

4. Examination:

The instruments should always be examined by the user prior to use in surgery.

Examination should be thorough, and in particular should take into account the presence of any visible defects, cracks, bending, or distortion, and that all components of the instruments are functional.

Never use instruments with obvious signs of excessive wear, damage, incomplete or which are otherwise non-functional.

5. Safety (active surgical instruments):

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K-wires are single use and must be replaced for each procedure.

The surgeon using the active surgical instrument is responsible for the proper operation of the instrument as well as any accessories or equipment, including the use of the K-wires.

Avoid using excessive force, twisting, or bending of the active surgical instrument in any unnatural or untrained way.

The active surgical instrument must be properly inserted and securely locked into the proper instrument before the instrument is turned on and/or engaged.

The surgical instrument may become heated from friction and the surgeon should take appropriate care to ensure that the patient is not harmed.

Monitor the heat output to prevent any possibility of tissue damage due to excessive heat effects.

The surgeon must make sure that the heat output is not used for any purpose other than its intended use in the orthopedic surgical procedure.

The surgical instrument must not be modified.

Resharpening of active surgical instrument should not be performed under any circumstances.

Contact with other metal objects could cause damage to the surgical instrument and may necessitate replacement.

DT MedTech, LLC informs the surgeon that repeated uses of the active surgical instrument can lead to incidents which would compromise the surgical technique or the results of the procedure.

6. Measuring instruments:

Some surgeries require instruments which incorporate a measuring function.

Ensure that these markings are not worn and any surface engravings are clearly visible.

Unless otherwise indicated on the instrument itself, the measurements provided by these instruments have the following characteristics:

Length Measurements Units: millimeter (mm) with accuracy in reading +/- 1mm.

Angle Measurement Units: degree (" with accuracy in reading +/- 1°.

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