

multiSTATION® LITA/RITA RetroSterno™ CURVED PADDLES

Technology Guide

READ THIS PRODUCT INSERT THOROUGHLY BEFORE USE




FIG. 1 – multiSTATION® LITA/RITA RetroSterno™ CURVED PADDLES

SOLUTIONS®


multiSTATION® LITA/RITA RetroSterno™ CURVED PADDLE DEVICE DESCRIPTION:

The RITA RetroSterno™ CURVED PADDLE and LITA RetroSterno™ CURVED PADDLE (FIG. 1) are reusable, sterilizable accessories to the SubX™ STERNAL ASCENDER™ SYSTEM. The LITA RetroSterno™ CURVED PADDLE 1 provides retraction of the sternum and sub-xiphoid access to the left internal thoracic artery. The RITA RetroSterno™ CURVED PADDLE 2 provides retraction of the sternum and sub-xiphoid access to the right internal thoracic artery. The textured zone 3 provides grip on the sub-sternal tissue during retraction. The RetroSterno™ CURVED PADDLES are offset laterally to maximize working space and concavely contoured to facilitate optimal visualization of the targeted vessel. The shape of the paddle post 4 interfaces with space lateral to the xiphoid process and opposite the targeted internal thoracic artery to create a low-profile anatomical insertion. The attachment studs 5 allow for quick connection to the accessory slot of the SubX™ STERNAL ASCENDER™ SYSTEM.

INDICATIONS FOR USE:


The multiSTATION® LITA RetroSterno™ CURVED PADDLE/multiSTATION® RITA RetroSterno™ CURVED PADDLE is indicated for use for retraction of bone and soft tissue.

1 SELECT



1. **SELECT** the LITA paddle for access to the left internal thoracic artery or the RITA paddle for access to the right internal thoracic artery.

2 ATTACH




2. **ATTACH** a paddle by inserting the paddle post into the accessory slot of the hand grip and pulling down to secure.

INSTRUCTIONS FOR USE

SETTING UP, ATTACHING, and POSITIONING the *multiSTATION*® Bridge

Note: Full setup instructions are found in the instructions for use for components 1-3 below.



multiSTATION® Bridge DESCRIPTION:

The *multiSTATION*® Bridge (FIG. 2) comprises two *RIGID miniARM*® Devices, two *RAIL CLAMPs*, a *SubX*™ *STERNAL ASCENDER*™ SYSTEM, and a compatible accessory such as a *LITA RetroSterno*™ *CURVED PADDLE*.


The bridge is assembled by attaching the rail clamps 1 to draped surgical table rails, connecting a *miniARM*® 2 to the rail clamp, and then connecting the ascender 3 to the *miniARM*® devices. This configuration provides a stable, rigid system with which compatible accessories may be securely raised and lowered. The *miniARM*® may be connected to either a side or center receiving DOC of the rail clamp to provide positional adjustability. The ascender, with a compatible accessory paddle 4 attached, is connected to the two *miniARM*® devices and adjusted into position. Once the desired position of the bridge is reached, the system is locked into place. The swivel bar of the ascender is used to raise the compatible accessory paddle, providing a stable opening to the surgical site. The accessory paddle can be lowered by rotating the swivel bar clockwise while pressing the ratcheting latch. The bridge can then be repositioned, or a different accessory paddle can be installed.

NOTE: The *multiSTATION*® *LITA RetroSterno*™ *CURVED PADDLE*, and *RITA RetroSterno*™ *CURVED PADDLE* must be cleaned and sterilized prior to use.

INSTRUCTIONS FOR USE

ATTACHING to the *SubX™ STERNAL ASCENDER™* and POSITIONING the *LITA/RITA RetroSterno™ CURVED PADDLE*


1 ATTACH



1. **ATTACH** a paddle to the sternal ascender by sliding the paddle attachment stud into the sternal ascender accessory slot and pulling down to secure.

NOTE: Reference the sternal ascender IFU for instructions on attaching the sternal ascender to the *miniARM®* devices.

2 POSITION




2. **POSITION** the paddle by adjusting the ascender and/or the *miniARM®* devices. The paddle post should interface with the space lateral to the xiphoid process, opposite the targeted internal thoracic artery.

NOTE: The paddle indicator on the top of the hand grip represents the location of the tip of the paddle. The paddle should always be kept parallel to the sternum.

Refer to the *SubX™ STERNAL ASCENDER™* instructions for use of the completed system.

PRECAUTION: Like all metal implements used in surgery, the components of this system conduct electricity and can be associated with electrical shock or arcing from cautery or other sources of current.


3 REMOVE






3. **REMOVE** the paddle by lifting the post up and out of the accessory slot of the hand grip.

multiSTATION® RetroSterno™ PADDLE REPROCESSING


- Disassemble and clean the device immediately after use. Do not allow a soiled device to dry.
- The RetroSterno™ PADDLES are not validated to be cleaned or sterilized while attached to the SubX™ Sternal Ascender.
- Cleaning agent used in validation: Steris Prolystica® 2X (enzymatic, neutral pH).
- Perform the final rinse using only freshly prepared purified water/highly purified water.
- Never use metal brushes or steel wool for cleaning.
- Prepare and reprocess other multiSTATION® devices according to each device's instructions for use.



multiSTATION® RetroSterno™ PADDLE REPROCESSING (continued)

 LUBRICATION	1. Apply instrument lubricant mixed to manufacturer's recommendations to prolong instrument life by submerging the entire device in the lubricant for a minimum of 30 seconds. If the hospital washer-disinfector has a lubrication cycle this can be used instead of manual lubrication. NOTE: LSI has validated the use of MicroLube™ C Instrument Lubricant on this device. Other instrument lubricant brands have not been tested and performance and results cannot be guaranteed.
 INSPECTION	1. Carefully inspect the device to assure that all visible soil has been removed. Generally, unmagnified visual inspection under good light conditions is sufficient. Repeat cleaning process if soil is detected. 2. Visually inspect the device for mild or excessive corrosion. If corrosion is present, discontinue use of the device in surgery, but complete reprocessing. 3. Visually inspect the device for any damage. If parts are damaged, discontinue use of the device in surgery, but complete reprocessing.
 PACKAGING	1. If using a sterilization pouch, place each instrument in its own individual pouch. If using a sterilization tray, read the sterilization tray instructions for use before proceeding. Ensure that the sterilization tray has been cleaned according to the sterilization tray instructions for use and load the tray bottom according to FIG. 4. 2. Package the device according to TABLE 1. The barrier system for sterilized re-usable instruments should meet the following requirements: • ISO 11607-1 • Suitable for pre-vacuum steam sterilization • Appropriate for medical use • Grade appropriate for weight of loaded tray per sterilization tray instructions for use and facility procedures 

multiSTATION® RetroSterno™ PADDLE REPROCESSING (continued)



STERILIZATION

1. The device must be properly cleaned prior to sterilization.
2. Perform sterilization cycle according to TABLE 1:

TABLE 1: multiSTATION® LITA/RITA RetroSterno™ CURVED PADDLE Sterile Packaging and Processing			
Method	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Moist heat (steam) sterilization according to ANSI/AAMI ST79	Immediate use steam sterilization according to ANSI/AAMI ST79
Container	multiSTATION® Sterilization Tray P/N 100034	No tray	No tray
Cycle	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)	Pre-vacuum (Pre-vac)
Packaging	2-layer polypropylene wrap	Pouch	No packaging
Temperature	132°C (270°F) (minimum)	132°C (270°F) (minimum)	132°C (270°F) (minimum)
Exposure Time	4 minutes (minimum)	4 minutes (minimum)	4 minutes (minimum)
Dry Time	65 minutes (minimum)	25 minutes (minimum)	N/A

Device(s) processed by immediate use sterilization should be transferred immediately, using aseptic technique, from the sterilizer to the point of use.

Refer to ANSI/AAMI ST79, Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.

The multiSTATION® LITA RetroSterno™ CURVED PADDLE, and RITA RetroSterno™ CURVED PADDLE have been validated to 20 reprocessing cycles. The useful lifespan of a surgical instrument is largely dependent on the care and handling of the instrument. Careful inspection and functional testing of the instrument should be used to determine the end of its serviceable life.



STORAGE

1. During storage, ensure the device remains in a sterile condition ready for reuse.
2. Shelf life is dependent on the sterile barrier employed, storage manner, and environmental and handling conditions.

CONTRAINDICATIONS

- Do not use with attachments other than accessories provided by LSI SOLUTIONS®.
- These devices are not intended for use except as indicated.

WARNINGS

- Federal law restricts this device to sale, distribution, and use by, or on, the order of a physician.
- Read and become familiar with all instructions, warnings, and cautions before using this product.
- The *LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE* shall be used in accordance with these instructions for use.
- Improper use of the *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE* may cause serious injury or death. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.
- When using the *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE*, patients must be immobilized or anesthetized.
- Discontinue use of the *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE* when patient is moving or being moved.
- Do not use this retraction system without adequate knowledge or experience regarding surgical retraction, including sternal and rib retraction. Avoid compressing, tensioning, displacing, or compromising tissue structures, such as the heart, that may be injured or lead to intraoperative or postoperative loss of function.
- Surgical or endoscopic procedures should be performed only by physicians having adequate training and familiarity with relevant techniques and anatomy. Medical literature relating to techniques, complications, and hazards should be consulted prior to use.


PRECAUTIONS

- **The *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE* are packaged as non-sterile.** Clean and sterilize prior to use.
- Proper cleaning, sterilization, packaging, storage and deployment are required to ensure sterility and safe use.
- If there are any variations between these instructions for use and the policies of your facility and/or the instructions of your cleaning/sterilizing equipment manufacturer, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your device.
- Use of the *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE* for a task other than what it is intended for can result in a damaged or broken device, and/or injury or death.
- Prior to use, inspect the *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE* to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.
- Ensure operating room table rails are of adequate strength to support use of this system to lift tissue structures as indicated for the individual patient.
- Ensure system connections, fittings and engagement features are secured appropriately, fully and accurately positioned to confer the required function.
- Ensure placement of system components to not compromise patient access during surgery or injure patient or tissue structures.
- Ensure retractors are placed in contact with tissue structures under adequate visualization. Do not force system components against non-targeted tissue.
- Ensure all system components are appropriately removed from the patient and operating room table prior to patient movement.
- Surgical instruments vary between manufacturers. Before instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility and ensure electrical isolation or grounding are not compromised.
- Like all metal implements used in surgery, the components of this system conduct electricity and can be associated with electrical shock or arcing from cautery or other sources of current.
- Avoid mechanical shock or overstressing the *multiSTATION® LITA RetroSterno™ CURVED PADDLE*, and *RITA RetroSterno™ CURVED PADDLE*.
- Only the cleaning and sterilization processes that are defined within these instructions for use have been validated.
- Store at room temperature.

ADVERSE REACTIONS

- No documented adverse reactions.

ORDERING INFORMATION

TABLE 2: <i>multiSTATION</i> ® LITA and RITA RetroSterno™ CURVED PADDLES PRODUCT ORDERING			
	REORDER	PRODUCT	DESCRIPTION
	REF 080900	<i>multiSTATION</i> ® SubX™ STERNAL ASCENDER™ SYSTEM	1 Shelf Box
	REF 081870	<i>multiSTATION</i> ® LITA RetroSterno™ CURVED PADDLE	1 Shelf Box
	REF 081880	<i>multiSTATION</i> ® RITA RetroSterno™ CURVED PADDLE	1 Shelf Box
	REF 100034	<i>multiSTATION</i> ® Sterilization Tray*	1 Shelf Box

LSI SOLUTIONS®

LSI SOLUTIONS, Inc.
7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Phone: +1 585.869.6600
Customer Service: +1 866.575.3493
Technical Support: +1 866.428.9092
Fax: +1 585.742.8086
www.lsisolutions.com

*The *multiSTATION*®
Sterilization Tray is
manufactured by Summit
Medical, 815 Vikings
Parkway, Suite 100, St. Paul,
MN 55121 U.S.A.

MADE IN THE USA

This Product Comes
with our LSI SOLUTIONS®
Perfect Performance Policy®
Call us at 866.575.3493 any time.

Patents: www.lsisolutions.com/patents

The LSI logo, LSI Solutions, Perfect Performance Policy, Miniarm, Multistation, Retrosternal, Sternal Ascender, and Subx are registered trademarks and trademarks of LSI SOLUTIONS®. Copyright © 2021, LSI SOLUTIONS®. All Rights Reserved.



www.lsisolutions.com/symbols

Prolystica is a registered trademark of Steris. MicroLube is a registered trademark of PurposeBuilt Brands. Use of these names, trademarks and brands does not imply endorsement. All company, product and service names used in this technology guide are for identification purposes only.