PURENEB® MANOSONIC

FR - Notice d'utilisation
EN - Instructions for use
DE - Bedienungsanleitung
NL - Gebruiksaanwijzing
IT - Manuale d'uso

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SUMMARY

1. INTENDED USE.....	14
a) Indications	14
b) Contra-indications	14
2. PRODUCT DESCRIPTION	14
a) Exploded view and description	14
b) Choice of nasal plug	15
3. CLEANING / DISINFECTION OF ACCESSORIES	16
a) Cleaning	16
b) Disinfection	16
4. USE	17
a) Accessories preparation	17
b) Compressor preparation	17
c) Use	17
d) After each use	18
5. AEROSOLTHERAPY.....	18
a) Performances of the MANOSONIC nebulizing system	18
6. DEVICE CLEANING	18
a) Device cleaning	18
7. IMPORTANT RECOMMENDATIONS.....	19
8. FAILURE REMEDIES.....	20
9. LIST OF ACCESSORIES AND SPARE PARTS	21
10. WARRANTY / RECYCLING	21
11. TECHNICAL SPECIFICATIONS	21
12. PICTOGRAMMES	22
13. COMPATIBILITE ELECTROMAGNETIQUE / ELECTROMAGNETIC COMPATIBILITY / ELEKTROMAGNETISCHE VERTRÄGLICHKEIT / ELEKTROMAGNETISCHE COMPATIBILITEIT / COMPATIBILITÀ ELETTRIMAGNETICA.....	56

1. INTENDED USE

This device is a manosonic nebulizer with tubo-tympanic target. That is to say, it will allow performing aerosoltherapy in the nasopharynx and eustachian tube.

What is an aerosol?

An aerosol is an air suspension of fine liquid particles. Their size is below 10 microns.

Manosonic aerosols = with 100 Hz sound vibrations and overpressure.

Aerosol diffusion is greatly improved by the periodic emission of sound vibrations. The application of overpressure facilitates diffusion in the eustachian tube.

a) Indications

The MANOSONIC nebulizing system is indicated for treatment of tubo-tympanic pathologies:

- secretory otitis media or seromucous otitis,
- tubotympanic sequelae (like retraction pocket...).


b) Contra-indications

- Child under 4 years old,
- Patient with a patulous Eustachian tube,
- Patient with cleft lip and palate.

2. PRODUCT DESCRIPTION

a) Exploded view and description

Nebulizing kit

1. Air inlet
2. Tank
3. Jet
4. Jet nozzle
5. Additional air inlet
6. Aerosol outlet
7. Vibrations and overpressure inlet

On the aerosol outlet (6), the user fixes the nasal plug, available in 3 sizes:

Big size


Medium size (standard)

Small size

Compressor

1. PRESSURE air outlet
2. Fuses
3. ON – OFF switch
4. “Overpressure ready” green indicator light
5. Adjustment knob for the overpressure
6. Power cord
7. VIBRATIONS and OVERPRESSURE outlet
8. Manometer indicating the value of the overpressure
9. Instructions for use
10. Nebulizing kit holder.

b) Choice of nasal plug

To make sure that the size is well suited to the nostrils of the patient and that airtightness is good:

Introduce the soft tips in the nostrils, obstruct the back of the nasal plug. The pressure rise should be felt in the nose without leakage to the edge of the nostrils. If not, humidify the cones of the nasal-plug and try again. Otherwise, try a bigger nasal plug.

3. CLEANING / DISINFECTION OF ACCESSORIES

a) Cleaning


Before the first use and after each session, for the elements on the picture beside:

Dismantle each element

Clean them in soapy hot water and rinse them

Let them dry on a clean cloth

b) Disinfection

Before the first use and once a day, for the elements present on the image beside:

Dismantle each element

Disinfect

- Place the elements in boiling water for 20 minutes,
- Or disinfect them cold according to the recommendations on the product box,
- Or in a microwave sterilizer according to the recommendations of the sterilizer,
- Or in the dishwasher at a maximum temperature of 70 ° C.

The nasal plug can only be **cold disinfected**, according to the recommendations on the product box.

Nota: cleaning and disinfection have to be made once a day.


Should the tubes be disinfected, let them perfectly dry before connecting them to the device.

Nota: after cleaning in boiling water the tubes colour can change (more white) but this would not change the good operating of your nebulizing kit.

4. USE

a) Accessories preparation

Respect the following steps:


Insert the jet nozzle on the bottom part.

Pour the solution into the tank, do not exceed the "full" level.

Screw the top part to the tank.

Fix the adapted nasal plug (see. § 2b) to the nebulizing kit.

Connect the pressure ("pression") narrow air tubing to the inlet (1) of the nebulizing kit and the "vibrations" large tube to the inlet (7) of the nebulizing kit.

b) Compressor preparation

Respect the following steps:

Connect the pressure tube to the pressure air outlet of the compressor.

Connect the vibrations tube to the vibrations outlet of the device.

Connect the device to the mains.

Start the device (switch on I).

c) Use

BLOW YOUR NOSE PROPERLY.

PHASE A: without overpressure – with sound vibrations (2 minutes)

As soon as the device is on, you should observe a mist coming out of the nebulizer.

Introduce the nasal plug air tightly into the nostrils.

Inhale calmly through the nose and exhale through the mouth for 2 minutes. Sound vibrations are emitted shortly before the green light (4) comes on.

PHASE B: with overpressure and periodic sound vibrations (around 6 min.)

Set the initial pressure to 2 or 3 (knob 5). When the green light is on, the overpressure is available.


Obstruct the additional air inlet of the nebulizing kit (5) with the thumb and swallow.

Swallowing releases the overpressure. The light is off. Release the additional air inlet of the nebulizing kit (5).

The patient must feel the release of the overpressure by a sensation of ear fullness.

If not, try again after having increased the overpressure with the adjustment knob (5). Breathe in through the nose and out through the mouth until the green light is on.

Before swallowing, it is recommended to gather one's saliva. If necessary, suck an acid drop or drink a sip of water with a straw.

The value of the overpressure reached can be read on the manometer. Adjust the overpressure between 10 and 50 mbar (0.15 and 0.75 PSI) according to the sensations of the patient and the indications of the doctor.

If the needle on the gauge does not move, check the airtightness of the nasal plug on the nose.

Start again the operation according to the doctor's indications.

Optimal duration: 6 min., that is to say 15 effective swallowings.

d) After each use

At the end of the treatment:

- Stop the device (switch on O).
- Disconnect the electric cord from the mains.
- Disconnect the tubes (pressure and vibrations) from the compressor.
- Follow the cleaning and disinfection instructions.


5. AEROSOLTHERAPY

a) Performances of the MANOSONIC nebulizing system

See granulometric curve opposite

MMAD: 4.5 μm

Cumulative particle size distribution of the 2.5% NaF aerosol produced by the MANOSONIC system

Information regarding **filling volumes**, **inhalable mass** and **air flow** are mentioned on the nebulizing kit packaging.

6. DEVICE CLEANING

a) Device cleaning

After each use, cleaning is limited to the external parts with a non-abrasive sponge soaked with alcohol or with a disinfecting product recommended by the manufacturer (Anios Surfa'Safe Premium).

7. IMPORTANT RECOMMENDATIONS

Use of electricity

- Before connecting the device to the mains, please check the device voltage and frequency (volts and Hz) and your distribution network.
- Warranty does not cover wrong connections.
- Connect the device to an easily accessible plug, to allow a rapid disconnection from the mains.
- PURENEB devices are designed to operate at altitudes less than or equal to 3 000 meters.
- Do not open the device when connected to the mains (electrical shock hazard).
- Caution: this medical device is directly connected to the electric power supply, do not use it in bath or under the shower, do not immerse.
- Keep this medical device away from any heat source.
- The mains plug is the disconnecting device and must always be accessible.

Use

- Keep this medical device out of the reach of children or animals (potential risks with the tubes and detachable part of the nebulizing kit).
- The device has to be installed and put into service according to the information supplied with the attached documents.
- Store properly the electric cord and tubes after each use.
- An animal can eat the electric cord. Check carefully the conditions of the equipment before use.
- This device is to be used in health care establishment or at home.
- Do not use the device outdoors, near a heat source or near a waterhole (bathtub, sink...).
- Make sure the device is not exposed to direct sunlight or excessive heat.
- After use, disconnect the device from the mains.
- If the compressor with the nebulizing kit is too noisy, do not hesitate to contact your local supplier.
- Fragile medical device, to be manipulated with care.
- Before connecting the device to the electric power supply, make sure that the nebulizing kit is connected to the unit in order to avoid an abnormal noise. (Compressor powered without its nebulizing kit is noisier).
- No special knowledge or training is required to use PURENEB devices.
- The casing of the device may reach 45 °C.
- Do not touch the device for more than 1 minute during use.
- Do not touch the tubes for more than 1 minute during the device use.
- Carefully keep these instructions for use booklet with the device as it provides necessary information for the good use of the equipment.
- The electric cord cannot be replaced by the maintenance staff.
- The patient is an intended operator. He can use all functions of the device safely.
- Contact the manufacturer if you need help to set up, use or maintain the medical device and to report any unexpected operation or event.
- A more detailed technical description is available on request

Remarks regarding aerosoltherapy

- Do not exceed 12 ml solution in the nebulizing kit.
- Do not exceed the mark "full" printed on the nebulizing kit.
- Clean if necessary the jets of the nebulizer with hot water, without using any metal wire.
- A modification of jets dimensions can have a bad effect on the particle size and on the aerosols therapeutic efficiency.
- Each nebulizing kit is intended for a single patient.
- Aerosol treatments can only be done on medical prescription and within the prescribed doses.
- Side effects: a risk of allergy of the body parts in touch with the nasal plug can appear.
- Unit not to be used in presence of a mixture of inflammable anaesthetic with air or oxygen or nitrous oxide.

CAUTION, in all cases:

- WARNING: Always disconnect the electrical cord from the mains supply before opening the body of the compressor (electrical shock hazard).
- WARNING: Only a professional repair shop is allowed to repair the device.
- WARNING: Modification forbidden.
- WARNING: Only use the accessories supplied and mentioned in the instructions for use.
- WARNING: Do not undertake service or maintenance operations while the device is in use.

8. FAILURE REMEDIES

Problem	Solution
No operating noise can be heard	Check that the electric cord is connected to the mains and to the device. Check that the switch is in position on. Check that main power supply is in accordance with the unit (Volt Hz). Check that electrical current arrives to the socket trying with another device.
No aerosol comes out of the nebulizing kit	Check that the nebulizer contains enough liquid: the minimum is indicated on the instructions for use of the nebulizing kit. Check that the tube is well connected to the compressed air inlet of the nebulizing kit and to the air outlet of the compressor. Check that the tube is not pinched or blocked. Check that the jet nozzle is completely inserted. Jet nozzle must be cleared (clean it with water).
Overpressure is not released	Check the airtightness of the set: nostrils – nasal plug – nebulizing kit.

9. LIST OF ACCESSORIES AND SPARE PARTS

Accessories

- Extraneb Manosonic adult nebulizing kit, with medium and big nasal plugs, ref. MS1A-GB
- Extraneb Manosonic paediatric nebulizing kit, with medium and small nasal plugs, ref. MS1E-GB

Spare parts

- Nasal plug, large size, ref. C28G
- Nasal plug, medium size, ref. C28E
- Nasal plug, small size, ref. C28F
- Tube 6x9, length 1 m, ref. C20P
- Tube 8,5x12, length 1 m, ref. T2
- Carrying case, ref. F64
- Fuses T 0,5 A H 250 V (according to IEC 60127-1)

10. WARRANTY / RECYCLING

Warranty: The compressor is warranted 3 years, parts and labour, against any manufacturing defect, except case of misuse. For after sales service under warranty, please contact your distributor.

Recycling

The device must be disposed of in accordance with local regulation on the disposal of electrical appliances. So you contribute to the human health and environmental protection.

Product on the market since August the 13th 2005.

11. TECHNICAL SPECIFICATIONS

Device class		IIa	Dimensions	L380xI250xh200 mm
Device class in electrical safety		2	Voltage	230 V ; 50 Hz
Type of applied parts		BF	Comsumption	105 VA
Expected shelf life of accessories		25 disinfection cycles	Using pressure	400 hPa
Expected shelf life of the device		5 years	Sound level	< 70 dB(A)
Weight (with carrying bag)		9.5 kg	Flow rate	5,15 L/min
Overpressure	Adjustable from 20 to 50 mbar		Running mode	Continuous
Electric supply		AC 230V (~ 230 V) 50 Hz		
Protection level		IP 30 - Protected against ingress of solid foreign objects ≥ 2.5 mm Not protected against vertical water drops falls		
Room temperature (T) – Transport and storage				-25°C < T < +40°C
Room temperature (T) – Operating				+5°C < T < +40°C
Relative humidity (H) – Transport, storage and operating				15% < H < 93%
Atmospheric pressure (P) – Transport, storage and operating				700 hPa < P < 1060 hPa

Applied part: nasal plug

Essential performance features: provide pressure to transform the solution into an aerosol.

Nota: CE certificate can be supplied on request.

Identification

- Each device has its own identification code (reference and unit) showing the manufacturing date, manufacturing order and personal code.

Example: 68 078 HE 03

68: manufacturing year (68 : 1950 + 68 = 2018)

078: day of the year

HE: personal code

03: number of order (in a series)

12. PICTOGRAMMES

Symbol on device	
	In accordance with Directive 93/42/CEE of 14th June 1993 regarding medical devices
	Manufacturer
	Device reference
	Device serial number
	Device class 2
IP 30	Level of protection
	Device type BF
	Alternative current (AC)
	See instructions for use
	Caution / Warning
	Dangerous voltage – Do not open the device
	Keep dry
	Fragile
	Keep away from light
○	Off
I	On
	Do not throw away. Product on the market since 13th August 2005

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Electro medical devices require special precautions against EMC and must be installed and set to operate in accordance with the EMC information provided in this manual.

WARNING: Avoid using this unit near or stacked onto other equipment as this may cause a malfunction. If this is necessary, this appliance and the other appliances should be monitored for normal operation.

WARNING: Use of accessories, transducers and cables other than those specified or supplied by the manufacturer of this device may cause an increase in electromagnetic emissions or a decrease in the immunity of this device and cause improper operation.

WARNING: Portable RF communication devices (including peripherals such as antenna cables and external antennas) should not be used closer than 12 inches (30 cm) to any part of the MANOSONIC device, including cables specified by the manufacturer. Otherwise, performance may be impaired.

WARNING: If the device malfunctions due to electromagnetic disturbances, the performance of the device may be impaired and result in delayed patient care.

Wireless communication devices such as wireless home networking devices, mobile phones, cordless phones and their bases, walkie-talkies can affect the device, and it is recommended to keep them at a distance of at least 0.5 m from the device.

Device:

Guidance and manufacturer's declaration – electromagnetic emissions		
Emissions test		Compliance
RF emissions - CISPR 11		Group 1
RF emissions - CISPR 11		Class B
Harmonic emissions - IEC 61000-3-2		Class A
Voltage fluctuations / flicker emissions - IEC 61000-3-3		Complies
Guidance and manufacturer's declaration – electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV ± 15 kV
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV ± 1 kV
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% U _T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
	0% U _T ; 1 cycle AND 70% U _T ; 25/30 cycles Single phase : at 0°	0% U _T ; 1 cycle AND 70% U _T ; 25/30 cycles Single phase : at 0°
Voltage interruptions IEC 61000-4-11	0% U _T : 250/300 cycles	0% U _T : 250/300 cycles
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m 80 MHz to 2,7 GHz
	385 MHz at 27 V/m 450 MHz at 28 V/m 710, 745, 780 MHz at 9 V/m 810, 870, 930 MHz at 28 V/m 1720, 1845, 1970 MHz at 28 V/m 2450 MHz at 28 V/m 5240, 5500, 5785 MHz at 9 V/m	385 MHz at 27 V/m 450 MHz at 28 V/m 710, 745, 780 MHz at 9 V/m 810, 870, 930 MHz at 28 V/m 1720, 1845, 1970 MHz at 28 V/m 2450 MHz at 28 V/m 5240, 5500, 5785 MHz at 9 V/m
Conducted RF IEC 61000-4-6	3V 0,15 MHz – 80 MHz 6V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% MA at 1kHz	3V 0,15 MHz – 80 MHz 6V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% MA at 1kHz