

Designed for the Workforce / Multiple Targeted Therapies in One Wrap

BENEFITS OF RED LIGHT THERAPY

Innovative, Patent-Pending Technology Our patent-pending therapeutic wavelengths deliver significant wellness benefits.

Safe LED-Based Therapy Safe and effective relief without the risks commonly associated with laser treatments.

Increased Circulation Stimulates blood flow to support healing and recovery.

Targeted Relief Targeted therapy helps soothe sore muscles and aching joints—right where you need it most.

BENEFITS OF HEAT/COMPRESSION/MASSAGE

Heat Therapy improves blood flow, relaxes muscles, eases pain and stiffness and reduces chronic pain

Compression Therapy improves circulation, reduces swelling, and promotes muscle recovery

Massage Therapy enhances circulation, reduces muscle tension and pain, and relaxes muscles and tissues

Combining these therapies can provide comprehensive benefits. This combined approach offers both pain relief and improved circulation to accelerate healing and enhance recovery.

SAFETY FIRST

Safety and Convenience Absolute RLT wraps are powered by convenient, rechargeable batteries. For added safety, these wraps feature a built-in controller with an auto-shut off function after a 20-minute therapy session. This automatic shut-off ensures a safe and controlled treatment.

USING YOUR ABSOLUTE RLT WRAP

Daily Use Use Absolute RLT for one 20-minute session per day, ideally 5 to 7 days per week.

Managing Discomfort If you're experiencing discomfort, you may increase usage to 2-3 sessions per day. As symptoms improve, gradually reduce to one session daily.

Discontinuing Use Once you're pain-free, you may stop using your Absolute RLT wrap. However, we recommend consulting your healthcare provider before making any changes.

Disclaimer: The information provided on this brochure about the Absolute RLT wrap is not a substitute for professional medical advice. Absolute DME does not claim to "cure" any specific type of musculoskeletal injury. A licensed physician should assess the nature of the patient's injury and determine if Absolute RLT is appropriate for their treatment plan. Always consult a licensed physician before using Absolute RLT or any light therapy device. Absolute RLT is only available with a valid Rx from a licensed physician.

NOTE: Available through Veterans First Contracting Program [VFCP – 38 U.S.C. § 8127(d)]. Accepts Government Purchasing Cards (GPCs).

Custom-Fit Comfort Our Ergonomically designed wraps contour to your body for maximum comfort and effectiveness.

Effortless Use One-touch operation. Just 20 minutes a day for noticeable results.

Therapy On the Go Lightweight and portable—take Absolute RLT with you for portable, effective therapy.


Non-Addictive Experience natural, drug-free relief without the risk of dependency.



KNEE WRAP
Model: 102351



SHOULDER WRAP
Model: 102370



LUMBAR WRAP
Model: 102371



Foot/Ankle Wrap
Model: 102463
NOTE: Coming Soon

Multiple Sessions If using Absolute RLT multiple times a day, allow at least 2 hours between sessions to ensure safe and effective treatment.

Clothing Guidelines Absolute RLT can be used over thin clothing, but for best results, apply wraps directly to bare skin.



660NM and 850NM
Infrared Light Therapy



Heat



Compression



Massage

LEARN MORE:

Contact Spartan Medical Inc.

(888) 240-8091

CS@SPARTANMEDICAL.COM



SPARTAN MEDICAL INC.

7537 Lindbergh Drive,
Gaithersburg, MD 20879

spartanmedical.com

UEI: LCJBD88KQ8G3

DUNS: 807438341 | CAGE:65VT5



Rev A 9/26