


User manual for the Health Care Provider (HCP)



A personalized treatment, with more insight and involvement in the care pathway.

App version : 2.0.0 and later

Manual version : 1.4

Data of release manual : 11/07/2025

Link to e-learning: <https://moveup.talentlms.com/index>

Table of Contents

Table of Contents.....	2
PRODUCT DESCRIPTION.....	3
<i>Intended use.....</i>	3
<i>Intended users.....</i>	3
<i>Target population / indication for use</i>	4
<i>Claims</i>	4
<i>moveUP app & devices</i>	4
<i>Warnings and precautions</i>	5
<i>Contra-indications and potential side effects</i>	5
<i>Notification of serious incidents.....</i>	5
<i>Human body contact</i>	6
<i>Accessories/products used in combination</i>	6
<i>Device lifetime</i>	6
<i>Company details</i>	7
PATIENT APP.....	8
MEDICAL MODULE	8
HEALTHCARE PLATFORM ACCOUNT CREATION.....	8
ONBOARDING OF PATIENTS	8
HEALTHCARE PLATFORM GENERAL FUNCTIONS	11
HEALTHCARE PLATFORM CARE PLAN FUNCTIONS	18
WELLNESS MODULE.....	21
WELLNESS MODULE ACTIVATION	21
WELLNESS MODULE OVERVIEW	22
Wellness module overview	23
Wellness report	23

PRODUCT DESCRIPTION

Intended use

For healthcare providers, moveUP provides efficient clinical management of pathologies and treatments. With the use of moveUP, valuable insights are given to the HCP/care team about the status of the patient condition and its progression. The HCP/care team can provide targeted therapy advice & adapt treatment (exercise) schemes based on the needs of the patient through the use of the moveUP app.

For patients, moveUP gives insights about their progress, provides instructions related to their disease or treatment and helps them manage their symptoms. Based on the patients' needs and timepoint in the patient pathway, patients can be followed-up intensively and receive active therapy advice through the moveUP app.

moveUP is intended for the clinical follow-up on adult patients who underwent (postoperative) or will undergo (preoperative) a Total Knee Arthroplasty (TKA), Total Hip Arthroplasty (THA), Chronic obstructive pulmonary disease (COPD) or Bariatric surgery.

Intended users

moveUP is intended to be used by patients and healthcare providers.

Patient

The main user of the mobile app and patient website is the patient. Patient users should be alert, mentally competent and capable of operating a tablet or smartphone.

Inclusion:

- Age:
 - minimum: 18 years
 - maximum: no limit
- Weight: NA
- Health & condition: capable of performing basic activities of daily living
- Language: understanding one of the available languages of the app (Dutch, French, German, English)

Exclusion:

- Patients who are mentally incompetent or having trouble expressing what they are feeling (for instance, mentally diseased people, people staying in elderly care centres.) are excluded.
- Patients who are not capable of operating a tablet or smartphone.
- Patients who can't understand one of the available languages of the app (Dutch, French, German, English)

Healthcare provider

The main user of the medical web interface is a healthcare practitioner (group) or clinical researcher (group), named the care team. The care team is able to operate a web interface via web browser on PC/tablet/smartphone. The healthcare practitioner needs to understand one of the available languages of the web interface (today only available in English).

Only qualified HCP with a diploma or certificate in the healthcare domain or qualified clinical researchers can access the platform and complete the moveUP e-learning trainings.

Target population / indication for use

moveUP is used by patients:

- Who underwent or will undergo a Knee (Resurfacing) arthroplasty
- Who underwent or will undergo a Hip (Resurfacing) arthroplasty
- Who underwent or will undergo a gastric sleeve or gastric bypass surgery
- Patients with COPD.

Claims

- Patients using moveUP are highly satisfied.
- Patients using moveUP show a high adherence/compliance level to their therapy.
- The use of moveUP is more cost-effective than standard of care for hip & knee arthroplasty
- The length of stay for Knee & hip arthroplasty can be reduced with the use of moveUP
- Patients can rehabilitate without the need to leave the home environment & without receiving F2F physical therapy session for TKA and THA
- Using moveUP is as effective as standard of care
- Using moveUP is as safe as standard of care

moveUP app & devices

The moveUP app runs on:

- iPhones/iPads (iOS 11.0 or newer)
- Android systems (Android 7.0.0 and newer).

If the patient has no compatible tablet/smartphone, moveUP can lease these devices to the patient.

- Lenovo Tab 3 - Android 6.0.1 or Higher
- Huawei MediaPad T3 - Android 7.0.0 or higher

The moveUP application needs an internet connection when in use.

The moveUP application is compatible with the following browsers:

- Chrome
- Safari
- Edge

- Firefox

We strongly recommend to use the latest version of the browser.

For your data safety please enable the 2-factor authentication(2FA) and use a safe network

Warnings and precautions

	• moveUP is not an emergency tool. In case of emergency please contact your doctor or call 112.
	• moveUP is not intended to monitor/treat vital parameters of critical diseases.
	• Patients need to be aware if questionnaires are not filled out on frequent basis or not filled out in a trustworthy, the care team has more difficulty to know their health situation & evolution.
	• moveUP is not a replacement of any treatment you need to follow, but is used as an addition.
	• The moveUP App offers personalized healthcare based on the patient's profile and therefore can not be shared, borrowed or traded between users. The information and advice of the moveUP application do not apply to a patient other than the one whose profile is stored in the application. The patient may not grant access to the moveUP App to third parties. Any other person can, even accidentally, send wrong information or change the profile of the patient, with the result that the moveUP-guided advice no longer corresponds to the needs of the patient.
	• The platform includes modules displaying wellness data. These data are not part of the medical device and are not intended to be used for clinical decision making

* This data originates from a non-medical device product & may not be validated as such and the data shown in the wellness module cannot be used for medical purposes. The wellness module information is not part of the CE certificate of moveUP.

Contra-indications and potential side effects

There are no contra-indications or known side effects

Notification of serious incidents

A serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or another person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health that resulted in any of the following:
 - life-threatening illness or injury,
 - permanent impairment of a body structure or a body function, hospitalisation
 - prolongation of patient hospitalisation, medical or surgical intervention to prevent life-threatening illness or injury
 - permanent impairment to a body structure or a body function,
 - chronic disease,
- a serious public health threat;

If you detect any of these cases please send an email to info@moveUP.care and your competent authority:

- Netherlands: meldpunt@igj.nl
- Belgium: meddev@fagg.be
- France: dedim.ugsv@ansm.sante.fr
- Germany: zlg@zlg.nrw.de

Human body contact

No human body contact with the patient or user, due to the nature of the product (software).

Accessories/products used in combination

There are no accessories.

Device lifetime

2-years, depending on the willingness of the user to update the app. We sent out a notification with the recommendation to update the app, to make sure the app will keep the performance & recently included features.

Company details

 moveUP NV/SA	Tel: +32 800 88 008	UDI: ++G166MOVEUPV200AU  (01)G166MOVEUPV200A (10)+\$+200MDR (11) /16D20251031 (21) HI
Head office Brussel: Cantersteen 47, B-1000 Brussels, Belgium Office Ghent: Oktrooiplein 1, B-9000 Ghent, Belgium	Email: info@moveUP.care www.moveup.care	

If you have any questions about the App or the tablet, you can contact moveUP via the following channels - in order of priority: Via the message system of the App. Via email to support@moveup.care - Via phone during office hours: +32 800 88 008



MoveUP guarantees delivery of the eIFU within seven (7) days upon request. Please contact info@moveUP.care if you want to obtain a free copy.

PATIENT APP

Please read the user manual for patients on www.moveup.care/appinfo to view the functionalities of the patient app.

MEDICAL MODULE

HEALTHCARE PLATFORM ACCOUNT CREATION

Step 1: The HCP registers on the healthcare platform on https://v1.medical.moveup.care/	
Step 2: moveUP will contact the HCP to add his/her profile to the correct organization.	
Step 3: Security measures: It is important to enable 2FA as you have access to sensitive health data.	
E-learning: For HCPs who will provide care through the use of the healthcare platform it is recommended to follow the e-learning courses for the indications listed on https://moveup.talentlms.com. We ask the HCPs to register for the e-learning course in their applicable language.	

ONBOARDING OF PATIENTS

Onboarding is the process of registering a new patient in the healthcare platform.

A crucial part is adding a care plan to the patient's profile. A care plan determines what will be visible for the patient in the application, e.g.: type of exercises, information modules, message functionality, etc. A care plan usually consists out of an intervention, a treatment and (a) care

team(s).

Onboarding patients

Patients Onboarding

The onboarding dashboard is easy retrievable in the health care platform in the upper left corner when clicking 'Onboarding'. It is also accessible through <https://patient-onboarding-admin-moveup.care> using the HCP credentials of the health care platform.

Click on 'new patient' to onboard a new patient and fill out the patient information.

After the patient information is filled out, press 'Save'.

New patient

Important note: Only one account can be coupled to an e-mail address. It is not possible to reuse an e-mail address that is already known within our database, e.g.: reusing the same e-mail address for the account of a partner or parent.

Patient's care plans

Next up, the care plan needs to be selected: press 'Start new care plan'. Drag and drop an intervention and treatment to the care plan and fill out the dates, name of the treating doctor and care team.

Drag and drop a study if applicable. Click Save.

Email variables

```

{{ firstname }}
{{ lastname }}
{{ code }}
{{ registrationlink }}
{{ organizationName }}
{{ surgeonName }}
{{ hosp_spec_par_start_NL }}
{{ organizationName }}
{{ bodyStructure }}




```






A mail template is shown. The HCP reviews the e-mail variables and clicks on 'send selected mails'.

An invitation mail will be sent to the patient with all the needed information to install and activate the application.

HEALTHCARE PLATFORM GENERAL FUNCTIONS

Below the HCP can retrieve all functions of the healthcare platform. The applicable functions will be shown based on the HCPs subscription. The extent of the available patient information is based on the patients subscription & his/her compliancy rate.

Patient overview		The health care platform gives a clear overview of the HCPs patients. In the top bar the relevant categories are shown that identify the patient and show relevant health parameters, unread messages, the lead clinician and patient compliance. By default, when opening the healthcare platform, the patient overview is sorted based on intervention date.
Tags		Select patients related to a specific pathology by using tags.
Filters		Select patients based on language, test patient, lead clinician or being activated or not. Activated can mean 2 things: either the patient registered him/herself with a code through the app or he has been marked as immediately activated in onboarding. Non activated: patients who haven't registered themselves through the app yet and haven't been marked as immediately activated.
Color codes		Color codes indicate which patients need closer attention.

Messaging		The number of unread messages are displayed in the overview.
Review / Flags		With the use of review flags it is known when the patient was checked or who needs to check the patient. See more info in Patient details – Review / Flags.
Notifications		Notifications can be sent directly to the patient's device home screen.
User management		Only limited HCPs of an organization have the permission to add users. In case a HCP has this permission, they can add users by clicking the 'New' button, fill out the credentials, select
		the applicable role and press 'Register'. A mail will be sent to the newly created user to follow the account creation steps, see above .
Patient details – Overview		In the patient details overview all relevant (health) information of the individual patient can be retrieved. With the use of the previous and next patient buttons the HCP can easily switch to another patient. In the left side bar, the HCP can navigate between status, log, care plan, profile, documents and more.

Patient details – Status

firstname yb/first lastname last/initials 26719

Gender: Female, Age: 76 years, Date of birth: , Language: Dutch

Start consultation, Exercises, Medical questionnaires (75%), Consultations (54.55%), Questionnaires, Reported events, Not verified

TKA General Therapy

Compliance OK for today or yesterday (Green dot)
Non compliant > 2 days (Yellow dot)
Non compliant > 7 days (Red dot)

The reported healthcare data are retrievable in the 'Status' section. From there it is possible to navigate to the recurring and medical questionnaire results, exercises and many more.

A **GREEN**, **YELLOW** or **RED** status dot shows how compliant the patient is to the care plan, for example filling out all the suggested recurring and medical questionnaires and performing the suggested exercises will result in a **GREEN** status dot.

Patient details – Medical questionnaires

firstname yb/first lastname last/initials 26719

Start consultation, Exercises, Medical questionnaires (30%), Consultations, Questionnaires, Reported events, Not verified

THA General Therapy

Via the 'medical questionnaire' button it is shown:

- ✓ Which medical questionnaires are active for a patient.
- ✓ How many of these are completed or not.
- ✓ What the reported results are on the individual questions within the medical questionnaire.

What the sub scores and total scores are to be able to see a quick indication on how the patient is progressing.

Patient details – Recurring questionnaires

firstname yb/first lastname last/initials 26719

Start consultation, Exercises, Medical questionnaires (100%), Consultations (28.57%), Questionnaires

Questionnaires

Answers

Day	Reported at	DD-Pain(1) Indicate how much pain you had in your affected joint during the night (0=not at all, 10=worst)	DD-Pain(2) Indicate how much pain you had in your affected joint during the day (0=not at all, 10=worst)	DD-Pain(3) Did you experience any pain other than in the affected joint?	DD-Pain(4) Indicate where and when you experienced pain.
30	29 Nov 2023, 16:30:59	Yes	Yes	Yes	
29	28 Nov 2023, 17:08:38	Yes	Yes	Yes	
28	27 Nov 2023, 16:07:33	Yes	Yes	Yes	
27	26 Nov 2023, 16:32:18	Yes	Yes	Yes	
26	25 Nov 2023, 16:30:46	Yes	Yes	Yes	

With these patient reported answers it is possible to follow the progress of important health parameters frequently.

Patient details - Reported events

firstname yb/first lastname last/initials 26719




Start consultation, Exercises, Medical questionnaires (100%), Consultations (28.57%), Questionnaires




Reported events

Answers

Day	Reported at	DD-Pain(1) Indicate how much pain you had in your affected joint during the night (0=not at all, 10=worst)	DD-Pain(2) Indicate how much pain you had in your affected joint during the day (0=not at all, 10=worst)	DD-Pain(3) Did you experience any pain other than in the affected joint?	DD-Pain(4) Indicate where and when you experienced pain.
242	Jan 25, 2023, 9:26:45 AM	0.52	0.81	Yes	Night: Back Exercises: Call No pain: Grains

The reported events are categorized per event. For each reported event you can view: (1) The date when the event was reported (+day in care plan) (2) The answers to all questions asked to the patient.

Patient details - Review / Flags		1st line review is for the HCP that is not the lead clinician (e.g. physical therapist), 2nd line review is for the lead clinician (e.g. surgeon)
		Click validate when a patient has been monitored. Monitoring can mean: ✓ Questions answered - Advice given ✓ Adjusted or checked exercises ✓ Adjusted or checked step goal ✓ Checked daily parameters: pain, swelling, activity,... ✓ Checked postoperative protocol: use of crutches/walker, medication,... Same principle for 2nd line review.
		Flag someone from first line ✓ Patient needs to be checked by someone from first line (e.g. a physical therapist) ✓ Give a reason for flag ✓ Select the discipline to flag ✓ Optional: add an e-mail notification to flag. Flagged discipline will receive an email about the flag. A flag will add a new line in the patient log with the reason for flagging(see further).


		Flag someone from second line ✓ Patient needs to be checked by someone from second line (e.g. a surgeon) ✓ Give a reason to flag. ✓ Select the discipline to flag. Optional: add an e-mail notification to flag. The flagged discipline will receive an email about that flag.
		Report a technical issue.
Patient details - Log		Flags and technical updates automatically appear in the log of the patient. In case of flagging: ✓ GREEN : solved flag ✓ RED : Unsolved flag for discipline of logged in user (e.g. PT) ✓ YELLOW : Unsolved flag for other discipline than logged in user (e.g. Surgeon) Other important information can be added to the log manually. The applicable tags are added to the log in view of easy log searches

Patient details – Report

The screenshot shows the 'Patient report' page with a sidebar on the left containing tabs: Home, Patient, Reports, and Documents. The main content area is divided into several sections:

- Personal information:** Includes patient name, age, gender, and surgery date.
- Surgery information:** Includes surgery type, region, physical therapist, and insurance.
- Patient-reported outcomes:** A table showing scores for various metrics like 'Outcome of Life', 'Pain', 'Quality of Life', and 'Patient-reported outcomes' across different time points (Pre-op, 6W, 2M, 6M, 1Y, 2Y).
- Clinical milestones:** A table showing milestones such as 'Drop (days of use after surgery)', 'Pain (days of use after surgery)', 'Activity (days of use after surgery)', 'Patient compliance', and 'Patient-reported outcomes'.

This report is used to visualize the main milestones and medical questionnaires of the patient, as well as the progression of the patient.



The report can be downloaded in pdf form by clicking the 'Download report' button. This way the HCP can send the report to the patient. It's also possible to share the link directly to the report, using the button 'Copy report link'.

Patient details – Profile

The screenshot shows the 'Patient profile' page with a sidebar on the left containing tabs: Basic info, Clinical data, and Care plan. The main content area is divided into several sections:

- Basic information:** Includes patient name, first name, surname, and last name.
- Clinical data:** Includes patient age, gender, and surgery date.
- Care plan:** Includes patient insurance, region, physical therapist, and insurance.






The HCP can view, edit or add all relevant basic, clinical, administrative or care plan specific data in the patient profile

Patient details – Messaging

The screenshot shows the 'Patient messaging' page with a sidebar on the left containing tabs: Home, Patient, Reports, and Documents. The main content area is divided into several sections:

- Messages:** A list of messages sent to the patient, including the date and time of the message.
- Message popup:** A modal window that appears when a message is clicked, showing the message content and a 'Push notification' button.

The message popup (which can be enlarged) is closed by default and the number of new messages is visible in the floating button. This functionality is used to give advice to patients and answer their questions. It can also be used to ask for pictures, for example to follow-up on wound healing, swelling, redness etc...


		When clicking 'messages', a set of prescribed personal messages is shown wherefrom the HCP can select an applicable message for that particular patient.
		Click on the "1st line symbol" to mark the question as read.
		Messages that are submitted by a HCP can be edited or removed. When hovering over a message, a button to edit or delete the message will become visible.
Patient details -- Documents		Documents can be uploaded in the 'Documents' section. Chat images and video recordings are automatically collected in this section in view of easy access to these materials.
Patient details - More		See which documents / contracts the patient has signed. See on which device the patient is using the application and when the application was first and last seen. Export data is only possible for limited HCPs due to security reasons and to protect data integrity.
Log out		In the upper right corner the HCP can log out. Know that the automatic log out time is 14 days, therefore it is recommended to log out more frequently to avoid unauthorized access.
Support		Reach out via support@moveup.care for questions or use the '?' button

		on the platform.
Care plan management		In care plan management it is possible to navigate between multiple functions that are applicable to the indication of the patient, such as exercises, step goal, recurring questionnaires, video analysis, automation, etc... For more information see further.

HEALTHCARE PLATFORM CARE PLAN FUNCTIONS

It is possible to provide care through the use of the healthcare platform. With the functions below the HCP can adapt the default treatment plan of their individual patient.

Only the applicable functions will be displayed conform the care plan.

Exercise management		The exercise management module includes the following possibilities of adaptation: ✓ Adapt description, frequency, and intensity of the exercises of the day. ✓ Adapt the daily message. ✓ Adapt the exercise scheme: ○ Blue check marks indicate active exercises, accompanied with frequency and intensity of these active exercises. ○ Schedule exercises over a longer period of time. Extend the standard exercise scheme with a set of exercises that are not part of the exercises defined in the care plan
---------------------	--	--

Recurring questionnaire

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis

Questionnaires

Using treatment scheme **Personalize**

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis Show in Patient language No Advanced Stop active follow-up Complete

Questionnaires

Name	Start day	End day	Frequency (days)	Active from	Reminders	Order	Preview	Remove
Discharge	3	10	1	10:00	No reminder set last	1	Preview	Remove
Postoperative rescue	0	3	1	10:00	No reminder set last	2	Preview	Remove
Activities	15	100	1	10:00	No reminder set last	3	Preview	Remove
Exercises	1	100	1	10:00	No reminder set last	4	Preview	Remove

Add questionnaire to personalized care plan

Questionnaire **Add**

Start date End date

Frequency Daily

Activation hour 15 Activation minutes 0

Save Cancel

Recurring questionnaires can be removed, reordered or a new recurring questionnaire can be added. The default list of recurring questionnaires, including their start day, end day, frequency and reminders is predetermined in the care plan.

Medical questionnaire

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** Video analysis Change Reset Stop active follow-up Complete

Medical scores scheme

Name	Start day	End day	Preview	Remove
HIP PROM collection month 6	170	210	Preview	Remove
HIP PROM collection Year 1	340	500	Preview	Remove
Discharge Questionnaire HIP	0	15	Preview	Remove

Add medical questionnaire to care plan

Medical score **Add**

Select a medical score

Activation type Day in treatment

Start Date End Date

0 0

Save Cancel

Medical questionnaires can also be removed or added to the default scheme. The default list of medical questionnaires, including their start and end day is predetermined in the care plan.

Video analysis

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** **Video analysis** Show in Patient language No Advanced Stop active follow-up Complete

Video analysis

Analysis enabled for this patient

Results NOT shown to this patient

Add request

Jan 12, 2024

Gait analysis

not recorded

Video requests

Name	Start day	End day	Preview	Remove
Gait analysis	-7	-1	Preview	Remove
Gait analysis	21	30	Preview	Remove
Gait analysis	40	45	Preview	Remove

Add video request to personalized care plan

Video request **Add**

Select a video type

Start date End date

0 0

Save Cancel





Request the patient to record a video for analysis

Video type **Add**

Select a Video type

Save Cancel

With the use of video requests, it is possible to gain more objective information, for example about the gait analysis or range of motion of an operated joint. The default list of video requests is predetermined in the care plan. A video request can be added directly by using the 'Add request' button or a video request can be scheduled for a longer period of time by using the 'Add' button and by adding a start and end day to put the video request active.

Goal manager		The HCP can enable and disable a behavioral change tool which enables patients to set their personal goal, define barriers and solutions in order to reach a healthy lifestyle. The HCP can see the progress on the patient's goals that are defined by the patient.
Change of care plan		The care plan can be adapted based on the patient's interventions, needs and wishes. This can be done via the onboarding tool or via the healthcare platform via the button 'Advanced'. Instructions how to change a care plan are clearly indicated within the tools.
Stop active follow-up		When the active follow-up is finished (what falls under active follow-up is predetermined in the care plan), press 'Stop active follow-up'. This will redirect the patient to a passive long term follow-up if applicable (also predetermined in the care plan if long term follow-up is needed).
		When no actions are required anymore from the HCP or the patient, the patient can be redirected to the 'Completed' status.

WELLNESS MODULE

The wellness module is an optional non-medical device module and therefore the data shown is not intended to be used for clinical decision making.

This data gives an overview of the general wellbeing of the patient.

The wellness module by default is blurred and needs to be activated to check the patient's wellness data.

WELLNESS MODULE ACTIVATION

To activate the wellness module you need to enter your profile settings and activate the option:

- Enable wellness data

moveUP Patients Onboarding Config Admin Laurens EN

Edit Laurens Baes

First name Laurens	Address Oktrooi plein 1
Surname Baes	City Gent
Username laurens@moveup.care	Postal code 9000
Email laurens@moveup.care	Country Belgium
Telephone number +32496123456	Title Title
Social security number RIZIV nr	Image Image
Authentication ☒ Two factor authentication	Language English
	Wellness module ☒ Enable wellness module

Save Set new password Manage consent/contracts

Upon activation, a mandatory acknowledgment pop-up is displayed. It informs you that the data from the Wellness Module must not be used to support clinical decision-making.

Platform notice

Wellness module terms of use

Dear User,
You are about to activate the Wellness Module.

What is the Wellness Module?

The wellness module provides you with additional data about the general wellbeing and activity levels of your patients. The data in the Wellness Module can come from:

- imported data from general-purpose devices such as e.g. activity trackers
- manual entry by patients

This data does not come from medical devices. This means that neither the devices from which the data is imported, nor the manual entries by the patients have been verified or validated by moveUP, and therefore, it should not be used for clinical decision-making.

What is the purpose of this data?

We provide the option to see wellness data to give you a broader idea about the general wellbeing of the patients and to give you the opportunity to discuss and give further guidance regarding wellbeing.

This module is an optional module and is not linked to the medical data presented in the "Medical data" tab of the dashboard.

The data in the wellness module is not intended to be used for clinical decisions.

How to recognize the wellness data?

Wellness data is displayed in a separate tab in the dashboard or in a separate report.

Whenever wellness data is displayed in a column, together with other data it has a grey background to visually indicate that this does not concern medical data. Furthermore, it has a warning sign next to it.

Enabling/Disabling the Wellness Module

At any time, a user can enable/disable the Wellness Module in the user settings.

Acknowledgement

By clicking on the button below, you acknowledge that you have read and understood that the Wellness Module is an optional module that provides you with additional data regarding the general wellbeing of the patient. This data is not validated by moveUP and is not intended for clinical decision making.

☐ I agree


Sign

Cancel

Once acknowledged, you gain access to the wellness data. This activation allows you to view the wellness data and export the wellness reports.


WELLNESS MODULE OVERVIEW

Wellness module overview



In this tab you'll be able to see the step data from the patient and see their evolution.

Wellness report



You can see the wellness reports when needed.