


# User manual for the Health Care Provider (HCP)



A personalized treatment, with more insight and involvement in the care pathway.

App version : 2.1.0 MDD and later

Manual version : 17.1

Data of release manual : 02/12/2025







Link to e-learning: <https://moveup.talentlms.com/index>

# Table of Contents

<b>PRODUCT DESCRIPTION.....</b>	<b>3</b>
Instructions .....	3
Intended use.....	4
Intended users .....	4
Target population / indication for use .....	5
Claims .....	5
moveUP companion/coach/therapy .....	5
moveUP therapy .....	5
moveUP app & devices .....	6
Warnings and precautions .....	6
Contra-indications and potential side effects .....	7
Notification of serious incidents .....	7
Human body contact .....	7
Accessories/products used in combination .....	8
Device lifetime .....	8
Company details .....	8
<b>PATIENT APP .....</b>	<b>9</b>
<b>MEDICAL MODULE .....</b>	<b>9</b>
HEALTHCARE PLATFORM ACCOUNT CREATION.....	9
ONBOARDING OF PATIENTS .....	9
HEALTHCARE PLATFORM GENERAL FUNCTIONS .....	12
HEALTHCARE PLATFORM CARE PLAN FUNCTIONS .....	19
<b>WELLNESS MODULE.....</b>	<b>22</b>
WELLNESS MODULE ACTIVATION .....	22
WELLNESS MODULE OVERVIEW .....	23
Wellness module overview .....	24
Wellness report .....	24

# PRODUCT DESCRIPTION

## Instructions

	<p>Read carefully this instruction for use. Make sure that you know how to use moveUP, but also understand the limitations of moveUP before you start using our service.</p> <p>If you do have any questions on the use, contact moveUP.</p>
	<p>In the final chapter of this document you find our contacts, which you can use during the use of our service, in case you would have technical problems with the application (moveUP dashboard or patient app) or when expertise is required about a care path.</p>
	<p>The Healthcare practitioner can't give access to third parties to access the moveUP dashboard. Every other person can, even by accident, send out wrong information or adapt the patient profile, with the result that the advice sent by moveUP no longer corresponds with the needs of the patient.</p>
	<p>In case the application is use while moveUP &amp; b.clinic is not asked to do the follow-up of patients, the individual healthcare provider or the health care facility is responsible for the follow-up of the patient and the action that are carried out.</p>
	<p>In compliance with the medical device legislation, you are required to report any complaints towards the application, any alteration regarding characteristics/specifications of the application and incidents/complications potentially linked to the use of the application to support@moveup.care. The AFMPS can also be contacted in case of incident via <a href="https://www.famhp.be/en/human_use/health_products/medical_devices_accessories/materiovigilance/how_notify">https://www.famhp.be/en/human_use/health_products/medical_devices_accessories/materiovigilance/how_notify</a></p>
	<p>In compliance with the GDPR legislation, please report eventual data leaks linked to the use of moveUP.</p>
	<p>Please, read the entire privacy policy: <a href="https://www.moveup.care/legal/privacy-policy">https://www.moveup.care/legal/privacy-policy</a></p>

## Intended use

For healthcare providers, moveUP provides efficient clinical management of pathologies and treatments. With the use of moveUP, valuable insights are given to the HCP/care team of the status of the patient and its evolution.

For patients, moveUP provides personalised information and instructions to help them managing their symptoms and progress in their rehabilitation. The intensity of follow-up is adapted based on patients' needs and timepoint in the patient pathway:

### **moveUP companion = monitoring and information, no active follow-up**

moveUP companion offers targeted information and evolution reports to patients. Their care team takes the evolution reports into account in their further recovery when appropriate.

### **moveUP coach = active follow-up by healthcare team**

moveUP coach offers targeted information and evolution reports to patients. Their care team is more actively involved and takes the evolution reports into account and can provide advice and exercise suggestions through the digital platform.

### **moveUP therapy = active follow-up by healthcare team, with data driven validated care protocol**

moveUP therapy offers a data-driven validated care protocol with certain category and level of exercises and activities, specifically targeted to the individual patient. Their care team can manually adapt the data driven validated protocol when needed. Patients can fully rehabilitation with moveUP without leaving their home environment.

\* Optional functionality that can be enabled: interoperability with Class IIa continues passive motion (CPM) medical device. For knee and hip patients who are using a Class IIa CPM medical device, moveUP can interoperate with the Class IIa CPM software. moveUP acts as a facilitator to easy assign a designated CPM exercise protocol chosen by a physician to a patient and to display the performed CPM exercises in the medical dashboard

## Intended users

moveUP is intended to be used by patients and healthcare providers. The main user of the mobile app and patient website is the patient.

Inclusion:

- Age: minimum 18 years / maximum no limit
- Health & condition: capable of performing basic activities of daily living
- Language: understanding one of the available languages of the app (Dutch, French, German, English)

## Exclusion:

- Patients who are mentally incompetent or having troubles to express what they are feeling (for instance, mentally diseased people, people staying in elderly care centres, ...) are excluded.
- Patients who are not capable of operating a tablet/smartphone and activity tracker.
- Patients who can't understand one of the available languages of the app (Dutch, French, German, English)

The main user of the medical web interface is a healthcare practitioner (group) or clinical researcher (group), named the care team. The care team is able to operate a web interface via web browser on PC/tablet/smartphone. The healthcare practitioner needs to understand one of the available languages of the web interface (today only available in English)

## Target population / indication for use

moveUP companion & coach is used by Musculo-skeletal, oncologic, respiratory, gastro-intestinal, cardiovascular, and neurologic patients, such as patients:

- who underwent or planning a gastric bypass or gastric sleeve operation
- who have or had a stroke, multiple myeloma, covid, familial hypercholesterolemia
- who have or had back or joint problems or operations
- moveUP therapy is used by knee & hip arthroplasty patients.
- moveUP is used for rehabilitation but is possible to include a patient before an intervention to know the previous state of health and establish goals for the rehabilitation.

## Claims

moveUP companion/coach/therapy

The intensity of follow-up is adapted based on the needs of the patient, via the moveUP Symptom & QoL monitoring tool.

moveUP enhances the clinical management of the patients, because early detection & management of complications is possible via the symptom & QoL monitoring tool

More efficient clinical management, such as the number of consultations, can be reduced

Enforces therapy compliance/adherence

The correct information is provided at the right time

moveUP therapy

With the use of moveUP therapy knee & hip arthroplasty patients can fully rehabilitate via the in-app care team without leaving their home environment

## moveUP app & devices

The moveUP app runs on:

- iPhones/iPads (iOS 11.0 or newer)
- Android systems (Android 7.0.0 and newer).

If the patient has no compatible tablet/smartphone, moveUP can lease these devices to the patient.

- Garmin Vivofit 3 and Vivofit 4
- Nokia Withings Go
- All the devices compatible with HealthKit (IOS)
- All the devices compatible with Google Fit (IOS and Android)

The moveUP application needs an internet connection when in use.





The moveUP application is compatible with the following browsers:


- Chrome
- Safari
- Edge
- Firefox

We strongly recommend to use the latest version of the browser.

For your data safety please enable the 2-factor authentication(2FA) and use a safe network

## Warnings and precautions

	<ul style="list-style-type: none"><li>• moveUP is not an emergency tool. In case of emergency please contact your doctor or call 112.</li></ul>
	<ul style="list-style-type: none"><li>• moveUP is not intended to monitor/treat vital parameters of critical diseases.</li></ul>
	<ul style="list-style-type: none"><li>• moveUP companion is not a replacement of any treatment you need to follow, but is used as an addition.</li></ul>
	<ul style="list-style-type: none"><li>• Patients need to be aware if questionnaires are not filled out on frequent basis or not filled out trustworthy, the care team has more difficulty to know their health situation &amp; evolution.</li><li>• The moveUP App offers personalized healthcare based on the patient's profile and therefore can not be shared, borrowed or traded between users. The information and advice of the moveUP application do not apply to a patient other than the one whose profile is stored in the application. The patient may not grant access to the moveUP App to third parties. Any</li></ul>

	<p>other person can, even accidentally, send wrong information or change the profile of the patient, with the result that the moveUP-guided advice no longer corresponds to the needs of the patient</p> <ul style="list-style-type: none"> <li>• With software, there is always a residual risk for bugs. If you have the feeling something is wrong or you don't get a response in the app, please reach out to us via support@ or call 0800 88 008.</li> </ul>
---	---

The moveUP IFU are only supplied in electronic format. If you need a full printed version please contact us at the e-mail address [info@moveup.care](mailto:info@moveup.care). This mail is available 24/7.

### Contra-indications and potential side effects

There are no contra-indications or known side effects

### Notification of serious incidents

A serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or another person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health that resulted in any of the following:
  - life-threatening illness or injury,
  - permanent impairment of a body structure or a body function, hospitalisation
  - prolongation of patient hospitalisation, medical or surgical intervention to prevent life-threatening illness or injury
  - permanent impairment to a body structure or a body function,
  - chronic disease,
- a serious public health threat;

If you detect any of these cases please send an email to [info@moveUP.care](mailto:info@moveUP.care) and your competent authority:

- Netherlands: [meldpunt@igj.nl](mailto:meldpunt@igj.nl)
- Belgium: [meddev@fagg.be](mailto:meddev@fagg.be)
- France: [dedim.ugsv@ansm.sante.fr](mailto:dedim.ugsv@ansm.sante.fr)
- Germany: [zlg@zlg.nrw.de](mailto:zlg@zlg.nrw.de)

### Human body contact

No human body contact with the patient or user, due to the nature of the product (software).

### Accessories/products used in combination

There are no accessories. If the patient has no compatible tablet/smartphone or activity tracker, moveUP can lease these devices to the patient. The leased devices are CE marked devices that meets the compatibility criteria outlined in the IFU for the app.

### Device lifetime

2-years, depending on the willingness of the user to update the app. We sent out a notification with the recommendation to update the app, to make sure the app will keep the performance & recently included features.

### Company details

 moveUP NV/SA	Tel: +32 800 88 008	UDI: ++G166MOVEUP210MDDDD0  (01) G166G166MOVEUP210MDDDD0 (10) ++\$+210MDD (11) /16D20251202 (21) Q
Head office Brussel: Cantersteen 47, B-1000 Brussels, Belgium  Secondary office: Kliniekstraat 27A 9050; Ghent, Belgium	Email: info@moveUP.care www.moveup.care	

If you have any questions about the App or the tablet, you can contact moveUP via the following channels - in order of priority: Via the message system of the App. Via email to support@moveup.care - Via phone during office hours: +32 800 88 008

MoveUP guarantees delivery of the eIFU within seven (7) days upon request. Please contact info@moveUP.care if you want to obtain a free copy.





# PATIENT APP

Please read the user manual for patients on [www.moveup.care/appinfo](http://www.moveup.care/appinfo) to view the functionalities of the patient app.

## MEDICAL MODULE

### HEALTHCARE PLATFORM ACCOUNT CREATION





<p>Step 1: The HCP registers on the healthcare platform on <a href="https://v1.medical.moveup.care/">https://v1.medical.moveup.care/</a></p>	
<p>Step 2: moveUP will contact the HCP to add his/her profile to the correct organization.</p>	
<p>Step 3: Security measures: It is important to enable 2FA as you have access to sensitive health data.</p>	
<p>E-learning: For HCPs who will provide care through the use of the healthcare platform it is recommended to follow the e-learning courses for the indications listed on <a href="https://moveup.talentlms.com">https://moveup.talentlms.com</a>. We ask the HCPs to register for the e-learning course in their applicable language.</p>	

### ONBOARDING OF PATIENTS

Onboarding is the process of registering a new patient in the healthcare platform.

A crucial part is adding a care plan to the patient's profile. A care plan determines what will be visible for the patient in the application, e.g.: type of exercises, information modules, message functionality, etc. A care plan usually consists out of an intervention, a treatment and (a) care





team(s).





<p><b>Onboarding patients</b></p>		<p>The onboarding dashboard is easy retrievable in the health care platform in the upper left corner when clicking 'Onboarding'. It is also accessible through <a href="https://patient-onboarding-admin-moveup.care">https://patient-onboarding-admin-moveup.care</a> using the HCP credentials of the health care platform.</p> <p>Click on 'new patient' to onboard a new patient and fill out the patient information.</p> <p>After the patient information is filled out, press 'Save'.</p>
	<p>New patient</p> 	<p>Important note: Only one account can be coupled to an e-mail address. It is not possible to reuse an e-mail address that is already known within our database, e.g.: reusing the same e-mail address for the account of a partner or parent.</p>
	<p>Patient's care plans</p> 	<p>Next up, the care plan needs to be selected: press 'Start new care plan'. Drag and drop an intervention and treatment to the care plan and fill out the dates, name of the treating doctor and care team.</p> <p>Drag and drop a study if applicable. Click Save.</p>
	<p>Email variables</p> <pre> {{ firstname }} {{ lastname }} {{ code }} {{ registrationlink }} {{ organizationName }} {{ surgeonName }} {{ hosp_spec_par_start_NL }} {{ organizationName }} {{ bodyStructure }} </pre> 	<p>A mail template is shown. The HCP reviews the e-mail variables and clicks on 'send selected mails'.</p> <p>An invitation mail will be sent to the patient with all the needed information to install and activate the application.</p>



## HEALTHCARE PLATFORM GENERAL FUNCTIONS

Below the HCP can retrieve all functions of the healthcare platform. The applicable functions will be shown based on the HCPs subscription. The extent of the available patient information is based on the patients subscription & his/her compliancy rate.

<p><b>Patient overview</b></p>		<p>The health care platform gives a clear overview of the HCPs patients.</p> <p>In the top bar the relevant categories are shown that identify the patient and show relevant health parameters, unread messages, the lead clinician and patient compliance. By default, when opening the healthcare platform, the patient overview is sorted based on intervention date.</p>
<p><b>Tags</b></p>		<p>Select patients related to a specific pathology by using tags.</p>
<p><b>Filters</b></p>		<p>Select patients based on language, test patient, lead clinician or being activated or not. Activated can mean 2 things: either the patient registered him/herself with a code through the app or he has been marked as immediately activated in onboarding. Non activated: patients who haven't registered themselves through the app yet and haven't been marked as immediately activated.</p>
<p><b>Color codes</b></p>		<p>Color codes indicate which patients need closer attention.</p>

Messaging		The number of unread messages are displayed in the overview.
Review / Flags		With the use of review flags it is known when the patient was checked or who needs to check the patient. See more info in Patient details – Review / Flags.
Notifications		Notifications can be sent directly to the patient's device home screen.
User management		Only limited HCPs of an organization have the permission to add users. In case a HCP has this permission, they can add users by clicking the 'New' button, fill out the credentials, select
		the applicable role and press 'Register'. A mail will be sent to the newly created user to follow the account creation steps, see <a href="#">above</a> .
Patient details – Overview		In the patient details overview all relevant (health) information of the individual patient can be retrieved. With the use of the previous and next patient buttons the HCP can easily switch to another patient. In the left side bar, the HCP can navigate between status, log, care plan, profile, documents and more.

## Patient details – Status

The reported healthcare data are retrievable in the 'Status' section. From there it is possible to navigate to the recurring and medical questionnaire results, exercises and many more.

A **GREEN**, **YELLOW** or **RED** status dot shows how compliant the patient is to the care plan, for example filling out all the suggested recurring and medical questionnaires and performing the suggested exercises will result in a **GREEN** status dot.

## Patient details – Medical questionnaires

Via the 'medical questionnaire' button it is shown:

- ✓ Which medical questionnaires are active for a patient.
- ✓ How many of these are completed or not.
- ✓ What the reported results are on the individual questions within the medical questionnaire.




What the sub scores and total scores are to be able to see a quick indication on how the patient is progressing.




## Patient details – Recurring questionnaires

With these patient reported answers it is possible to follow the progress of important health parameters frequently.

## Patient details - Reported events

The reported events are categorized per event. For each reported event you can view: (1) The date when the event was reported (+day in care plan) (2) The answers to all questions asked to the patient.


<b>Patient details</b> <b>- Review / Flags</b>		<p>1st line review is for the HCP that is not the lead clinician (e.g. physical therapist), 2<sup>nd</sup> line review is for the lead clinician (e.g. surgeon)</p>
		<p>Click validate when a patient has been monitored. Monitoring can mean:</p> <ul style="list-style-type: none"> <li>✓ Questions answered - Advice given</li> <li>✓ Adjusted or checked exercises</li> <li>✓ Adjusted or checked step goal</li> <li>✓ Checked daily parameters: pain, swelling, activity,...</li> <li>✓ Checked postoperative protocol: use of crutches/walker, medication,...</li> </ul> <p>Same principle for 2<sup>nd</sup> line review.</p>
		<p>Flag someone from first line</p> <ul style="list-style-type: none"> <li>✓ Patient needs to be checked by someone from first line (e.g. a physical therapist)</li> <li>✓ Give a reason for flag</li> <li>✓ Select the discipline to flag</li> <li>✓ Optional: add an e-mail notification to flag. Flagged discipline will receive an email about the flag.</li> </ul> <p>A flag will add a new line in the patient log with the reason for flagging(see <a href="#">further</a>).</p>

		<p>Flag someone from second line</p> <ul style="list-style-type: none"> <li>✓ Patient needs to be checked by someone from second line (e.g. a surgeon)</li> <li>✓ Give a reason to flag.</li> <li>✓ Select the discipline to flag.</li> </ul> <p>Optional: add an e-mail notification to flag. The flagged discipline will receive an email about that flag.</p>
		<p>Report a technical issue.</p>
<p>Patient details - Log</p>		<p>Flags and technical updates automatically appear in the log of the patient.</p> <p>In case of flagging:</p> <ul style="list-style-type: none"> <li>✓ <b>GREEN</b> : solved flag</li> <li>✓ <b>RED</b> : Unsolved flag for discipline of logged in user (e.g. PT)</li> <li>✓ <b>YELLOW</b> : Unsolved flag for other discipline than logged in user (e.g. Surgeon)</li> </ul> <p>Other important information can be added to the log manually. The applicable tags are added to the log in view of easy log searches</p>



## Patient details – Report

This report is used to visualize the main milestones and medical questionnaires of the patient, as well as the progression of the patient.










The report can be downloaded in pdf form by clicking the 'Download report' button. This way the HCP can send the report to the patient. It's also possible to share the link directly to the report, using the button 'Copy report link'.


## Patient details – Profile

The HCP can view, edit or add all relevant basic, clinical, administrative or care plan specific data in the patient profile

## Patient details – Messaging

The message popup (which can be enlarged) is closed by default and the number of new messages is visible in the floating button. This functionality is used to give advice to patients and answer their questions. It can also be used to ask for pictures, for example to follow-up on wound healing, swelling, redness etc...


		<p>When clicking 'messages', a set of prescribed personal messages is shown wherefrom the HCP can select an applicable message for that particular patient.</p>
		<p>Click on the "1st line symbol" to mark the question as read.</p>
		<p>Messages that are submitted by a HCP can be edited or removed. When hovering over a message, a button to edit or delete the message will become visible.</p>
<b>Patient details</b> <b>-- Documents</b>		<p>Documents can be uploaded in the 'Documents' section. Chat images and video recordings are automatically collected in this section in view of easy access to these materials.</p>
<b>Patient details</b> <b>- More</b>		<p>See which documents / contracts the patient has signed.</p> <p>See on which device the patient is using the application and when the application was first and last seen.</p> <p>Export data is only possible for limited HCPs due to security reasons and to protect data integrity.</p>
<b>Log out</b>		<p>In the upper right corner the HCP can log out. Know that the automatic log out time is 14 days, therefore it is recommended to log out more frequently to avoid unauthorized access.</p>
<b>Support</b>		<p>Reach out via <a href="mailto:support@moveup.care">support@moveup.care</a> for questions or use the '?' button</p>

		on the platform.
Care plan management		<p>In care plan management it is possible to navigate between multiple functions that are applicable to the indication of the patient, such as exercises, step goal, recurring questionnaires, video analysis, automation, etc... For more information see <a href="#">further</a>.</p>

## HEALTHCARE PLATFORM CARE PLAN FUNCTIONS

It is possible to provide care through the use of the healthcare platform. With the functions below the HCP can adapt the default treatment plan of their individual patient.

Only the applicable functions will be displayed conform the care plan.

Exercise management		<p>The exercise management module includes the following possibilities of adaptation:</p> <ul style="list-style-type: none"> <li>✓ Adapt description, frequency, and intensity of the exercises of the day.</li> <li>✓ Adapt the daily message.</li> <li>✓ Adapt the exercise scheme: <ul style="list-style-type: none"> <li>○ Blue check marks indicate active exercises, accompanied with frequency and intensity of these active exercises.</li> <li>○ Schedule exercises over a longer period of time.</li> </ul> </li> </ul> <p>Extend the standard exercise scheme with a set of exercises that are not part of the exercises defined in the care plan</p>
---------------------	--	--

## Recurring questionnaire

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis

Questionnaires

Using treatment scheme **Personalize**

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis Show in Patient language: No # Advanced Stop active follow-up Complete

Questionnaires

Name	Start day	End day	Frequency (times)	Active from	Reminders	Order	Preview	Remove
Discharge	3	10	1	10:00	No reminder set last		Preview	Remove
Postoperative reassess	0	3	1	10:00	No reminder set last		Preview	Remove
Activities	15	100	1	10:00	No reminder set last		Preview	Remove
Exercises	1	100	1	10:00	No reminder set last		Preview	Remove

Add questionnaire to personalized care plan

Questionnaire

Start date

End date

Frequency

Daily

Activation hour

15

Activation minutes

0

Save Cancel

Recurring questionnaires can be removed, reordered or a new recurring questionnaire can be added. The default list of recurring questionnaires, including their start day, end day, frequency and reminders is predetermined in the care plan.

## Medical questionnaire

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** Video analysis Change Reset Stop active follow-up Complete

Medical scores scheme

Name	Start day	End day	Preview	Remove
HIP PROM collection month 6	170	210	Preview	Remove
HIP PROM collection Year 1	340	500	Preview	Remove
Discharge Questionnaire HIP	0	15	Preview	Remove

Add medical questionnaire to care plan

Medical score

Select a medical score

Activation type

Day in treatment

Start Date

0

End Date

0

Save Cancel

Medical questionnaires can also be removed or added to the default scheme. The default list of medical questionnaires, including their start and end day is predetermined in the care plan.

## Video analysis

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** **Video analysis** Show in Patient language: No # Advanced Stop active follow-up Complete

Video analysis

Analysis enabled for this patient

Results NOT shown to this patient

Add request

Jan 12, 2024

Gait analysis

not recorded

Video requests

Name	Start day	End day	Preview	Remove
Gait analysis	-7	-1	Preview	Remove
Gait analysis	21	30	Preview	Remove
Gait analysis	40	45	Preview	Remove

Add video request to personalized care plan

Video request

Select a video type

Start date

0

End date

0

Save Cancel





Request the patient to record a video for analysis

Video type

Select a Video type

Save Cancel

With the use of video requests, it is possible to gain more objective information, for example about the gait analysis or range of motion of an operated joint. The default list of video requests is predetermined in the care plan. A video request can be added directly by using the 'Add request' button or a video request can be scheduled for a longer period of time by using the 'Add' button and by adding a start and end day to put the video request active.

<p>Goal manager</p>		<p>The HCP can enable and disable a behavioral change tool which enables patients to set their personal goal, define barriers and solutions in order to reach a healthy lifestyle. The HCP can see the progress on the patient's goals that are defined by the patient.</p>
<p>Change of care plan</p>		<p>The care plan can be adapted based on the patient's interventions, needs and wishes. This can be done via the onboarding tool or via the healthcare platform via the button 'Advanced'. Instructions how to change a care plan are clearly indicated within the tools.</p>
<p>Stop active follow-up</p>		<p>When the active follow-up is finished (what falls under active follow-up is predetermined in the care plan), press 'Stop active follow-up'. This will redirect the patient to a passive long term follow-up if applicable (also predetermined in the care plan if long term follow-up is needed).</p>
		<p>When no actions are required anymore from the HCP or the patient, the patient can be redirected to the 'Completed' status.</p>

# WELLNESS MODULE

The wellness module is an optional non-medical device module and therefore the data shown is not intended to be used for clinical decision making.

This data gives an overview of the general wellbeing of the patient.

The wellness module by default is blurred and needs to be activated to check the patient's wellness data.

## WELLNESS MODULE ACTIVATION

To activate the wellness module you need to enter your profile settings and activate the option:

- Enable wellness data

moveUP Patients Onboarding Config Admin Laurens EN

### Edit Laurens Baes

First name Laurens	Address Oktrooi plein 1
Surname Baes	City Gent
Username laurens@moveup.care	Postal code 9000
Email laurens@moveup.care	Country Belgium
Telephone number +32496123456	Title Title
Social security number Rizin nr	Image Image
Rizin nr Rizin nr	Language English
Authentication <input checked="" type="checkbox"/> Two factor authentication	Wellness module <input checked="" type="checkbox"/> Enable wellness module

Save Set new password Manage consent/contracts

Upon activation, a mandatory acknowledgment pop-up is displayed. It informs you that the data from the Wellness Module must not be used to support clinical decision-making.

## Platform notice

### Wellness module terms of use

Dear User,  
You are about to activate the Wellness Module.

### What is the Wellness Module?

The wellness module provides you with additional data about the general wellbeing and activity levels of your patients. The data in the Wellness Module can come from:

- imported data from general-purpose devices such as e.g. activity trackers
- manual entry by patients

**This data does not come from medical devices. This means that neither the devices from which the data is imported, nor the manual entries by the patients have been verified or validated by moveUP, and therefore, it should not be used for clinical decision-making.**

### What is the purpose of this data?

We provide the option to see wellness data to give you a broader idea about the general wellbeing of the patients and to give you the opportunity to discuss and give further guidance regarding wellbeing.

This module is an optional module and is not linked to the medical data presented in the "Medical data" tab of the dashboard.

**The data in the wellness module is not intended to be used for clinical decisions.**

### How to recognize the wellness data?

Wellness data is displayed in a separate tab in the dashboard or in a separate report.

Whenever wellness data is displayed in a column, together with other data it has a grey background to visually indicate that this does not concern medical data. Furthermore, it has a warning sign next to it.

### Enabling/Disabling the Wellness Module

At any time, a user can enable/disable the Wellness Module in the user settings.

### Acknowledgement

By clicking on the button below, you acknowledge that you have read and understood that the Wellness Module is an optional module that provides you with additional data regarding the general wellbeing of the patient. This data is not validated by moveUP and is not intended for clinical decision making.

☐ I agree


Sign

Cancel

Once acknowledged, you gain access to the wellness data. This activation allows you to view the wellness data and export the wellness reports.


## WELLNESS MODULE OVERVIEW

## Wellness module overview



In this tab you'll be able to see the step data from the patient and see their evolution.

## Wellness report



You can see the wellness reports when needed.