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






	<div><div><div><div><div><div>← Medical questionnaires</div><div>Questionnaires</div></div></div><div><div><div>Satisfaction hospital stay</div><div>Display: Satisfaction with your hospital stay</div><div>Preview</div></div><div><div><div>Forgotten Joint Score Heup</div><div>Display: Forgotten Joint Score - Hip</div><div>Preview</div></div><div><div><div>EQ-5D</div><div>Display: Health Questionnaire</div><div>Patient filed: 04/01/2024, 11:13</div><div>Preview</div></div><div><div><div>UCLA Activity Score</div><div>Display: UCLA Activity Score</div><div>Patient filed: 04/01/2024, 12:47</div><div>Preview</div></div></div></div></div><div><div>EQ-5D (NL-FR-EN-DE)</div><div><div>Scores</div><table><thead><tr><th>Name</th><th>Result</th></tr></thead><tbody><tr><td>Mobility</td><td>4</td></tr><tr><td>Self-care</td><td>3</td></tr><tr><td>Usual activities</td><td>4</td></tr><tr><td>Pain / discomfort</td><td>4</td></tr><tr><td>Anxiety / depression</td><td>1</td></tr><tr><td>VAS</td><td>52</td></tr></tbody></table><div>Answers</div><table><thead><tr><th>ID</th><th>Questions</th><th>Answers</th></tr></thead><tbody><tr><td></td><td>Day in care plan</td><td>-8</td></tr><tr><td></td><td>Started</td><td>4 Jan 2024, 11:13:22</td></tr><tr><td></td><td>Duration</td><td>0 hours 1 minutes 50 seconds</td></tr><tr><td>EQN00</td><td>MOBILITY</td><td>I have severe problems in walking about</td></tr><tr><td>EQI01</td><td>SELF-CARE</td><td>I have moderate problems washing or dressing myself</td></tr><tr><td>EQI02</td><td>USUAL ACTIVITIES (e.g. work, study, housework, family or leisure activities)</td><td>I have severe problems doing my usual activities</td></tr><tr><td>EQP03</td><td>PAIN / DISCOMFORT</td><td>I have severe pain or discomfort</td></tr></tbody></table></div></div></div></div></div></div>	Name	Result	Mobility	4	Self-care	3	Usual activities	4	Pain / discomfort	4	Anxiety / depression	1	VAS	52	ID	Questions	Answers		Day in care plan	-8		Started	4 Jan 2024, 11:13:22		Duration	0 hours 1 minutes 50 seconds	EQN00	MOBILITY	I have severe problems in walking about	EQI01	SELF-CARE	I have moderate problems washing or dressing myself	EQI02	USUAL ACTIVITIES (e.g. work, study, housework, family or leisure activities)	I have severe problems doing my usual activities	EQP03	PAIN / DISCOMFORT	I have severe pain or discomfort	<p>Gesamtpunktzahlen, um schnell erkennen zu können, wie der Patient vorankommt.</p>
Name	Result																																							
Mobility	4																																							
Self-care	3																																							
Usual activities	4																																							
Pain / discomfort	4																																							
Anxiety / depression	1																																							
VAS	52																																							
ID	Questions	Answers																																						
	Day in care plan	-8																																						
	Started	4 Jan 2024, 11:13:22																																						
	Duration	0 hours 1 minutes 50 seconds																																						
EQN00	MOBILITY	I have severe problems in walking about																																						
EQI01	SELF-CARE	I have moderate problems washing or dressing myself																																						
EQI02	USUAL ACTIVITIES (e.g. work, study, housework, family or leisure activities)	I have severe problems doing my usual activities																																						
EQP03	PAIN / DISCOMFORT	I have severe pain or discomfort																																						
<p>Patientendaten – Wiederkehrende Fragebögen</p>	<div><div><div><div><div>7</div><div>100%</div><div>28.57%</div></div><div><div>Start consultation</div><div>Exercises</div><div>Medical questionnaires</div><div>Consultations</div><div>Success</div></div></div><div><div>Questionnaires</div><table><thead><tr><th>Questionnaire</th><th>Answers</th><th>Status</th></tr></thead><tbody><tr><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td><td>Yes</td></tr><tr><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td><td>Yes</td></tr><tr><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td><td>Yes</td></tr><tr><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td><td>Yes</td></tr><tr><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td><td>Yes</td></tr><tr><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td><td>Yes</td></tr></tbody></table></div></div></div>	Questionnaire	Answers	Status	EQ-5D (NL-FR-EN-DE)	52	Yes	EQ-5D (NL-FR-EN-DE)	52	Yes	EQ-5D (NL-FR-EN-DE)	52	Yes	EQ-5D (NL-FR-EN-DE)	52	Yes	EQ-5D (NL-FR-EN-DE)	52	Yes	EQ-5D (NL-FR-EN-DE)	52	Yes	<p>Mit diesen Patientenantworten ist es möglich, den Verlauf wichtiger Gesundheitsparameter regelmäßig zu verfolgen.</p>																	
Questionnaire	Answers	Status																																						
EQ-5D (NL-FR-EN-DE)	52	Yes																																						
EQ-5D (NL-FR-EN-DE)	52	Yes																																						
EQ-5D (NL-FR-EN-DE)	52	Yes																																						
EQ-5D (NL-FR-EN-DE)	52	Yes																																						
EQ-5D (NL-FR-EN-DE)	52	Yes																																						
EQ-5D (NL-FR-EN-DE)	52	Yes																																						
<p>Patientendaten – Gemeldete Ereignisse</p>	<div><div><div><div><div>Reported events</div><div>Answers</div></div><table><thead><tr><th>Date</th><th>Event</th><th>Answer</th></tr></thead><tbody><tr><td>04 Jan 2024, 11:13:22</td><td>EQ-5D (NL-FR-EN-DE)</td><td>52</td></tr></tbody></table></div></div></div>	Date	Event	Answer	04 Jan 2024, 11:13:22	EQ-5D (NL-FR-EN-DE)	52	<p>Die gemeldeten Ereignisse werden nach Ereignis kategorisiert. Für jedes gemeldete Ereignis können Sie Folgendes anzeigen: (1) Das Datum der Meldung des Ereignisses (+ Behandlungsplan) (2) Die Antworten auf alle dem Patienten gestellten Fragen.</p>																																
Date	Event	Answer																																						
04 Jan 2024, 11:13:22	EQ-5D (NL-FR-EN-DE)	52																																						
<p>Patientendaten – Überprüfung/Markierungen</p>	<div><div><div><div><div>1st line review</div><div>2nd line review</div><div>No technical issue</div></div><div><div>Date of birth</div><div>Care path</div></div></div></div></div>	<p>Die Erstlinienprüfung ist für den Arzt bestimmt, der nicht der leitende Kliniker ist (z. B. Physiotherapeut), die Zweitlinienprüfung ist für den leitenden Kliniker (z. B. Chirurg).</p>																																						

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Patientendaten – Profil







The screenshot shows the 'Basic information' section of the patient profile. It includes fields for 'First Name' and 'Surname'. Below these are sections for 'National health insurance', 'Corona information', and 'Bariatrics information'. The 'First Name' field contains the text 'firstname.znetao9srk' and the 'Surname' field contains 'lastname.ix0nkv25fh'.

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
HEALTHCARE PLATFORM CARE PLAN FUNCTIONS

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The screenshot displays the 'Exercise Management' interface. At the top, there's a 'Manage patient scheme' table with columns for exercises (Icepacks Hip, Circulation, Heel raises, Glutes contraction) and days (Day 1, Day 2, Day 3). The table shows active exercises with blue checkmarks and durations (e.g., 6x 20 minutes). Below the table, there's a 'Plan exercise' modal for 'Glutes contraction' on 'Day 5'. The modal includes fields for 'Start day / End day' (5 / 5), 'Repeat' (Daily), 'Intensity start' (10), 'Intensity target' (10), 'Frequency start' (3), and 'Frequency target' (3). Buttons for 'Cancel', 'Apply', and 'Apply & save' are at the bottom.

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Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis

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Using treatment scheme **Personalize**

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis Show in Patient language: No Advanced + Stop active follow-up Complete

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Name	Start day	End day	Frequency (Week)	Active from	Reminders	Order	Preview	Remove
Discharge	9	10	1	10:00	Reminders not set		Preview	Remove
Postoperative recovery	0	3	1	10:00	Reminders not set		Preview	Remove
Activities	15	100	1	10:00	Reminders not set		Preview	Remove
Exercises	5	100	1	10:00	Reminders not set		Preview	Remove

Add questionnaire to personalized care plan

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Exercises Step goal **Recurring questionnaires** **Medical questionnaires** Video analysis Change Reset Stop active follow-up Complete

Medical scores scheme

Name	Start day	End day	Preview	Remove
HSP PROM collection month 6	170	250	Preview	Remove
HSP PROM collection Year 1	340	500	Preview	Remove
Discharge Questionnaire HSP	0	15	Preview	Remove

Add medical questionnaire to care plan

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Videoanalyse:

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** **Video analysis** Show in Patient language: No Advanced + Stop active follow-up Complete

Video analysis

Analysis enabled for this patient

Results NOT shown to this patient

Add request

Jan 12, 2024

Get analysis

not recorded

Video requests

Name	Start day	End day	Preview	Remove
Get analysis	-7	-1	Preview	Remove
Get analysis	21	30	Preview	Remove
Get analysis	40	45	Preview	Remove

Add video request to personalized care plan

Video request

Select a video type

Start date

0

End date

0

Save Cancel





Request the patient to record a video for analysis

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The screenshot shows the 'Edit Laurens Baes' page in the moveUP system. The page has a top navigation bar with 'moveUP', 'Patients', and 'Onboarding'. On the right, there are links for 'Config', 'Admin', and a user profile for 'Laurens'. The main content area is titled 'Edit Laurens Baes' and contains a form with two columns of fields. The left column includes 'First name' (Laurens), 'Surname' (Baes), 'Username' (laurens@moveup.care), 'Email' (laurens@moveup.care), 'Telephone number' (+32496123456), 'Social security number' (Social security number), and 'Ridv nr' (Ridv nr). The right column includes 'Address' (Oktrooi plein 1), 'City' (Gent), 'Postal code' (9000), 'Country' (Belgium), 'Title' (Title), 'Image' (Image), and 'Language' (English). At the bottom of the form, there are two toggle switches: 'Authentication' with 'Two factor authentication' (turned on) and 'Wellness module' with 'Enable wellness module' (turned off). The 'Enable wellness module' toggle is highlighted with a green circle. Below the form, there are three buttons: 'Save', 'Set new password', and 'Manage consent/contracts'.

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Platform notice

Wellness module terms of use

Dear User,
You are about to activate the Wellness Module.

What is the Wellness Module?

The wellness module provides you with additional data about the general wellbeing and activity levels of your patients. The data in the Wellness Module can come from:

- imported data from general-purpose devices such as e.g. activity trackers
- manual entry by patients

This data does not come from medical devices. This means that neither the devices from which the data is imported, nor the manual entries by the patients have been verified or validated by moveUP, and therefore, it should not be used for clinical decision-making.

What is the purpose of this data?

We provide the option to see wellness data to give you a broader idea about the general wellbeing of the patients and to give you the opportunity to discuss and give further guidance regarding wellbeing.

This module is an optional module and is not linked to the medical data presented in the "Medical data" tab of the dashboard.

The data in the wellness module is not intended to be used for clinical decisions.

How to recognize the wellness data?

Wellness data is displayed in a separate tab in the dashboard or in a separate report.

Whenever wellness data is displayed in a column, together with other data it has a grey background to visually indicate that this does not concern medical data. Furthermore, it has a warning sign next to it.

Enabling/Disabling the Wellness Module

At any time, a user can enable/disable the Wellness Module in the user settings.

Acknowledgement

By clicking on the button below, you acknowledge that you have read and understood that the Wellness Module is an optional module that provides you with additional data regarding the general wellbeing of the patient. This data is not validated by moveUP and is not intended for clinical decision making.

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