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Patientendetails – Bericht

The screenshot shows a web interface with a sidebar on the left containing icons for Home, Reports, LMS, Recent, Calendar, Profile, and Documents. The main content area is titled 'Patient's report - 17/09/2022'. It features a 'View report' button, which is highlighted with a red rectangle. Below this button is a link 'Print all (2 documents report)'. To the right of the button is a section titled 'Surf Information' containing details about the surgery type (Knee joint - total), surgeon (Olivier J. Serrano), therapist (Elise De Clerck), hospital (AZ Ziekenhuis Brussel), and treatment (TKA General Therapy). Below this is a section titled 'Patient-reported outcomes' with a table of scores.

Score	Pre op	Post op	Δ Pre
Overall knee scores (KSS)	50	70	20
Function (KSS)	53	75	22
Pain (KSS)	40	68	28
Activities of Daily Life (KSS)	39	61	22
Quality of Life (KSS)	31	55	24
Flexion/Extension (KSS)	4	4	0

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What are the scores meaning?
The mean score for the data using the oxford hip problem is 48, indicating full operation. **Red** supports for the knee 82%, and its average is 40 and for the patients with scores between 20, 40 and 60. **Blue** is for the knee between 60 and 80 and 82% of all patients for the patients above 80.

☆ Patient's scores over time

The radar chart displays scores for five metrics across six time points: Pre-op, 6 weeks, 3 months, 6 months, 1 year, and 2 years. The Test Group A (red) generally shows higher scores than the Other groups (blue) across most metrics and time points.

Metric	Pre-op	6 weeks	3 months	6 months	1 year	2 years
Oxford Hip Score	41	29	39			
Forgotten Joint Score	21		50			
Pain Score	63	48	70			
Symptoms	75	70	75			
Activities of Daily Life	72	68	70			

♣ Evolution of symptoms

The box plot shows the distribution of symptom scores over time. The Test Group A (red) shows a significant improvement in symptom scores over time, while the Other groups (blue) remain relatively stable.

Time Point	Test Group A (red)	Other groups (blue)
Pre-op	~41	~75
6 weeks	~29	~70
3 months	~39	~75
6 months	~40	~75
1 year	~40	~75
2 years	~40	~75

♣ Evolution of Activities of Daily Life

The box plot shows the distribution of Activities of Daily Life scores over time. The Test Group A (red) shows a significant improvement in Activities of Daily Life scores over time, while the Other groups (blue) remain relatively stable.

Time Point	Test Group A (red)	Other groups (blue)
Pre-op	~72	~72
6 weeks	~68	~72
3 months	~70	~72
6 months	~70	~72
1 year	~70	~72
2 years	~70	~72

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




Patientendaten – Profil


Basic information	Basic information
National health insurance	First Name *
Corona information	firstname.znetooa99
Bariatrics information	Surname *
	lastname.ix0nvk25f

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The screenshot displays the 'Exercise Management' interface. At the top, there's a 'Manage patient scheme' table with columns for exercises (Icepacks Hip, Circulation, Heel raises, Glutes contraction) and days (Day 1, Day 2, Day 3). The table shows active exercises with blue checkmarks and durations (e.g., 6x 20 minutes). Below the table, there's a 'Plan exercise' modal for 'Glutes contraction' on 'Day 5'. The modal includes fields for 'Start day / End day' (5 / 5), 'Repeat' (Daily), 'Intensity start' (10), 'Intensity target' (10), 'Frequency start' (3), and 'Frequency target' (3). Buttons for 'Cancel', 'Apply', and 'Apply & save' are at the bottom.

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Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis Show in Patient language: No Advanced + Stop active follow-up Complete

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Name	Start day	End day	Frequency (Week)	Active from	Reminders	Order	Preview	Remove
Discharge	9	10	1	10:00	Reminders not set		Preview	Remove
Postoperative recovery	0	3	1	10:00	Reminders not set		Preview	Remove
Activities	15	100	1	10:00	Reminders not set		Preview	Remove
Exercises	5	100	1	10:00	Reminders not set		Preview	Remove

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Videoanalyse:

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** **Video analysis** Show in Patient language: No Advanced + Stop active follow-up Complete

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



The screenshot shows the 'Edit Laurens Baes' page in the moveUP interface. The top navigation bar includes 'moveUP', 'Patients', 'Onboarding', 'Config', 'Admin', and a user profile for 'Laurens'. The main content area is titled 'Edit Laurens Baes' and contains two columns of form fields. The left column includes fields for First name (Laurens), Surname (Baes), Username (laurens@moveup.care), Email (laurens@moveup.care), Telephone number (+32496123456), Social security number, and a dropdown for Residence. The right column includes fields for Address (Oktrooi plein 1), City (Gent), Postal code (9000), Country (Belgium), Title, Image, and Language (English). At the bottom, there are two sections: 'Authentication' with a toggle for 'Two factor authentication' (currently off) and 'Wellness module' with a toggle for 'Enable wellness module' (currently on). At the very bottom, there are three buttons: 'Save', 'Set new password', and 'Manage consent/contracts'.

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Medical score: Select a medical score

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



	Stop treatment? Are you sure you want to stop the treatment? Stopping a treatment will mark the patient as being finished. Medical scores in the treatment will no longer be activated. Cancel Stop	patiënt worden doorgestuurd naar de status 'Voltooid'.
--	---	---

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- Wellnessgegevens inschakelen

The screenshot shows the 'Edit Laurens Baes' form in the moveUP system. The form is divided into two columns. The left column contains fields for First name (Laurens), Surname (Baes), Username (laurens@moveup.care), Email (laurens@moveup.care), Telephone number (+32496123456), Social security number (Social security number), and Rizin nr (Rizin nr). The right column contains fields for Address (Oktrooi plein 1), City (Gent), Postal code (9000), Country (Belgium), Title (Title), Image (Image), Language (English), and Wellness module (Enable wellness module). At the bottom, there are buttons for 'Save', 'Set new password', and 'Manage consent/contracts'.

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
User manual for the Health Care Provider (HCP)

A personalized treatment, with more insight and involvement in the care pathway.

App version : 2.2.0 MDD and later

Manual version : 17.2

Data of release manual : 27/01/2026

Link to e-learning: <https://moveup.talentlms.com/index>

Table of Contents

PRODUCT DESCRIPTION.....	3
Instructions	3
Intended use.....	4
Intended users	4
Target population / indication for use	5
Claims	5
moveUP companion/coach/therapy	5
moveUP therapy	5
moveUP app & devices	6
Warnings and precautions	6
Contra-indications and potential side effects	7
Notification of serious incidents	7
Human body contact	7
Accessories/products used in combination	8
Device lifetime	8
Company details	8
PATIENT APP	9
MEDICAL MODULE	9
HEALTHCARE PLATFORM ACCOUNT CREATION.....	9
ONBOARDING OF PATIENTS	9
HEALTHCARE PLATFORM GENERAL FUNCTIONS	12
HEALTHCARE PLATFORM CARE PLAN FUNCTIONS	19
WELLNESS MODULE.....	22
WELLNESS MODULE ACTIVATION	22
WELLNESS MODULE OVERVIEW	23
Wellness module overview	24
Wellness report	24

PRODUCT DESCRIPTION

Instructions

	Read carefully this instruction for use. Make sure that you know how to use moveUP, but also understand the limitations of moveUP before you start using our service. If you do have any questions on the use, contact moveUP.
	In the final chapter of this document you find our contacts, which you can use during the use of our service, in case you would have technical problems with the application (moveUP dashboard or patient app) or when expertise is required about a care path.
	The Healthcare practitioner can't give access to third parties to access the moveUP dashboard. Every other person can, even by accident, send out wrong information or adapt the patient profile, with the result that the advice sent by moveUP no longer corresponds with the needs of the patient.
	In case the application is use while moveUP & b.clinic is not asked to do the follow-up of patients, the individual healthcare provider or the health care facility is responsible for the follow-up of the patient and the action that are carried out.
	In compliance with the medical device legislation, you are required to report any complaints towards the application, any alteration regarding characteristics/specifications of the application and incidents/complications potentially linked to the use of the application to support@moveup.care. The AFMPS can also be contacted in case of incident via https://www.famhp.be/en/human_use/health_products/medical_devices_accessories/materiovigilance/how_notify
	In compliance with the GDPR legislation, please report eventual data leaks linked to the use of moveUP. Please, read the entire privacy policy: https://www.moveup.care/legal/privacy-policy

Intended use

For healthcare providers, moveUP provides efficient clinical management of pathologies and treatments. With the use of moveUP, valuable insights are given to the HCP/care team of the status of the patient and its evolution.

For patients, moveUP provides personalised information and instructions to help them managing their symptoms and progress in their rehabilitation. The intensity of follow-up is adapted based on patients' needs and timepoint in the patient pathway:

moveUP companion = monitoring and information, no active follow-up

moveUP companion offers targeted information and evolution reports to patients. Their care team takes the evolution reports into account in their further recovery when appropriate.

moveUP coach = active follow-up by healthcare team

moveUP coach offers targeted information and evolution reports to patients. Their care team is more actively involved and takes the evolution reports into account and can provide advice and exercise suggestions through the digital platform.

moveUP therapy = active follow-up by healthcare team, with data driven validated care protocol

moveUP therapy offers a data-driven validated care protocol with certain category and level of exercises and activities, specifically targeted to the individual patient. Their care team can manually adapt the data driven validated protocol when needed. Patients can fully rehabilitation with moveUP without leaving their home environment.

* Optional functionality that can be enabled: interoperability with Class IIa continuous passive motion (CPM) medical device. For knee and hip patients who are using a Class IIa CPM medical device, moveUP can interoperate with the Class IIa CPM software. moveUP acts as a facilitator to easily assign a designated CPM exercise protocol chosen by a physician to a patient and to display the performed CPM exercises in the medical dashboard

Intended users

moveUP is intended to be used by patients and healthcare providers. The main user of the mobile app and patient website is the patient.

Inclusion:

- Age: minimum 18 years / maximum no limit
- Health & condition: capable of performing basic activities of daily living
- Language: understanding one of the available languages of the app (Dutch, French, German, English)

Exclusion:

- Patients who are mentally incompetent or having troubles to express what they are feeling (for instance, mentally diseased people, people staying in elderly care centres, ...) are excluded.
- Patients who are not capable of operating a tablet/smartphone and activity tracker.
- Patients who can't understand one of the available languages of the app (Dutch, French, German, English)

The main user of the medical web interface is a healthcare practitioner (group) or clinical researcher (group), named the care team. The care team is able to operate a web interface via web browser on PC/tablet/smartphone. The healthcare practitioner needs to understand one of the available languages of the web interface (today only available in English)

Target population / indication for use

moveUP companion & coach is used by Musculo-skeletal, oncologic, respiratory, gastro-intestinal, cardiovascular, and neurologic patients, such as patients:

- who underwent or planning a gastric bypass or gastric sleeve operation
- who have or had a stroke, multiple myeloma, covid, familial hypercholesterolemia
- who have or had back or joint problems or operations
- moveUP therapy is used by knee & hip arthroplasty patients.
- moveUP is used for rehabilitation but is possible to include a patient before an intervention to know the previous state of health and establish goals for the rehabilitation.

Claims

moveUP companion/coach/therapy

The intensity of follow-up is adapted based on the needs of the patient, via the moveUP Symptom & QoL monitoring tool.

moveUP enhances the clinical management of the patients, because early detection & management of complications is possible via the symptom & QoL monitoring tool

More efficient clinical management, such as the number of consultations, can be reduced

Enforces therapy compliance/adherence

The correct information is provided at the right time

moveUP therapy

With the use of moveUP therapy knee & hip arthroplasty patients can fully rehabilitate via the in-app care team without leaving their home environment

moveUP app & devices

The moveUP app runs on:

- iPhones/iPads (iOS 11.0 or newer)
- Android systems (Android 7.0.0 and newer).

If the patient has no compatible tablet/smartphone, moveUP can lease these devices to the patient.

- Garmin Vivofit 3 and Vivofit 4
- Nokia Withings Go
- All the devices compatible with HealthKit (IOS)
- All the devices compatible with Google Fit (IOS and Android)

The moveUP application needs an internet connection when in use.

The moveUP application is compatible with the following browsers:

- Chrome
- Safari
- Edge
- Firefox

We strongly recommend to use the latest version of the browser.

For your data safety please enable the 2-factor authentication(2FA) and use a safe network

Warnings and precautions

	• moveUP is not an emergency tool. In case of emergency please contact your doctor or call 112.
	• moveUP is not intended to monitor/treat vital parameters of critical diseases.
	• moveUP companion is not a replacement of any treatment you need to follow, but is used as an addition.
	• Patients need to be aware if questionnaires are not filled out on frequent basis or not filled out trustworthy, the care team has more difficulty to know their health situation & evolution. • The moveUP App offers personalized healthcare based on the patient's profile and therefore can not be shared, borrowed or traded between users. The information and advice of the moveUP application do not apply to a patient other than the one whose profile is stored in the application. The patient may not grant access to the moveUP App to third parties. Any

	other person can, even accidentally, send wrong information or change the profile of the patient, with the result that the moveUP-guided advice no longer corresponds to the needs of the patient • With software, there is always a residual risk for bugs. If you have the feeling something is wrong or you don't get a response in the app, please reach out to us via support@ or call 0800 88 008.
---	---

The moveUP IFU are only supplied in electronic format. If you need a full printed version please contact us at the e-mail address info@moveup.care. This mail is available 24/7.

Contra-indications and potential side effects

There are no contra-indications or known side effects

Notification of serious incidents

A serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or another person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health that resulted in any of the following:
 - life-threatening illness or injury,
 - permanent impairment of a body structure or a body function, hospitalisation
 - prolongation of patient hospitalisation, medical or surgical intervention to prevent life-threatening illness or injury
 - permanent impairment to a body structure or a body function,
 - chronic disease,
- a serious public health threat;

If you detect any of these cases please send an email to info@moveUP.care and your competent authority:

- Netherlands: meldpunt@igj.nl
- Belgium: meddev@fagg.be
- France: dedim.ugsv@ansm.sante.fr
- Germany: zlg@zlg.nrw.de

Human body contact

No human body contact with the patient or user, due to the nature of the product (software).

Accessories/products used in combination

There are no accessories. If the patient has no compatible tablet/smartphone or activity tracker, moveUP can lease these devices to the patient. The leased devices are CE marked devices that meets the compatibility criteria outlined in the IFU for the app.

Device lifetime

2-years, depending on the willingness of the user to update the app. We sent out a notification with the recommendation to update the app, to make sure the app will keep the performance & recently included features.

Company details

 moveUP NV/SA	Tel: +32 800 88 008	UDI: ++G166MOVEUP220MDDDM  (01)G166MOVEUP220MDDDM (10)++\$7220 (11)/16D20260127 (21)\$Y
Head office Brussel: Cantersteen 47, B-1000 Brussels, Belgium Secondary office: Kliniekstraat 27A 9050; Ghent, Belgium	Email: info@moveUP.care www.moveup.care	

If you have any questions about the App or the tablet, you can contact moveUP via the following channels - in order of priority: Via the message system of the App. Via email to support@moveup.care - Via phone during office hours: +32 800 88 008



MoveUP guarantees delivery of the eIFU within seven (7) days upon request. Please contact info@moveUP.care if you want to obtain a free copy.

PATIENT APP

Please read the user manual for patients on www.moveup.care/appinfo to view the functionalities of the patient app.

MEDICAL MODULE

HEALTHCARE PLATFORM ACCOUNT CREATION



Step 1: The HCP registers on the healthcare platform on https://v1.medical.moveup.care/	
Step 2: moveUP will contact the HCP to add his/her profile to the correct organization.	
Step 3: Security measures: It is important to enable 2FA as you have access to sensitive health data.	
E-learning: For HCPs who will provide care through the use of the healthcare platform it is recommended to follow the e-learning courses for the indications listed on https://moveup.talentlms.com . We ask the HCPs to register for the e-learning course in their applicable language.	

ONBOARDING OF PATIENTS

Onboarding is the process of registering a new patient in the healthcare platform.




A crucial part is adding a care plan to the patient's profile. A care plan determines what will be visible for the patient in the application, e.g.: type of exercises, information modules, message functionality, etc. A care plan usually consists out of an intervention, a treatment and (a) care






team(s).

Onboarding patients		The onboarding dashboard is easy retrievable in the health care platform in the upper left corner when clicking 'Onboarding'. It is also accessible through https://patient-onboarding-admin-moveup.care using the HCP credentials of the health care platform. Click on 'new patient' to onboard a new patient and fill out the patient information. After the patient information is filled out, press 'Save'.
	New patient 	Important note: Only one account can be coupled to an e-mail address. It is not possible to reuse an e-mail address that is already known within our database, e.g.: reusing the same e-mail address for the account of a partner or parent.
	Patient's care plans 	Next up, the care plan needs to be selected: press 'Start new care plan'. Drag and drop an intervention and treatment to the care plan and fill out the dates, name of the treating doctor and care team. Drag and drop a study if applicable. Click Save.
	Email variables <pre> {{ firstname }} {{ lastname }} {{ code }} {{ registrationlink }} {{ organizationName }} {{ surgeonName }} {{ hosp_spec_par_start_NL }} {{ organizationName }} {{ bodyStructure }} </pre> 	A mail template is shown. The HCP reviews the e-mail variables and clicks on 'send selected mails'. An invitation mail will be sent to the patient with all the needed information to install and activate the application.

HEALTHCARE PLATFORM GENERAL FUNCTIONS

Below the HCP can retrieve all functions of the healthcare platform. The applicable functions will be shown based on the HCPs subscription. The extent of the available patient information is based on the patients subscription & his/her compliancy rate.

Patient overview		The health care platform gives a clear overview of the HCPs patients. In the top bar the relevant categories are shown that identify the patient and show relevant health parameters, unread messages, the lead clinician and patient compliance. By default, when opening the healthcare platform, the patient overview is sorted based on intervention date.
Tags		Select patients related to a specific pathology by using tags.
Filters		Select patients based on language, test patient, lead clinician or being activated or not. Activated can mean 2 things: either the patient registered him/herself with a code through the app or he has been marked as immediately activated in onboarding. Non activated: patients who haven't registered themselves through the app yet and haven't been marked as immediately activated.
Color codes		Color codes indicate which patients need closer attention.

Messaging		The number of unread messages are displayed in the overview.
Review / Flags		With the use of review flags it is known when the patient was checked or who needs to check the patient. See more info in Patient details – Review / Flags.
Notifications		Notifications can be sent directly to the patient's device home screen.
User management		Only limited HCPs of an organization have the permission to add users. In case a HCP has this permission, they can add users by clicking the 'New' button, fill out the credentials, select
		the applicable role and press 'Register'. A mail will be sent to the newly created user to follow the account creation steps, see above .
Patient details – Overview		In the patient details overview all relevant (health) information of the individual patient can be retrieved. With the use of the previous and next patient buttons the HCP can easily switch to another patient. In the left side bar, the HCP can navigate between status, log, care plan, profile, documents and more.

Patient details – Status

firstname yb/last first/lastname last/initials 26719

Gender: Female, Age: 76 years, Date of birth: , Language: Dutch

Start consultation, Exercises, Medical questionnaires (75%), Consultations (54.55%), Questionnaires, Reported events, Not verified

TKA General Therapy

Compliance OK for today or yesterday (Green dot)
Non compliant > 2 days (Yellow dot)
Non compliant > 7 days (Red dot)

The reported healthcare data are retrievable in the 'Status' section. From there it is possible to navigate to the recurring and medical questionnaire results, exercises and many more.

A **GREEN**, **YELLOW** or **RED** status dot shows how compliant the patient is to the care plan, for example filling out all the suggested recurring and medical questionnaires and performing the suggested exercises will result in a **GREEN** status dot.

Patient details – Medical questionnaires

Status, Start consultation, Exercises (7), Medical questionnaires (30%), Consultations, Questionnaires

THA General Therapy

Via the 'medical questionnaire' button it is shown:

- ✓ Which medical questionnaires are active for a patient.
- ✓ How many of these are completed or not.
- ✓ What the reported results are on the individual questions within the medical questionnaire.

What the sub scores and total scores are to be able to see a quick indication on how the patient is progressing.

Patient details – Recurring questionnaires

Status, Start consultation, Exercises (7), Medical questionnaires (100%), Consultations (28.57%), Questionnaires

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Day	Reported at	DD-Week(s) in the dressing dry?
30	29 Nov 2023, 16:30:59	Yes
29	28 Nov 2023, 17:08:38	Yes
28	27 Nov 2023, 16:07:33	Yes
27	26 Nov 2023, 16:32:18	Yes
26	25 Nov 2023, 16:30:46	Yes




With these patient reported answers it is possible to follow the progress of important health parameters frequently.




Patient details - Reported events

Reported events

Day	Reported at	DD-Pain(1) Indicate how much pain you had in your affected joint during the night	DD-Pain(2) Indicate how much pain you had in your affected joint during the day	DD-Pain(3) Did you experience any pain other than in the affected joint?	DD-Pain(4) Indicate where and when you experienced pain
242	Jan 25, 2023, 9:26:45 AM	0.52	0.81	Yes	Night: Back Exercises: Call Pain: Grains

The reported events are categorized per event. For each reported event you can view: (1) The date when the event was reported (+day in care plan) (2) The answers to all questions asked to the patient.

Patient details - Review / Flags		1st line review is for the HCP that is not the lead clinician (e.g. physical therapist), 2nd line review is for the lead clinician (e.g. surgeon)
		Click validate when a patient has been monitored. Monitoring can mean: ✓ Questions answered - Advice given ✓ Adjusted or checked exercises ✓ Adjusted or checked step goal ✓ Checked daily parameters: pain, swelling, activity,... ✓ Checked postoperative protocol: use of crutches/walker, medication,... Same principle for 2nd line review.
		Flag someone from first line ✓ Patient needs to be checked by someone from first line (e.g. a physical therapist) ✓ Give a reason for flag ✓ Select the discipline to flag ✓ Optional: add an e-mail notification to flag. Flagged discipline will receive an email about the flag. A flag will add a new line in the patient log with the reason for flagging(see further).

		Flag someone from second line ✓ Patient needs to be checked by someone from second line (e.g. a surgeon) ✓ Give a reason to flag. ✓ Select the discipline to flag. Optional: add an e-mail notification to flag. The flagged discipline will receive an email about that flag.
		Report a technical issue.
Patient details - Log		Flags and technical updates automatically appear in the log of the patient. In case of flagging: ✓ GREEN : solved flag ✓ RED : Unsolved flag for discipline of logged in user (e.g. PT) ✓ YELLOW : Unsolved flag for other discipline than logged in user (e.g. Surgeon) Other important information can be added to the log manually. The applicable tags are added to the log in view of easy log searches

Patient details – Report

Report

Personal information

Age: 22
Gender: Female
Height: 165 cm
Weight: 55 kg
Blood type: A Rh+

Surgery information

Surgery type: Laparoscopic cholecystectomy
Surgery date: 2023-10-10
Surgery location: General Surgery
Surgery status: Completed

Clinical milestones

Drop (days after surgery): 16
Discharge (days after surgery): 16
Pain (days after surgery): 16
Activity (days after surgery): 16
Solid food (days after surgery): 16
Bowel movement (days after surgery): 16
Catheter removal (days after surgery): 16
Wound dressing change (days after surgery): 16
Medication (days after surgery): 16
Vital signs (days after surgery): 16
Patient education (days after surgery): 16
Follow-up (days after surgery): 16

This report is used to visualize the main milestones and medical questionnaires of the patient, as well as the progression of the patient.

The report can be downloaded in pdf form by clicking the 'Download report' button. This way the HCP can send the report to the patient. It's also possible to share the link directly to the report, using the button 'Copy report link'.

Patient details – Profile

Basic info | Clinical data | Care plan

Basic information

National health insurance: [input field]
Corona information: [input field]
Bariatrics information: [input field]

Basic information

First Name: [input field]
Surname: [input field]

The HCP can view, edit or add all relevant basic, clinical, administrative or care plan specific data in the patient profile

Patient details – Messaging

Messages

11/01/2024






The Patient

11/01/2024

Messages 2

Push notification

The message popup (which can be enlarged) is closed by default and the number of new messages is visible in the floating button. This functionality is used to give advice to patients and answer their questions. It can also be used to ask for pictures, for example to follow-up on wound healing, swelling, redness etc...


		When clicking 'messages', a set of prescribed personal messages is shown wherefrom the HCP can select an applicable message for that particular patient.
		Click on the "1st line symbol" to mark the question as read.
		Messages that are submitted by a HCP can be edited or removed. When hovering over a message, a button to edit or delete the message will become visible.
Patient details -- Documents		Documents can be uploaded in the 'Documents' section. Chat images and video recordings are automatically collected in this section in view of easy access to these materials.
Patient details - More		See which documents / contracts the patient has signed. See on which device the patient is using the application and when the application was first and last seen. Export data is only possible for limited HCPs due to security reasons and to protect data integrity.
Log out		In the upper right corner the HCP can log out. Know that the automatic log out time is 14 days, therefore it is recommended to log out more frequently to avoid unauthorized access.
Support		Reach out via support@moveup.care for questions or use the '?' button

		on the platform.
Care plan management		In care plan management it is possible to navigate between multiple functions that are applicable to the indication of the patient, such as exercises, step goal, recurring questionnaires, video analysis, automation, etc... For more information see further.

HEALTHCARE PLATFORM CARE PLAN FUNCTIONS

It is possible to provide care through the use of the healthcare platform. With the functions below the HCP can adapt the default treatment plan of their individual patient.

Only the applicable functions will be displayed conform the care plan.

Exercise management		The exercise management module includes the following possibilities of adaptation: ✓ Adapt description, frequency, and intensity of the exercises of the day. ✓ Adapt the daily message. ✓ Adapt the exercise scheme: ○ Blue check marks indicate active exercises, accompanied with frequency and intensity of these active exercises. ○ Schedule exercises over a longer period of time. Extend the standard exercise scheme with a set of exercises that are not part of the exercises defined in the care plan
---------------------	--	--

Recurring questionnaire

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis

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Using treatment scheme **Personalize**

Exercises Step goal **Recurring questionnaires** Medical questionnaires Video analysis Show in Patient language No Advanced Stop active follow-up Complete

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Name	Start day	End day	Frequency (days)	Active from	Reminders	Order	Preview	Remove
Discharge	3	10	1	10:00	No reminder set last	1	Preview	Remove
Postoperative reassess	0	3	1	10:00	No reminder set last	2	Preview	Remove
Activities	15	100	1	10:00	No reminder set last	3	Preview	Remove
Exercises	1	100	1	10:00	No reminder set last	4	Preview	Remove

Add questionnaire to personalized care plan

Questionnaire **Add**

Start date End date

Frequency Daily

Activation hour 15 Activation minutes 0

Save Cancel

Recurring questionnaires can be removed, reordered or a new recurring questionnaire can be added. The default list of recurring questionnaires, including their start day, end day, frequency and reminders is predetermined in the care plan.

Medical questionnaire

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** Video analysis Change Reset Stop active follow-up Complete

Medical scores scheme

Name	Start day	End day	Preview	Remove
HIP PROM collection month 6	170	210	Preview	Remove
HIP PROM collection Year 1	340	500	Preview	Remove
Discharge Questionnaire HIP	0	15	Preview	Remove

Add medical questionnaire to care plan

Medical score **Add**

Select a medical score

Activation type Day in treatment

Start Date End Date

0 0

Save Cancel

Medical questionnaires can also be removed or added to the default scheme. The default list of medical questionnaires, including their start and end day is predetermined in the care plan.

Video analysis

Exercises Step goal **Recurring questionnaires** **Medical questionnaires** **Video analysis** Show in Patient language No Advanced Stop active follow-up Complete

Video analysis

Analysis enabled for this patient

Results NOT shown to this patient

Add request

Jan 12, 2024

Gait analysis

not recorded

Video requests

Name	Start day	End day	Preview	Remove
Gait analysis	-7	-1	Preview	Remove
Gait analysis	21	30	Preview	Remove
Gait analysis	40	45	Preview	Remove

Add video request to personalized care plan

Video request **Add**

Select a video type

Start date End date

0 0

Save Cancel





Request the patient to record a video for analysis

Video type **Add**

Select a Video type

Save Cancel

With the use of video requests, it is possible to gain more objective information, for example about the gait analysis or range of motion of an operated joint. The default list of video requests is predetermined in the care plan. A video request can be added directly by using the 'Add request' button or a video request can be scheduled for a longer period of time by using the 'Add' button and by adding a start and end day to put the video request active.

Goal manager		The HCP can enable and disable a behavioral change tool which enables patients to set their personal goal, define barriers and solutions in order to reach a healthy lifestyle. The HCP can see the progress on the patient's goals that are defined by the patient.
Change of care plan		The care plan can be adapted based on the patient's interventions, needs and wishes. This can be done via the onboarding tool or via the healthcare platform via the button 'Advanced'. Instructions how to change a care plan are clearly indicated within the tools.
Stop active follow-up		When the active follow-up is finished (what falls under active follow-up is predetermined in the care plan), press 'Stop active follow-up'. This will redirect the patient to a passive long term follow-up if applicable (also predetermined in the care plan if long term follow-up is needed).
		When no actions are required anymore from the HCP or the patient, the patient can be redirected to the 'Completed' status.

WELLNESS MODULE

The wellness module is an optional non-medical device module and therefore the data shown is not intended to be used for clinical decision making.

This data gives an overview of the general wellbeing of the patient.

The wellness module by default is blurred and needs to be activated to check the patient's wellness data.

WELLNESS MODULE ACTIVATION

To activate the wellness module you need to enter your profile settings and activate the option:

- Enable wellness data

The screenshot shows the 'Edit Laurens Baes' form in the moveUP system. The form is divided into two columns. The left column contains fields for First name (Laurens), Surname (Baes), Username (laurens@moveup.care), Email (laurens@moveup.care), Telephone number (+32496123456), Social security number (Rizin nr), and a section for Authentication with a toggle for 'Two factor authentication'. The right column contains fields for Address (Oktrooi plein 1), City (Gent), Postal code (9000), Country (Belgium), Title, Image, Language (English), and a section for the Wellness module with a toggle for 'Enable wellness module'. At the bottom, there are buttons for 'Save', 'Set new password', and 'Manage consent/contracts'.

Upon activation, a mandatory acknowledgment pop-up is displayed. It informs you that the data from the Wellness Module must not be used to support clinical decision-making.

Platform notice

Wellness module terms of use

Dear User,
You are about to activate the Wellness Module.

What is the Wellness Module?

The wellness module provides you with additional data about the general wellbeing and activity levels of your patients. The data in the Wellness Module can come from:

- imported data from general-purpose devices such as e.g. activity trackers
- manual entry by patients

This data does not come from medical devices. This means that neither the devices from which the data is imported, nor the manual entries by the patients have been verified or validated by moveUP, and therefore, it should not be used for clinical decision-making.

What is the purpose of this data?

We provide the option to see wellness data to give you a broader idea about the general wellbeing of the patients and to give you the opportunity to discuss and give further guidance regarding wellbeing.

This module is an optional module and is not linked to the medical data presented in the "Medical data" tab of the dashboard.

The data in the wellness module is not intended to be used for clinical decisions.

How to recognize the wellness data?

Wellness data is displayed in a separate tab in the dashboard or in a separate report.

Whenever wellness data is displayed in a column, together with other data it has a grey background to visually indicate that this does not concern medical data. Furthermore, it has a warning sign next to it.

Enabling/Disabling the Wellness Module

At any time, a user can enable/disable the Wellness Module in the user settings.

Acknowledgement

By clicking on the button below, you acknowledge that you have read and understood that the Wellness Module is an optional module that provides you with additional data regarding the general wellbeing of the patient. This data is not validated by moveUP and is not intended for clinical decision making.

☐ I agree


Sign

Cancel

Once acknowledged, you gain access to the wellness data. This activation allows you to view the wellness data and export the wellness reports.


WELLNESS MODULE OVERVIEW

Wellness module overview

In this tab you'll be able to see the step data from the patient and see their evolution.

Wellness report

You can see the wellness reports when needed.

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The screenshot displays the 'Exercices' (Exercises) module interface. It includes a sidebar with a list of exercises, a main area with a video player and a 'Circulation' section, and a 'Manage patient scheme' section. A 'Plan exercise' dialog box is open, showing settings for 'Start day / End day', 'Repeat', 'Intensity start (repetitions)', 'Intensity target (repetitions)', 'Frequency start', and 'Frequency target'. The 'Plan exercise' dialog box is highlighted with a red border.

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



The screenshot displays the 'Questionnaires' module interface. It includes a sidebar with a list of questionnaires, a main area with a 'Questionnaires' section, and a 'Personalize' button. A 'Add questionnaire to personalized care plan' dialog box is open, showing fields for 'Questionnaire', 'Start date', 'End date', 'Frequency', 'Activation hour', and 'Activation minutes'. The 'Add questionnaire to personalized care plan' dialog box is highlighted with a red border.

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


The screenshot displays the 'Medical scores scheme' module interface. It includes a sidebar with a list of medical scores, a main area with a 'Medical scores scheme' section, and an 'Add' button. A table lists medical scores with columns for 'Name', 'Start day', 'End day', 'Provider', and 'Remove'.

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The screenshot shows the 'Edit Laurens Baes' profile page. The page has a header with 'moveUP Patients Onboarding' and a user menu for 'Laurens'. The main content area is titled 'Edit Laurens Baes' and contains a form with the following fields:

- First name: Laurens
- Surname: Baes
- Username: laurens@moveup.care
- Email: laurens@moveup.care
- Telephone number: +32496123456
- Social security number: Social security number
- RIZIV nr: RIZIV nr
- Address: Oltroolplein 1
- City: Gent
- Postal code: 9000
- Country: Belgium
- Title: Title
- Image: Image
- Language: English
- Authentication: Two factor authentication (checked)
- Wellness module: Enable wellness module (unchecked)

At the bottom of the form, there are three buttons: 'Save', 'Set new password', and 'Manage consent/contracts'.

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