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


# User manual for the Health Care Provider(HCP)

# moveUP

A personalized treatment, with more insight and involvement  
in the care pathway.

App version: 2.3.0 MDD and later

Manual version: 18.0

Data of release manual: 04/03/2026

# Table of Contents

<b>Table of Contents</b> .....	<b>1</b>
<b>1. Product description</b> .....	<b>2</b>
1. Instructions .....	2
2. Intended Use .....	3
3. Intended Users .....	3
4. Target population / indication for use .....	4
5. Claims .....	4
6. moveUP app & devices .....	4
7. Warnings and precautions .....	5
8. Contra-indications and potential side effects .....	5
9. Notification of serious incidents .....	5
10. Human body contact .....	6
11. Accessories/products used in combination .....	6
12. Device lifetime .....	6
13. Company details.....	6
<b>2. Patient app</b> .....	<b>7</b>
<b>3. Medical module</b> .....	<b>7</b>
1. Healthcare platform account creation .....	7
2. Onboarding of patients.....	8
3. Healthcare platform general functions .....	9
3.1. Patient overview.....	9
3.2. Patient details.....	10
3.3. Log out.....	15
3.4. Support.....	15
3.5. Care plan management .....	15
4. Healthcare platform care plan functions.....	16
4.1. Exercise management.....	16
4.2. Recurring questionnaire .....	16
4.3. Medical questionnaire .....	17
4.4. Video analysis .....	17
4.5. Goal manager .....	17
4.6. Change of care plan .....	18
4.7. Stop active follow up .....	18
<b>4. WELLNESS MODULE</b> .....	<b>19</b>
1. WELLNESS MODULE ACTIVATION .....	19
2. WELLNESS MODULE OVERVIEW.....	20

# 1. Product description

## 1. Instructions

Read carefully this instruction for use. Make sure that you know how to use moveUP, but also understand the limitations of moveUP before you start using our service.

If you do have any questions on the use, contact moveUP.

In the final chapter of this document you find our contacts, which you can use during the use of our service, in case you would have technical problems with the application (moveUP dashboard or patient app) or when expertise is required about a care path.

The Healthcare practitioner can't give access to third parties to access the moveUP dashboard. Every other person can, even by accident, send out wrong information or adapt the patient profile, with the result that the advice sent by moveUP no longer corresponds with the needs of the patient.

In case the application is use while moveUP & b.clinic is not asked to do the follow-up of patients, the individual healthcare provider or the health care facility is responsible for the follow-up of the patient and the action that are carried out.

In compliance with the medical device legislation, you are required to report any complaints towards the application, any alteration regarding characteristics/specifications of the application and incidents/complications potentially linked to the use of the application to [support@moveup.care](mailto:support@moveup.care). The AFMPS can also be contacted in case of incident via this [link](#)

In compliance with the GDPR legislation, please report eventual data leaks linked to the use of moveUP.

Please, read the entire privacy policy: <https://www.moveup.care/legal/privacy-policy>

## 2. Intended Use

For healthcare providers, moveUP provides efficient clinical management of pathologies and treatments. With the use of moveUP, valuable insights are given to the HCP/care team of the status of the patient and its evolution.

For patients, moveUP provides personalised information and instructions to help them managing their symptoms and progress in their rehabilitation. The intensity of follow-up is adapted based on patients' needs and timepoint in the patient pathway:

- moveUP companion = monitoring and information, no active follow-up

moveUP companion offers targeted information and evolution reports to patients. Their care team takes the evolution reports into account in their further recovery when appropriate.

- moveUP coach = active follow-up by healthcare team

moveUP coach offers targeted information and evolution reports to patients. Their care team is more actively involved and takes the evolution reports into account and can provide advice and exercise suggestions through the digital platform.

- moveUP therapy = active follow-up by healthcare team, with data driven validated care protocol

moveUP therapy offers a data-driven validated care protocol with certain category and level of exercises and activities, specifically targeted to the individual patient. Their care team can manually adapt the data driven validated protocol when needed. Patients can fully rehabilitation with moveUP without leaving their home environment.

\* Optional functionality that can be enabled: interoperability with Class IIa continuous passive motion (CPM) medical device. For knee and hip patients who are using a Class IIa CPM medical device, moveUP can interoperate with the Class IIa CPM software. moveUP acts as a facilitator to easy assign a designated CPM exercise protocol chosen by a physician to a patient and to display the performed CPM exercises in the medical dashboard

## 3. Intended Users

moveUP is intended to be used by patients and healthcare providers. The main user of the mobile app and patient website is the patient.

Inclusion:

- Age: minimum 18 years / maximum no limit
- Health & condition: capable of performing basic activities of daily living
- Language: understanding one of the available languages of the app (Dutch, French, German, English)

Exclusion:

- Patients who are mentally incompetent or having troubles to express what they are feeling (for instance, mentally diseased people, people staying in elderly care centres, ...) are excluded.
- Patients who are not capable of operating a tablet/smartphone and activity tracker.
- Patients who can't understand one of the available languages of the app (Dutch, French, German, English)

The main user of the medical web interface is a healthcare practitioner (group) or clinical researcher (group), named the care team. The care team is able to operate a web interface via web browser on PC/tablet/smartphone. The healthcare practitioner needs to understand one of the available languages of the web interface (today only available in English)

## 4. Target population / indication for use

moveUP companion & coach is used by Musculo-skeletal, oncologic, respiratory, gastro-intestinal, cardiovascular, and neurologic patients, such as patients:

- who underwent or planning a gastric bypass or gastric sleeve operation
- who have or had a stroke, multiple myeloma, covid, familial hypercholesterolemia
- who have or had back or joint problems or operations
- moveUP therapy is used by knee & hip arthroplasty patients.
- moveUP is used for rehabilitation but is possible to include a patient before an intervention to know the previous state of health and establish goals for the rehabilitation.

## 5. Claims

### moveUP companion/coach/therapy

The intensity of follow-up is adapted based on the needs of the patient, via the moveUP Symptom & QoL monitoring tool.

moveUP enhances the clinical management of the patients, because early detection & management of complications is possible via the symptom & QoL monitoring tool

More efficient clinical management, such as the number of consultations, can be reduced

Enforces therapy compliance/adherence

The correct information is provided at the right time

### moveUP therapy

With the use of moveUP therapy knee & hip arthroplasty patients can fully rehabilitate via the in-app care team without leaving their home environment

## 6. moveUP app & devices

The moveUP app runs on:

- iPhones/iPads (iOS 11.0 or newer)
- Android systems (Android 7.0.0 and newer).

If the patient has no compatible tablet/smartphone, moveUP can lease these devices to the patient.

- Garmin Vivofit 3 and Vivofit 4
- Nokia Withings Go
- All the devices compatible with HealthKit (IOS)
- All the devices compatible with Google Fit (IOS and Android)

The moveUP application needs an internet connection when in use.

The moveUP application is compatible with the following browsers:

- Chrome
- Safari
- Edge
- Firefox

We strongly recommend using the latest version of the browser.

For your data safety, please enable the 2-factor authentication(2FA) and use a safe network

## 7. Warnings and precautions

moveUP is not an emergency tool. In case of emergency please contact your doctor or call 112.

moveUP is not intended to monitor/treat vital parameters of critical diseases.

moveUP companion is not a replacement of any treatment you need to follow, but is used as an addition.

Patients need to be aware if questionnaires are not filled out on frequent basis or not filled out trustworthy, the care team has more difficulty to know their health situation & evolution.

The moveUP App offers personalized healthcare based on the patient's profile and therefore can not be shared, borrowed or traded between users. The information and advice of the moveUP application do not apply to a patient other than the one whose profile is stored in the application. The patient may not grant access to the moveUP App to third parties. Any other person can, even accidentally, send wrong information or change the profile of the patient, with the result that the moveUP-guided advice no longer corresponds to the needs of the patient

With software, there is always a residual risk for bugs. If you have the feeling something is wrong or you don't get a response in the app, please reach out to us via support@ or call 0800 88 008.

The moveUP IFU are only supplied in electronic format. If you need a full printed version please contact us at the e-mail address info@moveup.care. This mail is available 24/7.

## 8. Contra-indications and potential side effects

There are no contra-indications or known side effects

## 9. Notification of serious incidents

A serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or another person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health that resulted in any of the following:
  - life-threatening illness or injury,
  - permanent impairment of a body structure or a body function, hospitalisation
  - prolongation of patient hospitalisation, medical or surgical intervention to prevent life-threatening illness or injury
  - permanent impairment to a body structure or a body function,
  - chronic disease,
- a serious public health threat;

If you detect any of these cases please send an email to info@moveUP.care and your competent authority:

- Netherlands: meldpunt@igj.nl
- Belgium: meddev@fagg.be
- France: dedim.ugsv@ansm.sante.fr
- Germany: zlg@zlg.nrw.de

## 10. Human body contact

No human body contact with the patient or user, due to the nature of the product (software).

## 11. Accessories/products used in combination

There are no accessories. If the patient has no compatible tablet/smartphone or activity tracker, moveUP can lease these devices to the patient. The leased devices are CE-marked devices that meet the compatibility criteria outlined in the IFU for the app.

## 12. Device lifetime

2-years, depending on the willingness of the user to update the app. We sent out a notification with the recommendation to update the app, to make sure the app will keep the performance & recently included features.

## 13. Company details

Tel: +32 800 88 008

UDI: + G166MOVEUP230MDD0

(01) G166MOVEUP230MDD0

(10) \$\$\$7230MDD

(11) /16D20260303

(21) 7\*

moveUP NV/SA

Head office Brussel:  
Cantersteen 47, B-1000  
Brussels, Belgium

Secondary office:  
Kliniekstraat 27A  
9050; Ghent, Belgium

Email:  
[info@moveUP.care](mailto:info@moveUP.care)  
[www.moveup.care](http://www.moveup.care)

If you have any questions about the App or the tablet, you can contact moveUP via the following channels  
- in order of priority: Via the message system of the App. Via email to [support@moveup.care](mailto:support@moveup.care) - Via phone during office hours: +32 800 88 008

moveUP guarantees delivery of the eIFU within seven (7) days upon request. Please contact [info@moveUP.care](mailto:info@moveUP.care) if you want to obtain a free copy.

## 2. Patient app

Please read the user manual for patients on [www.moveup.care/appinfo](http://www.moveup.care/appinfo) to view the functionalities of the patient app.

## 3. Medical module

### 1. Healthcare platform account creation

Step 1: The HCP registers on the healthcare platform on <https://v1.medical.moveup.care/>

**Login**

Username or Email  
Username

Password  
Password

**Login**

Request password reset

Login using an electronic identity

No account? [Register here](#)

Step 2: moveUP will contact the HCP to add his/her profile to the correct organization.

Step 3: Security measures: It is important to enable 2FA as you have access to sensitive health data.

**NAME** !

**! Enable 2FA**

[Click on course catalog in the top bar to see some example courses and get started. Subscribe for free.](#)

**Orthopaedics**

- Belgium**  
English (English)  
Français (French)
- The Netherlands**  
English (English)
- France**  
Français (French)
- Germany**  
Deutsch (German)

E-learning: For HCPs who will provide care through the use of the healthcare platform it is recommended to follow the e-learning courses for the indications listed on <https://moveup.talentlms.com>. We ask the HCPs to register for the e-learning course in their applicable language.

## 2. Onboarding of patients

Onboarding is the process of registering a new patient in the healthcare platform.

A crucial part is adding a care plan to the patient's profile. A care plan determines what will be visible for the patient in the application, e.g.: type of exercises, information modules, message functionality, etc. A care plan usually consists out of an intervention, a treatment and (a) care team(s).

The onboarding dashboard is easy retrievable in the health care platform in the upper left corner when clicking 'Onboarding'. It is also accessible through <https://patient-onboarding-admin-moveup.care> using the HCP credentials of the health care platform.

Click on 'new patient' to onboard a new patient and fill out the patient information.

After the patient information is filled out, press 'Save'.

Important note: Only one account can be coupled to an e-mail address. It is not possible to reuse an e-mail address that is already known within our database, e.g.: reusing the same e-mail address for the account of a partner or parent.


Next up, the care plan needs to be selected: press 'Start new care plan'. Drag and drop an intervention and treatment to the care plan and fill out the dates, name of the treating doctor and care team.

Drag and drop a study if applicable. Click Save.

A mail template is shown. The HCP reviews the e-mail variables and clicks on 'send selected mails'.

An invitation mail will be sent to the patient with all the needed information to install and activate the application.

### Patients Onboarding

#### New patient

#### Patient's care plans

#### Email variables

```

{{ firstname }}
{{ lastname }}
{{ code }}
{{ registrationlink }}
{{ organizationName }}
{{ surgeonName }}
{{ hosp_spec_par_start_HL }}
{{ organizationName }}
{{ bodyStructure }}

```

### 3. Healthcare platform general functions


Below the HCP can retrieve all functions of the healthcare platform. The applicable functions will be shown based on the HCPs subscription. The extent of the available patient information is based on the patients subscription & his/her compliancy rate.

#### 3.1. Patient overview

##### Patient overview


The health care platform gives a clear overview of the HCPs patients.

In the top bar the relevant categories are shown that identify the patient and show relevant health parameters, unread messages, the lead clinician and patient compliance. By default, when opening the healthcare platform, the patient overview is sorted based on intervention date


##### Tags

Select patients related to a specific pathology by using tags.

##### Filters

Select patients based on language, test patient, lead clinician or being activated or not. Activated can mean 2 things: either the patient registered him/herself with a code through the app or he has been marked as immediately activated in onboarding. Non activated: patients who haven't registered themselves through the app yet and haven't been marked as immediately activated.

##### Color codes

Color codes indicate which patients need closer attention.


Pain	Feeling ^	Pain v	Feeling
1.3	1	8.7	2
1.6	1.4	8.5	7.7
3	1.9	8.2	4.1

VS


## Messaging

The number of unread messages are displayed in the overview


## Review / Flags

With the use of review flags it is known when the patient was checked or who needs to check the patient. See more info in Patient details – Review / Flags.

## Notifications


Notifications can be sent directly to the patient's device home screen.

## User management

Only limited HCPs of an organization have the permission to add users. In case a HCP has this permission, they can add users by clicking the 'New' button, fill out the credentials, select


the applicable role and press 'Register'. A mail will be sent to the newly created user to follow the account creation steps, see above.

## 3.2. Patient details

### Overview

In the patient details overview all relevant (health) information of the individual patient can be retrieved. With the use of the previous and next patient buttons the HCP can easily switch to another patient. In the left side bar, the HCP can navigate between status, log, care plan, profile, documents and more.

## Status

The reported healthcare data are retrievable in the 'Status' section. From there it is possible to navigate to the recurring and medical questionnaire results, exercises and many more.

A GREEN, YELLOW or RED status dot shows how compliant the patient is to the care plan, for example filling out all the suggested recurring and medical questionnaires and performing the suggested exercises will result in a GREEN status dot.

Compliance OK for today or yesterday  
 Non compliant > 2 days  
 Non compliant > 7 days

## Medical questionnaires

Via the 'medical questionnaire' button it is shown:

Which medical questionnaires are active for a patient.

How many of these are completed or not.

What the reported results are on the individual questions within the medical questionnaire.

What the sub scores and total scores are to be able to see a quick indication on how the patient is progressing.

## Recurring questionnaires

With these patient reported answers it is possible to follow the progress of important health parameters frequently

Exercises	Day	Filled in at	Export
Crutches	30	29 Nov 2023, 16:30:59	(ID: WalkE) Is the dressing dry?
Activities	29	28 Nov 2023, 17:08:38	Yes
Postoperative heurta	28	27 Nov 2023, 16:07:33	Yes
Discharge	27	26 Nov 2023, 16:22:18	Yes
Activity tracker	26	25 Nov 2023, 16:30:46	Yes

## Reported events

The reported events are categorized per event. For each reported event you can view: (1) The date when the event was reported (+day in care plan) (2) The answers to all questions asked to the patient.

Day	Reported at	(ID: Pain-1) Indicate how much pain you had in your affected joint during the night (0=not at all - 10=worst)	(ID: Pain-2) Indicate how much pain you had in your affected joint during the day (0=not at all - 10=worst)	(ID: Pain-3) Did you experience any pain other than in the affected joint?	(ID: Pain-4) Indicate where and when you experienced pain.
242	Jan 23, 2023, 9:26:45 AM	0.52	0.81	Yes	Night: Back Exercises: Calf No pain: Green

## Review / Flags

1st line review is for the HCP that is not the lead clinician (e.g. physical therapist), 2nd line review is for the lead clinician (e.g. surgeon)

Click validate when a patient has been monitored. Monitoring can mean:

- Questions answered – Advice given
- Adjusted or checked exercises
- Adjusted or checked step goal
- Checked daily parameters: pain, swelling, activity,...
- Checked postoperative protocol: use of crutches/walker, medication,...

Same principle for 2nd line review.

Flag someone from first line

Patient needs to be checked by someone from first line (e.g. a physical therapist)

Give a reason for flag

Select the discipline to flag


Optional: add an e-mail notification to flag. Flagged discipline will receive an email about the flag.

A flag will add a new line in the patient log with the reason for flagging (see further).

Flag someone from second line

- Patient needs to be checked by someone from second line (e.g. a surgeon)
- Give a reason to flag.
- Select the discipline to flag.

Optional: add an e-mail notification to flag. The flagged discipline will receive an email about that flag.

Add first line flag

Enter flag reason

Disciplines to flag

Administrative / Other


PT

Notify users

Send mail

This will send an e-mail to all users in the selected disciplines with access to this patient. Keep in mind this might also include users outside of your own organization.

Add Cancel

Add second line flag

Enter flag reason

Disciplines to flag

PT 2nd line

Surgeon

Notify users

Send mail

This will send an e-mail to all users in the selected disciplines with access to this patient. Keep in mind this might also include users outside of your own organization.

Add Cancel

Report a technical issue.


Log

Flags and technical updates automatically appear in the log of the patient.

In case of flagging:


- GREEN : solved flag
- RED : Unsolved flag for discipline of logged in user (e.g. PT)
- YELLOW : Unsolved flag for other discipline than logged in user (e.g. Surgeon)

Other important information can be added to the log manually. The applicable tags are added to the log in view of easy log searches


Report

This report is used to visualize the main milestones and medical questionnaires of the patient, as well as the progression of the patient.

The report can be downloaded in pdf form by clicking the 'Download report' button. This way the HCP can send the report to the patient. It's also possible to share the link directly to the report, using the button 'Copy report link'.


## Profile

The HCP can view, edit or add all relevant basic, clinical, administrative or care plan specific data in the patient profile


## Messaging

The message popup (which can be enlarged) is closed by default and the number of new messages is visible in the floating button. This functionality is used to give advice to patients and answer their questions. It can also be used to ask for pictures, for example to follow-up on wound healing, swelling, redness etc...

When clicking 'messages', a set of prescribed personal messages is shown wherefrom the HCP can select an applicable message for that particular patient.

Click on the "1st line symbol" to mark the question as read.


Messages that are submitted by a HCP can be edited or removed. When hovering over a message, a button to edit or delete the message will become visible.

## Documents

Documents can be uploaded in the 'Documents' section. Chat images and video recordings are automatically collected in this section in view of easy access to these materials.


### More

See which documents / contracts the patient has signed.


See on which device the patient is using the application and when the application was first and last seen.

Export data is only possible for limited HCPs due to security reasons and to protect data integrity.

### 3.3. Log out

In the upper right corner the HCP can log out. Know that the automatic log out time is 14 days, therefore it is recommended to log out more frequently to avoid unauthorized access.


### 3.4. Support

Reach out via support@moveup.care for questions or use the '?' button on the platform.

### 3.5. Care plan management

In care plan management it is possible to navigate between multiple functions that are applicable to the indication of the patient, such as exercises, step goal, recurring questionnaires, video analysis, automation, etc... For more information see further.




## 4. Healthcare platform care plan functions

### 4.1. Exercise management

The exercise management module includes the following possibilities of adaptation:


- Adapt description, frequency, and intensity of the exercises of the day.
- Adapt the daily message.
- Adapt the exercise scheme:
  - Blue check marks indicate active exercises, accompanied with frequency and intensity of these active exercises.
  - Schedule exercises over a longer period of time.

Extend the standard exercise scheme with a set of exercises that are not part of the exercises defined in the care plan

### 4.2. Recurring questionnaire

Recurring questionnaires can be removed, reordered or a new recurring questionnaire can be added. The default list of recurring questionnaires, including their start day, end day, frequency and reminders is predetermined in the care plan.

### 4.3. Medical questionnaire

Medical questionnaires can also be removed or added to the default scheme. The default list of medical questionnaires, including their start and end day is predetermined in the care plan.

Name	Start day	End day	Preview	Remove
HIP PROM collection month 6	170	250	Preview	Remove
HIP PROM collection Year 1	340	500	Preview	Remove
Discharge Questionnaire HIP	0	15	Preview	Remove

**Add medical questionnaire to care plan**

Medical score: Select a medical score

Activation type: Day in treatment

Start Date: 0 End Date: 0

Save Cancel

### 4.4. Video analysis

With the use of video requests, it is possible to gain more objective information, for example about the gait analysis or range of motion of an operated joint. The default list of video requests is predetermined in the care plan. A video request can be added directly by using the 'Add request' button or a video request can be scheduled for a longer period of time by using the 'Add' button and by adding a start and end day to put the video request active.

**Video analysis**

Analysis enabled for this patient

Results NOT shown to this patient

Add request

Jan 12, 2024  
Gait analysis  
Not recorded

Name	Start day	End day	Remove	Preview	Message	Remove
Gait analysis	-7	-1	Remove			
Gait analysis	21	30	Remove			
Gait analysis	-0	45	Remove			

**Add video request to personalized care plan**

Video request: Select a video type

Start date: 0 End date: 0

Save Cancel

**Request the patient to record a video for analysis**

Video type: Select a Video type

Save Cancel

### 4.5. Goal manager

The HCP can enable and disable a behavioral change tool which enables patients to set their personal goal, define barriers and solutions in order to reach a healthy lifestyle. The HCP can see the progress on the patient's goals that are defined by the patient.

**Activities** Select all

Select activities that you want to display

- Jogging In the park
- Exercises (fitness, agility, etc.) at work In de polyvalente zaal
- Running In the forest

**Jogging In the park**

Selected month: 8/20

Activities done: 8/20

Hours recorded: 13h10m

Success rate: 53%

Selected week: 2/5

Activities done: 2/5

Hours recorded: 4h10m

Success rate: 100%

Planned: Mon Tue Wed Thu Fri

Patient's notifications: The day before at 11:39, The day it self at 11:40




See barriers and solutions, See report logs

Start date of this activity: 02-10-23



## 4.6. Change of care plan

The care plan can be adapted based on the patient's interventions, needs and wishes. This can be done via the onboarding tool or via the healthcare platform via the button 'Advanced'. Instructions how to change a care plan are clearly indicated within the tools.


## 4.7. Stop active follow up

When the active follow-up is finished (what falls under active follow-up is predetermined in the care plan), press 'Stop active follow-up'. This will redirect the patient to a passive long term follow-up if applicable (also predetermined in the care plan if long term follow-up is needed).

When no actions are required anymore from the HCP or the patient, the patient can be redirected to the 'Completed' status.

## 4. WELLNESS MODULE

The wellness module is an optional non-medical device module and therefore the data shown is not intended to be used for clinical decision making.

This data gives an overview of the general wellbeing of the patient.

The wellness module by default is blurred and needs to be activated to check the patient's wellness data.

### 1. WELLNESS MODULE ACTIVATION

To activate the wellness module you need to enter your profile settings and activate the option:

#### Enable wellness data

Upon activation, a mandatory acknowledgment pop-up is displayed. It informs you that the data from the Wellness Module must not be used to support clinical decision-making.


Once acknowledged, you gain access to the wellness data. This activation allows you to view the wellness data and export the wellness reports.

## 2. WELLNESS MODULE OVERVIEW


### Wellness module overview

In this tab you'll be able to see the step data from the patient and see their evolution.

### Wellness report

You can see the wellness reports when needed.

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Reported at		Answers			
Day	Reported at	00: Pain(1) Indicate how much pain you had in your affected joint during the night. (0=none; 1= Mild; 2= Moderate; 3= Severe)	01: Pain(2) Indicate how much pain you had in your affected joint during the day. (0=none; 1= Mild; 2= Moderate; 3= Severe)	02: Pain(3) Did you experience any pain other than in the affected joint?	03: Pain(4) Indicate where and when you experienced pain.
242	Jan 23, 2023, 9:26:45 AM	0.52	0.81	Yes	Night, Back Exercise, Call the pain Green

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
Care plan: ALL Type: MEDICAL

Showing 1 out of 2 total log messages.

Date	Type	Log	User
16/01/2024, 14:43	MEDICAL	This patient has an inflammation	

Patient log

Care plan: ALL Type: MEDICAL

Showing 2 out of 404 total log messages.

Date	Type	Log
25/09/2023, 02:11	FLAG (red line)	TRIGGER: The patient provided an improvement suggestion in medical scan: Satisfaction @ 1 week.
16/04/2023, 14:26	FLAG (red line)	TRIGGER: The patient contacted for issue of app index.

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Patient's report - 17/03/2023

Surgery Information

Surgeon type	Chronic Pain
Surgeon	Chronic Pain
Procedure	Chronic Pain
Procedure	Chronic Pain

Personal information

First name	John
Age	35
Gender	Female
Surgeon type	Chronic Pain

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

Clear row	Active exercises	Day	Icepacks Hip	Circulation	Heel raises	Glutes contraction
		-1				
		0	6x 20 minutes	6x 3 minutes		3x 15 repetitions
		1	6x 20 minutes	6x 3 minutes		3x 10 repetitions

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The screenshot shows the 'Recurring questionnaires' section of the application. At the top, there are tabs for 'Exercises', 'Step goal', 'Recurring questionnaires', 'Medical questionnaires', and 'Video analysis'. Below the tabs, there is a 'Questionnaires' section with a 'Personalize' button. A table lists existing questionnaires:

Name	Start day	End day	Frequency (times)	Active from	Reminders	Order	Preview	Remove
Discharge	3	30	5	14:00	No reminder set last		Preview	Remove
Postoperative names	0	3	5	14:00	No reminder set last		Preview	Remove
Activities	15	100	5	14:00	No reminder set last		Preview	Remove
Exercises	1	100	5	14:00	No reminder set last		Preview	Remove

Below the table is a modal titled 'Add questionnaire to personalized care plan'. It contains a dropdown menu for 'Questionnaire', input fields for 'Start date' and 'End date', a 'Frequency' dropdown set to 'Daily', and input fields for 'Activation hour' (set to 16) and 'Activation minutes' (set to 0). 'Save' and 'Cancel' buttons are at the bottom.

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The screenshot shows the 'Medical questionnaires' section. At the top, there are tabs for 'Exercises', 'Step goal', 'Recurring questionnaires', 'Medical questionnaires', and 'Video analysis'. Below the tabs, there is a 'Medical scores scheme' section with an 'Add' button. A table lists existing medical scores:

Name	Start day	End day	Preview	Remove
HIP PROM collection month 6	170	250	Preview	Remove
HIP PROM collection Year 1	340	500	Preview	Remove
Discharge Questionnaire HIP	0	15	Preview	Remove

Below the table is a modal titled 'Add medical questionnaire to care plan'. It contains a dropdown menu for 'Medical score', a dropdown for 'Activation type' set to 'Day in treatment', and input fields for 'Start Date' and 'End Date'. 'Save' and 'Cancel' buttons are at the bottom.

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The screenshot shows the 'Video analysis' section. At the top, there are tabs for 'Exercises', 'Step goal', 'Recurring questionnaires', 'Medical questionnaires', and 'Video analysis'. Below the tabs, there is a 'Video analysis' section with a 'Add request' button. A table lists existing video requests:

Name	Start day	End day	Remove	Messages	Remove
Gait analysis	-7	-1	Remove		Remove
Gait analysis	31	30	Remove	Messages	Remove
Gait analysis	-40	45	Remove		Remove

Below the table is a modal titled 'Add video request to personalized care plan'. It contains a dropdown menu for 'Video request', input fields for 'Start date' and 'End date', and 'Save' and 'Cancel' buttons. Below this modal is another modal titled 'Request the patient to record a video for analysis' with a dropdown for 'Video type' and 'Save' and 'Cancel' buttons.

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Add care plan?

Starting a new care plan will end the active follow-up of the current care plan.

Cancel End current, start new

Edit care plans for Patient

Patient's care plan

Intervention: Knee left - total

Name: 24/01/2024

Doctor: Joren Surgron

Status: Non-activated

Treatment: THA General Therapy

Start Date: 14/01/2024

Code: 163934

HCP: Advantia

Joren Surgron as surgeon is part of the care team

Change current care plan

Care plan: THA General Therapy

Note: Changing current Care plan type will cause the couple scheme to alter to reflect the default scheme of the New Care plan and also will reset all personalized MS, QS and Video Requests. If you also need to change the surgery date, please do this first. If it is changed after changing the care plan, there might be inconsistencies.

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Advanced Stop active follow up Complete

Stop active follow up

Stopping the active follow up will move the patient to the PROM stage. Personalisation of exercises after stop active follow up can be done in the tab 'care plan'. If exercises are not personalized, the patient will receive an overview with all exercises performed during his/her rehab.

Message

The message below will be sent to the patient using the messaging function of the app.

Type a message

Email

No emails configured

Stop Cancel

Advanced Back to active follow up Complete

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#### Stop treatment?

Are you sure you want to stop the treatment? Stopping a treatment will mark the patient as being finished. Medical scores in the treatment will no longer be activated.

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### Patients Onboarding

Dashboard Patients Support Account

Warning: Please verify all information provided here. There will be no email reminders until your information is verified. You will be able to edit your profile after you have been activated.

Patient Onboarding Dashboard  
Welcome

0 moveUP patients 0 not activated 0 activated patients

New patient

Managing Organization: [dropdown]  
Test Patient:

Given names: [input]  
Family name: [input]  
Gender:  Female  Male  Other  Unknown  
Date of Birth: [input]  
Language: [dropdown]

Address Lines: [input]  
Postal Code: [input]  
City: [input]  
Country: [dropdown]

Email: [input]  
Phone: [input]

Patient's care plans

Start a care plan to add intervention and treatment for this patient

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

Clear row	Active exercises	Day	Clear all	Clear all	Clear all	Clear all
		1				
		5		6x 20 minutes	6x 3 minutes	3x 10 repetitions
		5		6x 20 minutes	6x 3 minutes	3x 10 repetitions

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Name	Start day	End day	Frequency (times)	Active from	Reminders	Order
Discharge	3	30	5	16:00	No reminder set last	1
Postoperative status	0	3	5	16:00	No reminder set last	2
Activities	15	100	5	16:00	No reminder set last	3
Exercises	1	100	5	16:00	No reminder set last	4

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**Jogging In the park**

Selected month: 8/20  
Activities done: 8/20  
Hours recorded: 13h10m  
Success rate: 53%

Selected week: 2/5  
Activities done: 2/5  
Hours recorded: 4h10m  
Success rate: 100%

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**Advanced** **Stop active follow up** **Complete**

**Stop active follow up**

Stopping the active follow up will move the patient to the PROMM stage. Personalisation of exercises after stop active follow up can be done in the tab 'care plan'. If exercises are not personalized. The patient will receive an overview with all exercises performed during his/her rehab.

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No emails configured

**Stop** **Cancel**

**Advanced** **Back to active follow up** **Complete**

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**Platform notice**

**Wellness module terms of use**

Dear User,  
You are about to activate the Wellness Module.

**What is the Wellness Module?**  
The wellness module provides you with additional data about the general wellbeing and activity levels of your patients. The data in the Wellness Module can come from:

- imported data from general-purpose devices such as e.g. activity trackers
- manual entry by patients

This data does not come from medical devices. This means that neither the devices from which the data is imported, nor the manual entries by the patients have been verified or validated by moveUP, and therefore, it should not be used for clinical decision-making.

**What is the purpose of this data?**  
We provide the option to see wellness data to give you a broader idea about the general wellbeing of the patients and to give you the opportunity to discuss and give further guidance regarding wellbeing.

This module is an optional module and is not linked to the medical data presented in the "Medical data" tab of the dashboard.

The data in the wellness module is not intended to be used for clinical decisions.

**How to recognize the wellness data?**  
Wellness data is displayed in a separate tab in the dashboard or in a separate report.

Whenever wellness data is displayed in a column, together with other data it has a grey background to visually indicate that this does not concern medical data. Furthermore, it has a warning sign next to it.

**Enabling/Disabling the Wellness Module**  
At any time, a user can enable/disable the Wellness Module in the user settings.

**Acknowledgement**  
By clicking on the button below, you acknowledge that you have read and understood that the Wellness Module is an optional module that provides you with additional data regarding the general wellbeing of the patient. This data is not validated by moveUP and is not intended for clinical decision making.

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