

## WHITE HEROIN IDENTIFICATION

### Substance & Trace Test



#### Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the probe from test unit.	2	Sample suspected substance or surface. <b>Touch or Swipe</b>	3	Secure probe into test unit. Press firmly.	<p><b>Do not continue unless probe is secured properly in place.</b></p>

Identification	4	From both sides of the kit, squeeze on the circle of <b>Chamber 1</b> until the capsule breaks.	5	Gently shake test unit for 2-3 seconds.

Results		The appearance of a purple color indicates the presence of <b>Heroin</b> , if the substance is a white powder.
		If no purple color develops, Heroin is not present.
		<b>Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.</b>

Dilution	6	From both sides of the kit, squeeze on the circle of <b>Chamber 2</b> until the capsule breaks.	7	Gently shake test unit for 2-3 seconds.
			<p><b>Wait 5 minutes.</b> The test window will warm up. This is normal.</p>	

Dilution	8	From both sides of the kit, squeeze on the circle of <b>Chamber 3</b> until the capsule breaks.	9	Gently shake test unit for 2-3 seconds.
			<p>The test window will warm up. This is normal.</p>	



#### Disposal

After Dilution, the test kit is safe to dispose of in the trash.

01-19