

## BROWN HEROIN IDENTIFICATION

### Substance & Trace Test



#### Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	<b>1</b> Remove the probe from test unit. 	<b>2</b> Sample suspected substance or surface. <b>Touch or Swipe</b> Touch probe to the suspected substance.          Swipe either side of the probe across the suspect surface.	<b>3</b> Secure probe into test unit. Press firmly. 	<b>STOP</b> Do not continue unless probe is secured properly in place.

Chamber 1 Identification	<b>4</b> From both sides of the kit, squeeze on the circle of <b>Chamber 1</b> until the capsule breaks. 	<b>5</b> Gently shake test unit for 2-3 seconds. 

Chamber 1 Results	<p>If a brown powder is present, the appearance of a purple-red color strongly indicates the presence of <b>Brown Heroin</b>.</p> <p>If a brown powder is not present, the appearance of a purple-red color can indicate either the presence of Heroin or a basic alkaloid substance (possible drug).</p>
	<p>If a purple-red color develops, proceed to Disposal.</p> <p>If no purple-red color develops, Brown Heroin is not present. Proceed to Step 6 for further testing.</p>

Chamber 2 Identification	<b>6</b> From both sides of the kit, squeeze on the circle of <b>Chamber 2</b> until the capsule breaks. 	<b>7</b> Gently shake test unit for 2-3 seconds. 

Chamber 2 Results	<p>The appearance of a purple-red color indicates the presence of a <b>basic alkaloid substance (possible drug)</b>.</p>
	<p>If no purple-red color develops, Heroin is not present.</p>



#### Disposal

The test kit is safe to dispose of in the trash.

01-19