



Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the probe from test unit.	Sample suspected substance.				
	2		Touch	or	Swipe	or	Dig
			Touch probe to the suspected substance. 		Swipe either side of the probe across the suspected substance. 		Firmly press probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.

Test Preparation	3	Secure probe into test unit. Press firmly.	STOP Do not continue unless probe is secured properly in place.	Identification	4	From both sides of the kit, squeeze on the circle of Chamber 1 until the capsule breaks.	5	Gently shake test unit for 2-3 seconds.

Results		The appearance of a red-orange color <u>within one minute</u> indicates the presence of Amphetamine or Methamphetamine .	
		The appearance of a purple color <u>within one minute</u> indicates the presence of an opiate such as Heroin , Morphine or Codeine .	
		The appearance of a dark purple to black color indicates the presence of Ecstasy/MDMA or its analogues.	
		If none of the above colors develop, the above listed drugs are not present.	
		Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.	

Dilution	6	From both sides of the kit, squeeze on the circle of Chamber 2 until the capsule breaks.	7	Gently shake test unit for 2-3 seconds.
				Wait 5 minutes. The test window will warm up. This is normal.

Dilution	8	From both sides of the kit, squeeze on the circle of Chamber 3 until the capsule breaks.	9	Gently shake test unit for 2-3 seconds.
				The test window will warm up. This is normal.



Disposal

After Dilution, the test kit is safe to dispose of in the trash.