

KETAMINE IDENTIFICATION

Substance Test



Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the probe from test unit.	2 Sample suspected substance.		
			Touch or Dig or Dip		
			Touch probe to the suspected substance.	Firmly press probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.	Place the absorbent paper of the probe into the liquid. Make sure the paper is wet.

Test Preparation	3	Secure probe into test unit. Press firmly.	STOP Do not continue unless probe is secured properly in place.	Identification	4	From both sides of the kit, squeeze on the circle of Chamber 1 until the capsule breaks.	5	Gently shake test unit for 2-3 seconds.

Identification	6	From both sides of the kit, squeeze on the circle of Chamber 2 until the capsule breaks.	7	Gently shake test unit for 2-3 seconds.

Results		The appearance of a lavender precipitate within the liquid, in about one minute, indicates the presence of Ketamine .
	 	If no lavender precipitate develops within a minute, Ketamine is <u>not</u> present. * Also Note: If the liquid remains the dark red-purple of the reagent or a blue color appears on the paper, these do <u>not</u> indicate the presence of Ketamine.



Disposal

The test kit is safe to dispose of in the trash.

01-19