


EPHEDRINE IDENTIFICATION



Substance Test







Read Before Use




As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the probe from test unit.	2	Sample suspected substance.		
			Touch or Swipe or Dig	Touch probe to the suspected substance.	Swipe either side of the probe across the suspected substance.	Firmly press probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.

Identification	4	From both sides of the kit, squeeze on the circle of Chamber 1 until the capsule breaks.	5	Gently shake test unit for 2-3 seconds.
				

Identification	6	From both sides of the kit, squeeze on the circle of Chamber 2 until the capsule breaks.	7	Gently shake test unit for 2-3 seconds.
				A blue, blue-violet or purple color will develop.

Identification	8	From both sides of the kit, squeeze on the circle of Chamber 3 until the capsule breaks.	9	Vigorously shake test unit for 20 seconds.
				Hold test unit steady and keep upright, until 2 layers are visibly separated. This may take up to 2 minutes.

Results		The appearance of a purple color in the <u>upper</u> layer of liquid indicates the presence of Ephedrine or Pseudo-Ephedrine .
		If the <u>upper</u> layer of liquid does not turn a purple color, Ephedrine is not present.
		A blue color does not indicate a positive result.



Disposal

The test kit is safe to dispose of in the trash.

01-19