

## FENTANYL IDENTIFICATION

### Substance Test




#### Read Before Use

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

#### Test Preparation

1

Remove the probe from test unit.



2

Sample suspected substance.


#### Touch

Touch Smart-Tip™ probe to sample in different places on the sample. Make sure the tip is completely covered.



3

Secure probe into test unit. Press firmly.



**STOP**

**Do not continue unless probe is secured properly in place.**

#### Identification

4

From both sides of the kit, squeeze on the circle of **Chamber 1** until the capsule breaks.



5

Gently shake test unit for 2-3 seconds.



#### Results



The appearance of a red-purple (wine) color indicates the presence of **Fentanyl**.



The appearance of a purple color within one minute indicates the presence of an opiate such as **Heroin**, **Morphine** or **Codeine**.



The appearance of a red-orange color within one minute indicates the presence of **Amphetamine** or **Methamphetamine**.



If none of the above colors develop, the above mentioned substances are **not** present in the tested sample.



**Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.**



#### Dilution


6

From both sides of the kit, squeeze on the circle of **Chamber 2** until the capsule breaks.



7

Gently shake test unit for 2-3 seconds.



**Wait 5 minutes.**  
The test window will warm up. This is normal.

#### Dilution


8

From both sides of the kit, squeeze on the circle of **Chamber 3** until the capsule breaks.



9

Gently shake test unit for 2-3 seconds.



The test window will warm up. This is normal.



#### Disposal

After Dilution, the test kit is safe to dispose of in the trash.

0139