Electronic Instructions for Use (eIFU)

1. Introduction

The Byteflies Kit Model 2 is a wearable CE Class IIa medical device for continuous recording of non-invasive physiological and behavioral signals in healthcare and home settings. The Byteflies Kit Model 2 instructions for use (IFU) contains all the information required for safe and correct use of the Byteflies Kit. This document summarizes relevant safety and regulatory information.

- Byteflies Kit:** The collective name for all components of the Byteflies Kit Model 2 medical device.
- Sensor Dot:** A 2-channel biopotential and tri-axial inertial measurement device, indicated for recording electroencephalography (EEG), electrocardiography (ECG), heart rate (HR), respiratory rate (RR), electromyography (EMG), acceleration (ACC), angular velocity (GYR), and activity index (AI).
- Docking Station:** A component that can charge up to 5 Sensor Dots and provide a secure communication gateway to the Byteflies Cloud to upload recorded data.
- Byteflies Cloud:** A secure data storage and processing system for data recorded by Sensor Dots.

Intended Use

The Byteflies Kit Model 2 is a multisensor device intended to collect non-invasive physiological and behavioral signals to be used in the clinic and/or home for Users who need continuous monitoring.

Intended Users

- Trained Operator:** A healthcare professional trained in the usage of Byteflies Kits and its accessories.
- User:** A patient/subject with a medical need, as assessed by a healthcare professional, to monitor one or more of the supported signals; 3 years of age or older and without any contraindications (see Contraindications).

The Byteflies Kit is intended to be used by Trained Operators who configure the Byteflies Kit and its accessories in a medical environment (e.g. hospital). Users are wearers of the Sensor Dot, not operators of the Byteflies Kit, and are 3 years of age or older. Once the correct set up is configured, the User can take the Byteflies Kit home to continue recording data as instructed by the Trained Operator. The User may be required to return to the medical environment for removal or reconfiguration of Byteflies Kit components by the Trained Operator.

Clinical indications

The Byteflies Kit is intended for prescription use in the home or healthcare facility to acquire, record, and transmit physiological and behavioral signals from Users. It acquires, records, and transmits electroencephalogram (EEG), electrocardiogram (ECG), heart rate (HR), respiratory rate (RR), and/or electromyogram (EMG) signals, with optional accelerometer (ACC), gyroscope (GYR), and activity index (AI) signals. It only acquires these signals, no diagnostic claims are made about the User's condition.

Contraindications

The Byteflies Kit and its accessories should not be used on Users that

- 1) have been diagnosed with or are suspected to have life-threatening conditions that could result in immediate danger,
- 2) need to undergo an MRI scan or cardiac defibrillation,
- 3) are between 0 and 3 years of age.

Byteflies Kit can be used with extreme caution and under strict supervision on Users that have an implanted device such as a pacemaker, defibrillator, implantable loop recorder (ILR), vagal nerve stimulator (VNS), responsive neurostimulation device (RNS), or deep brain stimulator (DBS) on the condition that the prescribing physician can reasonably assert that the magnets in the Byteflies Kit components will not interfere with the function of the implanted device. The magnets inside Sensor Dot and Cradle are of type N50.

2. Safety Information

The following is important information for using the Byteflies Kit properly and safely. Carefully read this section before operating any component. The Byteflies Kit is a biopotential and motion recording device comprised of hardware and software that has been tested to assess its safety and effectiveness, and to establish substantial equivalence with predicate devices.

- The Byteflies Kit cannot provide diagnoses. It only records physiological and behavioral signals as configured by and under the authority of the Trained Operator.

- The Sensor Dot LED patterns only reflect battery status and operation mode. The LED does not reflect physiological status in any way and should never be interpreted as an indication of patient health.
- The Byteflies Kit is designed to be operated by a Trained Operator. The User should always receive proper instructions from a Trained Operator before using any component of the Byteflies Kit.
- All components of the Byteflies Kit are only intended to be used by the User as instructed by the Trained Operator and the User cannot configure any component of the Byteflies Kit.
- The data recorded by a Sensor Dot can only be used to assess heart, brain, muscle, or motion-related diseases when reviewed by a properly trained healthcare professional (e.g. a cardiologist, pulmonologist, neurologist, physiologist ...).
- Sensor Dots interface with the human body via Cradles. Cradles for monitoring of biopotential (ExG) signals are compatible with certain commercial ExG electrodes, such as but not limited to Byteflies ExG Adhesives. These commercial electrodes are not supplied as part of the Byteflies Kit. The use of any skin adhesives may cause mild skin irritation, such as redness, itching, or contact dermatitis in some Users. Consult the manufacturer of the ExG electrodes for more information.

Warnings

Do not use any component of the Byteflies Kit before reading these warnings, the IFU and - if applicable - a manual for the ExG electrodes or Byteflies compatible devices:

- Do not use the Byteflies Kit if any component is damaged.
- The Byteflies Kit is provided non-sterile.
- Please verify the LED status indicators to avoid misuse.
- If a Trained Operator or User experiences any adverse events while using a Byteflies Kit, discontinue the use and contact Byteflies immediately.
- Skin adhesives (either provided by Byteflies as a compatible device or from a compatible third party provider) should be replaced if they no longer stick firmly to the skin or as indicated in their respective manuals.
- Similar devices may cause signal interference during data transmission. Avoid operating the Byteflies Kit near interfering devices.
- Keep the Byteflies Kit away from children and pets. The Sensor Dot and the accessories may present a choking hazard, never place or let anyone place a component of the device in their mouth.
- If any component of the Byteflies Kit fails to operate after attempting all suggested troubleshooting steps, contact support as soon as possible.
- The battery used in a Sensor Dot may present a risk of fire, explosion, or chemical burn if mistreated. Do not expose Sensor Dots to excessive heat or fire. Do not crush, puncture, or incinerate as doing so can result in fire, explosion, or the release of toxic gasses. Do not use or charge if a Sensor Dot appears to be leaking, discolored, deformed, or in any way abnormal.
- All components of the Byteflies Kit should be returned to the healthcare provider or study staff at the conclusion of the prescribed period of use.
- Do not use a Sensor Dot on a User without first cleaning it according to the instructions.
- Do not use a Byteflies Kit without first charging the Sensor Dot(s) and completing proper set up as described in the IFU.
- Do not expose a Sensor Dot to strong sources of static electricity.
- Do not expose a Sensor Dot to strong electromagnetic fields.
- Do not leave a Sensor Dot on top of or next to other electrical equipment.
- Do not submerge any component of a Byteflies Kit in liquid.
- Do not clean a Sensor Dot with agents other than those listed in the cleaning instructions.
- Do not damage any component of the Byteflies Kit through drops, violent shaking or crushing.
- Do not allow Sensor Dots to come into contact with conductive materials, except when following instructions on how to recharge the battery.
- Sensor Dot is not a toy. Usage on children 3 years of age and older should be under strict supervision of an adult.
- Do not alter any component of the Byteflies Kit. Any modification is strictly prohibited.
- Do not use at temperatures lower than 15°C or greater than 35°C.
- Do not use where air humidity is lower than 10% or higher than 90%.
- Do not use in environments subjected to pressure lower than 700 hPa or greater than 1080 hPa.
- Do not use a Docking Station with a charger and charging cable other than the one supplied with the Byteflies Kit.



3. Device Information

LED notifications

A Sensor Dot and Docking Station each have a multicolor LED to notify the Trained Operator and User about normal operation, warning/status information, or errors. The meaning of the Sensor Dot LED changes depending on whether it is docked inside a slot on a Docking Station or undocked.

1) Sensor Dot in the Docking Station



When a Sensor Dot is docked, its LED can mean the following:

If the LED on one or more Sensor Dots is yellow (blinking or solid), please check the following:

- If you followed the Sensor Dot shut down procedure, the LED will be yellow until the Sensor Dot has shut down, which may take up to 15 seconds.
- If the Sensor Dot is receiving a software update, its LED will be yellow for up to 20 seconds. Please wait.
- Your Docking Station is connected to the original power adapter. If not, please use the original and see if that resolves the issue.
- The Sensor Dot may not be docked correctly. Take it off the Docking Station, gently wipe down the gold pins with a dry cloth, and place it back on the Docking Station. If the yellow LED persists, contact Byteflies.


2) Sensor Dot out of the Docking Station

- When the LED is red or will not turn on, do not use this particular Sensor Dot anymore and contact Byteflies.
- If one or more Sensor Dots have a yellow LED, check the troubleshooting steps above. If the problem persists, contact Byteflies.

3) Docking Station LED

The Docking Station can connect to the Byteflies Cloud over internet via a WiFi connection. An LED located near the micro-USB port communicates the state of the Docking Station:

- **No light:** Docking Station not connected to power.
- **Solid red:** Docking Station is connected to power but not to WiFi.
- **Blinking red:** Docking Station is in WiFi pairing mode.
- **Solid blue:** Docking Station is connected to WiFi.
- **Solid green:** Docking Station is connected to the Byteflies Cloud.
- **Blinking green:** Docking Station is uploading data to the Byteflies Cloud.
- **Solid purple:** Docking Station is in Configuration Mode
- **Solid yellow:** Docking Station is in Calibration Mode

Re-usability & Cleaning Procedures:

There is no need for routine cleaning or disinfection while the device is used for a single User. Before using any component on a new User, they should be disinfected by the Trained Operator.

In order to protect patients against the risk of cross-contamination, the Sensor Dot, Cradle, Docking Station and the Power Supply Unit should be cleaned and disinfected in preparation for a new patient (sterilization is not required):

To disinfect, use a fresh disposable Trionic™ D wipe and wipe the surface of the components until all visible soiling is removed. Ensure the entire surface remains wet for at least 1 minute (Trionic™ is bactericidal after 1 minute). Wipe off any remaining fluid with a clean disposable tissue. Do not use any abrasive cleaners such as acetone.

4. Contact

Email support: support@byteflies.com

EU office/manufacturer:
Borsbeeksebrug 22, 6th floor
2600 Antwerp
Belgium

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4. Kontakt

E-Mail Support: support@byteflies.com

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6. Regulatory Information

Declaration of Conformity

The Byteflies Kit Model 2 is in conformity with the essential requirements and provisions of the EU Medical Device Regulation (MDR)2017/45 concerning medical devices.

Electromagnetic Compatibility Information

1. Sensor Dot

Emissions test	Test level
Conducted and radiated RF emission (CISPR 11)	group 1, class B
Immunity tests	Test level
Electrostatic discharge (ESD) (IEC 61000-4-2)	8kV contact discharges & 15kV air discharges
Conducted disturbances induced by RF fields (IEC 61000-4-6)	3.0 Vrms on 0.15→ 80% AM, 1kHz
Radiated electromagnetic fields (IEC 61000-4-3)	10.0V/m 80MHz - 2.7GHz, 80% AM at 1kHz
Proximity field from wireless communications equipment (IEC 61000-4-3)	List of frequencies (Table 1) from 9V/m up to 28V/m PM (18Hz or 217Hz), FM 1kHz

2. Docking Station

Emissions test	Test level
Conducted and radiated RF emission (CISPR 11)	group 1, class B
Harmonic distortion (IEC 61000-3-2-0)	230 VAC
Voltage fluctuatoinis and flicker (IEC 61000-3-3)	230 VAC
Immunity tests	Test level
Electrostatic discharge (ESD) (IEC 61000-4-2)	8kV contact discharges & 15kV air discharges
Conducted disturbances induced by RF fields (IEC 61000-4-6)	3.0 Vrms on 0.15→ 80% AM, 1kHz
Radiated electromagnetic fields (IEC 61000-4-3)	10.0V/m 80MHz - 2.7GHz, 80% AM at 1kHz
Proximity field from wireless communications equipment (IEC 61000-4-3)	List of frequencies (Table 1) from 9V/m up to 28V/m PM (18Hz or 217Hz), FM 1kHz
Electrical Fast transient (EFT) (IEC 61000-4-4)	± 2kV on 230 VAC: Tr/Th – 5/50ns, 100kHz ±1kV on 230 VAC: Tr/Th – 5/50ns, 5kHz
Surge (IEC 61000-4-5)	±1.0 kV DMon 230 VAC: Tr/Th – 1.2/50 (8/20)µs
Power frequency magnetic field (IEC 61000-4-8)	30A/m, 50 & 60Hz
Voltage dips, short interruptions and voltage variations (IEC 61000-4-11)	230 & 100 VAC mains: 0% - 10ms; 0% - 20ms; 70% - 500ms; 0% - 5

Table 1 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Service	Modulation Immunity test level (V/m)	Immunity test level (V/m)
385	380 to 390	TETRA 400	Pulse modulation 18Hz	27
450	430 to 470	GMRS 460, FRS 460	FM +- 5kHz deviation 1kHz sin	28
710 745 780	704 to 787	LTE Band 13, 17	Pulse modulation 217Hz	9
810 870 930	800 to 960	GSM 800/900, TETRA 800 iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	28
1720 1845 1970	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	28
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	28
5240 5500 5785	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217Hz	9

Byteflies Kit Regulatory Standards

ISO 10993-1:2018	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
2017/745 EU MDR (2017)	EU Medical Device Regulation
2014/53/EU (2014)	Radio Equipment Directive
EN ISO 13485:2016	Medical Devices - Quality Management Systems - Requirements For Regulatory Purposes
FDA 21 CFR, part 820	Quality system regulation for medical device, federal register
EN ISO 14971:2019	Application of risk management to medical devices
EN 60601-1:2006 + A1:2013 (Ed.3.1)	Medical electrical equipment - Part 1: General requirements for safety
IEC 60601-1-2:2014 (Ed. 4.0)	Medical electrical equipment - part 1-2: General requirements for basic safety and essential performance - collateral standard: electromagnetic compatibility - requirements and tests
IEC 60601-1-11:2015	Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
IEC 60601-1-6:2010 + A1:2013	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
IEC 80601-2-26:2019	Medical electrical equipment - Part 2-26: Particular requirements for the basic safety and essential performance of electroencephalographs
IEC 60529:1989+ AMD1:1999+AMD2:2013	Degrees of protection provided by enclosures (IP code)
IEC 62366-1:2015	Medical devices - Part 1: Application of usability engineering to medical devices
IEC 62304:2006 + Amd 1:2015	Medical device software - Software life cycle processes + Amendment 1
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN ISO 15223-1: 2021	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
EN ISO 14155:2020	Clinical investigation of medical devices for human subjects – Good clinical practice
ISO 10993-1:2020	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Test for in vitro cytotoxicity
ISO 10993-10:2021	Biological evaluation of medical devices - Part 10: Tests for skin sensitization
ISO 10993-23:2021	Biological evaluation of medical devices - Part 23: Tests for skin irritation
MEDDEV 2.7/1 Rev. 4 (2016)	Clinical Evaluation: Guide for manufacturers and notified bodies

MEDDEV 2.12/1 rev 8 (2013)	Medical Devices Vigilance System
MEDDEV 2.12/2 rev 2 (2012)	Clinical Evaluation - Post Market Clinical Follow-up
MEDDEV 2.7/4 rev 12 (2010)	Guidelines on Clinical Investigation: A guide for manufacturers and notified bodies
MEDDEV 2.7/3 rev 5 (2015)	Clinical investigations: serious adverse event reporting
MDCG 2020-5 rev 1 (April 2020)	Clinical Evaluation - Equivalence A guide for manufacturers and notified bodies
ASTM D999: 2015	Test Methods for Vibration Testing of Shipping Containers
ASTM D4728: 2017	Standard Test Method for Random Vibration Testing of Shipping Containers
ASTM D5276: 2017	Standard Test Method for Drop Test of Loaded Containers by Free Fall
ASTM D4169: 2016	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM D4332: 2014	Standard Practice for Conditioning Containers, Packages or Packaging Components for Testing
ASTM F2213: 2017	Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment
ETSI EN 300-328: 2019	Electromagnetic Compatibility and radio spectrum matters (ERM) - Wideband Transmission Systems - Data transmission equipment operating In The 2.4 GHz ISM band and using wide band modulation techniques
IEC 62133-2 (2017)	Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems + Amendment 1
UN DOT 38.3	Safety of primary and secondary lithium cells and batteries during transport
HIPAA (1996)	Health Insurance Probability and Accountability Act
GDPR 2016/679 (2016)	General Data Protection Regulation
RoHS recast Directive 2011/65/EU (2011)	Restricting the use of hazardous substances in electrical and electronic equipment
WEEE Directive 2012/19/EU	Waste Electrical & Electronic Equipment (WEEE) Directive
ISO/IEC 27001:2013	Information technology - security techniques - information security management systems - requirements
ISO 27799:2016	Health informatics — Information security management in health using ISO/IEC 27002
FDA Guidance (2002)	General Principles of Software Validation; Final Guidance for Industry and FDA Staff
FDA Guidance (2016)	Applying Human Factors and Usability Engineering to Medical Devices

FDA Guidance (2005)	Guidance for the Content of Pre-Market Submissions for Software Contained in Medical Devices
FDA Guidance (2017)	Deciding When to Submit a 510(k) for a Change to an Existing Device
FDA Guidance (2016)	Use of International Standard ISO 10993-1, "Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process"
FDA Guidance (2019)	Safety and Performance Based Pathway
FDA Guidance (2013)	Radio Frequency Wireless Technology in Medical Devices
FDA Guidance (2018)	Content of Premarket Submissions for Management of Cybersecurity in Medical Devices (DRAFT)
FDA Guidance (2016)	Postmarket Management of Cybersecurity in Medical Devices
FDA Guidance (2016)	Unique Device Identification System: Form and Content of the Unique Device Identifier (UDI)
FDA Guidance (2017)	Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling
AAMI TIR 12 (2010)	Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
AAMI TIR 30 (2011 /(R)2016)	A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
AAMI TIR 69 (2017)	Risk Management of Radio-frequency Wireless Coexistence for Medical Devices and Systems
FCC part 15 Subpart C	FCC: Rules and regulations for Title 47 (Telecommunications)
EN 301 489-17 v3.2.4 (2020-09)	EMC Testing
FCC Part 15 Subpart B	EMC Testing
NBOG (2014)	Consensus paper for the Interpretation and Application of Annexes Z in EN ISO 14971 2012
NBOG BPG 2014-3 (2014)	Guidance for manufacturers and Notified Bodies on reporting of Design Changes and Changes of the Quality System