

PHANTOM XY SYSTEM

Technical Specifications & Clinical Performance Data

Dental Radiographic Quality Assurance | Powered by Science, Delivered by Subscription

87.63	98.71%	0.89%	29,367	53
Fleet-Wide Average Q-Score (out of 100)	Overall QA Test Pass Rate	Radiographic Retake Rate	QA Tests Performed	Dental Clinics Deployed

PRODUCT OVERVIEW

The Phantom XY System is a hardware-software quality assurance solution enabling dental clinics to perform standardized, objective digital quality assurance testing of digital panoramic, intraoral, and cephalometric radiographic equipment. The system comprises a precision-manufactured test phantom and an integrated cloud-based analytics platform (Phantom Cloud), together delivering continuous, manufacturer-agnostic radiographic quality assurance across all major digital dental X-ray platforms.

PHANTOM XY HARDWARE

Calibrated physical test phantom providing objective, reproducible image quality measurements across all digital X-ray modalities including intraoral sensors, panoramic OPG units, and cephalometric devices.

PHANTOM CLOUD PLATFORM

Cloud-based QA management delivering longitudinal performance trending, compliance documentation, regulatory reporting, and automated alert notifications when performance falls below diagnostic threshold.

Document Type	Prepared By	Date	Classification
Technical Specifications & Clinical Data	MVD Technologies Inc.	May 2026	Clinical / Sales

1. SYSTEM ARCHITECTURE

Hardware + Cloud Platform Overview

PHANTOM XY HARDWARE		PHANTOM CLOUD PLATFORM	
Physical Form Factor	Precision-manufactured calibrated test phantom	Platform Type	Cloud-based SaaS — no local installation required
Modality Compatibility	Intraoral sensors, panoramic OPG, cephalometric units	Access	Any device, any browser
Manufacturer Coverage	15+ manufacturers, 32+ device models validated	Data Logging	Real-time automated QA result logging
Setup Time	Under 2 minutes per device per test	Trending	Longitudinal performance trend analysis
Measurement Parameters	17 image quality parameters assessed per test	Alerts	Automated alerts when Q-Score falls below threshold
Q-Score Range	0–100 composite score; pass threshold ≥ 70	Reporting	Audit-ready regulatory compliance documentation
Hardware Cost	\$0 — included at no charge with 12-month subscription	Pricing Model	Subscription-based, predictable monthly cost
Operator Requirement	Any trained dental staff member	Data Coverage	Full audit trail; timestamped; exportable

Q-Score Methodology

The Q-Score is a composite metric (0–100) evaluating 17 image quality parameters from a standardized test phantom exposure. Parameters assessed include:

- Spatial resolution
- Density uniformity
- Artifact detection
- Electronic noise characterization
- Exposure calibration accuracy
- Contrast-to-noise ratio (CNR)
- Geometric accuracy
- Pixel dropout detection
- Software processing consistency
- Tube output consistency

A Q-Score ≥ 70 constitutes a passing result, indicating the radiographic device is performing within acceptable parameters for clinical diagnostic use. Failing devices are flagged immediately with detailed diagnostic recommendations delivered through the Phantom Cloud interface.

Standard QA Workflow — Under 2 Minutes Per Device

1	POSITION Place the Phantom XY test phantom in the radiographic field using standard positioning technique.
2	EXPOSE Acquire a standard exposure using the clinic's routine technique settings — no special exposure parameters required.
3	UPLOAD Upload the resultant image to the Phantom Cloud platform via any connected device — browser-based, no local software.

4

ANALYZE

Phantom Cloud automatically returns a quantitative Q-Score, pass/fail determination, trend data, and alerts within seconds.

2. DEVICE COMPATIBILITY

Panoramic & Intraoral Equipment — Validated Performance Data

The Phantom XY System has been validated across 15 distinct manufacturers and 32 unique radiographic device models, encompassing panoramic OPG units, intraoral sensors, and cephalometric systems. The following performance data were collected across 29,367 QA tests in 53 Ontario dental clinics over a 410-day observation period (January 2025 – March 2026).

Panoramic / OPG Device Performance by Manufacturer

Manufacturer	Model	Avg Q-Score	Pass Rate	Total Tests
Sirona	Orthophos S	94.87	98.76%	1,128
Sirona	Orthophos XG 3	90.97	99.14%	350
Sirona	Orthophos XG 5	92.14	99.07%	215
Sirona	Orthophos SL	95.99	95.89%	146
Sirona	Orthophos XG	94.86	100.00%	129
Carestream	CS 8100	80.50	99.74%	755
Carestream	CS 9600	94.69	99.18%	122
Planmeca	Pro Max	89.06	98.68%	606
Planmeca	Pro One	90.79	98.70%	307
Sorodex	Cranex	79.18	99.07%	430
Owandy	I-MAX	90.67	99.53%	213
Owandy	I-MAX Touch	88.72	99.37%	159
Kavo	OP 3D	91.05	100.00%	182
Instrumentarium	OP30	88.46	97.04%	203
Instrumentarium	OP300	75.83	100.00%	174
Acteon	X-Mind Prime	89.60	99.42%	172
Acteon	Trium	89.16	100.00%	118
Gendex	GXDP-300	80.60	100.00%	173
Gendex	GXDP-700	85.97	99.22%	128
Gendex	Orthoralix 9000	82.76	99.31%	144
Villa	Rotograph Evo	78.75	99.61%	255
Vatech	PaX-i3D Green	91.90	98.97%	97

Belmont	Bel-Cypher N	92.07	95.89%	73
Air Techniques	Pro Vecta S	92.51	98.48%	66

Table 1. QA performance metrics by panoramic/extraoral device. Q-Score range 0–100; pass threshold ≥ 70. Data from 29,367 tests across 53 Ontario clinics, Jan 2025–Mar 2026.

■ Q-Score ≥ 90: Excellent	■ Q-Score 80–89: Good	■ Q-Score 70–79: Acceptable (Monitor)	■ Q-Score < 70: FAIL — Immediate Action Required
---------------------------	-----------------------	---------------------------------------	--

Intraoral Sensor Market Distribution

Analysis of 96,861 radiographic exposures logged through Phantom Cloud reveals the intraoral sensor distribution across participating clinics:

Manufacturer	Models	Total Exposures	Market Share
Dexis	Titanium, Platinum, IXS	52,192	53.9%
Schick	33, 33 AE	8,593	8.9%
Carestream / Kodak	RVG 6200, RVG 5100	7,152	7.4%
Acteon	SOPIX	4,359	4.5%
Air Techniques	ScanX Plate	3,082	3.2%
Owandy	OPTEO	1,288	1.3%
Woodpecker	i-Sensor	929	1.0%
Planmeca	ProSensor	333	0.3%
Other / Unknown	Various	18,933	19.5%

Table 2. Intraoral sensor/imaging plate manufacturer distribution across 96,861 logged exposures.

The Phantom XY System's manufacturer-agnostic design enables consistent QA assessment across this highly fragmented sensor market, providing a uniform standard of image quality measurement regardless of the underlying hardware platform.

3. CLINICAL PERFORMANCE DATA

Retake Rate Analysis & Patient Safety Impact

Retake rate — the proportion of radiographic exposures repeated due to diagnostic inadequacy — is the primary measurable indicator of radiographic quality and its direct patient safety impact. Each retake represents an additional, unnecessary radiation exposure. Published literature reports typical retake rates of 5–10% in facilities without systematic QA programs, with some environments reporting rates as high as 15%.

0.89
%

System-Wide Retake Rate

96,861 total exposures | 866 retakes | 53 clinics

5–11× lower than industry-reported baseline of 5–10%

~3,977 to 8,820 unnecessary radiation exposures prevented during study period

Retake Reason Analysis

Categorized analysis of 866 retake events provides actionable intelligence for clinic quality improvement:

Retake Reason	Count	Percentage	Category
Patient moved during imaging	233	26.91%	Operator-dependent
Tooth apex not in frame	176	20.32%	Operator-dependent
Wrong positioning	129	14.90%	Operator-dependent
Overlapping	110	12.70%	Operator-dependent
Other	89	10.28%	Mixed
Cone cut	63	7.27%	Operator-dependent
Poor image quality	39	4.50%	Equipment-related
Software error during acquisition	20	2.31%	Equipment-related
Over/under exposed	7	0.81%	Equipment-related

Table 3. Categorized retake reason analysis, 866 retake events.

Equipment-related retake causes (poor image quality, over/underexposure, software errors) total only 7.62% of all retakes — directly validating the Phantom XY System's effectiveness in maintaining radiographic equipment within optimal operating parameters.

4. DEPLOYMENT & MARKET ADOPTION

Ontario Clinical Network — 410-Day Observation Period

53	24	171	29,367
Participating Clinics	Ontario Municipalities	Registered Clinical Users	QA Tests Performed
96,861	71.63	15	32
Radiographic Exposures	Avg QA Tests / Day	Manufacturers Tested	Unique Device Models

Quarterly Growth Trajectory

Quarter	New Clinics Enrolled	Cumulative Total	Q-o-Q Growth
Q1 2025 (Jan–Mar)	10	10	—
Q2 2025 (Apr–Jun)	8	18	-20%
Q3 2025 (Jul–Sep)	8	26	0%
Q4 2025 (Oct–Dec)	15	41	+87.5%
Q1 2026 (Jan–Mar)	12	53	-20%
Q1 2026 Additional*	+25	78+	—

*25 additional Ontario clinics confirmed enrolled in March 2026. Table 4. Quarterly clinic enrollment trajectory.

High-Compliance Clinic Performance — Top 15 by Weekly Testing Frequency

Clinic	Avg Tests / Week	Total Tests	Weeks Active
Monarch Fennell	49.88	1,696	34
Monarch Dentistry Bunting	41.26	1,444	35
Rowntree Gate Dental	34.60	173	5
Monarch Erin Mills	34.41	1,273	37
Monarch Dentistry Tremont	33.33	1,200	36
Monarch Kitchener	29.06	1,046	36
Monarch Dentistry Colborne	27.31	983	36
Bogart Dental	26.54	345	13
Royal City Dental	23.84	1,311	55
Monarch Niagara Kalar	21.62	800	37
Monarch Dentistry Centennial	21.42	707	33
Sun Dent Family Dental Etobicoke	20.68	1,158	56
Don River Dental	20.69	269	13
The Children's Dental Centre	20.06	1,023	51
Wexford Dental	19.35	600	31

Table 5. Top-performing clinics by weekly QA testing frequency.

5. REGULATORY COMPLIANCE

Health Canada Safety Code 30 Alignment

Health Canada Safety Code 30 (Radiation Protection in Dentistry) establishes comprehensive requirements for quality assurance in dental radiography. The Phantom XY System directly addresses all key Safety Code 30 requirements:

Safety Code 30 Requirement	Phantom XY Capability	Study Evidence
Routine QA testing of X-ray equipment	Automated QA testing with standardized phantom and quantitative Q-Score analysis	29,367 QA tests; 71.63 tests/day avg. across 53 clinics
Documentation and record-keeping of QA results	Cloud-based audit trail with timestamped results, device tracking, and exportable reports	410 days of continuous data collection with full traceability
Monitoring of radiographic retake rates	Automated retake tracking with categorized reason codes and trend analysis	0.89% retake rate across 96,861 exposures with detailed reason analysis
Equipment performance monitoring	Quantitative Q-Score trending with automated alerts for performance degradation	32 device models from 15 manufacturers continuously monitored
Radiation dose optimization (ALARA)	Identification of underperforming equipment that may cause overexposure or retakes	Est. 3,977+ unnecessary exposures prevented during study period
Staff training and quality awareness	Real-time feedback on equipment performance; data-driven quality culture	171 registered clinical users actively engaged in QA workflows

Table 6. Safety Code 30 compliance mapping for the Phantom XY System.

Additional Regulatory Frameworks Supported

- Health Canada Radiation Protection Bureau — periodic performance testing of dental X-ray equipment
- IAEA Safety Reports Series — quality control testing protocols including image quality assessment using test phantoms
- Canadian Dental Association (CDA) — evidence-based radiographic protocols including equipment performance monitoring
- Provincial Dental Regulatory Colleges — documented QA protocols as condition of dental facility licensure
- American Dental Association (ADA) / AAOMR — regular performance testing and image quality monitoring programs
- Ontario Occupational Health and Safety Act — radiation workplace requirements for occupational exposure reduction

6. SUBSCRIPTION MODEL & PRICING

Phantom XY System + Phantom Cloud

INCLUDED WITH SUBSCRIPTION

Phantom XY Hardware	Calibrated test phantom — \$0 equipment cost
Phantom Cloud Platform	Full cloud QA management, trending, alerts, compliance reporting
Unlimited QA Tests	No per-test fees — test as frequently as required

Multi-Device Coverage All intraoral sensors and OPG units at the facility

Regulatory Documentation Audit-ready reports for provincial regulatory inspections

The subscription model eliminates capital expenditure on QA hardware. The Phantom XY physical test phantom is provided at no charge with a 12-month subscription — removing all financial barriers to implementing evidence-based radiographic quality assurance in every dental practice.

7. SUPPORTING CLINICAL EVIDENCE

Peer-Reviewed Literature Summary

The clinical case for dental radiographic quality assurance is established by a robust body of peer-reviewed literature. The following publications form the evidentiary foundation for the Phantom XY System's clinical positioning:

PMC 11447568	Quality Assessment of Periapical Radiographs	<i>FGDP Guidelines Study, 2024</i>
Finding:	Only ~2/3 of intraoral radiographs meet diagnostic acceptability standards. Inadequate contrast and geometric distortion are the primary failure modes.	
PMC 11486910	Australian Dentists' Knowledge of Radiographic Interpretive Errors	<i>Journal of Dental Research, 2024</i>
Finding:	Interpretation errors cause direct, documented patient harm including unnecessary extractions, missed abscesses, false-positive diagnoses, and medicolegal liability.	
PMC 10891473	Image Quality Assessment — Hand-Held vs. Wall-Mounted Devices	<i>Dentomaxillofacial Radiology, 2024</i>
Finding:	Equipment drift causes image quality degradation that is not perceptible during routine clinical use without objective phantom-based testing.	
PMC 9974235	Factors Affecting Interpretation of Dental Radiographs	<i>BMC Oral Health, 2023</i>
Finding:	Radiographic artifacts misidentified as pathology drive unnecessary invasive treatment — a documented risk amplified by suboptimal image quality.	

Your patients deserve diagnostically reliable imaging. Every. Single. Time.

Contact us: info@mvdxraytech.com | 1-855-MVD-XRAY (683-9729) | www.mvdxraytech.com