


Case Study: μ Mass™ identified a human-relevant signal that sharpened the mechanism of action (MoA) before IND

Objective: suppress a functional human biomarker with free-uptake siRNA delivery in 3D


Step 1: 2D transfection optimisation



Transfection reagents are toxic at effective concentrations in primary cells

siRNA uptake in 2D primary human cells confirmed by fluorescent labelling. Standard reagents (RNAiMax, TransIT) showed dose-dependent toxicity as per expectations


Step 2: Free-uptake delivery in 3D



Free-floating siRNA integrated into primary, human 3D aggregates with μ Mass™

Free-uptake delivery confirmed by fluorescence across 2D-to-3D transition; reducing cellular toxicity vs. standard transfection. μ Mass™ retains siRNA in the aggregate environment throughout culture

Step 3: Functional disease readout



A trigger-induced functional biomarker, read across candidates and doses

Disease phenotype induced by a defined trigger, generating a 4-log rise in a functional biomarker. Biomarker suppression confirmed functional disease switching

Set up: 1 week

Uptake optimization: 2 weeks


Results assays: 1 week

*siRNA: small interfering RNA

*Free uptake: siRNA delivered without transfection reagent

*qPCR: quantitative polymerase chain reaction

Outcome: assumption challenged. Human 3D aggregates and μ Mass™ flagged the candidate the mouse data didn't



*MoA: mechanism of action, *siRNA: small interfering RNA, *Tr.: disease phenotype trigger

How it was tested

siRNA delivered by free uptake into human, primary 3D aggregates in μ Mass™
Disease phenotype was triggered across three concentrations, and functional human biomarker measured at Day 5

What the data showed


$>80\%$ phenotypic reduction of functional biomarker across all tested concentrations
siRNA#1 had no effect, unambiguous lead selection on a functional human endpoint.
The MoA observed in μ Mass™ differed from what the mouse model predicted

Next steps for the program

Characterize the MoA divergence — Verify potency in multiple donors — Test durability — Add mechanical load (Re-plate™)


Assay results within four weeks — GLP-like, SOP-based workflows — ISO-certified facilities

Case Study: μ Mass™ identified a human-relevant signal that sharpened the mechanism of action (MoA) before IND



Gain access to μ Mass™ for your program

Whether you need study execution, in-house platform access, or a tailor-made solution—chiron offers multiple ways to engage




Device Sales

Device sales via technology transfer: bring μ Mass™ into your own laboratory as a consumable platform; compatible with standard workflows and equipment



Contract Research Services

Outsource studies to chiron for rapid, expert execution, generating human-relevant, decision-ready data with no internal burden



Model Co-development

Co-develop a bespoke model with chiron that fits your specific therapeutic, disease, or translational application

Why work with chiron?

Human-relevant insights, for more informed decision making, resulting in lower a financial burden for drug development

Difficult-to-access insights, made measurable

- Direct comparison of multiple treatments and donors on one platform
- Cell-therapy and advanced-modality testing where animal models are poorly predictive
- 3D solid tumour models that better capture tumour structure

More learning, at a lower cost per data point

- Generate large, comparable datasets that can be fed to predictive platforms
- Enable leaner clinical trials, reducing the development timeline
- Rapid application development and assay processing

Thank You

