



PORTFOLIO
2025

WWW.SPHERE-ENERGY.EU

NEXT GENERATION BATTERY CELL TESTING



20
25

TABLE OF CONTENTS




FOCUS AREAS IN SSB RESEARCH	01	ASC-P: POUCH CELL BATTERY TESTING	12
PRECISION IN RESEARCH	02	ASC-P+: ADVANCED POUCH CELL BATTERY TESTING	13
SOLID-STATE & POUCH CELL SOLUTIONS	03	ASC PRESSURE FRAME	18
EXPLORING OUR SETUPS	04	ASC-P MAPPING SETUP	19
ASC-T+: ADVANCED TESTS	06	PRESSURE & TEMPERATURE SOLUTIONS	20
ASC-X: SCREENING TESTS	07	THICKNESS MONITORING SOLUTIONS	21
ASC-AD: PRESSURE BALANCING	08	ASC PACKAGES	22
ASC-AD+: LOW-PRESSURE RANGES	09	ASC PUBLICATIONS	26
ASC-G: GAS ANALYSIS	10		





FOCUS AREAS IN SOLID-STATE BATTERY RESEARCH

The field of solid-state battery research often faces the challenge of conducting electrochemical tests on complex samples under specific mechanical conditions. Sphere Energy has responded to this challenge by designing a range of testing equipment that is customized to meet the requirements of experiments in this area of research. Our equipment is designed to provide support in various aspects, including sample preparation, pressure, and thickness monitoring, making testing more convenient and streamlined.




Key aspects for reproducibility and data reliability


Pressure tracking from the initial state of sample preparation until the end of the electrochemical tests.

Direct access to relevant data


We have designed our setups to offer you convenient access to essential data while conducting tests on solid-state battery materials. Through our testing cells, pressure can be monitored not only during the sample preparation protocol but also along the entire electrochemical test phase, generating important insights about the electro-chemo-mechanical performance of your studied systems.



Pressure monitoring data of electrochemical tests.



Sample size (thickness) monitoring during electrochemical tests.



EIS measurements of a sample over different temperatures.

Pressure, thickness, and temperature are some of the parameters covered by our setup. You will find many other options to enrich your R&D activities in the following pages of this catalog. Enjoy it!



SIMPLIFYING RESEARCH WITH PRECISION AND PERFORMANCE

Compact & Lab-Optimized Design

- Product dimensions designed for optimal use under lab conditions.
- Perfectly suitable for glove-box handlings and transfers.
- User-friendly assembly.

Simple usability under difficult conditions

- Reproducible results by monitoring all important parameters during testing and preparation steps.
- Prepare battery materials at their measurement location to avoid cross-contamination and pellet damage.
- A broad range of features covering: airtight cells, fixed distance pressurization, dynamic pressure management, gas analysis, temperature control, resistance measurement probe, reference-electrode, pressure and thickness monitoring.
- Simplified sample material loading thanks to the symmetric cell design.
- Long-lasting design thanks to high-end materials and replaceable wear parts.
- Time-saving mechanisms paving the way to high throughput battery testing

Modular approach

- The ASC family is made from intercompatible products - so you can easily enlarge your testing possibilities whether it's by integrating new modules, expanding sample sizes or incorporating real-time environmental monitoring systems.





COMPLETE SOLUTIONS FOR SOLID-STATE & POUCH CELL TESTING

01 ELECTROCHEMICAL TEST UNDER PRECISE SAMPLE GEOMETRY AND PRESSURE CONTROL

PRESSURE BEHAVIOUR OF SAMPLES DURING TEST

All ASC setups facilitate solid-state battery sample preparation by enabling in-situ pelletization of complex material layers. Integrated pressure monitoring allows continuous tracking during sample formation and testing.

Adapted setups: **ASC-T+ / **ASC-AD+** / **ASC-X** / **ASC-AD** / **ASC-G** / **ASC-C** / **ASC-P+***

**In combination with ASC pressure frame*

02 SWELLING TESTS UNDER FIXED PRESSURE PROTOCOLS

VOLUME EXPANSION AND CONTRACTION DURING ELECTROCHEMICAL TESTS

In order to get insights into breathing and swelling behavior of battery cells during electrochemical tests, our setups can be combined with thickness monitoring options and the pressure can be applied via springs for fixed pressure testing protocols.

Adapted setups: **ASC-T+ / **ASC-AD+** / **ASC-X** / **ASC-P+***

03 IONIC CONDUCTIVITY MEASUREMENTS

EIS MEASUREMENTS WITH PRECISE PRESSURE CONTROL AND THICKNESS MONITORING

Our systems enable accurate monitoring of sample dimensions under applied pressure, allowing for precise correlation with impedance values and enabling exact estimations of ionic conductivity in your experiments.

Adapted setups: **ASC-T+ (practical for Arrhenius plots) / **ASC-X***

04 ADVANCED ELECTROCHEMICAL TESTS WITH REFERENCE-ELECTRODE

ANODE AND CATHODE DECONVOLUTED VOLTAGE INFORMATION

The ASC-T+ is equipped with metallic-based reference-electrode embedded in the internal insulation sleeve and placed at the center of the sample.

Adapted setups: **ASC-T+*

05 ADVANCED CHARACTERIZATION VIA GAS ANALYSIS

REALTIME GAS ANALYSIS OF DECOMPOSITION REACTIONS DURING ELECTROCHEMICAL TESTS

The ASC-G's built-in gas flow channels enable direct and efficient interfacing with gas analysis equipment for comprehensive gas evaluation in electrochemical tests.

Adapted setups: **ASC-G*



COVERING KEY ASPECTS OF YOUR RESEARCH


OVERVIEW OF PARAMETERS COVERED BY SPHERE'S SETUPS

Type	Technical purpose	Technical purpose	Pressure application	Airtightness	Max pressure/force	Pressure monitoring	Swelling monitoring
ASC-X		Screening tests	Fixed thickness	Yes	400 MPa	Optional	Optional
ASC-AD		Long term testing	Fixed thickness	Yes	400 MPa	No	No
ASC-G		Mass spectroscopy	Fixed thickness	Yes	400 MPa	Optional	Optional
ASC-T+		Advanced tests and high temperature	Fixed thickness Or fixed pressure	Optional	400 MPa	Yes	Optional
ASC-AD+		Low pressure range	Fixed thickness Or fixed pressure	Yes	100 MPa	Yes	Optional
ASC-P		Pouch format	Fixed thickness fixed pressure or hybrid	No	20 kN	Optional	Optional
ASC-P+		Pouch format	Fixed thickness fixed pressure or hybrid	No	20 kN	Optional	Optional



FOR SOLID-STATE BATTERY TESTING EQUIPMENT

**Precision Testing for
Innovative Energy Materials**



Dimensions: 210 x 120 x 100 mm
Weight: 5.5 kg
Pressure range: 0-400 MPa
Temperature range: 25 to 200°C

The ASC-T+

Advanced tests on solid-state batteries and battery materials

The ASC-T+ is a versatile setup designed specifically for tests on solid-state battery research. It provides great flexibility to meet various research requirements, which is crucial in this fast-evolving battery industry.

The setup comes with a built-in temperature control system, ensuring accurate experiments up to 200°C. Pressure application and monitoring are achieved through the use of a mechanical press (pressure frame) also provided with the setup.

Additionally, the system is capable of accommodating multiple sample sizes, ranging from 8 mm to 14.5 mm, resulting in a 3.3 times increase in the active surface area. This adaptability is essential for the development of new materials and advanced battery systems.

MAIN FEATURES


- Multiple sample sizes
- Pressure Monitoring
- Temperature control
- Exchangeable int. sleeves
- Reference electrode
- Compact dimensions
- Fixed pressure mode


ADD-ON OPTIONS

- Fixed thickness mode
- Swelling monitoring
- Multiple spring constants
- 4-point probe resistivity

ASC-T+ setup contains:

- 1x ASC-T+ test cell equipped with built-in temperature control system
- 1x Set of pistons and internal sleeves for multiple sample diameters (8, 10, 12 and 14.5mm)
- 1x Ceramic sleeve for samples of 8 mm of diameter
- 1x Internal sleeve with reference-electrode based on Cu
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Pressure frame equipped with a pressure sensor for material pressing and pressure monitoring
- 1x Control unit equipped with a temperature control system, pressure display and analogical output for pressure data transmission to potentiostat, cycler or computer
- 1x Glove box feedthrough adaptor, set of cables and connection leads





Dimensions: 65 x 65 x 57 mm

Weight: 750 g

Pressure range: 0-400 MPa

Temperature range: -20 to 120°C

The ASC-X

Precision Airtight Testing with Interchangeable Pistons

The ASC-X is a new generation electrochemical cell developed to meet the rigorous demands of solid-state battery R&D. Designed for precision under extreme conditions, ASC-X enables researchers to apply and control pressures of up to 400 MPa within an airtight environment—crucial for probing the mechanical and electrochemical interplay in solid-state systems. Its flexible sample holder accommodates diameters from 8 to 12 mm, supporting both material-scarce fundamental studies and signal-intensive advanced testing. Whether working with minimal material quantities or exploring larger interfaces, ASC-X adapts to the evolving challenges of solid-state innovation.

What truly sets ASC-X apart is its redesigned sealing mechanism that decouples the pressure applied to the sample from that applied to the ferrules. This novel architecture ensures highly accurate and consistent pressure control, essential for reproducibility and in-depth analysis. Its small form factor and simplified assembly protocol reduce setup time, enabling more consistent results under controlled conditions. With ASC-X, researchers can now focus on what matters most: unlocking the next breakthroughs in solid-state battery performance.

ASC-X setup contains:



- 1x ASC-X test cell
- 1x Set of pistons and internal sleeves for multiple sample diameters (8, 10, 12 mm)
- 1x Torque wrench for manual pressure setting
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Set of tools, sealing ferrules and cables


MAIN FEATURES

- Airtight Precision Design
- Rapid Interchangeable Pistons (8, 10, 12 mm)
- Fixed thickness mode
- Exchangeable int. sleeves
- Enhanced Reproducibility

ADD-ON OPTIONS

- Pressure monitoring
- Swelling monitoring
- Fixed pressure mode
- Multiple spring constants





Dimensions: 65 x 65 x 60.5 mm
Weight: 770 g
Sample preparation: 0-400 MPa
Test with constant P: 5-100 MPa
Temperature range: -20 to 120°C

The ASC-AD

Airtight electrochemical cell for tests under constant pressure

The ASC-AD is a practical setup for electrochemical tests under stable conditions of applied pressure. This setup offers an active pressure balancing system, ensuring that the initially set stacking pressure is kept constant during the expansion or contraction of the tested materials. This is key for optimal results in long-term electrochemical cycling.

The dynamic pressure applied during the tests can be perfectly customized from 5 to 100 MPa via a variety of plungers and a selection of springs having different spring coefficients. During the sample preparation phase, the setup can apply up to 400 MPa using the pressure tool delivered with the setup. This setup is also airtight based on Swagelok ferrules and can be used for testing air-sensitive materials outside the glove box.

The ASC-AD can be used in combination with the pressure frame for sample preparation only, for example during the stage of sample densification inside the glove box. Once the cell is assembled with the sealing ferrules, it becomes a stand-alone testing cell for electrochemical tests.

ASC-AD setup contains:


- 1x ASC-AD electrochemical test cell
- 1x Compression tool
- 3x Dynamic plungers with different lengths
- 2x Internal sleeves in PEEK for samples of 8 mm of diameter
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Set of tools, springs, sealing ferrules and cables


MAIN FEATURES

- Airtight
- Large pressure range
- Fixed pressure mode
- Sample diameter 8 mm
- Exchangeable int. sleeves

ADD-ON OPTIONS

- Multiple spring constants





Dimensions: 90 x 90 x 145 mm

Weight: 2.2 kg

Pressure range: 0-100 MPa

Temperature range: -20 to 75°C

The ASC-AD+

Airtight electrochemical testing cell tailored for low-pressure conditions and under multiple cell configurations

The ASC-AD+ is our latest product for solid-state battery testing, designed to effortlessly explore multiple sample dimensions while providing crucial data on pressure and swelling through a streamlined, user-friendly interface.

This airtight setup allows for increasing the surface area of the samples up to 3.3-fold. The built-in pressure sensor generates real-time pressure information during the tests. All this with the possibility of tuning the setup for experiments under various pressure values, from a few kPa to 100 MPa, via the fixed thickness mode (pistons blocked) or via the constant pressure mode (piston connected with spring).

The ASC-AD+ setup can also be equipped with a thickness monitoring system for precisely tracking the swelling of battery cells, providing valuable insights into the mechanical properties of battery materials.

MAIN FEATURES


- Multiple sample sizes
- Pressure monitoring
- Swelling monitoring
- Built-in pressure sensor
- Fixed thickness mode
- Fixed pressure mode
- Exchangeable int. sleeves
- Multiple spring constants


ADD-ON OPTIONS

- Swelling monitoring

ASC-AD+ setup contains:

- 1x ASC-AD+ electrochemical test cell
- 1x Control unit equipped with pressure display and analogical output for pressure data transmission to potentiostat, cycler or computer
- 1x Set of pistons and internal sleeves for multiple sample diameters (8, 10, 12 and 14.5mm)
- 1x Glove box feedthrough adaptor
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Set of loading tools and adaptors for sample preparation
- 1x Set of tools, springs, sealing o-rings and cables





Dimensions: 65 x 65 x 56 mm
Weight: 740 g
Pressure range: 0-400 MPa
Temperature range: -20 to 80°C

The ASC-G

Airtight electrochemical test cell equipped with gas flow channels for combined gas analysis

The ASC-G is designed to combine the functionalities of our setups with an external gas analysis characterization technique, enabling real-time analysis of decomposition reactions and gasses released during your electrochemical measurements.

The setup is equipped with four gas channels that can be used in multiple configurations, for example for separating the gasses from anode and cathode reactions. The gas flow circulating in the channels is distributed inside the cell through a small gap placed between the pistons and the internal wall of the cell.

The setup comes with a calibrated torque wrench and a table that correlates torque values with pressure, which enables precise adjustment of the internal pressure at the sample level.

ASC-G setup contains:


- 1x ASC-G electrochemical test cell
- 1x Torque wrench for manual pressure setting
- 4x Inlet/outlet tubes in PEEK with 50 µm of internal diameter
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Set of tools, sealing ferrules, channel stoppers, cables and accessories

MAIN FEATURES

- Gas flow channels
- Airtight
- Large pressure range
- Fixed thickness mode
- Sample diameter 8 mm

ADD-ON OPTIONS

- Pressure monitoring
- Swelling monitoring
- Fixed pressure mode
- Multiple spring constants





FOR POUCH CELL TESTING EQUIPMENT

**Optimized Solutions for Pouch
Cell Evaluation**



The ASC-P

Battery pouch cell testing setup


The ASC-P is a specialized setup for battery researchers focused on prototyping, enabling precise control over the pressure applied to pouch cells. Its design allows for fine-tuning of mechanical stress simulations through adjustable pressure screws and springs, crucial for testing battery durability and performance. This system simplifies the electrical connection to testing devices with its intuitive tab contacts, ensuring reliable data collection.

In its versatile operation, the ASC-P features two modes: fixed thickness and fixed pressure, each catering to different experimental needs. The fixed thickness mode offers insights into the cell's behavior under constant dimensions, while the fixed pressure mode, with adjustable spring coefficients, accommodates variations in cell thickness, ensuring consistent pressure. This adaptability makes the ASC-P an indispensable tool for researchers aiming to push the boundaries of battery innovation.

Additionally, the ASC-P offers the possibility for precise pressure control and monitoring when coupled with the ASC pressure frame, or via a built-in pressure sensor. The setup can also be equipped with a thickness monitoring system for tracking the swelling of the cells assessing valuable information about volume expansion.

ASC-P setup contains:

- 1x ASC-P testing setup
- 1x Torque wrench for manual pressure application
- 1x Set of medium-force springs (spring rate 216N/mm)
- 1x Set of low-force springs (spring rate 98N/mm)
- 1x Set of extra low-force springs (spring rate 32N/mm)
- 1x Set of tools, screws and accessories



Dimensions: 157 x 112 x 70.5 mm

Weight: 3.8 kg

Force range: 0-20 kN


Temperature range: -20 to 80°C


MAIN FEATURES

- Fixed thickness mode
- Fixed pressure mode
- Multiple spring options
- Facilitated tab connection
- High force (up to 20 kN)
- Max cell size 108 x 75 x 8
- Compatible with ASC pressure frame

ADD-ON OPTIONS

- Pressure monitoring
- Swelling monitoring
- Built-in pressure sensor





Dimensions: 157 x 112 x 100 mm

Weight: 4.8 kg

Force range: 0-20 kN

Temperature range: -20 to 80°C

The ASC-P+

Precision Testing with Real-Time Pressure Control

Meet the ASC-P+, our compression jig dedicated to pouch and prismatic cell testing. Designed for precision, this upgraded setup comes equipped with a built-in pressure monitoring system that provides real-time data – giving you unparalleled control over pressure dynamics throughout your cell testing. With this system, you'll gain immediate insights, helping you fine-tune your experiments and optimize battery cell performance.

The ASC-P+ adapts to your research needs with fixed pressure and fixed thickness modes, along with adjustable spring options that help you maintain consistent pressure during cell swelling. Whether you're studying cell degradation, mechanical stress, or performance under varied conditions, the ASC-P+ ensures accurate and reproducible results.

The ASC-P+ is a compact system engineered for efficiency. It combines robust performance with a user-friendly design, providing cutting-edge additional functionalities such as swelling monitoring via thickness, pressure and temperature mapping, and control over cell heat extraction. For researchers seeking to demonstrate breakthroughs in battery cell technology over prototype cell sizes, the ASC-P+ offers the precision and adaptability required to advance your work.

ASC-P+ setup contains:


- 1x ASC-P+ testing setup
- 1x Built-in pressure sensor
- 1x Analogical data output
- 1x Set of medium-force springs (spring rate 216 N/mm)
- 1x Set of low-force springs (spring rate 98N/mm)
- 1x Set of extra low-force springs (spring rate 32 N/mm)
- 1x Set of tools, screws and accessories

MAIN FEATURES

- Pressure monitoring
- Fixed thickness mode
- Fixed pressure mode
- Multiple spring options
- Facilitated tab connection
- High force (up to 20 kN)
- Max cell size of 108 x 75 x 8 mm

ADD-ON OPTIONS

- Swelling monitoring
- Pressure mapping
- Temperature mapping





Pressure applied via a fixed thickness mode

The ASC-P+ offers users the flexibility to select different methods of applying pressure to cells during testing. The "fixed thickness mode" is designed to maintain a constant distance between the compression plates, secured by pressure screws. In this setup, the cell's swelling is restricted by the fixed position of the plates, leading to significant pressure variations throughout charge and discharge cycles. The plot below illustrates these pressure fluctuations in a Li-ion battery cell.




Fig. 1 - Setup assembled for tests under the "fixed thickness" protocol (left), where the plates are locked in place after applying the desired pressurizing force. Charge and discharge profiles, along with the force response to cell swelling, are shown on the right.

Pressure applied via a fixed force mode

In contrast, the "fixed force mode" is a testing protocol in which pressure is stabilized by incorporating springs into the testing setup. The springs act as a mechanical buffer, absorbing changes in thickness caused by cell swelling with minimal variation in applied pressure. This simple and affordable yet powerful feature enables pressure stabilization without the need for complex automated setups. Additionally, the use of springs allows for coupling the tests with in-operando thickness monitoring, providing valuable insights into battery cell behavior.




Fig. 2 - Setup assembled for tests under the "fixed force" protocol (left), with the jig equipped with pressurizing springs. Charge and discharge profiles, along with the stabilized force response that actively compensates for cell swelling, are shown on the right.



Multiple options of pressurizing springs

Pressure stabilization has a significant impact on electrochemical performance and battery pack design. The key factor behind pressure stabilization is the spring coefficient of the springs used in the setup. This critical parameter can be adjusted by selecting different springs from the wide range of spring coefficients offered by Sphere Energy.



Precise characterization of volume changes at the cell level




Battery performance is shaped by the electrochemical evolution with impact on cell geometry and mechanics – swelling, breathing, and structural changes that directly impact efficiency, safety, and lifespan. The ASC-P+ provides high-resolution, in-operando thickness monitoring to reveal how cells evolve under real-world conditions.

Tracking Long-Term Swelling Trends

Battery aging and capacity fade are often linked to structural transformations at the cell level. Battery cells gradually swell over time due to SEI formation, gas generation, metal plating and electrochemical grinding. ASC-P+ enables real-time tracking of cell thickness, helping researchers correlate swelling with electrochemical performance and develop strategies to mitigate long-term degradation.


Fig. 3 - Cell thickness evolution during cycling under fixed force testing protocol.



Electrochemical Breathing: The Hidden Dynamics of Charge-Discharge Cycles

During cycling, cells undergo repeated expansion and contraction due to the reversible volume changes at the electrode level. This breathing effect plays a key role in mechanical fatigue and interfacial stability. By quantifying these subtle deformations in real time, ASC-P+ provides key insights for cell design, pressure management, and reliability improvements.


Fig. 4 - Voltage and cell thickness during cycling under fixed force testing protocol.



Sub-Micron Sensitivity for In-Operando Characterization


Mechanical changes at the micro-scale often dictate macroscopic battery behavior. The ASC-P+ detects thickness variations as small as 0.1 μm , revealing fine mechanical fluctuations often undetectable. This level of precision allows researchers to correlate mechanical responses with electrochemical performance, uncovering stability shifts and early signs of degradation.

Fig. 5 - Zoom into the thickness data collected during cell testing showing the capability of the system to detect minimal changes in the cell volume.





Detecting homogeneity and evolution for improved cell performance and safety




Advanced characterization for pouch and prismatic cells

Understanding how pressure and temperature distribute across a battery cell is essential for optimizing performance, enhancing durability, and ensuring safety. Even small variations in mechanical or thermal conditions can influence ion transport, electrolyte behavior, and electrode stability, ultimately shaping a cell's efficiency and lifespan. By mapping these critical parameters, researchers can move beyond static testing and gain a real-time, spatially resolved view of how a cell responds under different conditions. The ASC-P+ can be equipped with up to **56 surface sensors** enabling a complete characterization of pressure and temperature distribution at the cell level.


Attesting pressure distribution on ASC-P+

Not all cells distribute pressure evenly — cell geometry, electrode stacking, and internal structure all contribute to variations in mechanical stress. Ensuring that your test equipment is capable of providing homogenous pressure distribution is key for characterizing the performance of battery cells.

Here below you can find the pressure distribution of the ASC-P+ measured assessed via prescale pressure films (Fujifilm), revealing how force is applied across: i. metallic plates, ii. pressure pads and a iii. pouch cell samples.



Pressure distribution on bare metallic plates



Pressure distribution on ASC-P+ equipped with pressure pads



Pressure distribution on a pouch cell sample using pressure pads



Sphere

ADD-ONS & ACCESSORIES

**Precision Control for Pressure
and Temperature**



MAX FLEXIBILITY THROUGH THE ASC PRESSURE FRAME



In solid-state battery R&D, precise control over a wide range of parameters is critical for consistent innovation. Recognizing the importance of accurate pressure management, our test cells are designed for seamless integration with our mechanical press, the pressure frame. This tool is capable of delivering and monitoring up to 20 kN of force, all within a compact design that easily fits inside the small antechamber of a glove box.





Compactability of ASC Pressure frame with the main cells

The pressure frame plays a crucial role in enhancing the flexibility of our ASC setups. It smoothly integrates all the different variations of cell tests under the same pressure monitoring platform, allowing standard and accurate measurements of pressure under different tests applications.



ASC-P – PRESSURE AND TEMPERATURE MAPPING SETUP



✓ Mapping sensors directly in contact with the cell body



✓ 56 detection points



Advanced characterization for pouch and prismatic cells

The distribution of pressure and temperature in pouch and prismatic cells significantly affects their performance and longevity. Characterizing these parameters is crucial for optimizing the operating conditions of battery cells. This optimization leads to higher energy storage capacities, better swelling management, and ensures that the cells operate safely.

Our ASC-P and ASC-P+ setups, designed for pouch and prismatic cell testing, can be equipped with pressure and temperature mapping systems. These mapping systems utilize flexible circuits that can accommodate a large number of sensors to monitor pressure and temperature distribution during battery cell testing. Developed in collaboration with Flexoo (Germany), the mapping systems consist of thin pads that are incorporated into Sphere's testing jigs, ensuring direct contact with the surface of the cells.



Pressure: 1kPa to 5 MPa / 20° to 100°C
Temperature: <1°C accuracy on T_{meas} . / <3% accuracy on P_{meas} .

A complete turnkey solution for advanced insights

The pressure and temperature mapping setup provides a comprehensive solution for assessing the uniformity of pouch and prismatic cells. Signals collected from multiple sensors are processed using specialized hardware and software, enabling the measurement of up to three sensor maps simultaneously.



PRESSURE MONITORING AND TEMPERATURE CONTROL

The control box plays a vital role in powering and processing all the information collected from the pressure and temperature sensors. It displays the value of the pressure applied to the sample in real-time. The data can be transmitted to external devices via an analog output, ensuring full synchronization of pressure information with the electrochemical data from a potentiostat/cycler. This integration is important for achieving optimal performance and accuracy. The control box is also tasked with monitoring and regulating the temperature of the ASC-T+ setup, which is crucial for exploring advanced tests in the field of solid-state batteries.

Two different options of control boxes are available depending on your application:



ASC CONTROL BOX with pressure frame and ASC-T+ cell

1

1. ASC CONTROL BOX FOR TEMPERATURE

This control box is designed to monitor the mechanical pressure applied to the samples inside the different cells and to control the ASC-T+ temperature. It includes an analogical output for external communication of mechanical pressure values.

PROGRAMMABLE PID OPTION

The ASC CONTROL BOX can also be equipped with a programmable PID that receives temperature control protocols directly from a computer connection.



ASC CONTROL BOX FOR PRESSURE with pressure frame and ASC-T+ cell

2

2. ASC CONTROL BOX FOR PRESSURE

This simplified control box is designed solely for monitoring the mechanical pressure applied to samples within various cell types. It features an analog output for the external communication of mechanical pressure values.



THICKNESS MONITORING SOLUTIONS

Sample thickness plays a critical role in the research and development of solid-state batteries, particularly when studying solid-state electrolytes or swelling phenomena. Our options for thickness measurement provide a range of solutions that enable you to precisely assess this parameter. You can measure sample thickness in-situ and in-operando within most ASC cells using our pressure frame.



1

1. IN-OPERANDO THICKNESS MONITORING SYSTEM

The thickness measurement relies on a displacement sensor attached to the pressure frame through an articulated arm. This setup is designed for in-operando thickness measurements, featuring data display and automated data recording via software or analogical output. It offers a precision of 1 μm and 0.1 μm of resolution.



2

2. IN-SITU THICKNESS MONITORING SYSTEM

The thickness measurement is based on a micrometer attached to the pressure frame through an articulated arm. This setup is specifically designed for in-situ thickness measurements, with a data display integrated into the probe. It achieves a precision of 3 μm and a resolution of 1 μm .



3

3. WIRELESS CONNECTION FOR THICKNESS OPTION

Wireless connection between micrometer and computer for punctual data recording. Each time the triggering button is pressed, a data point is sent directly to an Excel table. The product package includes a transmitter, receiver, and necessary cables.



ASC PACKAGE

ADVANCED ELECTROCHEMICAL TESTS WITH HIGH-THROUGHPUT



Screening tests / Long term testing / Advanced tests and high temperature / Fixed thickness / fixed pressure / Airtightness / High pressure / Pressure monitoring

The ASC Package is an essential toolkit for conducting comprehensive research related to solid-state batteries. It comprises the ASC-T+, designed for temperature- and pressure-controlled electrochemical tests; the ASC-X, a next-generation airtight cell offering precise pressure control up to 400 MPa with interchangeable pistons (8, 10, and 12 mm); and the ASC-AD, utilized for dynamic pressure control. Moreover, it is equipped with a mechanical press (pressure frame), a control box (control unit), a calibrated torque wrench, and all the necessary accessories required for smooth handling. With the ability to run three electrochemical cells in parallel, it enables detailed analysis of battery materials while ensuring optimal throughput.

ASC PACKAGE contains:

- 1x ASC-T+ electrochemical test cell with temperature control and pressure monitoring
- 1x ASC-X airtight electrochemical test cell
- 1x ASC-AD airtight electrochemical test cell with dynamic pressure control system
- 1x Pressure frame (mechanical press)
- 1x ASC control box (control unit)
- 1x Calibrated torque wrench
- 1x Glove box feedthrough adaptor, set cables and connection leads
- 1x Set of tools, cables and accessories



ASC+ PACKAGE

ADVANCED ELECTROCHEMICAL TESTS WITHIN LARGE PRESSURE RANGES



Advanced tests and high temperature / Low pressure range / Fixed thickness / Fixed pressure / High & low pressure / Pressure monitoring / Swelling monitoring (Optional)

The ASC+ package represents the pinnacle of our testing solutions, amalgamating the most sophisticated setups in our portfolio: the ASC-T+ and ASC-AD+. This package is best suited for deep analysis of key parameters of solid-state battery R&D. The ASC-T+ is specifically designed to conduct temperature-controlled tests and can achieve a maximum temperature of 200°C and 400 MPa pressure. Conversely, the ASC-AD+ is optimized for low-pressure tests in an airtight environment. Both setups provide real-time pressure monitoring data and facilitate test upscaling with multiple sample sizes.

ASC-T+ setup contains:

- 1x ASC-T+ test cell equipped with built-in temperature control system
- 1x Set of pistons and internal sleeves for multiple sample diameters (8, 10, 12 and 14.5mm)
- 1x Ceramic sleeve for samples of 8 mm of diameter
- 1x Internal sleeve with reference-electrode based on Cu
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Pressure frame equipped with a pressure sensor for material pressing and pressure monitoring
- 1x Control unit equipped with a temperature control system, pressure display and analogical output for pressure data transmission to potentiostat, cycler or computer
- 1x Glove box feedthrough adaptor, set of cables and connection leads*

ASC-AD+ setup contains:

- 1x ASC-AD+ electrochemical test cell
- 1x Control unit equipped with pressure display and analogical output for pressure data transmission to potentiostat, cycler or computer
- 1x Set of pistons and internal sleeves for multiple sample diameters (8, 10, 12 and 14.5mm)
- 1x Glove box feedthrough adaptor
- 1x Sample removing tool
- 1x Holder for weighing and loading materials in the cell
- 1x Set of loading tools and adaptors for sample preparation
- 1x Set of tools, springs, sealing o-rings and cables



ASC-X PACKAGE

PRESSURE MONITORED SOLID-STATE BATTERY TESTS



Screening tests / Fixed thickness / Airtight / High Pressure / Swelling monitoring (Optional)

The ASC-X package is specifically designed for testing air-sensitive materials outside the glovebox, combining mechanical compression with precise, reproducible pressure control in an airtight environment. This package includes the ASC-X electrochemical test cell—featuring interchangeable pistons (8, 10, and 12 mm), enhanced sealing architecture, and compatibility with add-ons like pressure and swelling monitoring—a pressure frame, and a control box. It offers a robust and flexible solution for advanced solid-state battery research, whether working with limited sample quantities or exploring larger interfaces.


ASC-X PACKAGE contains:

- 1x ASC-X airtight electrochemical test cell
- 1x Pressure frame (mechanical press)
- 1x ASC control box (control box for pressure)
- 1x Calibrated torque wrench
- 1x Glove box feedthrough adaptor, set cables and connection leads
- 1x Set of tools, cables and accessories



ASC-P PACKAGE

ADVANCED TESTS ON POUCH CELL BATTERIES



Pouch format / Fixed thickness / Fixed pressure or hybrid / Swelling monitoring (Optional) / Pressure mapping (Optional) / Temperature mapping (Optional)

The ASC-P package offers essential tools for battery researchers focused on prototyping. The ASC-P enables precise control over mechanical stress simulations and pressure adjustments, critical for investigating battery durability and optimizing performance under varying conditions. With its dual operational modes (fixed thickness / fixed pressure), it adapts to different experimental needs. The integration with the Pressure Frame and ASC Control Box facilitates consistent pressure monitoring and enhanced data reliability. This comprehensive approach supports researchers in advancing battery technology through rigorous, accurate tests under controlled key parameters.

ASC-P PACKAGE contains:

- 1x ASC-P testing setup
- 1x Pressure frame (mechanical press)
- 1x ASC control box (control box for pressure)
- 1x Torque wrench for manual pressure application
- 1x Glove box feedthrough adaptor, set cables and connection leads
- 1x Set of medium-force springs (spring rate 216N/mm)
- 1x Set of low-force springs (spring rate 98N/mm)
- 1x Set of extra low-force springs (spring rate 32N/mm)
- 1x Set of tools, cables and accessories



PUBLICATIONS FEATURING ASC SETUPS

Neumann, T., Cerny, A. T., Kalthoff, S., Eisele, L., Krossing, I., & Dsoke, S. (2025). Chemical Prelithiation of Silicon Powder and its Role as Anode Material for All-Solid-State Batteries. *Batteries & Supercaps*, 2500332.

Yang, J., Schaller, M., Cherkashinin, G., Zhang, R., Indris, S., Corte, D. A. D., ... & Strauss, F. (2025). Glass-Ceramic Lithium Thiophosphate Electrolytes with Enhanced Conductivity and (Chemo) mechanical Properties for All-Solid-State Batteries. *Chemistry of Materials*.

Sivavec, J., Baumgärtner, J. F., Stoian, D. C., Klimpel, M., van Beek, W., Kovalenko, M. V., & Kravchyk, K. V. (2025). On the Feasibility of Pairing Pyrochlore Iron (III) Hydroxy Fluoride Cathode with Argyrodite Li₆PS₅Cl Solid-State Electrolyte for Low-Cost All-Solid-State Batteries. *Batteries & Supercaps*, e202400731.

Shen, F., McGahan, M., Pietras, J. D., Doeff, M. M., Battaglia, V. S., & Tucker, M. C. (2025). Optimization of catholyte for halide-based all-solid-state batteries. *Journal of Power Sources*, 640, 236709.

Artal, R., Andersen, H. L., Del Olmo, R., Villaluenga, I., Sobrados, I., Diez-Gómez, V., ... & Aguadero, A. (2025). Cation disorder and lithium conductivity in mechanochemically synthesized chloride solid electrolytes. *Solid State Ionics*, 428, 116952.

Neumann, T., Alexander Dold, L., Thomas Cerny, A., Tröster, E., Günthel, M., Fischer, A., ... & Biro, D. (2025). Towards Sustainable Sulfide-Based All-Solid-State-Batteries: An Experimental Investigation of the Challenges and Opportunities Using Solid Electrolyte Free Silicon Anodes. *Batteries & Supercaps*, 8(2), e202400412.

Samanta, S., Bera, S., Biswas, R. K., Mondal, S., Mandal, L., & Banerjee, A. (2024). Ionocovalency of the Central Metal Halide Bond-Dependent Chemical Compatibility of Halide Solid Electrolytes with Li₆PS₅Cl. *ACS Energy Letters*, 9(8), 3683-3693.

Mandal, L., Biswas, R. K., Bera, S., Ogale, S. B., & Banerjee, A. (2024). Evolution of Interfacial Electro-Chemo-Mechanics between Lithium Metal and Halide Solid Electrolyte. *Chemistry of Materials*, 36(20), 10336-10350.

Shen, F., McGahan, M., Pietras, J. D., Lau, G. Y., Doeff, M. M., Battaglia, V. S., & Tucker, M. C. (2023). Tape casting of thin electrolyte and thick cathode for halide-based all-solid-state batteries. *Journal of the Electrochemical Society*, 170(10), 100505.

Karkera, G., Soans, M., Akbaş, A., Witter, R., Euchner, H., Diemant, T., ... & Fichtner, M. (2023). A structurally flexible halide solid electrolyte with high ionic conductivity and air processability. *Advanced Energy Materials*, 13(30), 2300982.

Brighi, M., Murgia, F., & Černý, R. (2022). Mechanical behavior and dendrite resistance of closo-hydroborate solid electrolyte. *Advanced Materials Interfaces*, 9(3), 2101254.


Dugas, R., Dupraz, Y., Quemin, E., Koç, T., & Tarascon, J. M. (2021). Engineered three-electrode cells for improving solid state batteries. *Journal of The Electrochemical Society*, 168(9), 090508.

Marchini, F., Saha, S., Alves Dalla Corte, D., & Tarascon, J. M. (2020). Li-rich layered sulfide as cathode active materials in all-solid-state Li-metal batteries. *ACS applied materials & interfaces*, 12(13), 15145-15154.



ABOUT

SPHERE ENERGY



The world is shifting towards sustainable energy solutions with a strong focus on reducing CO2 emissions. At the core of this revolution is global R&D, crucial for discovering new materials and technologies. However, a significant challenge remains: the lack of specialized testing equipment and advanced data solutions, hindering the exploration of materials that require novel testing methods and deeper data analysis.

Sphere Energy addresses this challenge by providing testing solutions and cutting-edge data analytics. As a team of battery experts, we support the entire battery value chain, from cell manufacturers to OEMs, offering **state-of-the-art hardware**, **customized battery testing**, and **AI-driven insights**. Our Battery Technology Centre serves as the hub for testing and validating advanced battery technologies, integrating AI and data analytics to enable unmatched innovation and optimization in the pursuit of sustainable, next-generation energy solutions.



Our products can be fast delivered worldwide when in stock. If not immediately available, we commit to delivering within 30 to 45 days.



Guaranteed responses within 24 hours, ensuring your queries are addressed promptly and efficiently.



All products are backed by an official 24-month warranty.



Our products are crafted with the highest standards of quality, durability and performance.



Info@sphere-energy.eu

+33 7 69 07-62-70

Sphere Energy SAS

250bis Boulevard Saint Germain

75007 Paris

France