

Comparison of Whole Blood Preservation Methods for Cytometric Analysis

Carly Lancaster, McKenna Coughlin, Gage Black, Justin Jarrell, Ramji Srinivasan, Li-Chun Cheng
Teiko Bio, Salt Lake City, UT, United States

Introduction

Clinical trials rely on high-dimensional immune profiling to evaluate treatment responses, but inconsistencies in blood collection, shipping, and processing result in a 10–15% sample loss. Approximately 80% of samples are processed live, yet transport under variable conditions causes a 49% reduction in total cell count with disproportionate loss of monocytes and neutrophils. These challenges highlight the need for reliable, scalable whole blood stabilization methods that preserve immune composition and support both surface and intracellular marker detection.

We compared three blood stabilization methods: Cyto-Chex, TransFix, and CellSave, to Teiko's two-step blood preservation kit (TokuKit). Cyto-Chex supported convenient single-tube collection but resulted in reduced B cell frequencies and increased CD4⁺CD8⁺ T cell artifacts after 7 days. TransFix preserved total population frequencies but showed poor resolution of monocyte subsets (CD14/CD16) and failed to detect intracellular markers such as T-bet. CellSave maintained overall immune frequencies and achieved the highest correlation to frozen controls at 7 days ($r = 0.97$ with BD Lyse) but exhibited slightly elevated granulocyte background.

TokuKit samples can be stored at -80°C for more than 6 months and showed strong correlation to fresh samples ($R > 0.97$). No measurable skewing of immune population frequencies or signal loss was observed over time, and both surface and intracellular markers were reliably detected. These results demonstrate that TokuKit provides robust long-term whole blood stabilization, maintains accurate immune population frequencies, and supports detection of both surface and intracellular markers. This approach enables more flexible sample handling for multi-site clinical trials and retrospective immune monitoring studies.

Teiko's TokuKit Workflow






Figure 1. Schematic of the whole blood collection process using the TokuKit

Immune Cell Frequencies are Comparable between TokuKit-Fixed and Live Whole Blood Samples



Validation Parameter	Metric	Result	Acceptance Criteria
Intra-run precision	Average %CV across all measurable populations	2.16%	≥95% of measurements CV ≤ 25%
Inter-run precision	Average %CV across all measurable populations	6.83%	≥95% of measurements CV ≤ 30%
Post-collection stability (frozen at -80°C)	Percent difference between $t=0$ and $t=x$ month	up to 15-month	≤25% change between timepoints

Figure 2. (A) Schematic of gating strategy for Live (left) and TokuKit-Fixed (right) whole blood samples. (B) Correlation of population frequencies (left) and functional marker subset frequencies (right) between live and TokuKit-fixed samples stored for 14 days at -80°C . (C) Chart of intra- and inter-run precision across all measurable populations as well as long-term stability measures for TokuKit fixed whole blood samples.

Method and Panel for Comparison of Stabilization Methods




Figure 3. (A) Schematic of collection and acquisition methods. Whole blood from four healthy donors was collected into Cyto-Chex, Transfix, CellSave, or TokuKit (contains buffers "Stable Lyse" and "Stable Store," (SLST) buffers manufactured by Smart Tube®). TokuKit samples were frozen at -80°C for 14 days and served as the gold standard reference control. Cyto-Chex, Transfix, and CellSave samples were stored at 4°C or room temperature (RT) for 1, 2, 3, or 7 days, lysed with ACK or BD Lyse/Fix, and stained with (B) a 41-marker mass cytometry immune profiling panel, and acquired on a Helios mass cytometer.

Transfix Preserves Total Event Counts During Acquisition




Figure 4. Comparison of total events acquired by CyTOF after blood collection in Cyto-Chex (blue), Transfix (red), and CellSave (green) compared to TokuKit (black).

Comparison of Whole Blood Preservation Methods to TokuKit




Figure 5. (A) Graphs of immune cell frequencies represented as percent of total leukocytes and corresponding dot plots comparing blood preservation methods. (B) Heat maps of r values for each experimental condition compared to TokuKit-Fixed whole blood.

Conclusions

CytoChex	Transfix	CellSave	TokuKit (-80°C storage)
✓ Single tube collection	✓ Preserves cell population frequencies similar to TokuKit	✓ Single tube collection	✓ Can store for long term (6+ months)
✓ Preservation at ambient temperature up to 14 days	✗ Struggle to separate monocyte subpopulations (CD14 ⁺ /- and CD16 ⁺ /-) in gating	✓ Preservation at ambient temperature up to 4 days	✓ Highly correlated with fresh, live cell staining
✗ Does not preserve B cell population frequencies	✗ Failed to detect intracellular markers (e.g. T-bet)	✗ Does not preserve B cell population frequencies	✓ Preserves immune cell frequencies
✗ Increase double positive T (CD4 ⁺ CD8 ⁺) that may be artifacts when stored for >7 days		✗ Reduced total event count by CyTOF	

Explore our dashboard: app.teiko.bio | teiko.bio | lichun@teiko.bio | SMS at (385) 355-9108