



Novel CellCage™ technology integrates image-based phenotyping and single-cell transcriptomics to study dynamic behaviors of living cells

Frank Charbonier¹, Yunmin Li¹, Nick Elder¹, Olaia Vila¹, Sina Farahvashi¹, Shreya Dashmukh¹, Mehdi Mohseni¹, Camilla Valente², Justin Poelma¹, Nirvan Rouzbeh¹, Elaheh Farjami¹, Meng Taing¹, Mika Maenaga¹, Teresa Al¹, Shan Sabri¹, Liz Wu¹, Pier Federico Gherardini¹, Sina Moeinzhadeh¹, Matthew Spitzer², Tarun Khurana¹, Gary Schroth¹, Mostafa Ronaghi¹

¹Cellanome, Foster City, CA ²Dept. of Otolaryngology-Head and Neck Surgery, UCSF School of Medicine, San Francisco, CA

ABSTRACT Imaging-based technologies are widely used to study cellular phenotypes in vitro because they provide multi-modal readouts of morphological, molecular, and functional cell states. However, integrated technologies which can directly link dynamic behaviors of living cells with high quality transcriptomic data at single-cell resolution are lacking. To address these limitations, we developed a novel CellCage™ Enclosure (CCE) technology: biocompatible and degradable hydrogel microenvironments formed by light-guided polymerization around individual cells or defined groups. CCEs enable long term culture of tens of thousands of cells, longitudinal imaging of dynamic cell behaviors, and direct mapping of image-based phenotypes to transcriptomic profiles on Cellanome's R3200 platform. In contrast to standard well plate, nanowell, or droplet-based methods, the properties of CCEs (i.e., size, shape, location, porosity, substrate) can be configured depending on the cell and assay type. CCEs can be formed on the order of seconds and selectively degraded on the order of minutes using cell-compatible chemistry. Permeable CCE walls allow diffusion of nutrients, waste, and fluorescence reagents, enabling the culture of suspension or adherent cells for weeks. We demonstrate CCE versatility across multiple applications. First, we correlated phagocytosis activity with transcriptomes in microglial cells, linking functional phenotypes to gene expression states. Second, we integrated image-based morphological embeddings of adherent cells with single-cell transcriptomic data to reveal functionally relevant cellular heterogeneity. Third, we co-enclosed dendritic cells with T cells to capture time-resolved activation dynamics during cell-cell interactions. Beyond these established applications, we also highlight emerging capabilities enabled by CCE technology. These include selective enrichment or selective retention of target cells from mixed populations for functional profiling, e.g., retention of CD56+ NK cells from an unlabeled PBMC cell mix. Additionally, we demonstrated integrated calcium imaging and transcriptome profiling from single astrocytes. Together, these results demonstrate that CCE technology enables novel multi-modal measurements linking dynamic cellular behaviors, morphological phenotypes, and single-cell transcriptomics across diverse experimental contexts.

