

INTRODUCTION

Single-cell technologies have revolutionized our understanding of cellular heterogeneity, yet current high-throughput methods face critical limitations. Conventional approaches rely on enzymatic dissociation and harsh microfluidic processing that disrupt functional networks and artificially alter transcriptomic states. These constraints are especially problematic for brain cells relevant to neuro-immunological and neurodegenerative diseases, which are large, fragile, and dependent on intact functional units. We developed a high-throughput, multi-modal platform enabling longitudinal studies (days to weeks) of individual cells, co-cultures, and neuronal networks, followed by endpoint transcriptomics without cell dissociation. The platform captures thousands of cells in CellCage™ Enclosures (CCE), generating linked imaging, functional, and transcriptome data. Flexible CCE designs support diverse strategies—from isolating single cells to capturing somas within intact neuronal networks or configurations supporting axonal outgrowth. Deep-learning morphology quantification, cell tracking, functional assays, and mRNA sequencing are linked across thousands of individual CCE per experiment.

Using this platform, we establish assays combining image-based phenotypic analysis with transcriptomics, including microglia phagocytosis and fluorescent imaging in neurons and astrocytes. Neurons remain functionally intact within networks until lysis, preserving cell-cell interactions lost during dissociation. We also establish glial-neuronal co-cultures where neurons form functional networks first, then microglial cells are introduced to the network to permit glial-neuron interactions such as phagocytic activity and surveillance behaviors while capturing transcriptomic changes in neurons. By enabling transcriptomic study of clinically relevant functions while preserving cell integrity, this platform provides a powerful tool for pathway analysis and target discovery in neuro-immunological research.

A NOVEL MICRO-3D PRINTING TECHNOLOGY ENABLES HIGH-THROUGHPUT EVALUATION OF CELL PHENOTYPES AND TRANSCRIPTOMES




Fig 1. Schematic of workflow for novel Cellanome technology enabling the measurement of multiple phenotypic and functional assays from the same cells in CellCage™ enclosures. Tens of thousands of cells in suspension are mixed with hydrogel precursor and loaded on a 4 or 8-lane flow cell. Positions of cells are identified and CellCage™ enclosures are generated around cells with light-guided polymerization in an automated fashion. CellCage™ enclosures are permeable to reagents (e.g., small molecules, antibodies) enabling long-term culturing and a variety of imaging-based, longitudinal phenotypic and functional assays to be performed on the same cells. Cells can be lysed within CellCage™ enclosures to generate barcoded cDNA for downstream library prep and sequencing off the instrument.

CAPTURING INTACT NEURON SOMA IN CELLCAGE™ ENCLOSURES






Fig 2. Using CellCage™ Enclosures to study connected neurons. (A) bit.bio ioGlutamatergic neurons were thawed and loaded into the flow cell. Cells were cultured for 7 days to allow for network formation. On the instrument, cells were automatically detected and CellCage™ Enclosures (CCE) created around individual neuron soma. After they were enclosed, neurons were imaged and transcriptome collected without cell dissociation. (B) Representative image showing 5 neuron soma captured inside of CCE. (C) Quality control of transcriptome data showing high number of molecular barcodes (MBC) collected for the majority of CCE. (D) Merged data from two lanes of neuron data. Top plot shows clusters identified by transcriptomic differences. Bottom plot shows the data points split for each lane. (E) Summary statistics of two lanes of neuron data showing the number of CCE, molecular barcodes, and genes recovered.

NEURON-MICROGLIA CO-CULTURES WITH SINGLE NEURON TRANSCRIPTOME



CHANGES IN ASTROCYTE TRANSCRIPTOME AND MORPHOLOGY WITH ACTIVATION



MORPHOLOGY PREDICTS PHAGOCYTOSIS IN MICROGLIA

